

StarLine®

*People From Your Community
Who Serve Your Community*

This catalog is brought to you courtesy of your authorized StarLine healthcare distributor. StarLine Dealers Association is a national association comprised of independently-owned medical supply distributors committed to high levels of service and satisfaction. StarLine distributors have provided physicians' offices, hospitals, nursing homes, laboratories and clinics with our extensive product offering for 20 years.

As an authorized member of StarLine, your distributor is able to negotiate volume pricing and value programs on medical products and equipment. Their affiliation with StarLine provides you, the healthcare provider, the benefit of the **StarLine Advantage** – high quality, competitively priced medical supplies delivered with reliable service.

StarLine Dealers Association was founded in 1986 to provide member distributors with quality products and competitive pricing through cooperative purchasing. StarLine is a long-recognized leader within medical distribution, with a focus giving local independent distributors the same access to national brand and private label products as other large regional competitors.

Throughout this catalog you will see **StarLine®** brand products, along with our proprietary brands **Courier®**, **ProAdvantage®**, **RediScreen™**, **Omnicide®** and **Dermodine™**, as well as other brand name products from reputable, well-known manufacturers. Although most, if not all, of the national-brand products in this catalog are available through your StarLine distributor, manufacturers do reserve the right to sell through their respective distributors. Your StarLine distributor will advise you if there are any products they cannot supply and will be able to show you similar products that are available to fit your needs.

*Thank you for supporting your
local independent healthcare distributor.*

Please see the next page for more catalog information.

StarLine®

The Quality Line in Healthcare

Look for the StarLine® brand throughout this catalog wherever you see the black outlined box. →

StarLine brand products, exclusive to StarLine distributors, are field tested and carefully screened to meet our high standards before entering the marketplace. This kind of diligent scrutiny, plus thousands of satisfied customers over the years, assures you this brand is the best value.

Quality... Value... Availability...

StarLine®

Gloves • Disposable Paper Products
Elastic Bandages • Gauze Sponges
Gauze Bandages • Cotton Balls
Adhesive Bandages • Applicators
Tongue Depressors • Penlights
Non-Adherent Pads • Lancets
Alcohol Prep Pads • Stethoscopes
Sphygmomanometers • Thermometry
Disposable Electrode Dermal Tips
Speculum • Detergents & Solutions
Microscope Slides • Cover Glass
Furniture & Equipment



Sterile Towel Drapes
Elastic Bandages
Sterile Gauze Sponges
Stretch Gauze Bandages
Stockinettes • Personal Care Items
Stainless Steel Instruments
Sterilization Pouches & Tubing
Wheelchairs

Pro ADVANTAGE®

Tongue Depressors • Applicators
Medical Spoons • First Aid Kit
Adult Wipes • Microscope Slides

RediScreen™

Rapid Diagnostics

Omnicide®

Glutaraldehyde Solutions

Dermodine™

Prep and Scrub Solutions

StarLine®

People From Your Community Who Serve Your Community

HEALTHCARE PRODUCT CATALOG INFORMATION

Please note the following when using this catalog.

1. While most standard items are featured, some available items may not be listed. Please contact your StarLine distributor representative to assist you in finding exactly what you need.
2. During the two year shelf-life of this catalog, certain products and lines will be deleted and new merchandise will be added. Please feel free to contact your StarLine distributors representative about the availability of any item.
3. Authorized distributors vary geographically. Therefore, in some instances, your StarLine distributor may not represent a specific manufacturer that appears in this catalog.
4. Your StarLine distributor may carry products, lines and manufacturers that do not appear in this catalog. Your distributor can inform you of their complete manufacturer offerings, specialty products and additional private label products if available.
5. An effort was made to best match the catalog numbers and respective photos together; however, some photos shown may represent a similar item or an item with accessories or additional parts. Some products will vary slightly from the photo. **Please consult your StarLine representative about specific products before purchasing based on a photograph.**

TERMS

Regular terms of sales are extended by StarLine distributors to most customers with favorable credit. Long-term financing can be arranged. Please ask your StarLine distributor representative for financing details and terms and conditions unique to your distributor.

GUARANTEE

The products we have presented in this catalog have been carefully chosen from reputable, well-known manufacturers. In addition to their standard guarantee, the StarLine name should assure your complete satisfaction.

SHIPPING

Purchases are normally shipped the most economical way possible unless your order specifies the type of carrier desired (parcel post, express, freight or other). Should damage or shortages occur, please notify your StarLine distributor representative immediately so they can initiate the proper adjustment procedures.

PRICING

This catalog does not provide pricing information due to the constantly changing prices in the healthcare products industry which are affected by individual pricing policies, periodic cost adjustments and competitive factors existing in this distributor's market. Please contact your StarLine distributor representative for a quote with their current pricing.

We look forward to servicing your needs with a winning combination of

- ★ **Quality Assurance**
- ★ **Local Service**
- ★ **Prompt Delivery**

- ★ **Competitive Pricing**
- ★ **In-Stock Warehousing**
- ★ **Personal Attention**

© 2006 National Distribution and Contracting, Inc. All Rights Reserved.

Not responsible for typographical or photographic errors. Under no circumstances shall the language, graphic content or other information contained within this catalog be reproduced or otherwise distributed in any form without the express written consent of National Distribution and Contracting, Inc.

StarLine is a subsidiary of National Distribution and Contracting, Inc. (NDC).

Printed and produced in the USA.

www.ndc-inc.com



STARLINE POWDER-FREE LATEX EXAM GLOVES

- 145054 X-Small (5 - 5 1/2) Gloves, 100/bx, 10 bx/cs
- 145052 Small (6 - 6 1/2) Gloves, 100/bx, 10 bx/cs
- 145050 Medium (7 - 7 1/2) Gloves, 100/bx, 10 bx/cs
- 145051 Large (8 - 8 1/2) Gloves, 100/bx, 10 bx/cs
- 145053 X-Large (9 - 9 1/2) Gloves, 100/bx, 10 bx/cs



STARLINE LOW MODULUS LATEX EXAM GLOVES

- 145002 Small Gloves, Non-Sterile,, 100/bx, 10 bx/cs
- 145001 Large Gloves, Non-Sterile,, 100/bx, 10 bx/cs
- 145000 Unisize Gloves, Non-Sterile, 100/bx, 10 bx/cs



STARTEX BRAND POWDERED LATEX EXAM GLOVES

- 145034 X-Small Exam Gloves, 100/bx, 10 bx/cs
- 145032 Small Exam Gloves, 100/bx, 10 bx/cs
- 145030 Medium Exam Gloves, 100/bx, 10 bx/cs
- 145031 Large Exam Gloves, 100/bx, 10 bx/cs
- 145033 X-Large Exam Gloves, 100/bx, 10 bx/cs



STARLINE SYNTHETIC VINYL POWDER-FREE EXAM GLOVES

- 145060 Medium Vinyl Gloves, 100/bx, 10 bx/cs
- 145061 Large Vinyl Gloves, 100/bx, 10 bx/cs
- 145063 X-Large Vinyl Gloves, 100/bx, 10 bx/cs



STARLINE SYNTHETIC VINYL POWDERED EXAM GLOVES

Lightly powdered satin stretch synthetic gloves. Manufactured in a strict quality control environment to exceed the ASTM and FDA standards for watertight testing and pinholes. Made of a unique formulation for a softer, more tactile feel. Non-sterile. Medical grade.

- 129502 Small Gloves, 6 - 6 1/2, 100/bx, 10 bx/cs
- 129504 Medium Gloves, 7-7 1/2, 100/bx, 10 bx/cs
- 129506 Large Gloves, 8 - 8 1/2, 100/bx, 10 bx/cs
- 129508 X-Large Gloves, 100/bx, 10 bx/cs



COURIER® STERILE TOWEL DRAPE

Disposable. Poly-interlined. Sterile.

- N207100 Plain, 18" x 26", Folded to 5" x 7 1/2", 50/bx, 6 bx/cs
- N207105 Fenestrated, 2 3/4" Opening, 18" x 26", Folded to 5" x 7 1/2", 50/bx, 6 bx/cs



STARLINE ECONOMY EXAMINATION TABLE PAPER

Soft crepe finish is pleasant to touch, thicker, more absorbent and protects tender skin areas. Smooth finish is strong, clean and easy to use. Clings to table surface. All types have substantial wet strength.

- 122063 18" x 200 ft Economy Table Paper, Smooth Finish, White, 12/cs
- 122065 21" x 200 ft Economy Table Paper, Smooth Finish, White, 12/cs
- 122062 18" x 225 ft Economy Table Paper, Smooth Finish, White, 12/cs
- 122064 21" x 225 ft Economy Table Paper, Smooth Finish, White, 12/cs
- 122072 18" x 125 ft Economy Table Paper, Crepe Finish, White, 12/cs
- 122074 21" x 125 ft Economy Table Paper, Crepe Finish, White, 12/cs



STARLINE EXAMINATION TABLE PAPER

Soft crepe finish is pleasant to touch, thicker, more absorbent and protects tender skin areas. Smooth finish is strong, clean and easy to use. Clings to table surface. All types have substantial wet strength.

- 151003 18" x 225 ft Table Paper, Smooth Finish, White, 12/cs
- 151007 21" x 225 ft Table Paper, Smooth Finish, White, 12/cs
- 151013 18" x 125 ft Table Paper, Crepe Finish, White, 12/cs
- 151017 21" x 125 ft Table Paper, Crepe Finish, White, 12/cs



STARLINE EXAMINATION DRAPE SHEETS

Made of soft tissue embossed for strength and durability. Lengths for both partial draping or more complete draping procedures.

- 1510302 2-Ply Drape Sheet, White, 40" x 48, 100/cs



STARLINE TISSUE/POLY PILLOWCASES

Constructed of soft, absorbent tissue embossed to extra thick poly for comfort and protection against moisture. Each case is vented on the bottom seam to allow pillow insertion without the usual air build up.

- 1510360 Tissue/Poly Pillowcase, 21" x 30", White, 100/cs



STARLINE DISPOSABLE TOWELS

Embossed towels that offer the combination of softness and absorbency for any clinical setting.

- 1510101 Fully Embossed Towels, 3-Ply Tissue, 13.5" x 18", White, 500/cs



1510315



1510515

STARLINE EXAMINATION CAPES

3-ply tissue capes feature embossed patterns for softness and durability. TPT capes offer comfort and strength and accommodate office procedures involving fluids, extended exams and patient perspiration. Front/back opening.

- 1510315 3-Ply Exam Cape, 30" x 21", White, 100/cs
- 1510515 TPT Exam Cape, 30" x 21", Blue, 100/cs



1510320



1510520

STARLINE EXAMINATION GOWNS

3-ply tissue gown ideal for shorter use procedures. TPT gowns provide added strength for extended exams. Gowns feature front/back openings.

- 1510320 3-Ply Exam Gown, 30" x 42", White, 50/cs
- 1510520 TPT Exam Gown, 30" x 42", Blue, 50/cs



STARLINE ELASTIC BANDAGES WITH DELUXE CLIPS

Non-sterile knitted high compression latex free. Manufactured with synthetic fiber and polyester. Offers consistent controlled pressure for sports injuries and other applications. Bandages secure with elastic clips. 5 yard lengths.

- 207369 2" Tan Elastic Bandage, 10/bx, 6 bx/cs
- 207376 3" Tan Elastic Bandage, 10/bx, 6 bx/cs
- 207377 4" Tan Elastic Bandage, 10/bx, 6 bx/cs
- 207378 6" Tan Elastic Bandage, 10/bx, 6 bx/cs



STARLINE ELASTIC BANDAGE WITH STANDARD CLIPS

Non-sterile knitted high compression latex free. Manufactured with synthetic fiber and polyester. Offers consistent controlled pressure for sports injuries and other applications. Bandages secure with standard clips. 5 yard lengths.

- 122562 2" Natural Elastic Bandage, 10/bx, 6 bx/cs
- 122563 3" Natural Elastic Bandage, 10/bx, 6 bx/cs
- 122564 4" Natural Elastic Bandage, 10/bx, 6 bx/cs
- 122566 6" Natural Elastic Bandage, 10/bx, 6 bx/cs



COURIER® ELASTIC BANDAGES WITH SELF-CLOSURE

Non-sterile knitted premium quality and latex free. Knitted with a blend of synthetic and 100% cotton yarns and increased absorbency. They feature convenient self-closure at each end for easy application. High level of compression for treatment of a variety of conditions.

- N877002 2" x 5 yds Self Closure Elastic Bandage, 10/bx, 6 bx/cs
- N877003 3" x 5 yds Self Closure Elastic Bandage, 10/bx, 6 bx/cs
- N877004 4" x 5 yds Self Closure Elastic Bandage, 10/bx, 6 bx/cs
- N877006 6" x 5 yds Self Closure Elastic Bandage, 10/bx, 6 bx/cs



STARLINE GAUZE SPONGES - NON-STERILE

- 155228 2" x 2" Gauze Sponge, 8-Ply, 200/slv, 25 slv/cs
- 155222 2" x 2" Gauze Sponge, 12-Ply, 200/slv, 40 slv/cs
- 155338 3" x 3" Gauze Sponge, 8-Ply, 200/slv, 20 slv/cs
- 155332 3" x 3" Gauze Sponge, 12-Ply, 200/slv, 20 slv/cs
- 155448 4" x 4" Gauze Sponge, 8-Ply, 200/slv, 20 slv/cs
- 155442 4" x 4" Gauze Sponge, 12-Ply, 200/slv, 10 slv/cs
- 155446 4" x 4" Gauze Sponge, 16-Ply, 200/slv, 10 slv/cs



COURIER® GAUZE SPONGES – STERILE

- N157002 1s, 2" x 2" Gauze Sponge, 12-Ply, 100/bx, 24 bx/cs
- N157003 1s, 3" x 3" Gauze Sponge, 12-Ply, 100/bx, 24 bx/cs
- N157004 1s, 4" x 4" Gauze Sponge, 12-Ply, 100/bx, 12 bx/cs
- N157022 2s, 2" x 2" Gauze Sponge, 8-Ply, 50/tr, 30 tr/cs
- N157023 2s, 3" x 3" Gauze Sponge, 12-Ply, 40/tr, 30 tr/cs
- N157024 2s, 4" x 4" Gauze Sponge, 8-Ply, 25/tr, 24 tr/cs
- N157025 2s, 4" x 4" Gauze Sponge, 12-Ply, 25/tr, 24 tr/cs



COURIER® STRETCH GAUZE BANDAGES

Non-sterile and white. Ideal for use as secondary dressings to hold primary dressings in place. It can also be used with IV's and splints or for mild compression or support. Clings to itself eliminating use of tape. Woven edges minimize lint and unraveling. High absorbency, softness and conformability.

- N157042 2" x 75" Stretch Gauze, 12/bx, 8 bx/cs
- N157043 3" x 75" Stretch Gauze, 12/bx, 8 bx/cs
- N157044 4" x 75" Stretch Gauze, 12/bx, 8 bx/cs
- N157046 6" x 75" Stretch Gauze, 6/bx, 8 bx/cs



STARLINE CONFORMING GAUZE BANDAGES

Non sterile. Clings to itself and stretches to conform. 60% cotton for breathability, absorbency and comfort blended with 40% polyester for strength and flexibility. Low linting and highly absorbent, helps reduce lint from entering the wound. Ease of application - no need to tuck or twist. Used in the operating, emergency rooms, as well as surgical floors, orthopedics and dermatology. Used to secure dressings and excellent for general bandaging, secondary burn dressings. 4.1 yards stretched.

- 1295412 2" Gauze Bandage, 12/bg, 8 bg/cs
- 1295413 3" Gauze Bandage, 12/bg, 8 bg/cs
- 1295414 4" Gauze Bandage, 12/bg, 8 bg/cs
- 1295416 6" Gauze Bandage, 12/bg, 8 bg/cs



STARLINE TONGUE DEPRESSORS

Non-Sterile

- 256500 5 1/2" Junior Tongue Depressor, 500/bx, 10 bx/cs
- 256520 6" x 11/16" Standard Tongue Depressor, 500/bx, 10 bx/cs



PROADVANTAGE® TONGUE DEPRESSORS

Non-Sterile

- 70000 5 1/2" Junior Tongue Depressor, 500/bx, 10 bx/cs
- 72000 6" x 11/16" Standard Tongue Depressor, 500/bx, 10 bx/cs
- 72200 6" Super Blade Tongue Depressor, 500/bx, 10 bx/cs

Sterile

- 76900 6" Standard Blade Tongue Depressor, Sterile, 1/pk, 100 pk/bx, 10 bx/cs
- 77000 5 1/2" Junior Tongue Depressor, Sterile, 1/pk, 100 pk/bx, 10 bx/cs



STARLINE COTTON TIPPED APPLICATORS

Non-Sterile

- 256561 3" x 1/12", Wood Stick, 100/bg, 10 bg/bx, 10 bx/cs
- 256562 6" x 1/12", Wood Stick, 100/bg, 10 bg/bx, 10 bx/cs

Sterile

- 256568 6" x 1/12", Wood Stick, Sterile, 2/pk, 100 pk/bx, 10 bx/cs



76100



77600

PROADVANTAGE® APPLICATORS

Non-Sterile

- 75900 Plain Applicator, 6" x 1/12", 1000/bx, 20 bx/cs
- 76100 Cotton-Tipped Applicator, 3" x 1/12", Wood Stick, 100/bg, 10 bg/bx, 10 bx/cs
- 76200 Cotton-Tipped Applicator, 6" x 1/12", Wood Stick, 100/bg, 10 bg/bx, 10 bx/cs
- 77500 16" Proctoscopic Polypropylene Shaft, Rayon-Tipped Applicator, 50/bx, 10 bx/cs
- 77600 8" OB/GYN Polypropylene Shaft, Rayon-Tipped Applicator, 50/bx, 10 bx/cs

Sterile

- 76700 Cotton-Tipped Applicator, Sterile, 6" x 1/12", Wood Stick, 1/pk, 100 pk/bx, 10 bx/cs
- 76800 Cotton-Tipped Applicator, Sterile, 6" x 1/12", Wood Stick, 2/pk, 100 pk/bx, 10 bx/cs
- 76801 Cotton-Tipped Applicator, Sterile, 6" x 1/10", Plastic Stick, 2/pk, 100 pk/bx, 10 bx/cs
- 77100 Cotton-Tipped Applicator, Sterile, 3" x 1/12", Wood Stick, 2/pk, 100 pk/bx, 10 bx/cs



PROADVANTAGE® MEDICAL SPOON

Non-Sterile

- 71510 3 5/8" Wooden Medical Spoon, 1000/bx, 10 bx/cs



COURIER® COTTON STOCKINETTES

Unbleached cotton provides maximum absorbency and patient comfort. For use under and over casts and other dressing applications. Non-sterile.

- N0781220 2" x 25 yds Cotton Stockinette
- N0781320 3" x 25 yds Cotton Stockinette
- N0781420 4" x 25 yds Cotton Stockinette
- N0781620 6" x 25 yds Cotton Stockinette

COURIER® SYNTHETIC STOCKINETTES

- N0786212 2" x 25 yds Synthetic Stockinette
- N0786312 3" x 25 yds Synthetic Stockinette
- N0786412 4" x 25 yds Synthetic Stockinette
- N0786612 6" x 25 yds Synthetic Stockinette



STARLINE NON-ADHERENT PADS

- 10707565 2" x 3" Pad, 100/bx, 12 bx/cs
- 10707575 3" x 4" Pad, 100/bx, 12 bx/cs



STARLINE FABRIC ADHESIVE BANDAGES

- Sterile**
- 10701580 3/4" x 3" Bandages, 100/bx, 12 bx/cs
 - 10701595 1" x 3" Bandages, 100/bx, 12 bx/cs
 - 10701570 2" x 4" Bandages, 50/bx, 24 bx/cs
 - 10701602 Knuckle Bandages, 100/bx, 12 bx/cs
 - 10701607 Fingertip Bandages, 100/bx, 12 bx/cs



STARLINE SHEER ADHESIVE BANDAGES

- Sterile**
- 10701275 3/4" x 3" Adhesive Bandages, 100/bx, 12 bx/cs
 - 10701290 1" x 3" Adhesive Bandages, 100/bx, 12 bx/cs
 - 10701303 7/8" Spots, 100/bx, 12 bx/cs



STARLINE COTTON BALLS

Pure white cotton balls made from absorbent 100% bleached cotton. Soft, compact and consistent in size and weight. They will not come apart when saturated. Non-sterile.

- 10705510 Medium Cotton Balls, 2M/bg, 2 bg/cs
- 10705530 Large Cotton Balls, 1M/bg, 2 bg/cs



STARLINE ALCOHOL PREP PADS

Prep pad impregnated with 70% USP Grade isopropyl alcohol. Used as a skin antiseptic. Packaged individually in a hermetically sealed foil pouch.

- 178002 Medium Prep Pads, 1 1/8" x 2 3/8", 200/bx, 10 bx/cs



PROADVANTAGE® FIRST AID KIT

This 25-person industrial/emergency first aid kit meets federal OSHA standard §1910.151b and ANSI Z308.1-2003. Kit contains 158 pieces: 1- plastic case, 100 - 1" x 3" adhesive plastic bandages, 3 - knuckle fabric bandages, 3 - fingertip fabric bandages, 1 set - triangular sling/bandage with safety pins, 4 - 3" x 3" gauze pads, 5 - 2: x 3" non-stick pads, 1 - 5" x 9" trauma pad, 1 - 3" conforming gauze roll bandage, 1 - 2" conforming gauze roll bandage, 2 - sterile eye pads, 15 - antiseptic cleansing wipes, 3 - ammonia inhalant ampoules, 6 - 0.9 gm first aid / burn cream packs, 1 - 4" x 5" instant cold compress, 1 - 1/2" x 10 yd. first aid tape roll, 1 - 4-1/2" scissors, 1 - 4" plastic tweezers, 4 - vinyl exam gloves, 1 - 1/2 oz. eye wash, 1 - first aid guide. Made in the USA.

- P440025 25-Person Industrial/Emergency First Aid Kit, 10/cs



**PROADVANTAGE® ADULT WIPES
PREMIUM SPUNLACE CLEANSING WASHCLOTHS**

Pre-moistened washcloths of spunlace with aloe and lanolin. Hypoallergenic, alcohol free, latex free. Soft to the skin, easy dispensing. Mild scent.

- P90001 Flip Top Soft Pack, Adult Wipes, 9" x 13", 64/pk, 8 pk/cs
- P90002 Easy Open Hard Tub, Adult Wipes, 9" x 13", 64/tub, 8 tub/cs



COURIER® BATHING SYSTEM

Premium pre-moistened needle punch cleansing and bathing washcloths with aloe. Hypoallergenic, alcohol free, latex free. Soft to the skin, accordion folded for easy dispensing. Mild scent.

- N90003 Bathing System Needle Punch Material, 8 1/2" x 8", 8/pk, 24 pk/cs



COURIER® ANTISEPTIC PERINEUM CLEANSER

A gentle, spray-on cleanser designed specifically for incontinence care. Gently cleans and deodorizes with a light, fresh scent. It maintains a neutral pH, is low sudsing for easy rinsing, and contains an antimicrobial agent to prevent perineal irritation from bacteria.

- N909150 No-Rinse Antiseptic Perineal Wash, 8 oz, 12/cs



COURIER® PERINEUM CLEANSER

Freshly scented no-rinse incontinence and stoma care products. Soothing and gentle spray wash. Pleasantly scented with water lily papaya. Its gentle soothing formula quickly cleans, refreshes and eliminates odors. Rich in lanolin and aloe, it replenishes natural skin oils.

- N909110 No-Rinse Perineal Wash w/aloe & lanolin, 8 oz, 48/cs
- N909120 No-Rinse Perineal Wash w/aloe & lanolin, 128 oz, 4/cs



COURIER® PERSONAL CARE SAMPLE KIT

This cake box sampler contains (1) 8 oz bottle each of the Courier No-Rinse Perineal Wash with aloe and lanolin, Courier No-Rinse Antiseptic Perineal Wash, Courier Body Wash and Shampoo with aloe vera, and Courier Conditioning Deodorizing Shampoo and Body Wash.

- N90CCSP Cake Box Sampler Contains: (1) each of N909110, N909150, N909210, N909250



COURIER® BODY WASH & SHAMPOO GEL

An ultra mild, aloe enriched all-in-one body bath and shampoo product. Leaves hair clean and skin soft, moist, and smooth. Mild enough for everyday use. The green pear fragrance is light and pleasant, never over-powering. The flip top cap means dispensing is never a problem.

- N909210 Body Wash & Shampoo w/aloe vera, 8 oz, 48/cs
- N909220 Body Wash & Shampoo w/aloe vera, 128 oz, 4/cs



COURIER® PREMIUM BODY WASH & SHAMPOO GEL

An effective odor reducing, non-alkaline, free-rinsing shampoo, body wash and skin cleanser. Its special combination of conditioners and emollients, in a pearlescent lotion base, provide this gentle product with its unique rich lather. Contains aloe vera to smooth and condition skin.

- N909250 Conditioning, Deodorizing Shampoo & Body Wash, 8 oz, 48/cs
- N909260 Conditioning, Deodorizing Shampoo & Body Wash, 128 oz, 4/cs



DERMIDINE PVP PREP SOLUTION

Formulated to USP 10% available iodine. For use as an antiseptic for preparation of the skin prior to surgery. Also used as a first aid antiseptic to help prevent infection in minor cuts, scrapes, and burns.

- 81763900930 4 oz PVP Prep Solution, 48/cs
- 81763900931 8 oz PVP Prep Solution, 24/cs
- 81763900932 16 oz PVP Prep Solution, 24/cs
- 81763900934 Gallon PVP Prep Solution, 4/cs



DERMIDINE PVP SURGICAL SCRUB SOLUTION

Formulated to USP 10% available iodine. For washing whenever a germicidal soap is required.

- 81763900940 4 oz PVP Surgical Scrub, 48/cs
- 81763900941 8 oz PVP Surgical Scrub, 24/cs
- 81763900942 16 oz PVP Surgical Scrub, 24/cs
- 81763900944 Gallon PVP Surgical Scrub, 4/cs



STARLINE SINGLE HEAD NURSE STETHOSCOPES

Diaphragm only aluminum chestpiece. Adjustable chrome-plated binaurals. Flexible 22" PVC Y tubing. White plastic eartips. Overall length 31". Latex free.

- | | |
|----------------------------------|-----------------------------------|
| 104266030 Black Stethoscope | 104266049 White Stethoscope |
| 104266031 Gray Stethoscope | 104266052 Neon Orange Stethoscope |
| 104266033 Light Blue Stethoscope | 104266053 Neon Pink Stethoscope |
| 104266035 Royal Blue Stethoscope | 104266054 Neon Yellow Stethoscope |
| 104266037 Green Stethoscope | 104266056 Burgundy Stethoscope |
| 104266039 Red Stethoscope | 104266057 Dark Green Stethoscope |
| 104266040 Pink Stethoscope | 104266059 Teal Stethoscope |
| 104266042 Gold Stethoscope | 104266060 Purple Stethoscope |
| 104266044 Lavender Stethoscope | 104266061 Navy Stethoscope |
| 104266045 Orchid Stethoscope | |
| 104266046 Magenta Stethoscope | |
| 104266047 Yellow Stethoscope | |



STARLINE DUAL-HEAD STETHOSCOPES

Combination diaphragm/bell chestpiece for wide frequency response. Bell and diaphragm equipped with non-chill ring for patient comfort. Lightweight, durable construction. Adjustable metal binaurals. Flexible 22" PVC tubing. Plastic eartips. Overall length 31 1/2". Latex free.

- | | |
|----------------------------------|-----------------------------------|
| 104267030 Black Stethoscope | 104267047 Yellow Stethoscope |
| 104267031 Gray Stethoscope | 104267049 White Stethoscope |
| 104267033 Light Blue Stethoscope | 104267052 Neon Orange Stethoscope |
| 104267035 Royal Blue Stethoscope | 104267053 Neon Pink Stethoscope |
| 104267037 Green Stethoscope | 104267054 Neon Yellow Stethoscope |
| 104267039 Red Stethoscope | 104267056 Burgundy Stethoscope |
| 104267040 Pink Stethoscope | 104267057 Dark Green Stethoscope |
| 104267042 Gold Stethoscope | 104267059 Teal Stethoscope |
| 104267044 Lavender Stethoscope | 104267060 Purple Stethoscope |
| 104267045 Orchid Stethoscope | 104267061 Navy Stethoscope |
| 104267046 Magenta Stethoscope | |



STARLINE INFANT STETHOSCOPES

Infant sized, combination diaphragm/bell chestpiece for wide frequency response. 25mm diaphragm/bell equipped with non-chill ring for patient comfort. Chrome-plated brass chestpiece. Adjustable chrome-plated binaurals. Flexible 22" PVC Y tubing. White plastic eartips. Overall length 31". Latex free.

104267630	Black Stethoscope	104267647	Yellow Stethoscope
104267631	Gray Stethoscope	104267649	White Stethoscope
104267633	Light Blue Stethoscope	104267652	Neon Orange Stethoscope
104267635	Royal Blue Stethoscope	104267653	Neon Pink Stethoscope
104267637	Green Stethoscope	104267654	Neon Yellow Stethoscope
104267639	Red Stethoscope	104267656	Burgundy Stethoscope
104267640	Pink Stethoscope	104267657	Dark Green Stethoscope
104267642	Gold Stethoscope	104267659	Teal Stethoscope
104267644	Lavender Stethoscope	104267660	Purple Stethoscope
104267645	Orchid Stethoscope	104267661	Navy Stethoscope
104267646	Magenta Stethoscope		



STARLINE PEDIATRIC STETHOSCOPES

Pediatric sized, combination diaphragm/bell chestpiece for wide frequency response. 1 1/4" diaphragm/bell equipped with non-chill ring for patient comfort. Lightweight aluminum construction. Adjustable chrome-plated binaurals. Flexible 22" PVC Y tubing. White plastic eartips. Overall length 31 1/2". Latex free.

104267530	Black Stethoscope	104267547	Yellow Stethoscope
104267531	Gray Stethoscope	104267549	White Stethoscope
104267533	Light Blue Stethoscope	104267552	Neon Orange Stethoscope
104267535	Royal Blue Stethoscope	104267553	Neon Pink Stethoscope
104267537	Green Stethoscope	104267554	Neon Yellow Stethoscope
104267539	Red Stethoscope	104267556	Burgundy Stethoscope
104267540	Pink Stethoscope	104267557	Dark Green Stethoscope
104267542	Gold Stethoscope	104267559	Teal Stethoscope
104267544	Lavender Stethoscope	104267560	Purple Stethoscope
104267545	Orchid Stethoscope	104267561	Navy Stethoscope
104267546	Magenta Stethoscope		



STARLINE SPRAGUE RAPPAPORT STETHOSCOPES

Threaded chestpiece drum of chrome plated zinc permits the use of five interchangeable chestpiece fittings. 22" tubing length; overall length 30". Latex free, double tube configuration for improved acoustic separation. Proprietary valve mechanism eliminates acoustic leakage. Adjustable chrome plated brass binaurals. Complete with an accessory pouch containing the bell chestpieces, chestpiece key, 2 spare diaphragms, one pair of mushroom eartips and one pair of silicon eartips. Weighs 7.5 oz. Inspected and packaged in the U.S.A. Five year warranty.

104264130	Black Stethoscope	104264150	Neon Blue Stethoscope
104264131	Gray Stethoscope	104264151	Neon Green Stethoscope
104264133	Light Blue Stethoscope	104264152	Neon Orange Stethoscope
104264135	Royal Blue Stethoscope	104264153	Neon Pink Stethoscope
104264137	Green Stethoscope	104264154	Neon Yellow Stethoscope
104264139	Red Stethoscope	104264155	Boysenberry Stethoscope
104264140	Pink Stethoscope	104264156	Burgundy Stethoscope
104264142	Gold Stethoscope	104264157	Dark Green Stethoscope
104264144	Lavender Stethoscope	104264159	Teal Stethoscope
104264145	Orchid Stethoscope	104264160	Purple Stethoscope
104264146	Magenta Stethoscope	104264161	Navy Stethoscope



STARLINE SPHYG/SPRAGUE KIT

Black enamel 300mmHg no-pin stop manometer. Color matched sprague-type stethoscope. Calibrated nylon cuff. Latex bladder and bulb. Chrome plated valve. Color matched nylon carrying case that is extra large with one large zippered compartment, outside Velcro® pocket, end strap and handy key ring. Case measures 10" x 6" x 3". Assembled, inspected and packaged in the U.S.A. Two year inflation system warranty. 20 year calibration warranty.

30761	Black Sphyg/Sprague Kit
30762	Burgundy Sphyg/Sprague Kit
30763	Dark Green Sphyg/Sprague Kit
30764	Magenta Sphyg/Sprague Kit
30765	Silver Sphyg/Sprague Kit
30766	Teal Sphyg/Sprague Kit
30767	Purple Sphyg/Sprague Kit



STARLINE PREMIUM ANEROID SPHYGMOMANOMETER

This is our premier sphygmomanometer. Features #10428000 lifetime warranted gauge and top quality inflation system, including a nylon Velcro® cuff and highest quality bulb and valve. Complete with leatherette case.

- 104272311 Adult Sphygmomanometer, Gray
- 104272351 Adult Sphygmomanometer, Royal Blue
- 104272371 Adult Sphygmomanometer, Green
- 104272391 Adult Sphygmomanometer, Red
- 104272431 Adult Sphygmomanometer, Orange
- 104272461 Adult Sphygmomanometer, Magenta
- 104272561 Adult Sphygmomanometer, Burgundy
- 104272571 Adult Sphygmomanometer, Dark Green
- 104272591 Adult Sphygmomanometer, Teal
- 104272601 Adult Sphygmomanometer, Purple
- 104272611 Adult Sphygmomanometer, Black
- 30740 Adult Sphygmomanometer, Navy
- 104272615 Large Adult Sphygmomanometer, Black
- 104272612 Newborn Sphygmomanometer, Black
- 104272613 Infant Sphygmomanometer, Black
- 104272614 Child Sphygmomanometer, Black
- 104272616 Thigh Sphygmomanometer, Black

Latex-Free

- 30741 Adult Sphygmomanometer, Black, Latex Free
- 30742 Child Sphygmomanometer, Black, Latex Free
- 30743 Large Adult Sphygmomanometer, Black, Latex Free

Accessories

- 10428000 Gauge Only, Latex Free



STARLINE ANEROID SPHYGMOMANOMETER

Black enamel 300mmHg no-pin stop manometer. Range markings to facilitate selection of correct cuff size. Artery label and gauge holder. Latex inflation bladder and bulb. Chrome plated deflation valve. Complete with leatherette zippered carrying case and operating instructions. Assembled, inspected and packaged in the U.S.A. Ten year calibration warranty.

- 104277311 Adult Aneroid Sphygmomanometer, Gray
- 104277551 Adult Aneroid Sphygmomanometer, Blue



STARLINE POCKET ANEROID SPHYGMOMANOMETER

This is our "standard" sphygmomanometer. Features are the same as the premier aneroid series, except for a standard bulb/valve and bladder combination.

- 104246571 Adult Sphygmomanometer, Dark Green, Latex
- 104276311 Adult Sphygmomanometer, Gray, Latex
- 104276351 Adult Sphygmomanometer, Royal Blue, Latex
- 104276371 Adult Sphygmomanometer, Green, Latex
- 104276391 Adult Sphygmomanometer, Red, Latex
- 104276431 Adult Sphygmomanometer, Orange, Latex
- 104276461 Adult Sphygmomanometer, Magenta, Latex
- 104276561 Adult Sphygmomanometer, Burgundy, Latex
- 104276591 Adult Sphygmomanometer, Teal, Latex
- 104276601 Adult Sphygmomanometer, Purple, Latex
- 104276611 Adult Sphygmomanometer, Black, Latex
- 30940 Adult Sphygmomanometer, Navy, Latex
- 104276615 Large Adult Sphygmomanometer, Black, Latex
- 104276612 Newborn Sphygmomanometer, Black, Latex
- 104276613 Infant Sphygmomanometer, Black, Latex
- 104276614 Child Sphygmomanometer, Black, Latex
- 104276616 Thigh Sphygmomanometer, Black, Latex

Latex-Free

- 104230941 Adult Sphygmomanometer, Black, Latex Free
- 30942 Child Sphygmomanometer, Black, Latex Free
- 30943 Large Adult Sphygmomanometer, Black, Latex Free

STARLINE 2-TUBE BLADDERS

Dip molded latex or latex free. 55cm tubing length. Size and model number printed on bladder. For use with most aneroid and mercurial blood pressure instruments. Individually boxed. Two year warranty.

- 10428010 Two-Tube Bladder, Latex, Adult

STARLINE CALIBRATED NYLON CUFF

Featuring the marking system. For use with most aneroid or mercurial blood pressure instruments. Two year warranty.

- 104283610 Adult Nylon Cuff, Black, Latex

STARLINE BLOOD PRESSURE ACCESSORIES

For use with most aneroid or mercurial blood pressure instruments. Individually boxed.

- 104284610 Adult Cuff, 2-Tube, Black
- 104285610 Adult Cuff, 1-Tube, Black
- 30840 Adult Cuff, 2-Tube, Black, Latex Free
- 30850 Adult Cuff, 1-Tube, Black, Latex Free
- 10428075 Standard Bulb & Valve, Latex

STARLINE INFLATION SYSTEM

Includes nylon cuff, deluxe 2-tube bladder, bulb and valve. For use with most aneroid or mercurial blood pressure instruments. Individually boxed.

- 104286610 Adult Inflation System, Black, Latex
- 30860 Adult Inflation System, Black, Latex Free



STARLINE DIGITAL THERMOMETER

Accurately records body temperature in 60 seconds. Liquid crystal display read-out displays temperature within 1/10th of a degree. Audible beep signals peak temperature. A replaceable battery provides 200 hours of use. Each package comes complete with thermometer, 5 probe covers, battery, storage case, instructions, and one year limited warranty. Can be used orally, rectally or under the arm. Latex free.

10420410 Digital Thermometer w/beeper, °F



STARLINE PHYSICIAN'S PENLIGHT

Penlight has high-intensity prefocused bulb powered by two AAA batteries. Aluminum contact strip assures good contact for brighter light. Extruded case has a unique V-clip. 1-year shelf life. Useful life is 3- to 6-months with intermittent use.

6666S Dr. Pack Disposable Penlight, 6/pk



STARLINE DISPOSABLE ACTIVE ELECTRODES — DERMAL TIPS

Can be used with most electrosurgical generators that use 3/32" electrodes.

- 1010A804 Sharp Tip Electrodes, Non-Sterile, 100/bx
- 1010A805 Sharp Tip Electrodes, Sterile, 50/bx
- 1010A806 Blunt Tip Electrodes, Non-Sterile, 100/bx
- 1010A807 Blunt Tip Electrodes, Sterile, 50/bx



STARLINE LANCETS

5925001 General Purpose Lancets, Orange, Standard Cannula 23ga, 200/bx

COURIER® INSTRUMENTS

Manufactured to German standards with exceptional value. Lifetime guarantee against defects in material or workmanship assures years of serviceability from each instrument. Selection of the most popular patterns of surgical instruments.



COURIER® ADSON DRESSING FORCEPS

Manufactured to German standards with exceptional value. Lifetime guarantee against defects in material or workmanship assures years of serviceability from each instrument.

N407120 4 3/4" Adson Dressing Forceps



COURIER® ADSON TISSUE FORCEPS

Manufactured to German standards with exceptional value. Lifetime guarantee against defects in material or workmanship assures years of serviceability from each instrument.

N407125 4 3/4" Adson Tissue Forceps, 1x2 Teeth



N407100



N407105

COURIER® DRESSING FORCEPS

Manufactured to German standards with exceptional value. Lifetime guarantee against defects in material or workmanship assures years of serviceability from each instrument.

N407100 5" Dressing Forceps

N407105 5 1/2" Dressing Forceps



COURIER® FOERSTER FORCEPS

Manufactured to German standards with exceptional value. Lifetime guarantee against defects in material or workmanship assures years of serviceability from each instrument.

N407225 9 1/2" Straight Serrated Foerster Forceps



N407200



N407205

COURIER® HALSTED MOSQUITO FORCEPS

Manufactured to German standards with exceptional value. Lifetime guarantee against defects in material or workmanship assures years of serviceability from each instrument.

- N407200 5" Straight Halsted Mosquito Forceps
- N407205 5" Curved Halsted Mosquito Forceps



COURIER® BROWN-ADSON FORCEPS

Manufactured to German standards with exceptional value. Lifetime guarantee against defects in material or workmanship assures years of serviceability from each instrument.

- N407130 4 3/4" Brown-Adson Forceps, 7x7 Teeth



N407210



N407215

COURIER® KELLY FORCEPS

Manufactured to German standards with exceptional value. Lifetime guarantee against defects in material or workmanship assures years of serviceability from each instrument.

- N407210 5 1/2" Straight Kelly Forceps
- N407215 5 1/2" Curved Kelly Forceps



N407135



N407140



N407145

COURIER® SPLINTER FORCEPS

Manufactured to German standards with exceptional value. Lifetime guarantee against defects in material or workmanship assures years of serviceability from each instrument.

- N407135 3 1/2" Splinter Forceps
- N407140 4" Splinter Forceps
- N407145 4 1/2" Splinter Forceps



N407110



N407115

COURIER® TISSUE FORCEPS

Manufactured to German standards with exceptional value. Lifetime guarantee against defects in material or workmanship assures years of serviceability from each instrument.

- N407110 5" Tissue Forceps 1x2 Teeth
- N407115 5 1/2" Tissue Forceps 1x2 Teeth



COURIER® KEVORKIAN-YOUNGE BIOPSY FORCEPS

Manufactured to German standards with exceptional value. Lifetime guarantee against defects in material or workmanship assures years of serviceability from each instrument.

- N407430 Kevorkian-Youngue Biopsy Forceps



COURIER® BACKHAUS TOWEL CLAMPS

Manufactured to German standards with exceptional value. Lifetime guarantee against defects in material or workmanship assures years of serviceability from each instrument.

N407220 3 1/2" Backhaus Towel Clamp



COURIER® BANDAGE & UTILITY SCISSORS

Manufactured to German standards with exceptional value. Lifetime guarantee against defects in material or workmanship assures years of serviceability from each instrument.

N407035 7 1/2" Bandage & Utility Scissors



COURIER® LITTAUER STITCH SCISSORS

Manufactured to German standards with exceptional value. Lifetime guarantee against defects in material or workmanship assures years of serviceability from each instrument.

N407045 4 1/2" Littauer Stitch Scissors



N407005



N407010



N407015



N407020

COURIER® IRIS SCISSORS

Manufactured to German standards with exceptional value. Lifetime guarantee against defects in material or workmanship assures years of serviceability from each instrument.

- N407005 4" Straight Iris Scissors
- N407010 4" Curved Iris Scissors
- N407015 4 1/2" Straight Iris Scissors
- N407020 4 1/2" Curved Iris Scissors



N407025

N407030

COURIER® LISTER BANDAGE SCISSORS

Manufactured to German standards with exceptional value. Lifetime guarantee against defects in material or workmanship assures years of serviceability from each instrument. Selection of the most popular patterns of surgical instruments.

- N407025 5 1/2" SS Lister Bandage Scissors
- N407030 7 1/4" SS Lister Bandage Scissors



COURIER® O.R. SCISSORS

Manufactured to German standards with exceptional value. Lifetime guarantee against defects in material or workmanship assures years of serviceability from each instrument.

N407000 5 1/2" Straight Sharp/Blunt O.R. Scissors



COURIER® SPENCER STITCH SCISSORS

Manufactured to German standards with exceptional value. Lifetime guarantee against defects in material or workmanship assures years of serviceability from each instrument.

N407040 3 1/2" Spencer Stitch Scissors

The Courier® Guarantee

All Courier instruments are guaranteed to be free from defects in workmanship and material. Any Courier instrument which proves defective in workmanship and/or material will either be repaired or replaced at our discretion without charge.



N407300



N407305



N407310



N407315



N407320

COURIER® NEEDLE HOLDERS

Manufactured to German standards with exceptional value. Lifetime guarantee against defects in material or workmanship assures years of serviceability from each instrument.

- N407300 Webster 5" Needle Holder
- N407305 Halsey 5" Needle Holder
- N407310 Baumgartner 5 1/2" Needle Holder
- N407315 Mayo-Hegar 5" Needle Holder
- N407320 Mayo-Hegar 6" Needle Holder



STARLINE VAGINAL SPECULA

Recommended for routine gynecological exams. Disposable. Allows the examination of the vaginal mucosa at 360° angle without rotation. Wide range of adjustability. Ease of operation and manipulation with one hand. Self-holding in various opening positions. Made from hard, clear and non-toxic plastic resin of high strength.

- 16108309 Small Vaginal Specula, 10/bg, 10 bg/cs
- 16108310 Medium Vaginal Specula, 10/bg, 10 bg/cs
- 16108311 Large Vaginal Specula, 10/bg, 10 bg/cs



COURIER® PEDERSON VAGINAL SPECULA

Manufactured to German standards with exceptional value. Lifetime guarantee against defects in material or workmanship assures years of serviceability from each instrument. Selection of the most popular patterns of surgical instruments.

- N407415 Small Vaginal Specula
- N407420 Medium Vaginal Specula
- N407425 Large Vaginal Specula



COURIER® GRAVES VAGINAL SPECULA

Manufactured to German standards with exceptional value. Lifetime guarantee against defects in material or workmanship assures years of serviceability from each instrument. Selection of the most popular patterns of surgical instruments.

- N407400 Small Vaginal Specula
- N407405 Medium Vaginal Specula
- N407410 Large Vaginal Specula



N407435



N407440

COURIER® TISCHLER BIOPSY PUNCH

Manufactured to German standards with exceptional value. Lifetime guarantee against defects in material or workmanship assures years of serviceability from each instrument. Selection of the most popular patterns of surgical instruments.

- N407435 Baby-Tischler Biopsy Punch
- N407440 Tischler Biopsy Punch 3x7mm



COURIER® SELF-SEAL STERILIZATION POUCHES

- N012790 2.75" x 9" Pouch, 200/bx, 6 bx/cs
- N013554 3.5" x 5.25" Pouch, 200/bx, 6 bx/cs
- N013590 3.5" x 9" Pouch, 200/bx, 6 bx/cs
- N015015 5" x 15" Pouch, 200/bx, 6 bx/cs
- N011016 10" x 16" Pouch, 200/ctn
- N015410 5.25" x 10" Pouch, 200/bx, 6 bx/cs
- N017513 7.5" x 13" Pouch, 200/bx, 6 bx/cs
- N018024 8" x 24" Pouch, 200/ctn
- N011318 Cassette Pouch, 13" x 18", 200/ctn
- N011217 Cassette Pouch, 12" x 17", 200/ctn
- N013522 Cassette Pouch, 3.5" x 22.5", 200/ctn



COURIER® STERILIZATION TUBING

- N01102 2" x 100 ft Tubing, 16 rl/cs
- N01103 3" x 100 ft Tubing, 8 rl/cs
- N01104 4" x 100 ft Tubing, 16 rl/cs
- N01106 6" x 100 ft Tubing, 8 rl/cs
- N01108 8" x 100 ft Tubing, 4 rl/cs
- N01110 10" x 100 ft Tubing, 4 rl/cs



STARCLENZ 3000 INSTRUMENT DETERGENT/LUBRICANT

Multipurpose detergent for cleaning and lubricating stainless steel surgical instruments in one step. Incorporates two nonsilicone, water soluble lubricants, that clean and lubricate hinged and metal to metal surfaces. Eliminates a messy and time consuming lubrication step. Instrument life is prolonged by preventing corrosion, pitting and rusting that can destroy valuable instruments. May be used in manual, ultrasonic and washer sterilizer/decontaminators.

12283000 StarClenz 3000 Instrument Detergent, 1 Gallon w/1 oz pump, 4/cs



STARZYMATIC 2000 DUAL ENZYMATIC DETERGENT

Low foaming for cleaning instruments. Fast and effective against complex soils. Effectively attacks both protein and starch based soils. Reduces unpleasant odors with a built-in odor suppressant. Ultrasonic solution. Concentrated. Neutral pH formula.

12282000 StarZymatic 2000 Detergent, 1 Gallon w/1oz pump, 4/cs



STARZYMATIC SINGLE ENZYMATIC DETERGENT

Low foaming for cleaning instruments. Designed to attack the most common protein based instrument soils. pH neutral and non-corrosive. Ultrasonic solution. Concentrated.

12282001 StarZymatic Detergent, 1 Gallon w/1oz pump, 4/cs



STARLINE HIGH SUDS INSTRUMENT DETERGENT

Concentrated powder and liquid are biodegradable with neutral pH. Contain rust inhibitors and rinsing agents. Hemolytic to remove blood.

B2S High Suds Instrument Detergent Liquid For Manual Cleaning, Gallon, 4/cs
 B3S High Suds Instrument Detergent Powder, 5 lb Tub, 6/cs



STARLINE LOW SUDS INSTRUMENT DETERGENT

Biodegradable with neutral pH. Hemolytic to remove blood. Contains rust inhibitor and rinsing agent. Highly concentrated.

B1S Low Suds Instrument Detergent Liquid For Ultrasonic Cleaners, Gallon, 4/cs





STARLINE ENZYMATICAL CLEANING SOLUTION

Disintegrates blood products. Neutral pH, bacteriostatic, double enzyme. Non-corrosive, safe for scopes. Cleans hard to reach surfaces. Can be used in ultrasonic cleaners.

B9S Enzymatic Cleaner, Gallon, 4/cs



STARLINE INSTRUMENT LUBRICANT

Non-silicone and non-toxic solution that eases stiffness. Steam permeable and water soluble. Prevents electrolysis on points and edges. Helps prevent rust.

T2S Instrument Lubricant, Gal, 4/cs



OMNICIDE® 14-DAY HIGH LEVEL DISINFECTANT/STERILANT

A 2.6% buffered glutaraldehyde that can be reused as a sterilant and high-level disinfectant on immersible lensed and heat sensitive instruments for up to 14 days. Safe for use on lensed instruments and is compatible with plastic, rubber, metal and other material used in medical and dental instrumentation and equipment requiring disinfection or sterilization. No dilution is required.

N099001 Omnicide® Gallon, 2.6% Buffered Glutaraldehyde, 4/cs



OMNICIDE® 28-DAY HIGH LEVEL DISINFECTANT/STERILANT

A long-life, buffered glutaraldehyde that can be reused as a sterilant and high-level disinfectant for up to 28 days. Recommended for immersible and heat sensitive instruments. Contains a surfactant to enhance the disinfection process and contains a rust-inhibitor to reduce corrosion. No dilution is required.

N099002 Omnicide® Gallon, 2.8% Buffered Glutaraldehyde, 4/cs



STARLINE MICROSCOPE SLIDES

Made from finest sheet glass to assure smooth, flat surface. No striations, bubbles or ripples. Ground, smooth edges. Pre-cleaned and ready to use. 1.00 mm thick. 3" x 1". Packaging: 10 gross/case.

1506420 Plain Microscope Slides
1506421 Frosted Microscope Slides

STARLINE COVER GLASS

Superior corrosion resistance. Uniform flatness and thinness with no warping. Crush-resistant plastic box has semi-soft foam interior for lint free protection. Packaging: 1 ounce per plastic carton, 10 cartons/box

150612222 22 mm x 22 mm Cover Glass, Thickness No. 1 (0.13 to 0.17 mm)
150622222 22 mm x 22 mm Cover Glass, Thickness No. 2 (0.17 to 0.25 mm)

PROADVANTAGE MICROSCOPE SLIDES

Made from finest sheet glass to assure smooth, flat surface. No striations, bubbles or ripples. Ground, smooth edges. Pre-cleaned and ready to use. 1.00 mm thick. 3" x 1". Packaging: 10 gross/case.

N467300 Plain Microscope Slides
N467311 Frosted Microscope Slides

REDISCREEN™ DRUGS OF ABUSE CASSETTE

The CLIA Waived, 6 panel test includes individually wrapped test cassette, dropper and test for AMP/COC/THC/M-AMP/OPI/PCP in urine.

2060 Drug of Abuse, CLIA Waived 6 Test Cassette, 25/bx



REDISCREEN™ REDI CUP™ DRUGS OF ABUSE - DOA

Cut-off levels meet SAMSHA standards. Performance comparable to laboratory screening methods. Clear and easy-to-read results. Self-contained cup minimizes collector exposure to urine. Results in five minutes. Flexible drug test combinations. Collect specimen and secure lid. Invert cup and place on a flat surface for 30 seconds. Place cup upright on flat surface. Read results in five minutes.

2054 Drugs of Abuse, 6 Test Includes M-AMP, COC, THC, OPI BZO & PCP, 25/bx
2055 Drugs of Abuse, 10 Test Includes: M-AMP, COC, THC, OPI BZO, PCP, TCA, MTD, AMP & BAR, 25/bx



REDISCREEN™ URINALYSIS REAGENT STRIPS - URS 10

Qualitative and semi-quantitative results. Screen for carbohydrate, metabolism, kidney and liver function and bacteria. Clear and easy-to-read results. Simultaneous determination of up to 10 parameters. One step, no added reagents. Dip test strip in urine sample for one second. Draw strip across the rim of the container to remove excess urine. Wait 30-60 seconds. Compare test strip to color scale.

2001 Urinalysis Reagent Strips, 10 Test Includes: Blood, Urobilinogen, Bilirubin, Protein, Nitrite, Ketones, Glucose, Specific Gravity, pH & Leukocytes, CLIA Waived, 100/bx



REDISCREEN™ FECAL OCCULT BLOOD TEST

2006 Fecal Occult Blood Test, Immunoassay Cassette Device, CLIA Waived, 20/bx



2080

2081

REDISCREEN™ ULTRA GUAIAIC OCCULT BLOOD TEST & ACCESSORIES

2080 Patient Pack, comes with 100 Triple Slide Tests, 6 x 10mL Developers, Applicators, Foil Lined Envelopes, CLIA Waived, 50/pk, 2 pk/bx
2081 Lab Pack, comes with 100 Single Slide Tests, 2 x 10mL Developers, Applicators, CLIA Waived, 100 test/bx
2084 Office Pack, comes with 34 Triple Slide Tests, 2 x 10mL Developers, Applicators, 34/bx

Accessories

2082 8 x 10mL Developers and Applicators, CLIA Waived, 1/bx
2083 2 x 100 Test Tape, 2 x 10mL Developers, Applicators, CLIA Waived, 200/bx

REDISCREEN™ H. PYLORI DEVICE

Individually packaged devices with disposable sample droppers, 20 disposable heparinized capillary tubes and dispensing bulb, 1 each Positive & Negative Control, 1 sample buffer.

2059 CLIA Waived, 20 tests/bx



REDISCREEN™ STREP A TEST THROAT SWAB

97% sensitivity. Equivalent results to a traditional swab culture. Add sample, wait five minutes and read results. Kit includes reagents, controls, tubes, swabs, and test strips. Remove strip from foil pouch. Add reagents to tube. Swab throat and add to reagent solution. Place test strip into solution. Read results in five minutes.

2033 Strep A Test Throat Swab Dipstick, CLIA Waived, 25/bx



REDISCREEN™ URINE HCG PREGNANCY CASSETTE DEVICE

99% clinical accuracy. 25 mIU/ml sensitivity. Some positive results can be seen in as soon as one minute. One step, no reagents. Remove cassette from foil pouch and lay on flat surface. Add three drops of urine to the sample well. Read results in three minutes for urine.

2036 Urine hCG Pregnancy Cassette Device, CLIA Waived, 25/bx
2036-100 Urine hCG Pregnancy Cassette Device, 25/bx, 4bx/cs



REDISCREEN™ URINE HCG PREGNANCY TEST STRIPS

99% clinical accuracy. 25 mIU/ml sensitivity. Some positive results can be seen as soon as one minute. One step, no reagents. Remove strip device from canister. Immerse the strip in urine sample. Read results in three minutes.

2032 Urine hCG Pregnancy Test Strips, CLIA Waived, 100/bx



REDISCREEN™ URINE/SERUM HCG PREGNANCY CASSETTE DEVICE

99% clinical accuracy. 25 mIU/ml sensitivity. Some positive results can be seen in as soon as one minute. One step, no reagents. Remove cassette from foil pouch and lay on flat surface. Add three drops of urine or serum to the sample well. Read results in three minutes for urine and five minutes for serum.

2012 Urine/Serum hCG Pregnancy Cassette Device, CLIA Waived, 25/bx

StarLine®

*Quality Products
Manufactured by
Quality People,
Available through
Quality Distribution*



STARLINE 4-LEG IV STAND

Posts and extensions are constructed of chrome plated steel tubing. Extensions are secured with a sturdy friction knob and allow height adjustment of 47" to 84 1/2". Epoxy coated steel base and twin wheel casters provide ease of transport. No loss knob.

11353001 IV Stand, 2-Rams Horn Hooks



STARLINE 5-LEG IV STAND

Epoxy coated steel base with twin wheel casters for ease of transport. Posts are constructed of 1" diameter chrome plated steel tubing. Extensions are made with 7/8" OD chrome plated steel tubing and are secured with a sturdy friction knob. No loss knob. Height adjustment is 52 1/2" to 93 3/4".

11353501 IV Stand, 2-Rams Horn Hooks

11353601 IV Stand, 4-Rams Horn Hooks



STARLINE I.V. STAND

Heavy bar steel base to achieve the lowest possible center of gravity. Base is epoxy coated to resist corrosion. Ram's horn type hooks prevent accidental spilling. Height adjusts 51.75"-94". No loss knob

11350403 Two-Hook Model I.V. Stand



STARLINE 18" ROUND HAMPER STAND

Made of chrome-plated steel with satin finish on the lid. The foot-pedal mechanism is activated by a spring for easy, hands-free operation. The M-shaped base design makes this hamper very stable and allows the four ball-bearing swivel casters to roll easily over different floor surfaces.

11353334 Hamper Stand, 18" Round



Cart shown with optional basket

STARLINE MULTIPURPOSE CART

Suited for almost any task. Carries two 6" wire baskets or one 12" basket (not included). Folds flat for storage. Chrome plated steel frame. No assembly required

11353400 Cart
11353020 6" Basket, 2/ctn



STARLINE OVERBED TABLE

Table features walnut laminate top with vinyl trim. Spring loaded mechanism allows quick height adjustment from 28"-43". Sturdy chrome four point base features dual wheel nylon casters for easy positioning. Top measures 15" x 30".

11350600 Walnut Laminate Overbed Table



STARLINE UTILITY CART

Plastic, all purpose and sturdy to accommodate almost all sensitive equipment and disposable supplies. 200 lb capacity. Three ABS plastic shelves and 3" swivel casters.

11353530 Utility Cart, 27 1/4" X 17" Shelf Size, 35 1/2"H, 11" Between Cart Shelves #65350



STARLINE STEP ON METAL WASTE CAN

Steel outer can and lid provide secure waste disposal. Features removable plastic inner liners for easy cleaning. Step-on foot pedal opens lid hygienically. White enamel finish complements any office decor.

11350266 32 Qt Waste Can



11350220

11352002

STARLINE STEP STOOL

9" high single step stool is ideal for examination rooms, x-ray, and other areas where sturdy support is required. Stool legs are constructed of 18 gauge tubular steel with chrome plated finish. Topped with a skid resistant mat.

11350220 Step Stool w/handrail
11352002 Step Stool, 2/ctn



STARLINE GOOSENECK EXAM LAMP

Durable chrome finish. Adjustable height 50"-74". 12" flexible gooseneck arm allows easy, one-handed 360 degree adjustment. Uses one 10-40 incandescent watt bulb (not included).

11350500 Adjustable 50"-74" Exam Lamp



11350001

11350002

STARLINE PNEUMATIC STOOL

Features a 5-leg, die brushed aluminum 23" diameter base. The base has a satin finish top with complementary black sides and black/chrome hooded casters. Equipped with a smooth operating air cylinder lift mechanism, operated by a convenient semi-ring under seat. Seat cushion is 4" thick, 16" diameter. 40 standard colors available, plus a variety of special order colors.

11350001 5-Leg Pneumatic Stool
11350002 5-Leg Pneumatic Stool w/backrest



STARLINE VALUE PLUS PNEUMATIC STOOL

Pneumatic height adjustment of 17 1/2" - 22 1/2". Features 5-leg, 21" diameter composite base. Poly foam seat cushion is 4" thick, 15" diameter. Available in upholstery colors to match any exam room.

11350800 5-Leg Pneumatic Stool



STARLINE THREADLESS SCREW ADJUSTABLE STOOL

Thread height adjustment. 15" diameter. 4" thick poly-foam seat cushion. 21" diameter black composite base.

11350214 5-Leg Stool



11350501

11350511



11350521

STARLINE SCREW ADJUSTABLE REVOLVING STOOLS

Features 4" thick polyfoam cushions upholstered in durable vinyl for long lasting use. 14" cushion diameter. 19.5" diameter aluminum base with solid steel, machine thread screw for height adjustments of 20"-28". Dual wheel casters. 40 standard colors available, plus a variety of special order colors.

- 11350501 4-Leg Stool
- 11350511 4-Leg Stool, w/round footing
- 11350521 4-Leg Stool, w/square footing



COURIER® WHEELCHAIRS STANDARD

The perfect combination of both quality and economy. Meets RESNA standards for manual wheelchairs in product performance. Low maintenance, dual sealed precision bearings. Heavy-duty inner liners keep seat and back from stretching. 16-gauge crossbraces add strength and durability. 250 lbs weight capacity. Features 8" solid rubber casters, 24" urethane rear wheels, chrome handrim, mid-night blue vinyl upholstery, padded upholstered arms, carry pocket and heavy duty inner liner upholstery. Limited lifetime warranty. Front riggings sold separately on the swingaway adult frames (N51050 or N51052).

- N51001 16" Wide Wheelchair w/fixed frame, full permanent arm & aluminum footplates
- N51002 18" Wide Wheelchair w/fixed frame, full permanent arm & aluminum footplates
- N51003 16" Wide Wheelchair w/swingaway adult frame & full permanent arms
- N51004 18" Wide Wheelchair w/swingaway adult frame & full permanent arms
- N51005 16" Wide Wheelchair w/swingaway adult frame & full removable arms
- N51006 18" Wide Wheelchair w/swingaway adult frame & full removable arms
- N51007 16" Wide Wheelchair w/swingaway adult frame & desk removable arms
- N51008 18" Wide Wheelchair w/swingaway adult frame & desk removable arms
- N51009 20" Wide Wheelchair w/swingaway adult frame & full removable arm
- N51010 20" Wide Wheelchair w/swingaway adult frame & desk removable arm



COURIER® WHEELCHAIRS HEAVY DUTY

Features superior durability, rollability and streamlined looks. It is designed for the active individual who needs the strength of a reinforced heavy-duty frame. Meets RESNA standards for manual wheelchairs in product performance. Heavy-duty inner liners keep seat and back from stretching. Features swingaway adult frame, 8" solid casters, 24" urethane rear wheels, chrome handrim, blue vinyl upholstery, padded upholstered armrests and heavy duty inner liner. 350 lbs weight capacity. Dual-axle frame allows for easy adjustment for adult (19 1/2") or hemi (17 1/2") seat-to-floor height. Lifetime limited warranty. Front riggings sold separately (N51050 or N51052).

- N52027 22" Wide Wheelchair w/full removable arm
- N52028 22" Wide Wheelchair w/desk removable arm



COURIER® WHEELCHAIRS HEMI

A lower seat to floor height of 17 1/2". Meets RESNA standards for manual wheelchairs in product performance. Low maintenance, dual sealed precision bearings. Flame retardant upholstery is mildew and bacteria resistant. 250 lbs weight capacity. Features swingaway adult frame, 8" solid rubber casters, 24" urethane rear wheels, chrome handrim, midnight blue vinyl upholstery, padded upholstered arms, carry pocket and heavy duty inner liner upholstery. Lifetime limited warranty. Front riggings sold separately (N51050 or N51052).

- N51023 16" Wide Wheelchair w/full removable arm
- N51024 18" Wide Wheelchair w/full removable arm
- N51025 16" Wide Wheelchair w/desk removable arm
- N51026 18" Wide Wheelchair w/desk removable arm



COURIER® WHEELCHAIRS LIGHTWEIGHT

The perfect wheelchair for both short and long term use. This versatile, lightweight frame meets RESNA standards for manual wheelchairs in product performance. Low maintenance, dual sealed precision bearings. Urethane upholstery is mildew and bacteria resistant. 250 lbs weight capacity. Features swingaway adult frame, 8" solid rubber casters, 24" urethane rear wheels, composite handrim, lightweight black nylon upholstery, padded upholstered arms and heavy duty inner liner upholstery. Limited lifetime warranty. Front riggings sold separately (N51050 or N51052).

- N51019 16" Wide Wheelchair w/full removable arm
- N51020 18" Wide Wheelchair w/full flip back arms
- N51021 16" Wide Wheelchair w/desk removable arm
- N51022 18" Wide Wheelchair w/desk flip back arms

COURIER® WHEELCHAIR ACCESSORIES

- N51050 Footrests w/aluminum footplate (For standard chairs)
- N51052 Elevating Legrests w/aluminum footplate & padded calf pads w/black upholstery & hemi spacing (For standard, lightweight & hemi chairs)

