



STARLINE DIGITAL THERMOMETER

Accurately records body temperature in 60 seconds. Liquid crystal display readout displays temperature within 1/10th of a degree. Audible beep signals peak temperature. A replaceable battery provides 200 hours of use. Each package comes complete with thermometer, 5 probe covers, battery, storage case, instructions, and one year limited warranty. Can be used orally, rectally or under the arm. Latex free.

10420410 Digital Thermometer w/beeper, °F

ADC® ADTEMP™ I DIGITAL THERMOMETER

An electronic thermometer designed and priced to be a viable replacement for mercury filled glass thermometers, which are gradually being phased out of healthcare. This new model features a striking translucent housing with integral case. It is dual-scale (C/F) switchable, water resistant and offers a last reading memory. It does not have a beeper, but it will lock in the reading once measurement is complete. This makes an ideal item as a true mercury thermometer replacement. Features a storage case, long life battery and auto-off feature.

412-00 Single Patient Use Digital Thermometer, Bulk, 20/pk
 412K Single Patient Use Digital Thermometer SPU Kit

ADC® ADTEMP™ II DIGITAL THERMOMETER

Features a water resistant design, last reading memory and is dual scale (C/F) switchable. Replaces the 410 and 6001 series thermometers.

413B Oral Digital Thermometer
 413B-00 Oral Digital Thermometer, 25/pkg
 413BK Oral Digital Thermometer SPU Kit
 413R Rectal Digital Thermometer
 413R-00 Rectal Digital Thermometer, 25/pkg
 413RK Rectal Digital Thermometer SPU Kit

ADC® ADTEMP™ III DIGITAL THERMOMETER

An electronic thermometer with a flexible probe designed and priced to be a viable replacement for mercury filled glass thermometers, which are gradually being phased out of healthcare. The flexible probe provides greater patient safety and is ideal for pediatric applications. This new model features a striking translucent housing with integral case. It is dual-scale (C/F) switchable, water resistant and offers a last reading memory. Beeps when measurement is complete. Features a storage case, long life battery and auto-off feature.

414 Digital Thermometer Kit, Flex Tip, Retail Boxed
 414K Digital Thermometer Kit, Flex Tip, Boxed w/5 sheaths


ADC® ADTEMP™ IV DIGITAL THERMOMETER

415 Digital Thermometer
 415K Digital Thermometer SPU Kit

ADC® ADTEMP™ V DIGITAL THERMOMETER

Represents state-of-the-art in micro digital fever thermometers. This new unit features a stylish European design, extra large display and flexible tip. It is water resistant, dual scale (C/F) switchable and has a last reading memory. The response time is 10 seconds rectal and 18 seconds oral. Beeps when measurement is complete. Features a storage case, long life battery and auto-off feature.

418 Super Fast Read Digital Thermometer, Flex Tip, Boxed w/5 sheaths



ADC® ADTEMP™ THERMOMETER SHEATHS

Disposable design eliminates cross contamination. Cost effective. Easy to apply. Suitable for oral or rectal use. Digital model designed for use with ADC® and most comparably sized digital thermometers. Glass model designed for use on all glass thermometers. Made in the U.S.A.

416-100 ADTEMP™ Digital Sheaths, 100/bx
 416-1000 ADTEMP™ Digital Sheaths, 1000/bx
 416-50 ADTEMP™ Digital Sheaths, 50/bx

BD THERMOMETERS

522862 Basic Digital Thermometer, 12/bx

Accessories

522870 Probe Covers For Oral Thermometer, 50/sp, 20 sp/cs
 522871 Probe Covers For Rectal Thermometer, 50/sp, 20 sp/cs



OMRON COMPACT DIGITAL THERMOMETER

No fear of broken glass or mercury poisoning. Fast, accurate reading in 60 seconds. Premium long life battery (up to 400 hours of continuous use). Waterproof for easy cleaning.

MC-102 Digital Thermometer (Spanish/Portuguese) Centigrade
 MC-110B Digital Rectal Thermometer, 20/bx



MC-130MO



MC-130MR

OMRON DIGITAL FEVER THERMOMETER

Oral or rectal available. Quick accurate temperature in 60 seconds. Beeps at peak temperature. Memory recalls last temperature. Meets ASTM (E1112) standard. Waterproof for easy cleaning.

MC-130MO Digital Oral Thermometer
 MC-130MR Digital Rectal Thermometer



OMRON DIGITAL FLEXIBLE THERMOMETER

Flexible tip for added comfort and safety. Beeps at peak temperature. Memory recalls last temperature. 20 second reading. Meets ASTM (E1112) standard.

MC-206 Flexible Digital Thermometer

OMRON DIGITAL RIGID THERMOMETER

3 times faster than traditional digital thermometers. Features include jumbo display for easy reading, memory and long-life 3-volt battery. For use oral, underarm and rectal locations. Safe, no fear of glass breakage or mercury. Includes 5 probe covers.

MC-106N Rigid Digital Thermometer

OMRON DIGITAL THERMOMETER & PROBE COVERS

Individually wrapped disposable probe covers fit securely over the sensor portion on any Omron or Marshall thermometer.

0007PROB	Probe Covers, 50/bx
0009PROB	Probe Covers, 100/bx
600RPROBE	Probe Covers For MC-600, 25/bx
601MPROB500	Probe Covers For MC-601, 500/bx

ADC® ADTEMP™ 422 VETERINARY DIGITAL THERMOMETER

Auto-off function conserves battery life. Range 90°F ± 2°F or 32°C - 43.9°C ± 1°C depending upon scale selection. Replaceable 1.55v (LR41) type battery provides up to 1,500 measurements. Specially designed housing with lanyard attachment for enhanced removal. Integral carry case.

422 Adtemp™ 422 Veterinary Digital Thermometer

ADC® ADTEMP™ HYPOTHERMIA THERMOMETER

The first digital extended range thermometer. Auto off function conserves battery life. Replaceable 1.55v type battery provides up to 1,500 measurements. Clear two-piece plastic storage case.

419 Digital Hypothermia Thermometer



OMRON BASAL DIGITAL THERMOMETER

Identify ovulation timing. Includes temperature tracking chart. Beeps at peak temperature. Waterproof for easy cleaning. Measures in 60 seconds and is three times faster than basal glass thermometers. Meets ASTM (E1112) standard.

MC-301N Basal Digital Thermometer



TIDI SANITHERM THERMOMETER SHEATHS

20521	Mercury Oral Thermometer Sheaths, 100/bx, 50 bx/cs
20524	Mercury Oral Thermometer Sheaths, 500/bx, 10 bx/cs
20525	Mercury Oral Thermometer Sheaths, 5000/cs
20538	Mercury Oral Thermometer Sheaths, Bulk, 11000/cs
20631	White Digital Oral Thermometer Sheaths, Bulk, 11000/cs
20633	Digital Oral Thermometer Sheaths, 100/bx, 50 bx/cs
20634	Digital Oral Thermometer Sheaths, 500/bx, 10 bx/cs
20636	Digital Oral Thermometer Sheaths, 5000/cs
20657	Digital Oral Thermometer Sheaths, Packed in Sheets of 10, Indexed Every 50, Bulk, 10000/cs
20710	Prelubricated Mercury Rectal Thermometer Sheaths, 100/bx, 50 bx/cs
20713	Prelubricated Mercury Rectal Thermometer Sheaths, 500/bx, 10 bx/cs
20714	Prelubricated Mercury Rectal Thermometer Sheaths, 5000/cs
20775	Digital Rectal Thermometer Sheaths, 100/bx, 50 bx/cs
20777	Digital Rectal Thermometer Sheaths, 500/bx, 10 bx/cs



8884221000



8884309800

KENDALL FILAC F-1500 ELECTRONIC THERMOMETER

Large, clear LCD temperature display. Automatic display. Auto Off feature. 30-35 second temperature reading. Audible signal. In-house two point calibration. Two lithium batteries. Lightweight.

8884309800	Thermometer w/oral probe
8884103650	Oral Probe Only
8884103660	Rectal Probe Only

Accessories

8884221000	Probe Covers, 20/pkg, 25 pkg/bx, 10 bx/cs
------------	---



KENDALL FILAC FASTEMP ELECTRONIC THERMOMETER

4-10 second reading. Isolation chambers. Global iconography. Oversized backlit LCD display. Fahrenheit/Celsius button. Last temperature recall. 30 second pulse timer. Temperature range of 85°F-110°F/30C-43C. Electronic anti-theft system. ASTM and CEN compliance.

202000	Complete Oral/Ax System Including Electric Thermometer, 4 Lang Ops Man, Blue Oral/Ax Iso Chamber, Blue Probe w/4' cord, Wall Mountable Base, Mounting Hardware
202018	Fastemp Monitor, Including 4 Language Operating Manuals
202026	Fastemp Blue Oral/Axillary Probe w/blue 4' cord
202027	Fastemp Blue Oral/Axillary Probe w/blue 9' cord
202028	Fastemp Blue Oral/Axillary Isolation Chamber
202034	Fastemp Red Rectal Probe w/red 4' cord
202036	Fastemp Red Rectal Isolation Chamber
202037	Fastemp Red Rectal Probe w/red 9' cord
202038	Fastemp Red Rectal Isolation Chamber, Red Rectal Probe w/red 4' cord, Wall Mountable BS, Hardware
202099	Fastemp Calibration Plug
2020PC	Primary Care Kit Includes: FastTemp Thermometer, Blue Oral/Axillary Isolation Chamber, Blue Oral/Axillary Probe w/4' cord, Wall Mountable Base, Mounting Hardware, (500) Probe Covers




KENDALL FASTEMP PROBE COVERS

Disposable, single-use. Smooth, rigid plastic. Dispenser package fits in the isolation chamber. Contains no latex or natural rubber.

202020	Probe Covers for FASTEMP Thermometer, 20/pk, 24 pk/bx, 10 bx/cs
--------	---



01690-200



01690-400

**WELCH ALLYN SURETEMP® PLUS
ELECTRONIC THERMOMETER**

- 01690-200 Model 690 Electronic Thermometer, Oral Probe w/oral probe well, 2-Year Limited Warranty
- 01690-201 Model 690 Electronic Thermometer, Rectal Probe w/rectal probe well, 2-Year Limited Warranty
- 01690-300 Model 690 Electronic Thermometer, One Per Room, 9 ft Oral Probe w/oral probe well, 2-Year Limited Warranty (must be locked to the wall)
- 01690-301 Model 690 Electronic Thermometer, One Per Room, 9 ft Rectal Probe
- 01690-400 Model 690 Electronic Thermometer, Wall Mount, Oral Probe w/oral probe well, 2-Year Limited Warranty
- 01690-700 Model 690 Electronic Thermometer, Wall Mount, 9 ft Oral Probe w/oral probe well, rolling stand, spare probe & well holde, 2-Year Limited Warranty
- 01692-200 Model 692 Electronic Thermometer, Wall Mount, Security System w/ID location field, Oral Probe w/oral probe well, 3-Year Limited Warranty
- 01692-201 Model 692 Electronic Thermometer, Wall Mount, Security System w/ID location field, Rectal Probe w/rectal probe well, 3-Year Limited Warranty
- 01692-300 Model 692 Electronic Thermometer, One Per Room, Wall Mount, 9 ft Oral Probe w/oral probe well, 3-Year Limited Warranty (must be locked to the wall)
- 01692-700 Model 692 Electronic Thermometer, Wall Mount, 9 ft Oral Probe w/oral probe well, rolling stand, spare probe & well holder, 3-Year Limited Warranty



WELCH ALLYN SURETEMP® THERMOMETER

Oral response time: approx. 4 seconds. Axillary response time: approx. 10 seconds. Rectal response time: approx. 15 seconds. Continuous temperature reading. Measurement range: 84.0°F to 108.0°F. Unit dimensions: 2.25" x 3" x 7". Weight: 11.5 oz. Three AA alkaline batteries included. Three-year limited warranty on instrument. One-year limited warranty on probe.

- 01692-301 Model 692 SureTemp Plus Thermometer, One Per Room, Wall Mount, 9 ft Rectal Probe w/rectal probe well (must be locked to the wall)

**Contact your
representative if you
need more information.**



02678-000



05031-101

WELCH ALLYN SURETEMP® THERMOMETER ACCESSORIES

- 02678-000 4 ft Oral/Axillary Probe
- 02678-100 9 ft Oral/Axillary Probe
- 02679-000 4 ft Rectal Probe
- 02679-100 9 ft Rectal Probe
- 05031-101 Disposable Probe Covers, 1000/cs
- 06133-000 Calibration Key For M600
- 06137-000 Calibration Key For 767T, M678, M679
- 08043-200 Spare Probe Holder (does not include probe)
- 08070-000 Stand For SureTemp Units
- 08071-000 Infared Stand For SureTemp Units
- 08072-000 Thermometer Stand For Model M678
- 21004-0000 Additional Wall Mount
- 21004-1000 SureTemp Wall Mount w/security magnet
- 21004-2000 Wall Holder Assembly
- 21012-0000 Neck & Wrist Straps
- 21013-0000 Wall Holder For 250 Probe Covers

ADC® ADTEMP™ TYMPANIC THERMOMETER

Measures core temperature instantly from the ear reliably and accurately. Patented technology scans 512 times per second. Probe cover free design is less costly and more convenient to use. Large 6mm display is easy to read. Ultra compact ergonomic design fits comfortably in operator's hand and patient's ear. Fever alarm at 100.4°F. Complete with lithium battery, zippered nylon storage/travel case, and operating instructions.

- 421 Tympanic Thermometer



8884810055



8884830004

KENDALL GENIUS MODEL 3000A TYMPANIC THERMOMETER

Takes an accurate, non-invasive temperature. Peak Select System records and analyzes measurements. Scan button instantly activates thermometer. Release button ejects disposable probe. LCD menu prompt window. Incorporates 15-second pulse timer.

- 8884830004 Genius Model 3000A Thermometer
- Accessories**
- 8884810055 Cassette of 21 Probe Covers, 5 cassettes/pkg, 10 pkg/bx, 2 bx/cs
- 8884813711 Side Utility Basket for Mobile Thermometer Stand, 1/cs
- 8884814065 Wall Mount
- 8884830301 Mobile Thermometer Stand, 1/cs

OMRON MC-514 EAR THERMOMETER

Instant, one second reading (takes 250 scans per second). Convenient storage case keeps thermometer and lens covers together. Engineered button placement ensures, safe comfortable use. Sleek design fits comfortably in the ear and in the palm of any size hand. Meets accuracy standards of ASTM (E 1965-98). Memory feature recalls last temperature reading. Measuring mode reads in Fahrenheit and Celsius. Lens cover detection for hygienic use. Five-year warranty.

MC-514 Ear Thermometer



OMRON INSTANT EAR THERMOMETER LENS COVER

MC-509NLC Thermometer Lens Covers, 40/bx

OMRON LENS COVER FOR GENTLE TEMP® INSTANT EAR THERMOMETER

MC-505LC Lens Covers For Ear Thermometer, 40/bx



WELCH ALLYN PRO4000 EAR THERMOMETER

- 04000-200 Braun ThermoScan Pro4000 Ear Thermometer w/3 yr limited warranty
- 04000-300 Braun ThermoScan Pro4000 Ear Thermometer w/security base station and 3 yr limited warranty
- 04000-400 Braun ThermoScan Pro4000 Ear Thermometer w/roll stand and 3 yr limited warranty

Accessories

- 01802-110 Braun Black Body Calibration Device
- 24001-0000 Rolling Thermometer Stand and Rechargeable Base Station
- 24001-1000 Recharging Base Station
- 53009-000 Replacement Alkaline Battery
- 53020-0000 Braun ThermoScan Pro4000 Rechargeable Battery Pack
- 80199-0000 Braun ThermoScan Pro4000 Replacement Power Cord for Base Station
- 85378-0000 Rolling Stand Mounting Bracket (Mounting Bracket for WA Universal Stand)


WELCH ALLYN EAR THERMOMETER ACCESSORIES

- 05075-005 Probe Covers, 5000/cs
- 05075-800 Probe Covers, 800/cs
- 08075-000 Pro3000 Security Rolling Stand w/security base
- 23001-0000 Wall Holder & Probe
- 23001-1000 Wall Holder
- 80195-0000 Security Tether

WELCH ALLYN 9020 TYMPANIC THERMOMETER ACCESSORIES

While the 9020 Tympanic Thermometer is no longer available, certain accessories for it are.


05051-001 Probe Covers, 1000/bx



3M™ TEMPA-DOT™ AX THERMOMETERS

Single-use and disposable with an adhesive strip to hold itself in place in the underarm area when measuring axillary body temperature. Individually wrapped.


5001 Fahrenheit Thermometer, 80/bx, 20 bx/cs
5002 Celsius Thermometer, 80/bx, 20 bx/cs



3M™ TEMPA-DOT™ PLUS THERMOMETERS

Single-use and disposable pre-assembled in a sheath for use in measuring rectal body temperature.

5192 Fahrenheit Thermometer, 100/bx, 20 bx/cs
5502 Celsius Thermometer, 100/bx, 20 bx/cs



3M™ TEMPA-DOT™ SINGLE-USE THERMOMETERS

Single-use and disposable thermometer that can measure oral or axillary body temperature.

5122 Fahrenheit Thermometer, Sterile, Individually Wrapped, 100/bx, 20 bx/cs
5125 Fahrenheit Thermometer, Sterile, Individually Wrapped, 100/bx, 5 bx/cs
5126 Fahrenheit Thermometer, Sterile, Individually Wrapped, 24/evn, 84 env/cs
5147 Fahrenheit Thermometer, Non-Sterile, Unwrapped, 250/bx, 24 bx/cs
5532 Celsius Thermometer, Sterile, Individually Wrapped, 100/bx, 20 bx/cs
5557 Celsius Thermometer, Non-Sterile, Unwrapped, 250/bx, 24 bx/cs



STARLINE ANEROID SPHYGMOMANOMETER

Black enamel 300mmHg no-pin stop manometer. Range markings to facilitate selection of correct cuff size. Artery label and gauge holder. Latex inflation bladder and bulb. Chrome plated deflation valve. Complete with leatherette zippered carrying case and operating instructions. Assembled, inspected and packaged in the U.S.A. Ten year calibration warranty.

- 104277311 Adult Aneroid Sphygmomanometer, Gray
- 104277551 Adult Aneroid Sphygmomanometer, Blue



STARLINE POCKET ANEROID SPHYGMOMANOMETER

This is our "standard" sphygmomanometer. Features are the same as the premier aneroid series, except for a standard bulb/valve and bladder combination.

- 104230941 Adult Sphygmomanometer, Black, Latex Free
- 104246571 Adult Sphygmomanometer, Dark Green, Latex
- 104276311 Adult Sphygmomanometer, Gray, Latex
- 104276351 Adult Sphygmomanometer, Royal Blue, Latex
- 104276371 Adult Sphygmomanometer, Green, Latex
- 104276391 Adult Sphygmomanometer, Red, Latex
- 104276431 Adult Sphygmomanometer, Orange, Latex
- 104276461 Adult Sphygmomanometer, Magenta, Latex
- 104276561 Adult Sphygmomanometer, Burgundy, Latex
- 104276591 Adult Sphygmomanometer, Teal, Latex
- 104276601 Adult Sphygmomanometer, Purple, Latex
- 104276611 Adult Sphygmomanometer, Black, Latex
- 104276612 Newborn Sphygmomanometer, Black, Latex
- 104276613 Infant Sphygmomanometer, Black, Latex
- 104276614 Child Sphygmomanometer, Black, Latex
- 104276615 Large Adult Sphygmomanometer, Black, Latex
- 104276616 Thigh Sphygmomanometer, Black, Latex
- 30940 Adult Sphygmomanometer, Navy, Latex
- 30942 Child Sphygmomanometer, Black, Latex Free
- 30943 Large Adult Sphygmomanometer, Black, Latex Free



STARLINE PREMIUM ANEROID SPHYGMOMANOMETER

This is our premier sphygmomanometer. Features #10428000 lifetime warranted gauge and top quality inflation system, including a nylon Velcro cuff and highest quality bulb and valve. Complete with leatherette case.

- 104272311 Adult Sphygmomanometer, Gray
- 104272351 Adult Sphygmomanometer, Royal Blue
- 104272371 Adult Sphygmomanometer, Green
- 104272391 Adult Sphygmomanometer, Red
- 104272431 Adult Sphygmomanometer, Orange
- 104272461 Adult Sphygmomanometer, Magenta
- 104272561 Adult Sphygmomanometer, Burgundy
- 104272571 Adult Sphygmomanometer, Dark Green
- 104272591 Adult Sphygmomanometer, Teal
- 104272601 Adult Sphygmomanometer, Purple
- 104272611 Adult Sphygmomanometer, Black
- 104272612 Newborn Sphygmomanometer, Black
- 104272613 Infant Sphygmomanometer, Black
- 104272614 Child Sphygmomanometer, Black
- 104272615 Large Adult Sphygmomanometer, Black
- 104272616 Thigh Sphygmomanometer, Black
- 10428000 Gauge Only, Latex Free
- 30740 Adult Sphygmomanometer, Navy
- 30741 Adult Sphygmomanometer, Black, Latex Free
- 30742 Child Sphygmomanometer, Black, Latex Free
- 30743 Large Adult Sphygmomanometer, Black, Latex Free



STARLINE SPHYG/SPRAGUE KIT

Black enamel 300mmHg no-pin stop manometer. Color matched sprague-type stethoscope. Calibrated nylon cuff. Latex bladder and bulb. Chrome plated valve. Color matched nylon carrying case that is extra large with one large zippered compartment, outside Velcro® pocket, end strap and handy key ring. Case measures 10" x 6" x 3". Assembled, inspected and packaged in the U.S.A. Two year inflation system warranty. 20 year calibration warranty.

- 30761 Black Sphyg/Sprague Kit
- 30762 Burgundy Sphyg/Sprague Kit
- 30763 Dark Green Sphyg/Sprague Kit
- 30764 Magenta Sphyg/Sprague Kit
- 30765 Silver Sphyg/Sprague Kit
- 30766 Teal Sphyg/Sprague Kit
- 30767 Purple Sphyg/Sprague Kit



ADC® E-SPHYG DIGITAL ANEROID SPYGMOMANOMETER

Virtually indestructible design can withstand repeated drops of 5 feet or more. Extra large backlit LCD with 7/8" digits, is easy to read in virtually any light. Splash resistant housing for use in harsh environments. Swivel cuff clip allows the manometer to be re-positioned for easy viewing. Auto off function conserves battery strength. Dual low battery warnings. Recalibration reminder after 10,000 cycles (recalibration is free for the life of the instrument). Systolic assist alerts practitioner to the onset of Korotkoff sounds - makes a great training aid. Deflation rate display promotes proper technique in accordance with AHA guidelines. Pulse rate display at completion of measurement. Manometer works with virtually every cuff on the market. Requires familiar and trusted manual auscultation to obtain blood pressure readings (use of stethoscope). 2 "AAA" batteries (included). May be ordered as a complete sphygmomanometer or gauge only for use with existing inflation systems. Lifetime calibration warranty.

7002	Digital Aneroid, Adult, Black
7002C	Digital Aneroid, Child, Black
7002CLF	Digital Aneroid, Child, Latex Free, Black
7002DB	Digital Aneroid, Adult, Navy
7002DBLF	Digital Aneroid, Adult, Latex Free, Navy
7002I	Digital Aneroid, Infant, Black
7002LF	Digital Aneroid, Adult, Latex Free, Black
7002N	Digital Aneroid, Newborn, Black
7002OR	Digital Aneroid, Adult, Orange
7002ORLF	Digital Aneroid, Adult, Latex Free, Orange
7002RB	Digital Aneroid, Adult, Royal Blue
7002RBLF	Digital Aneroid, Adult, Latex Free, Royal Blue
7002V	Digital Aneroid, Adult, Purple
7002VLF	Digital Aneroid, Adult, Latex Free, Purple
7002X	Digital Aneroid, Large Adult, Black
7002XLF	Digital Aneroid, Large Adult, Latex Free, Black
8002	Digital Aneroid, Gauge Only



ADC® GENERAL PRACTICE MULTICUFF KIT

Adult, child, and large adult ADCUFF™ calibrated nylon cuffs with one tube latex inflation bladders and palm style manometer in a compact storage tray. Rugged gauge features precision crafted German mechanism housed in a patented double shell ABS impact polymer casing. Ambidextrous trigger release air valve. Quick release luer slip connectors. Four section smoke gray plexiglass storage tray features a small 14" x 3" footprint to save on counter space or may be wall mounted with included hardware. Assembled, inspected, and packed in the U.S.A. from foreign components. Two year inflation system warranty. Lifetime calibration warranty.

705GPK	Multicuff Kit, Black, Latex
--------	-----------------------------



ADC® SYSTEM 5™ BLOOD PRESSURE SYSTEM

A multipocketed, heavy duty 1000D nylon attaché style case that keeps each component securely in place yet instantly accessible. Dual zipper closure permits rapid access to the contents, protects against the elements. Measures a compact 10" x 9" x 4" when closed, weighs just 1 1/2 lbs. Compact enough to fit in most emergency bags or trauma boxes. Reinforced carry handle doubles as a hang tab for wall mounting in the open position. Includes infant, child, adult, large adult and thigh ADCUFF™ nylon cuffs featuring the Size Guide™ marking system with a one tube latex inflation bladder. Cuffs equipped with quick release luer slip connector to permit rapid and secure selection of the appropriate cuff. Both systems utilize the revolutionary DIAGNOSTIX™ 804 palm style manometer with precision crafted German mechanism housed in a patented double shell ABS impact polymer casing, Ambidextrous trigger release air valve. Case and cuffs available in color coordinated navy or orange. Assembled, inspected and packaged in the U.S.A. Two year inflation system warranty. Lifetime calibration warranty.

740	5 Cuff Blood Pressure System, Navy
740LF	5 Cuff Blood Pressure System, Navy, Latex Free
746	5 Cuff Blood Pressure System, Orange
746LF	5 Cuff Blood Pressure System, Orange, Latex Free



ADC® DIAGNOSTIX™ 703 SERIES SPHYGMOMANOMETER

Revolutionary palm style manometer (#804) with integral latex free bulb and trigger style release valve for one hand convenience. Precision crafted 300mmHg no-pin stop German mechanism tested to 200,000 cycles. Mechanism housed in a patented double shell ABS impact polymer casing. Ambidextrous infinite control trigger release air valve. Quick release luer slip connector. Six clinical cuff sizes available. ADCUFF™ nylon cuff with Size Guide™ marking system ensures proper cuff positioning and use of correct cuff size. Dip molded inflation bladder withstands heavy use. Complete with leatherette zippered carrying case with nylon zipper. Assembled, inspected, and packaged in the U.S.A. from foreign components. Two year inflation system warranty. Lifetime calibration warranty. Aneroids completely latex free have a LF in their number.

703	Adult Aneroid, Black, Latex
-----	-----------------------------



760GARCQ

ADC® PROSPHYG™ 760 SERIES ANEROID

- 760BOOCQ Child, Aneroid, Betty Boop
- 760GARCQ Child, Aneroid, Garfield
- 760PEACQ Child, Aneroid, Peanuts
- 760BOOQ Adult, Aneroid, Betty Boop
- 760GARQ Adult, Aneroid, Garfield
- 760PEAQ Child, Aneroid, Peanuts
- 760BOOXQ Large Adult, Aneroid, Betty Boop
- 760GARXQ Large Adult, Aneroid, Garfield
- 760PEAXQ Large Adult, Aneroid, Peanuts
- 768-641BOO Adult, Aneroid, Betty Boop
- 768-641GAR Adult, Aneroid, Garfield
- 768-641PEA Adult, Aneroid, Peanuts

ADC® PRO'S COMBO II ANEROID

- 768-641BOO Adult, Aneroid, Betty Boop
- 768-641GAR Adult, Aneroid, Garfield
- 768-641PEA Adult, Aneroid, Peanuts



OMRON ECONOMY ANEROIDS

Calibrated nylon or cotton cuff with artery label and manometer tab. Bold, black gradations on white faceplate. Gray or blue enamel, no stop pin, 300 mm Hg gauge. Latex inflation system. Standard air release valve. Zippered carrying case. 10 year calibration warranty > 3 mm Hg.

- 11-200 Adult Aneroid w/nylon cuff, Blue, Blue Enamel
- 11-200C Child Aneroid w/nylon cuff, Blue, Blue Enamel
- 11-200X Large Adult Aneroid w/nylon cuff, Blue, Blue Enamel



OMRON MARSHALL® SPHYGMOMANOMETER WITH COTTON CUFF

Long-life cotton cuff with artery label and manometer tab. Bold, black gradations on white faceplate. Gray enamel, no stop pin, 300 mm Hg gauge. Latex inflation system. Standard air release valve. Zippered carrying case. Lifetime calibration warranty > 3 mm Hg.

- 0115M Adult Sphygmomanometer, Gray
- 0115MC Child Sphygmomanometer, Gray
- 0115MO Large Adult Sphygmomanometer, Gray



OMRON MARSHALL® SPHYGMOMANOMETER WITH NYLON CUFF

Long-life, calibrated nylon cuff, with artery label and manometer tab. Bold, black gradations on white faceplate. Gray enamel, no stop pin, 300 mm Hg gauge. Latex inflation system. Deluxe air release valve. Zippered, leatherette carrying case. Lifetime calibration warranty > 3 mm Hg.

- 0108M Adult Sphygmomanometer, Black
- 0108MC Child Sphygmomanometer, Black
- 0108MI Infant Sphygmomanometer, Black
- 0108ML Thigh Sphygmomanometer, Black
- 0108MNL Adult Sphygmomanometer, Black, Latex Free
- 0108MO Large Adult Sphygmomanometer, Black



OMRON SELF-TAKING BLOOD PRESSURE KIT

Long-lasting cotton D-ring cuff. One-hand operation. Easy-to-read gauge. Professional stethoscope attached to cuff. Includes vinyl carrying case, illustrated instructions, and blood pressure chart. One year warranty.

- 0104 Adult Blood Pressure Kit, Black
- 0104MAJ Large Adult Blood Pressure Kit, Blue
- 0116 Adult Blood Pressure Kit, Two-Party, Black
- HEM-18 Self Taking Manual Blood Pressure
- HEM-18S Self Taking Manual Blood Pressure Kit w/vinyl bladder



5098-02



5098-18



5098-20

WELCH ALLYN CLASSIC HAND ANEROID

Unsurpassed accuracy and durability in the traditional classic design. For more than 30 years, the Tycos Classic Hand Aneroid has delivered performance and reliability in worldwide clinical use. Precision, jewelled bearings resist wear for consistently accurate readings. Every dial face is individually laser-engraved for exact calibration. Crisp white graphics on large charcoal-gray scale for better readability. Unique phosphor-bronze aneroid bellows ensure consistent pressure measurement. Recessed dial brings needle and dial markings close together and eliminates parallax. Single-tube leuc connector allows for easy interchange of cuff sizes. Lifetime calibration and manufacturing warranty. All latex free.

- 5098-02 Aneroid w/Adult Cuff, Latex Free
- 5098-05 Aneroid w/Child Cuff, Latex Free
- 5098-06 Aneroid w/Infant Cuff, Latex Free
- 5098-11 Aneroid w/Adult Cuff, w/o Case, Latex Free
- 5098-18 Pediatric Blood Pressure Kit w/case, Includes Newborn, Infant & Child Size Cuffs, Latex Free
- 5098-20 Family Practice Kit w/case, Includes Child, Adult & Large Adult Cuffs, Latex Free
- 5098-23 Family Practice Kit w/case, Includes Child Print Cuff, Adult & Large Adult, Latex Free

Accessories

- 5098-03 Gauge & Bulb Only, Latex Free



WELCH ALLYN CLASSIC POCKET ANEROID

Provides protection to the manometer during the rigors of clinical use. Precision, jewelled bearings resist wear for consistently accurate readings. Gauge attached to cuff protects instrument and prevents cuff and gauge from being easily separated. Laser-engraved dial for precise calibration. Unique phosphor-bronze bellows ensure consistent pressure measurement. Crisp, white graphics on large charcoal-gray scale for better readability. Lifetime warranty on calibration and manufacturing quality. All latex free.

- 5090-01 Aneroid w/Adult Hook Cuff, Latex Free
- 5090-02 Aneroid w/Adult Cuff, Latex Free
- 5090-06 Aneroid w/Child Cuff, Latex Free
- 5090-14 Aneroid w/Adult Cuff, w/o case, Latex Free
- 5090-22 Aneroid w/adult conductive cuff, Latex Free
- 5090-37 Aneroid w/Thigh Cuff w/o case, Latex Free
- 5090-41 Aneroid w/Large Adult Cuff, Latex Free

Accessories

- 5090-03 Gauge Only, Latex Free



WELCH ALLYN TYCOS® DS66 HAND ANEROID

More accurate blood pressure readings & lower lifetime ownership cost due to award-winning design that is gear-free & shock-resistant. Welch Allyn's line of aneroid sphygmomanometers featuring gear-free, shock-resistant DuraShock technology offer practical benefits compared to traditional sphygmomanometers w/gears. This is especially true in busy medical environments where instruments may be handled roughly or dropped. With DuraShock technology, Welch Allyn sphygmomanometers will remain in calibration longer, providing more accurate readings than traditional sphygmomanometers w/gears. They won't need to be replaced as often as traditional gauges w/gears, providing lower lifetime ownership cost. Only Welch Allyn aneroids feature gear-free, shock resistant DuraShock technology to deliver unprecedented durability. The only gauge capable of passing AAMI Shock Resistant standard ensures your sphygmomanometer will remain accurately calibrated longer. Lower lifetime cost of ownership since DuraShock gauges don't need replacing as often due to drops & rough handling. Only Welch Allyn offers a jeweled movement w/an industrial sapphire at the key friction point for long product life. Only Welch Allyn calibrates by computer & laser engraves a unique, custom-matched dial face for unprecedented accuracy over the instrument life. Lighter weight than traditional sphygmomanometers so it won't weigh down your lab coat. 1-piece cuff design is more durable than traditional cuff & bladder systems & it eliminates the need to remove or replace bladders. Latex-free cuff is easy to clean, disinfect, or sterilize for safety. 15-year calibration warranty for overall value.

- 5098-27 Aneroid w/Adult Cuff in Zippered Case, Latex Free Cuff
- 5098-28 Aneroid w/Large Adult Cuff in Zippered Case, Latex Free Cuff
- 5098-29 Aneroid w/Child Cuff in Zippered Case, Latex Free Cuff
- 5098-30 Multi-Cuff Kit, Includes Adult, Large Adult & Child-Print Inflation Systems in Zippered Case, Latex Free Cuff



5086-05



5098-70



5098-73

WELCH ALLYN TYCOS® TR-2 HAND ANEROID

Lightweight and compact, the Tycos TR-2 hand aneroid provides portability in an affordably priced, hand-held design. The trigger-release air valve is positioned for optional left or right-handed use and allows for easy fingertip control of cuff deflation. Large easy-to-read dial face with bright graphics. Balanced palm design allows for comfortable use over repeated measurements. Meets AAMI accuracy standard of ± 3 mmHg. Comes in a zippered case. Five year warranty. All latex free.

- 5098-70 Aneroid w/Adult Cuff, Latex Free
- 5098-71 Aneroid w/Large Adult Cuff, Latex Free
- 5098-72 Aneroid w/Child Print Cuff, Latex Free
- 5098-73 Aneroid Multi-Cuff Kit, Includes Adult, Large Adult & Child Print Cuff, Latex Free

Accessories

- 5082-63 Child Print Cuff, Latex Free
- 5086-05 Bulb Assembly, Latex Free
- 5098-79 Crystal, Latex Free




WELCH ALLYN TYCOS® TR-2 HOME ANEROID

Designed for simple and accurate personal blood pressure monitoring at home or work. Built-in stethoscope chestpiece allows for clear detection of blood pressure sounds with minimal training. Easy-to-use trigger deflation valve provides smooth release of air pressure. D-ring cuff for convenient self-application. Large & bright dial face for easy viewing of numbers. Instructional audio cassette and blood pressure recording forms. Meets AAMI accuracy standard of ± 3 mmHg. Five year warranty. All latex free.

- 5098-75 Aneroid w/Adult Cuff, Latex Free
- 5098-76 Aneroid w/Large Adult Cuff, Latex Free

Accessories

- 5082-67 Adult D-Ring Cuff & Bag, Latex Free
- 5082-68 Large Adult D-Ring Cuff & Bag, Latex Free
- 5098-74 Gauge & Bulb Only, Latex Free



750CW



752CM



752M

ADC® DIAGNOSTIX™ 750 SERIES ANEROID

Large, easy to read 6" dial with white scale on contrasting black face for high visibility. Luminescent dial for visibility in low light conditions. Lifetime calibration warranty. Contemporary Euro-style black casing complements any decor. Rear mounting bracket doubles as an inflation system storage compartment. ADCUFF™ inflation system with black nylon cuff, Size Guide™ marking system, premium latex bladder and bulb and 8' (when extended) coiled tubing with luer connectors. Includes all necessary hardware and instructions. Wall bracket swivels through 110° arc for enhanced visibility. Includes an ABS cuff storage compartment in main unit.

- 750CW Child Wall-Mounted Aneroid, Black, Latex
- 750CWLF Child Wall-Mounted Aneroid, Black, Latex Free
- 750IW Infant Wall-Mounted Aneroid, Black, Latex
- 750IWLF Infant Wall-Mounted Aneroid, Black, Latex Free
- 750NW Newborn Wall-Mounted Aneroid, Black, Latex
- 750TW Thigh Wall-Mounted Aneroid, Black, Latex
- 750TWLF Thigh Wall-Mounted Aneroid, Black, Latex Free
- 750W Adult Wall-Mounted Aneroid, Black, Latex
- 750WLF Adult Wall-Mounted Aneroid, Black, Latex Free
- 750XW Large Adult Wall-Mounted Aneroid, Black, Latex
- 750XWLF Large Adult Wall-Mounted Aneroid, Black, Latex Free
- 752CM Child Mobile Aneroid, Black, Latex
- 752CMLF Child Mobile Aneroid, Black, Latex Free
- 752IM Infant Mobile Aneroid, Black, Latex
- 752IMLF Infant Mobile Aneroid, Black, Latex Free
- 752M Adult Mobile Aneroid, Black, Latex
- 752MLF Adult Mobile Aneroid, Black, Latex Free
- 752NM Newborn Mobile Aneroid, Black, Latex
- 752TM Thigh Mobile Aneroid, Black, Latex
- 752TMLF Thigh Mobile Aneroid, Black, Latex Free
- 752XM Large Adult Mobile Aneroid, Black, Latex
- 752XMLF Large Adult Mobile Aneroid, Black, Latex Free

**Can't find what
you are looking for?
Give us a call for fast,
personal service.**



7670-01



7670-03

WELCH ALLYN 767 SERIES WALL & MOBILE ANEROIDS

The Tycos 767 Wall and Mobile Aneroids offer all the dependability and accuracy of our classic aneroids, along with an updated look. Precision, jeweled-movement manometer that is certified to be accurate within one percent of scale ($\pm 3\text{mmHg}$). Laser-engraved dial face for precision calibration. Built-in basket for securely storing the inflation system. High-impact housing and bezel to protect the movement and crystal. Large dial face (6 1/4") and bright orange pointer for easy reading. 40° side-to-side swivel for clear viewing at all angles. The Mobile aneroid has adjustable height stand and locking collar, 4' coiled tubing and 360° manometer swivel for ease of reading. Lifetime calibration and manufacturing warranty.

- 7670-01 Wall Aneroid w/Adult Cuff, Latex Free
- 7670-02 Wall Aneroid Manometer Only w/8 ft tubing, Latex Free
- 7670-03 Mobile Aneroid w/Adult Cuff, Latex Free
- 7670-04 Mobile Aneroid w/weighted base & adult cuff, Latex Free

Accessories

- 7670-06 Mobile Stand & Basket Assembly Only (base weight not included), Latex Free
- 7670-07 Inflation System Basket, White, Latex Free
- 7670-08 Crystal, Latex Free



5091-13



5091-38



5091-41



5091-47

WELCH ALLYN TYCOS® 509 SERIES WALL & MOBILE ANEROIDS

Tycos 509 Series aneroids provide all the accuracy and reliability you expect from Tycos - at an affordable price. Lifetime calibration and manufacturing warranty. Precision, jeweled bearings resist wear for consistently accurate readings. Unique phosphor-bronze aneroid bellows ensure consistent pressure measurement. Visual check of accuracy on the dial. Strong die-cast zinc internal housing protects mechanism. The Wall style has an 8' coiled tubing with 90° swivel for flexible reading, 10° forward tilt to reduce glare. The Mobile style has 360° swivel with adjustable height stand, and 4' coiled tubing.

- 5091-23 Burnt Orange Wall Aneroid Manometer Only w/8 ft tubing, Latex Free
- 5091-38 Charcoal Wall Aneroid w/Adult Cuff, Latex Free
- 5091-39 Charcoal Wall Aneroid Manometer Only w/8 ft tubing, Latex Free
- 5091-41 Charcoal Mobile Aneroid w/Adult Cuff, Latex Free

Accessories

- 5091-13 Inflation System Basket, Small, Stainless Steel, Latex Free
- 5091-47 Inflation System Basket, Large, Black, Latex Free
- 5091-48 Mobile Aneroid Stand & Basket Assembly Only, Latex Free
- 7670-05 Mobile Stand Base Weight Only, Latex Free

ADC® ADVANTAGE™ 6012 SEMI-AUTOMATIC BP MONITOR

Compact and easy-to-use. Measures systolic, diastolic pressure and pulse. Manual inflation. Self adjusting latex free D-bar cuff. Automatic deflation valve. Auto off.

- 6012 Advantage™ 6012 BP Monitor

ADC® ADVANTAGE™ 6013 AUTOMATIC BP MONITOR

Smart Logic™ Technology determines ideal inflation level. Measures systolic, diastolic pressure, and pulse. Automatic inflation and deflation control. Self adjusting latex free D-bar cuff. Recalls last 14 readings with date and time stamp. Auto off. Optional AC adapter (sold separately).

- 6013 Advantage™ 6013 BP Monitor, 9" to 13" Limb Range
- 6013SA Advantage™ 6013 BP Monitor, 7.5" to 9.4" Limb Range
- 6013X Advantage™ 6013 BP Monitor, 13" to 17" Limb Range

Accessories

- 850-6013 Cuff & Bladder, Adult For 6013 BP Monitor, Navy, Latex Free
- 850-6013SA Cuff and Bladder, Adult for 6013SA BP Monitor, Navy, Latex Free
- 850-6013X Cuff & Bladder, Large Adult For 6013 BP Monitor, Navy, Latex Free

ADC® ADVANTAGE™ 6014 ADVANCED BP MONITOR

Smart Logic™ Technology determines ideal inflation level. Measures systolic, diastolic pressure, and pulse. Automatic inflation and deflation control. Average mode automatically averages three readings for more meaningful results. Self adjusting latex free D-bar cuff. Advanced two zone memory with date and time stamp recalls last 30 readings for two different people. Auto off. Optional thermal printer (sold separately). Optional AC adapter (sold separately).

- 6014 Advantage™ 6014 BP Monitor

Accessories

- 6013ZAC AC Adapter For 6013, 6014, 6014P
- 6014P Printer For 6014 & 6016 Monitors
- 6014PAPER Paper For 6014P Printer, 5/pkg

ADC® ADVANTAGE™ 6016 DIGITAL WRIST BP MONITOR

Smart Logic™ Technology determines ideal inflation level. Measures systolic, diastolic pressure, and pulse. New 8 bit microprocessor and advanced algorithm for more accurate readings. Average mode automatically averages three readings for more meaningful results. Latex free wrist cuff fits 5" to 7" wrists. Advanced two zone memory with date and time stamp recalls last 30 readings for two different people. Auto off. Compact storage case. Optional thermal printer (sold separately).

- 6016 Advantage™ 6016 BP Monitor

OMRON AUTOMATIC BLOOD PRESSURE MONITOR WITH COMFIT™ CUFF

Comfit™ cuff is pre-formed for quick and proper fit (9" to 17" arms). Clinically proven accuracy. Technology offers personalized inflation for maximum comfort. 90 memory recall with date and time stamp. Plastic case for convenient storage included. Blood pressure and pulse rate displayed on a large digital panel. Large, clear digital display. AC adapter and 4 AA batteries included. Five-year warranty.

- HEM-780 Automatic Blood Pressure Monitor w/Comfit™ Cuff



OMRON DIGITAL BP MONITOR WITH EZ WRAP CUFF

In the clinical setting, the EZ Wrap Cuff (9"-17") allows the unit to be used on a wider range of arm sizes, often negating the need for separate, larger sized cuffs. For home use, the EZ Wrap Cuff makes applying the cuff much easier for an individual without assistance. The IntelliSense technology allows for comfortable, personalized inflation and accurate results. Two person, 21 memory recall with date and time stamped for each person. AC adapter and 4 "AA" batteries included. Five-year warranty.

- HEM-773AC Digital BP Monitor w/EZ Wrap Cuff



OMRON DIGITAL WRIST BP MONITOR WITH APS

APS (Advanced Positioning Sensor) technology means that unless the wrist is held in the proper position, slightly above the heart, it will not provide a reading, eliminating mixed results due to user technique. 90 memory recall. IntelliSense technology offers personalized inflation for maximum comfort, graphing capabilities with date and time stamp. Five-year warranty.

HEM-637 Digital Wrist BP Monitor w/APS



OMRON INTELLISENSE™ AUTOMATIC PLUS BLOOD PRESSURE MONITOR

IntelleSense™ technology offers automated level setting and electronic deflation control. Completely automatic, operates with the push of a button. Large, clear digital display. 14 memory storage with AC adapter and illustrated instruction manual included. Quickly and accurately measures blood pressure and pulse. Tailored arm cuff is durable and fits for arms 9"-13" in circumference (large cuff, 13" to 17" available). Five year warranty.

HEM-711AC Automatic Inflation Monitor w/fuzzy logic & AC adapter



HEM-422CRLC




HEM-432CN

OMRON REGENCY SERIES MANUAL INFLATION BLOOD PRESSURE MONITOR

Large digital panel with 14 reading memory. Easy squeeze inflation bulb quickly fills cuff with air. Complete with zippered carrying case. Operates on 4 AA batteries (not included).

HEM-422CRLC Manual Inflation Monitor w/large adult cuff
HEM-432CN Manual Inflation Monitor



ADC® DIAGNOSTIX™ 922 SERIES SPHYGMOMANOMETER

Rugged diecast aluminum housing with textured black baked enamel finish complements any decor. Fold down cover with spring lock mechanism. Graduated unbreakable plastic cartridge tube features 4.00mm inside bore for unrestricted, friction-free mercury flow. High visibility red graduations. Filter system keeps mercury clean for reliable performance. Locking reservoir secures mercury during storage, transport, or maintenance. ADCUFF™ black nylon cuff with ADC's proprietary Size Guide™ marking system ensures proper cuff positioning and use of correct cuff size. Dip molded latex inflation bladder withstands heavy use. Latex inflation bulb fitted with filter screen end valve reduces dust build-up. ADFLOW™ chrome plated brass air release valve with filter screen protection permits precise deflation control and years of reliable service. Comprehensive owner's manual. Assembled, inspected, and packaged in the U.S.A. from foreign components. Two year inflation system warranty. Lifetime warranty on cartridge tube.

- 922 Adult Sphygmomanometer, Black, Latex
- 922C Child Sphygmomanometer, Black, Latex
- 922CLF Child Sphygmomanometer, Black, Latex Free
- 922I Infant Sphygmomanometer, Black, Latex
- 922ILF Infant Sphygmomanometer, Black, Latex Free
- 922LF Adult Sphygmomanometer, Black, Latex Free
- 922N Newborn Sphygmomanometer, Black, Latex
- 922T Thigh Sphygmomanometer, Black, Latex
- 922X Large Adult Sphygmomanometer, Black, Latex
- 922XLF Large Adult Sphygmomanometer, Black, Latex Free
- 922TLF Thigh Sphygmomanometer, Black, Latex Free

Accessories

- 922-00 Manometer Only For 922, Latex Free
- 922-126 Plastic Cartridge Tube For 922, Latex Free

ADC® DIAGNOSTIX™ 932 SERIES SPHYGMOMANOMETER

Rugged housing constructed from a combination of virtually indestructible ABS and diecast aluminum. Striking black finish complements any decor. Contemporary upright design with weighted base and integral wire inflation system basket. Graduated unbreakable plastic cartridge tube features a massive 5.00mm inside bore for unrestricted, friction-free mercury flow. White oversized numerals on contrasting black scale combined with high visibility red graduations for easy reading. EZ-Tube™ system allows safe, simplified maintenance of the cartridge tube or top cap filters. Additional safety provided by the locking reservoir which secures mercury for transport or maintenance. ADCUFF™ black nylon cuff with ADC's proprietary Size Guide™ marking system ensures proper cuff positioning and use of correct cuff size. Dip molded latex inflation bladder withstands heavy use. Latex inflation bulb fitted with filter screen end valve reduces dust build-up. ADFLOW™ chrome plated brass air release valve with filter screen protection permits precise deflation control. Comprehensive owner's manual. Assembled, inspected, and packaged in the U.S.A. from foreign components. Two year inflation system warranty. Lifetime warranty on cartridge tube.

- 932 Adult Sphygmomanometer, Black, Latex
- 932C Child Sphygmomanometer, Black, Latex
- 932CLF Child Sphygmomanometer, Black, Latex Free
- 932I Infant Sphygmomanometer, Black, Latex
- 932ILF Infant Sphygmomanometer, Black, Latex Free
- 932LF Adult Sphygmomanometer, Black, Latex Free
- 932N Newborn Sphygmomanometer, Black, Latex
- 932T Thigh Sphygmomanometer, Black, Latex
- 932TLF Thigh Sphygmomanometer, Black, Latex Free
- 932X Large Adult Sphygmomanometer, Black, Latex
- 932XLF Large Adult Sphygmomanometer, Black, Latex Free



ADC® DIAGNOSTIX™ 952B SERIES SPHYGMOMANOMETER

Rugged diecast aluminum housing. Striking black finish complements any decor. Graduated unbreakable plastic cartridge tube (recessed for safety) features 5.00mm inside bore for unrestricted, friction-free mercury flow. Wall bracket swivels through 110° arc for enhanced visibility and viewer comfort. EZ-Tube™ system allows safe, simplified maintenance of the cartridge tube or top cap filters. Locking reservoir secures mercury during transport or maintenance. ADCUFF™ nylon cuff with ADC®'s proprietary Size Guide™ marking system ensures proper cuff positioning and use of correct cuff size. Complete with vinyl coated wire cuff basket, 8 ft length (when extended) of coiled tubing, all mounting hardware and comprehensive owner's manual. Assembled, inspected and packaged in the U.S.A. Two year inflation system warranty. Lifetime warranty on cartridge tube.

- 952B Adult Sphygmomanometer, Black, Latex
- 952BC Child Sphygmomanometer, Black, Latex
- 952BCLF Child Sphygmomanometer, Black, Latex Free
- 952BI Infant Sphygmomanometer, Black, Latex
- 952BILF Infant Sphygmomanometer, Black, Latex Free
- 952BLF Adult Sphygmomanometer, Black, Latex Free
- 952BN Newborn Sphygmomanometer, Black, Latex
- 952BT Thigh Sphygmomanometer, Black, Latex
- 952BTLF Thigh Sphygmomanometer, Black, Latex Free
- 952BX Large Adult Sphygmomanometer, Black, Latex
- 952BXLF Large Adult Sphygmomanometer, Black, Latex Free

Accessories

- 952-025 Cuff Storage Basket For 952
- 952-129 Plastic Cartridge Tube For 952
- 952B-00 Manometer Only

ADC® ADVANCED DIGITAL BP MONITOR WITH PC LINK

- 6017 Advanced Digital BP w/PC link (software not included), Navy
- 6017PC PC Link Software & USB Cable

**Call for
current pricing!**



H-003D



HEM-907-STAND

OMRON DIGITAL BLOOD PRESSURE PARTS & ACCESSORIES

- 004IC Interface Cable RS232 For Model HEM-7471C
- 0080AGRAY D-Ring Contour Cuff For HEM-412C, HEM-712C & HEM-703CP, Gray
- 0090TRP Replacement Roll of Thermal Paper For HEM-705CP, 5 rl/bx
- 0094AGRAY D-Ring Curved Contour Cuff, Gray
- 0097IR Ink Roller For HEM-703CP & 97
- 0097RP Replacement Roll Paper For HEM-703CP & 97, 5/bx
- 0146GRAY Gray Bulb & Valve For HEM-412C & HEM-425C
- H-003D Large Adult Cuff, Gray
- H-003DS Small Size Cuff For IntelliSense Monitors
- HEM-907-BLDRL Bladder Large Cuff For HEM-907 & 907XL
- HEM-907-BLDRM Bladder Medium Cuff For HEM-907 & 907XL
- HEM-907-BLDRS Bladders Small Cuff For HEM-907 & 907XL
- HEM-907-CL19 Cuff & Bladder Set, Large 32-42cm For HEM-907 & 907XL
- HEM-907-CR19 Cuff & Bladder Set, Medium 22-32cm For HEM-907 & 907XL
- HEM-907-CS19 Cuff & Bladder Set, Small 17-22cm For HEM-907 & 907XL
- HEM-907-CUFFL Cuff For Large Bladder 32-42cm For HEM-907 & 907XL
- HEM-907-CUFFM Cuff For Medium Bladder 22-32cm For HEM-907 & 907XL
- HEM-907-CUFFS Cuff For Small Bladder 17-22cm For HEM-907 & 907XL
- HEM-907-PBAT Battery Pack For HEM-907 & 907XL BP Unit
- HEM-907-PKIT Pole Mount For HEM-907 & 907XL BP Unit
- HEM-907-STAND Stand For HEM-907 & 907XL BP Unit
- HEM-907-WKIT Wall Mount Kit For HEM-907 & 907XL BP Unit
- HEM-AC-J HEM-907 & 907XL Adpater
- HEM-ADPT1 AC Adapter For Omron & Marshall Auto Inflate Monitors
- HEM-TUBE-100 Air Tube 1.0m WL For HEM-907 & 907XL
- HEM-TUBE-130 Air Tube 1.3m WL For HEM-907 & 907XL

ADC® ADCUFF™ 1-TUBE BLADDERS

Dip molded latex or latex free. 55cm tubing length. Size and model number printed on bladder. For use with most palm style aneroids. Individually boxed. Two year warranty.

- 820 One-Tube Bladder, Latex, Adult
- 820C One-Tube Bladder, Latex, Child
- 820CLF One-Tube Bladder, Latex Free, Child
- 820I One-Tube Bladder, Latex, Infant
- 820ILF One-Tube Bladder, Latex Free, Infant
- 820LF One-Tube Bladder, Latex Free, Adult
- 820N One-Tube Bladder, Latex, Newborn
- 820T One-Tube Bladder, Latex, Thigh
- 820TLF One-Tube Bladder, Latex Free, Thigh
- 820X One-Tube Bladder, Latex, Large Adult
- 820XLF One-Tube Bladder, Latex Free, Large Adult

ADC® ADCUFF™ 2-TUBE BLADDERS

Dip molded latex or latex free. 55cm tubing length. Size and model number printed on bladder. For use with most aneroid and mercurial blood pressure instruments. Individually boxed. Two year warranty.

810	Two-Tube Bladder, Latex, Adult
810C	Two-Tube Bladder, Latex, Child
810CH	Two-Tube Bladder, Latex, Hospital Pack, Child, 10/pkg
810CHLF	Two-Tube Bladder, Latex Free, Hospital Pack, Child, 10/pkg
810CLF	Two-Tube Bladder, Latex Free, Child
810E	Two-Tube Bladder, Latex, Adult, Standard
810H	Two-Tube Bladder, Latex, Hospital Pack, Adult, 10/pkg
810HLF	Two-Tube Bladder, Latex Free, Hospital Pack, Adult, 10/pkg
810I	Two-Tube Bladder, Latex, Infant
810IH	Two-Tube Bladder, Latex, Hospital Pack, Infant, 5/pkg
810IHLF	Two-Tube Bladder, Latex Free, Hospital Pack, Infant, 5/pkg
810ILF	Two-Tube Bladder, Latex Free, Infant
810LF	Two-Tube Bladder, Latex Free, Adult
810N	Two-Tube Bladder, Latex, Newborn
810NH	Two-Tube Bladder, Latex, Hospital Pack, Newborn, 5/pkg
810T	Two-Tube Bladder, Latex, Thigh
810TH	Two-Tube Bladder, Latex, Hospital Pack, Thigh, 5/pkg
810THLF	Two-Tube Bladder, Latex Free, Thigh, 5/pkg
810TLF	Two-Tube Bladder, Latex Free, Thigh
810X	Two-Tube Bladder, Latex, Large Adult
810XH	Two-Tube Bladder, Latex, Hospital Pack, Large Adult, 10/pkg
810XHLF	Two-Tube Bladder, Latex Free, Hospital Pack, Large Adult, 10/pkg
810XLF	Two-Tube Bladder, Latex Free, Large Adult

STARLINE INFLATION BULBS & VALVES

Dip molded latex or latex free. For use with most aneroid and mercurial blood pressure instruments. End valve included. Individually boxed. Two year warranty.

10428075	Standard Bulb & Valve, Latex
----------	------------------------------

STARLINE 2-TUBE BLADDERS

Dip molded latex or latex free. 55cm tubing length. Size and model number printed on bladder. For use with most aneroid and mercurial blood pressure instruments. Individually boxed. Two year warranty.

10428010	Two-Tube Bladder, Latex, Adult
----------	--------------------------------

STARLINE CALIBRATED NYLON CUFF

Featuring the marking system. For use with most aneroid or mercurial blood pressure instruments. Two year warranty.

104283610	Adult Nylon Cuff, Black, Latex
-----------	--------------------------------

STARLINE CUFFS & BLADDERS


For use with most aneroid or mercurial blood pressure instruments. Individually boxed.

104284610	Adult Cuff, 2-Tube, Black
104285610	Adult Cuff, 1-Tube, Black
30840	Adult Cuff, 2-Tube, Black, Latex Free
30850	Adult Cuff, 1-Tube, Black, Latex Free

STARLINE INFLATION SYSTEM

Includes nylon cuff, deluxe 2-tube bladder, bulb and valve. For use with most aneroid or mercurial blood pressure instruments. Individually boxed.

104286610	Adult Inflation System, Black, Latex
30860	Adult Inflation System, Black, Latex Free



ADC® ADCUFF™ & BLADDERS, 1-TUBE

Durable nylon construction. Exclusive Size Guide™ marking system. Latex or latex free inflation bladder with one 55cm tubes. For use with most palm style aneroid instruments. Individually boxed. Two-year warranty.

850	ADCUFF™ & Bladder, 1-Tube, Latex, Black, Adult
850BD	ADCUFF™ & Bladder, 1-Tube, Latex, Burgundy, Adult
850BDLF	ADCUFF™ & Bladder, 1-Tube, Latex Free, Burgundy, Adult
850C	ADCUFF™ & Bladder, 1-Tube, Latex, Black, Child
850CLF	ADCUFF™ & Bladder, 1-Tube, Latex Free, Black, Child
850DB	ADCUFF™ & Bladder, 1-Tube, Latex, Navy, Adult
850DBC	ADCUFF™ & Bladder, 1-Tube, Latex, Navy, Child
850DBCLF	ADCUFF™ & Bladder, 1-Tube, Latex Free, Navy, Child
850DBI	ADCUFF™ & Bladder, 1-Tube, Latex, Navy, Infant
850DBILF	ADCUFF™ & Bladder, 1-Tube, Latex Free, Navy, Infant
850DBLF	ADCUFF™ & Bladder, 1-Tube, Latex Free, Navy, Adult
850DBN	ADCUFF™ & Bladder, 1-Tube, Latex, Navy, Newborn
850DBT	ADCUFF™ & Bladder, 1-Tube, Latex, Navy, Thigh
850DBTLF	ADCUFF™ & Bladder, 1-Tube, Latex, Navy, Thigh, Latex-Free
850DBX	ADCUFF™ & Bladder, 1-Tube, Latex, Navy, Large Adult
850DBXLF	ADCUFF™ & Bladder, 1-Tube, Latex Free, Navy, Large Adult
850DG	ADCUFF™ & Bladder, 1-Tube, Latex, Dark Green, Adult
850DGLF	ADCUFF™ & Bladder, 1-Tube, Latex Free, Dark Green, Adult
850I	ADCUFF™ & Bladder, 1-Tube, Latex, Black, Infant
850ILF	ADCUFF™ & Bladder, 1-Tube, Latex Free, Black, Infant
850LF	ADCUFF™ & Bladder, 1-Tube, Latex Free, Black, Adult
850M	ADCUFF™ & Bladder, 1-Tube, Latex, Magenta, Adult
850MLF	ADCUFF™ & Bladder, 1-Tube, Latex Free, Magenta, Adult
850N	ADCUFF™ & Bladder, 1-Tube, Latex, Black, Newborn
850RB	ADCUFF™ & Bladder, 1-Tube, Latex, Royal Blue, Adult
850RBLF	ADCUFF™ & Bladder, 1-Tube, Latex Free, Royal Blue, Adult
850T	ADCUFF™ & Bladder, 1-Tube, Latex, Black, Thigh
850TL	ADCUFF™ & Bladder, 1-Tube, Latex, Teal, Adult
850TLF	ADCUFF™ & Bladder, 1-Tube, Latex Free, Black, Thigh
850TLF	ADCUFF™ & Bladder, 1-Tube, Latex Free, Teal, Adult
850V	ADCUFF™ & Bladder, 1-Tube, Latex, Purple, Adult
850VLF	ADCUFF™ & Bladder, 1-Tube, Latex Free, Purple, Adult
850X	ADCUFF™ & Bladder, 1-Tube, Latex, Purple, Large Adult
850XLF	ADCUFF™ & Bladder, 1-Tube, Latex Free, Purple, Large Adult
851	ADCUFF™ & Bladder, 1-Tube, Latex, Gray, Adult
851LF	ADCUFF™ & Bladder, 1-Tube, Latex Free, Gray, Adult
854	ADCUFF™ & Bladder, 1-Tube, Latex, Green, Adult
854LF	ADCUFF™ & Bladder, 1-Tube, Latex Free, Green, Adult
855	ADCUFF™ & Bladder, 1-Tube, Latex, Red, Adult
855LF	ADCUFF™ & Bladder, 1-Tube, Latex Free, Red, Adult
856	ADCUFF™ & Bladder, 1-Tube, Latex, Orange, Adult
856C	ADCUFF™ & Bladder, 1-Tube, Latex, Orange, Child
856CLF	ADCUFF™ & Bladder, 1-Tube, Latex Free, Orange, Child
856I	ADCUFF™ & Bladder, 1-Tube, Latex, Orange, Infant
856ILF	ADCUFF™ & Bladder, 1-Tube, Latex Free, Orange, Infant
856LF	ADCUFF™ & Bladder, 1-Tube, Latex Free, Orange, Adult
856T	ADCUFF™ & Bladder, 1-Tube, Latex, Orange, Thigh
856TLF	ADCUFF™ & Bladder, 1-Tube, Latex Free, Orange, Thigh
856X	ADCUFF™ & Bladder, 1-Tube, Latex, Orange, Large Adult
856XLF	ADCUFF™ & Bladder, 1-Tube, Latex Free, Orange, Large Adult



ADC® ADCUFF™ & BLADDERS, 2-TUBE

Durable nylon construction. Exclusive Size Guide™ marking system. Latex or latex free inflation bladder with two 55cm tubes. For use with most aneroid or mercurial blood pressure instruments. Individually boxed. Two-year warranty.

- 840 ADCUFF™ & Bladder, 2-Tube, Latex, Black, Adult
- 840BD ADCUFF™ & Bladder, 2-Tube, Latex, Burgundy, Adult
- 840BDLF ADCUFF™ & Bladder, 2-Tube, Latex Free, Burgundy, Adult
- 840C ADCUFF™ & Bladder, 2-Tube, Latex, Black, Child
- 840CH ADCUFF™ & Bladder, 2-Tube, Latex, Black, Child, 10/pkg
- 840CHLF ADCUFF™ & Bladder, 2-Tube, Latex Free, Black, Child, 10/pkg
- 840CLF ADCUFF™ & Bladder, 2-Tube, Latex Free, Black, Child
- 840DB ADCUFF™ & Bladder, 2-Tube, Latex, Navy, Adult
- 840DBC ADCUFF™ & Bladder, 2-Tube, Latex, Navy, Child
- 840DBCLF ADCUFF™ & Bladder, 2-Tube, Latex Free, Navy, Child
- 840DBI ADCUFF™ & Bladder, 2-Tube, Latex, Navy, Infant
- 840DBILF ADCUFF™ & Bladder, 2-Tube, Latex Free, Navy, Infant
- 840DBLF ADCUFF™ & Bladder, 2-Tube, Latex Free, Navy, Adult
- 840DBN ADCUFF™ & Bladder, 2-Tube, Latex, Navy, Newborn
- 840DBT ADCUFF™ & Bladder, 2-Tube, Latex, Navy, Thigh
- 840DBTLF ADCUFF™ & Bladder, 2-Tube, Latex Free, Navy, Thigh
- 840DBX ADCUFF™ & Bladder, 2-Tube, Latex, Navy, Large Adult
- 840DBXLF ADCUFF™ & Bladder, 2-Tube, Latex Free, Navy, Large Adult
- 840DG ADCUFF™ & Bladder, 2-Tube, Latex, Dark Green, Adult
- 840DGLF ADCUFF™ & Bladder, 2-Tube, Latex Free, Dark Green, Adult
- 840H ADCUFF™ & Bladder, 2-Tube, Latex, Black, Adult, 10/pkg
- 840HLF ADCUFF™ & Bladder, 2-Tube, Latex Free, Black, Adult, 10/pkg
- 840I ADCUFF™ & Bladder, 2-Tube, Latex, Black, Infant
- 840IH ADCUFF™ & Bladder, 2-Tube, Latex, Black, Infant, 5/pkg
- 840IHLF ADCUFF™ & Bladder, 2-Tube, Latex Free, Black, Infant, 5/pkg
- 840ILF ADCUFF™ & Bladder, 2-Tube, Latex Free, Black, Infant
- 840LF ADCUFF™ & Bladder, 2-Tube, Latex Free, Black, Adult
- 840M ADCUFF™ & Bladder, 2-Tube, Latex, Magenta, Adult
- 840MLF ADCUFF™ & Bladder, 2-Tube, Latex Free, Magenta, Adult
- 840N ADCUFF™ & Bladder, 2-Tube, Latex, Black, Newborn
- 840NH ADCUFF™ & Bladder, 2-Tube, Latex, Black, Newborn, 5/pkg
- 840RB ADCUFF™ & Bladder, 2-Tube, Latex, Royal Blue, Adult
- 840RBLF ADCUFF™ & Bladder, 2-Tube, Latex Free, Royal Blue, Adult
- 840T ADCUFF™ & Bladder, 2-Tube, Latex, Black, Thigh
- 840TH ADCUFF™ & Bladder, 2-Tube, Latex, Turquoise, Adult, 10/pkg
- 840THLF ADCUFF™ & Bladder, 2-Tube, Latex Free, Black, Thigh, 5/pkg
- 840TL ADCUFF™ & Bladder, 2-Tube, Latex, Teal, Adult
- 840TLF ADCUFF™ & Bladder, 2-Tube, Latex Free, Black, Thigh
- 840TLLF ADCUFF™ & Bladder, 2-Tube, Latex Free, Teal, Adult
- 840V ADCUFF™ & Bladder, 2-Tube, Latex, Purple, Adult
- 840VLF ADCUFF™ & Bladder, 2-Tube, Latex Free, Purple, Adult
- 840X ADCUFF™ & Bladder, 2-Tube, Latex, Black, Large Adult
- 840XH ADCUFF™ & Bladder, 2-Tube, Latex, Black, Large Adult, 10/pkg
- 840XHLF ADCUFF™ & Bladder, 2-Tube, Latex Free, Black, Large Adult, 10/pkg
- 840XLF ADCUFF™ & Bladder, 2-Tube, Latex Free, Black, Large Adult
- 841 ADCUFF™ & Bladder, 2-Tube, Latex, Gray, Adult
- 841LF ADCUFF™ & Bladder, 2-Tube, Latex Free, Gray, Adult
- 844 ADCUFF™ & Bladder, 2-Tube, Latex, Green, Adult
- 844LF ADCUFF™ & Bladder, 2-Tube, Latex Free, Green, Adult
- 845 ADCUFF™ & Bladder, 2-Tube, Latex, Red, Adult
- 845LF ADCUFF™ & Bladder, 2-Tube, Latex Free, Red, Adult
- 846 ADCUFF™ & Bladder, 2-Tube, Latex, Orange, Adult
- 846C ADCUFF™ & Bladder, 2-Tube, Latex, Orange, Child

ADC® ADCUFF™ & BLADDERS, 2-TUBE (Continued)


- 846CLF ADCUFF™ & Bladder, 2-Tube, Latex Free, Orange, Child
- 846I ADCUFF™ & Bladder, 2-Tube, Latex, Orange, Infant
- 846ILF ADCUFF™ & Bladder, 2-Tube, Latex Free, Orange, Infant
- 846LF ADCUFF™ & Bladder, 2-Tube, Latex Free, Orange, Adult
- 846T ADCUFF™ & Bladder, 2-Tube, Latex, Orange, Thigh
- 846TLF ADCUFF™ & Bladder, 2-Tube, Latex Free, Orange, Thigh
- 846X ADCUFF™ & Bladder, 2-Tube, Latex, Orange, Large Adult
- 846XLF ADCUFF™ & Bladder, 2-Tube, Latex Free, Orange, Large Adult



875PLF-10



8971-10



8972-10

ADC® ADCUFF™ BLOOD PRESSURE CUFF ACCESSORIES

- 611-20-10 Y Adaptor, No Fittings For Use w/two tube cuffs, 10/pkg
- 875PLF-10 Bulb & Valve, S.P.U., 10/pkg
- 891 Luer Connector Set, Metal, Includes 891F & 891M
- 891F Luer Connector, Female, Aluminum
- 891F-10 Female Luer Slip Connector w/barbed end for use w/ADC & Baum (metal), 10/pkg
- 891M Luer Connector, Male Metal
- 891M-10 Male Luer Slip Connector w/barbed end for use w/ADC & Baum (metal), 10/pkg
- 8971-10 Screw Adaptor w/barbed end for use w/Dynamap, 10/pkg
- 8972-10 Male luer slip connector w/barbed end (plastic), 10/pkg
- 8973-10 Female Locking Luer Connector w/barbed end (plastic), 10/pkg
- 8974-10 Bayonet Connector For Use w/Hewlett Packard & Siemens Monitors, 10/pkg
- 8975-10 Neonatal Cuff Connector, Male Luer Slip Fitting For 3/32" ID Tubing, 10/pkg
- 920-011 Double Male Connector
- 920-011-10 Double Male Connector For 5/32" ID Tubing (plastic), 10/pkg
- 920-039 Filter/Diaphragm/Gasket Set



860



860BD

ADC® ADCUFF™ INFLATION SYSTEMS

Includes ADCUFF™ nylon cuff featuring the Size Guide™ marking system. Latex or latex free inflation bladder with two 55cm tubes. ADFLOW™ bulb and valve with filter screen protection. For use with most aneroid or mercurial blood pressure instruments. Individually boxed. Two-year warranty.

- 860 Inflation System, Latex, Black, Adult
- 860BD Inflation System, Latex, Burgundy, Adult
- 860BDF Inflation System, Latex Free, Burgundy, Adult
- 860C Inflation System, Latex, Black, Child
- 860CLF Inflation System, Latex Free, Black, Child
- 860DB Inflation System, Latex, Navy, Adult
- 860DBC Inflation System, Latex, Navy, Child
- 860DBCLF Inflation System, Latex Free, Navy, Child
- 860DBI Inflation System, Latex, Navy, Infant
- 860DBILF Inflation System, Latex Free, Navy, Infant
- 860DBLF Inflation System, Latex Free, Navy, Adult
- 860DBN Inflation System, Latex, Navy, Newborn
- 860DBT Inflation System, Latex, Navy, Thigh
- 860DBTLF Inflation System, Latex Free, Navy, Thigh
- 860DBX Inflation System, Latex, Navy, Large Adult
- 860DBXLF Inflation System, Latex Free, Navy, Large Adult
- 860DG Inflation System, Latex, Dark Green, Adult
- 860DGLF Inflation System, Latex Free, Dark Green, Adult
- 860I Inflation System, Latex, Black, Infant
- 860ILF Inflation System, Latex Free, Black, Infant
- 860LF Inflation System, Latex Free, Black, Adult
- 860M Inflation System, Latex, Magenta, Adult
- 860MLF Inflation System, Latex Free, Magenta, Adult
- 860N Inflation System, Latex, Black, Newborn
- 860RB Inflation System, Latex, Royal Blue, Adult
- 860RBLF Inflation System, Latex Free, Royal Blue, Adult
- 860T Inflation System, Latex, Black, Thigh
- 860TL Inflation System, Latex, Teal, Adult
- 860TLF Inflation System, Latex Free, Black, Thigh
- 860TLLF Inflation System, Latex Free, Teal, Adult
- 860V Inflation System, Latex, Purple, Adult
- 860VLF Inflation System, Latex Free, Purple, Adult
- 860X Inflation System, Latex, Black, Large Adult
- 860XLF Inflation System, Latex Free, Black, Large Adult
- 861 Inflation System, Latex, Gray, Adult
- 861LF Inflation System, Latex Free, Gray, Adult
- 864 Inflation System, Latex, Green, Adult
- 864LF Inflation System, Latex Free, Green, Adult
- 865 Inflation System, Latex, Red, Adult
- 865LF Inflation System, Latex Free, Red, Adult
- 866 Inflation System, Latex, Orange, Adult
- 866C Inflation System, Latex, Orange, Child
- 866CLF Inflation System, Latex Free, Orange, Child
- 866I Inflation System, Latex, Orange, Infant
- 866ILF Inflation System, Latex Free, Orange, Infant
- 866LF Inflation System, Latex Free, Orange, Adult
- 866T Inflation System, Latex, Orange, Thigh
- 866TLF Inflation System, Latex Free, Orange, Thigh
- 866X Inflation System, Latex, Orange, Large Adult
- 866XLF Inflation System, Latex Free, Orange, Large Adult

ADC® ADCUFF™ NYLON CUFF

Featuring the Size Guide™ marking system. For use with most aneroid or mercurial blood pressure instruments. Two year warranty.

- 830 Nylon Cuff, Black, Adult
- 830BD Nylon Cuff, Burgundy, Adult
- 830C Nylon Cuff, Burgundy, Child
- 830CH Nylon Cuff, Black, Child, 10/pkg
- 830DB Nylon Cuff, Navy, Adult
- 830DBC Nylon Cuff, Navy, Child
- 830DBI Nylon Cuff, Navy, Infant
- 830DBN Nylon Cuff, Navy, Newborn
- 830DBT Nylon Cuff, Navy, Thigh
- 830DBX Nylon Cuff, Navy, Large Adult
- 830DG Nylon Cuff, Dark Green, Adult
- 830H Nylon Cuff, Black, Adult, 10/pkg
- 830I Nylon Cuff, Black, Infant
- 830IH Nylon Cuff, Black, Infant, 5/pkg
- 830M Nylon Cuff, Magenta, Adult
- 830N Nylon Cuff, Black, Newborn
- 830NH Nylon Cuff, Black, Newborn, 5/pkg
- 830RB Nylon Cuff, Royal Blue, Adult
- 830T Nylon Cuff, Black, Thigh
- 830TH Nylon Cuff, Black, Thigh, 5/pkg
- 830TL Nylon Cuff, Teal, Adult
- 830V Nylon Cuff, Purple, Adult
- 830X Nylon Cuff, Black, Large Adult
- 830XH Nylon Cuff, Black, Large Adult, 10/pkg
- 831 Nylon Cuff, Gray, Adult
- 834 Nylon Cuff, Green, Adult
- 835 Nylon Cuff, Red, Adult
- 836 Nylon Cuff, Orange, Adult
- 836C Nylon Cuff, Orange, Child
- 836I Nylon Cuff, Orange, Infant
- 836T Nylon Cuff, Orange, Thigh
- 836X Nylon Cuff, Orange, Large Adult

ADC® INFLATION BULBS & VALVES

Dip molded latex or latex free. For use with most aneroid or mercurial blood pressure instruments. End valve included. Individually boxed. Two-year warranty.

- 870 ADFLOW™ Bulb w/filter screen end valve, Latex
- 870-1 Filter Screen End Valve Only
- 870H ADFLOW™ Bulb w/filter screen end valve, Latex, 10/pkg
- 870HLF ADFLOW™ Bulb w/filter screen end valve, Latex Free, 10/pkg
- 870L Large ADFLOW™ Bulb w/filter screen end valve, Latex
- 870LF ADFLOW™ Bulb w/filter screen end valve, Latex Free
- 870LH Large ADFLOW™ Bulb w/filter screen end valve, Latex, 10/pkg
- 872 ADFLOW™ Bulb & Valve (870 & 871), Latex
- 872H ADFLOW™ Bulb & Valve (870 & 871), Latex, 10/pkg
- 872HLF ADFLOW™ Bulb & Valve (870 & 871), Latex Free, 10/pkg
- 872L ADFLOW™ Bulb & Valve (870 & 871), Large, Latex
- 872LF ADFLOW™ Bulb & Valve (870 & 871), Latex Free
- 872LH ADFLOW™ Bulb & Valve (870 & 871), Large, Latex, 10/pkg
- 873 Standard Bulb w/regular end valve, Latex
- 873-1 Standard End Valve
- 873LF Standard Bulb w/regular end valve, Latex Free
- 875 Standard Bulb & Valve, Latex
- 875LF Standard Bulb & Valve, Latex Free

OMRON CUFF AND BLADDER REPLACEMENTS

21-071	Calibrated Nylon Cuff, Blue, Adult
21-063	Cotton Cuff, Adult, Gray
21-063C	Cotton Cuff, Child, Gray
21-063X	Cotton Cuff, Large Adult, Gray
21-071C	Calibrated Nylon Cuff, Blue, Child
21-071I	Calibrated Nylon Cuff, Blue, Infant
21-071L	Calibrated Nylon Cuff, Blue, Thigh
21-071X	Calibrated Nylon Cuff, Blue, Large Adult
21-073	Calibrated Nylon Cuff, Gray, Adult
21-073C	Calibrated Nylon Cuff, Gray, Child
21-073I	Calibrated Nylon Cuff, Gray, Infant
21-073L	Calibrated Nylon Cuff, Gray, Thigh
21-073X	Calibrated Nylon Cuff, Gray, Large Adult
21-081	Double Tube Bladder, Adult
21-081NL	Double Tube Bladder, Adult, Latex Free (Neoprene)
21-082	Double Tube Bladder, Infant
21-082NL	Double Tube Bladder, Infant, Latex Free (Neoprene)
21-083	Double Tube Bladder, Child
21-083NL	Double Tube Bladder, Child, Latex Free (Neoprene)
21-084	Double Tube Bladder, Large Adult
21-084NL	Double Tube Bladder, Large Adult, Latex Free (Neoprene)
21-085	Double Tube Bladder, Thigh
21-085NL	Double Tube Bladder, Thigh, Latex Free (Neoprene)
21-091	Single Tube Bladder, Adult
21-092	Single Tube Bladder, Infant
21-093	Single Tube Bladder, Child
21-094	Single Tube Bladder, Large Adult
21-095	Single Tube Bladder, Thigh

OMRON CUFF AND BLADDER SETS

21-031	Calibrated Blue Nylon Cuff w/double-tube bladder, Adult
21-031C	Calibrated Blue Nylon Cuff w/double-tube bladder, Child
21-031CNL	Calibrated Blue Nylon Cuff w/double-tube bladder, Child, Latex Free (Neoprene)
21-031I	Calibrated Blue Nylon Cuff w/double-tube bladder, Infant
21-031INL	Calibrated Blue Nylon Cuff w/double-tube bladder, Infant, Latex Free (Neoprene)
21-031L	Calibrated Blue Nylon Cuff w/double-tube bladder, Thigh
21-031LNL	Calibrated Blue Nylon Cuff w/double-tube bladder, Thigh, Latex Free (Neoprene)
21-031NB	Calibrated Blue Nylon Cuff w/double-tube bladder, Newborn
21-031NL	Calibrated Blue Nylon Cuff w/double-tube bladder, Adult, Latex Free (Neoprene)
21-031X	Calibrated Blue Nylon Cuff w/double-tube bladder, Large Adult
21-031XNL	Calibrated Blue Nylon Cuff w/double-tube bladder, Large Adult, Latex Free (Neoprene)
21-033	Calibrated Gray Nylon Cuff w/double-tube bladder, Adult
21-033C	Calibrated Gray Nylon Cuff w/double-tube bladder, Child
21-033I	Calibrated Gray Nylon Cuff w/double-tube bladder, Infant
21-033X	Calibrated Gray Nylon Cuff w/double-tube bladder, Large Adult
21-037	Calibrated Black Nylon Cuff w/double-tube bladder, Adult
21-037C	Calibrated Black Nylon Cuff w/double-tube bladder, Child
21-037X	Calibrated Black Nylon Cuff w/double-tube bladder, Large Adult
21-053	Calibrated Gray Nylon Cuff w/single-tube bladder, Adult
21-053C	Calibrated Gray Nylon Cuff w/single-tube bladder, Child
21-053I	Calibrated Gray Nylon Cuff w/single-tube bladder, Infant
21-053L	Calibrated Gray Nylon Cuff w/single-tube bladder, Thigh
21-053X	Calibrated Gray Nylon Cuff w/single-tube bladder, Large Adult
21-201	Calibrated Blue Nylon Cuff w/double-tube bladder, inflation bulb, air release valve, Adult
21-201C	Calibrated Blue Nylon Cuff w/double-tube bladder, inflation bulb, air release valve, Child
21-201I	Calibrated Blue Nylon Cuff w/double-tube bladder, inflation bulb, air release valve, Infant
21-201L	Calibrated Blue Nylon Cuff w/double-tube bladder, inflation bulb, air release valve, Thigh
21-201X	Calibrated Blue Nylon Cuff w/double-tube bladder, inflation bulb, air release valve, Large Adult

OMRON INFLATION BULBS

21-122R	Inflation Bulb For 107M Palm Aneroid & R-1077
21-125	Inflation Bulb w/end valve only
21-125NL	Inflation Bulb, Latex Free w/end valve only
21-140	Inflation Bulb w/standard air release valve
21-141	Inflation Bulb w/deluxe air release valve

OMRON MANOMETER PARTS & ACCESSORIES

0179	Calibration Tool For Standard Aneroids
0180	Gauge Guard, Blue
148	Standard Aneroid Gauge, Gray
156	Marshall® Professional Gauge For 108M
21-001	Exactus II Gauge, Blue
21-150	Crystal For Standard Aneroids

OMRON MERCURIAL SPHYGMOMANOMETER PARTS

21-321	8 ft Latex Coiled Tubing
21-321NL	Coiled Tubing 8 ft, Latex Free (Neoprene)
21-347	BMS Wall Basket

TRIMLINE BLADDERLESS SOFT CUFFS, LATEX FREE

Economical, soft fabric blood pressure cuffs are comfortable for patients, even during lengthy procedures. Each size is printed in a different color print for easy identification. Cuffs feature the exclusive Pre-Gaged® BALANCED™ design with an exterior sizing guide and a "Use Larger Cuff" area to help prevent miscuffing, and a longer bladder length to facilitate right or left arm use. Cuffs and tubing are completely latex free. Can be used with all manual and NIBP monitors; connectors/adapters available. Conforms to American Heart Association and AAMI recommendations. Made in the USA.


8602XL	Two Tubes Cuff w/bulb, valve & clamp, Large Adult, Extra Long, Red Print, 5/bx
8900	Two Tubes cuff w/no connectors (Add "S" to # For Single Tube), Thigh, Brown Print, 5/bx



5080-01



5082-01H



5082-43



5086-01



5087-01



5090-65

WELCH ALLYN ANEROID ACCESSORIES & PARTS

5080-01	Adult, Non-Conductive Hook Cuff
5082-01	Adult Non-Conductive Cuff
5082-01H	Adult Cuff, Latex Free, 12/bx
5082-02	Child Non-Conductive Cuff
5082-03	Infant Cuff w/1-tube bag
5082-06	Adult Non-Conductive Cuff, Blue
5082-07	Newborn Cuff w/1-tube bag
5082-08	Infant Cuff w/2-tube bag
5082-11	Thigh Non-Conductive Cuff

WELCH ALLYN ANEROID ACCESSORIES & PARTS

(Continued)

- 5082-16 Large Adult Non-Conductive Cuff
- 5082-16H Large Adult Cuff, Latex Free, 6/bx
- 5082-21 Child Inflation System w/2 tube bag, Latex Free
- 5082-21H Inflation System, Child, Latex Free, 6/bx
- 5082-22 Adult Inflation System w/2 tube bag, Latex Free
- 5082-22H Inflation System, Adult, Latex Free, 6/bx
- 5082-23 Large Adult Inflation System w/2 tube bag, Latex Free
- 5082-23H Inflation System, Large Adult, Latex Free, 6/bx
- 5082-24 Child Cuff w/2-tube bag
- 5082-25 Adult Cuff w/2-tube bag
- 5082-257 Integrated Gauge Adapter Only
- 5082-26 Large Adult Cuff w/2-tube bag
- 5082-42 Child Cuff w/1-tube bag
- 5082-43 Adult Cuff w/1-tube bag
- 5082-44 Large Adult Cuff w/1-tube bag
- 5082-45 Thigh Inflation System
- 5085-07 Sphyg Carrying Case, Polyester
- 5085-08 Two-Sided Sphyg Carrying Case, Polyester
- 5086-01 Latex Inflation Bulb
- 5086-01H Latex Inflation Bulb, Latex Free, 12/bx
- 5086-03H Large Inflation Bulb, Latex Free, 6/bx
- 5087-01 Air Release Valve
- 5087-01H Air Release Valve, Latex Free, 12/bx
- 5087-03 Valve Thumb Screw For Hand Model (pre-1979)
- 5087-04 Valve Assembly For Hand Model (pre-1979)
- 5087-08 Valve Assembly For Hand Model (after Jan. 1, 1979)
- 5087-09 Valve Thumb Screw For Hand Model (after Jan. 1, 1979)
- 5088-01 Air Release Valve & Bulb
- 5088-01H Bulb & Valve Combination, Latex Free, 12/bx
- 5089-01 Adult Inflation Bags, Non-Conductive, 2 Tubes
- 5089-01H Inflation Bag, Adult, 2 Tubes, Latex Free, 12/bx
- 5089-02 Child Inflation Bags, Non-Conductive, 2 Tubes
- 5089-02H Inflation Bag, Child, 2 Tubes, Latex Free, 12/bx
- 5089-03 Adult Inflation Bags, Non-Conductive, 1 Tube
- 5089-04 Child Inflation Bags, Non-Conductive, 1 Tube
- 5089-09 Metal Y Adaptor
- 5089-11 Metal Twist Fit Connector Set
- 5089-12 Tubing 27" Length
- 5089-13H Coil Tubing 8 ft Length, Latex Free, 12/bx
- 5089-18 Thigh Inflation Bags, Non-Conductive, 1 Tube
- 5089-19 Thigh Inflation Bags, Non-Conductive, 2 Tubes
- 5089-19H Inflation Bag, Thigh, 2 Tubes, Latex Free, 6/bx
- 5089-21 Large Adult Inflation Bags, Non-Conductive, 1 Tube
- 5089-22 Large Adult Inflation Bags, Non-Conductive, 2 Tubes
- 5089-22H Inflation Bag, Large Adult, 2 Tubes, Latex Free, 6/bx
- 5089-25 Infant Inflation Bags, Non-Conductive, 1 Tube
- 5089-26 Infant Inflation Bags, Non-Conductive, 2 Tubes
- 5090-29 Bulb End Valve
- 5090-30 Screen Assembly
- 5090-44 Front Ring Bezel For Aneroid
- 5090-45 Shatterproof Crystal For Aneroid
- 5090-64 Clip Assembly For Pocket Aneroid
- 5090-65 Sphyg Repair Kit
- 5091-16 Crystal For Wall Gauge, Square
- 5091-33 Cuff Clip
- 5097-15 Filter & Gasket Set For Mercurial Gauge
- 5097-16 Tube Cap, Wall Mercurial
- 5097-25 Incline Tube, Mobile Mercurial
- 5097-32 Vertical Tube, Wall Mercurial, Marked 4.2 mm
- 5097-33 Incline Tube, Mobile Mercurial, Marked 4.2 mm
- 5098-07 Luer Lock (set)
- 5098-17 Housing Screws For Hand Gauge (Packet of 12)
- 7052-25 Printer Paper For 7052-23 Only, 5 rolls/pack
- DS44A Model DS44 Gauge w/integrated gauge adapter only
- DS45A Model DS45 Gauge w/integrated gauge adapter only



WELCH ALLYN TYCOS® DS66 HAND ANEROID ACCESSORIES

- 5082-18 Child Cuff, Print, Latex Free
- 5082-56 Adult Cuff & Black Inflation Bag, Latex Free
- 5082-58 Large Adult Cuff & Black Inflation Bag, Latex Free
- 5082-59 Adult Cuff, Black, Latex Free
- 5082-61 Large Adult Cuff, Black, Latex Free
- 5082-62 Child Print Cuff & Black Inflation Bag, Latex Free
- 5086-04 Bulb Assembly, Latex Free
- 5098-31 Gauge, Latex Free
- 5098-40 Crystal, Latex Free

BCI MINI-TORR PLUS® NON-INVASIVE BLOOD PRESSURE MONITOR

Provides a base to build on for your monitoring needs. Offers cost-effective, non-invasive blood pressure monitoring. Options which can be added include pulse oximetry, a built-in printer, oral or rectal temperature, and/or a lightweight roll stand or IV pole mount. The ability to combine the options needed in your clinical environment at an affordable price, makes this monitor an important addition to your equipment needs.

- 6004-000 Mini-Torr Plus® NIBP Monitor
- 6004-001 Mini-Torr Plus® NIBP/Oximeter
- 6004-004 Mini-Torr Plus® NIBP/Electronic Thermometry
- 6004-005 Mini-Torr Plus® NIBP/Printer
- 6004-006 Mini-Torr Plus® NIBP/Oximeter/Printer
- 6004-009 Mini-Torr Plus® NIBP/Electronic Thermometry/Printer
- 6004-011 Mini-Torr Plus® NIBP/SpO2/Electronic Thermometry
- 6004-013 Mini-Torr Plus® NIBP/SpO2/Thermometry/Printer



143265



143268



90121-1

BURDICK AMBULATORY BLOOD PRESSURE SYSTEMS

Provides unmatched accuracy for identifying and managing your hypertensive patients.

- 143265 Ultralite ABP Monitor w/software
- 143266 Ultralite ABP Monitor w/report generator
- 143268 ABP Monitor w/report generator
- 90121-1 ABP Report Management System Base Station Interface (for Windows), Features Remote Communication Capability For Data Retrieval & Monitor Programming via Modem
- 90239A ABP Report Generator/Printer
- 143267 ABP Monitor w/software



BURDICK AMBULATORY BLOOD PRESSURE SYSTEM ACCESSORIES

- 000-0027-02 Patient Diary For ABP
- 006-0210-02 Thermal Paper for 90239A ABP Report Generator, 57mm x 32, 6 rl/bx
- 012-0034-03 90239-80 Interface Cable (90239-80 to PC Serial Port)
- 012-0096-01 Modem Interface Cable (ABP Monitor to Modem)
- 012-0097-02 ABP Interface Cable (ABP Monitor to PC)
- 012-0105-02 Base Station Interface Cable (PC to Modem)
- 015-0070-00 Cuff Support Harness For All Cuffs
- 016-0077-01 Large Adult Cuff, 32-42cm
- 016-0080-00 Waist Belt for ABP
- 016-0109-01 XL Adult Cuff, 38-50cm (includes cuff support harness)
- 016-0262-00 Strap For Use w/vinyl pouch for model 90207 ABP
- 016-0340-00 Cloth Pouch for Model 90217 ABP
- 016-0410-00 Cloth Pouch for Model 90207 ABP (SL logo)
- 040-0549-00 Data Management Option Kit, Permits Storage of up to 40 Reports on the Report Generation Data Card
- 146-0008-00 Lithium Battery, 3-Volt, For Use w/90207 ABP
- 672-0167-00 Report Generation Data Card



90207-1Q



90217-1Q

BURDICK AMBULATORY BLOOD PRESSURE SYSTEMS WITH QUICK CONNECT

Our ABP monitors are among the smallest and lightest weight on the market today, with advanced reporting software that allows rapid analysis of patient information. Extensive testing by independent clinicians has validated ABPM accuracy time and time again, making them widely recognized and the "gold standard" for ABP. Both the 90217 Ultralite and 90207 ABP monitors feature a compact, lightweight design that optimizes patient comfort—a key factor in successful ABP monitoring. Comfort helps ensure patient compliance and, consequently, may decrease the number of artifactual readings. Collects measurements for a period of 24 hours or up to 8 days. Lightweight recorder operates quietly for patient comfort. Oscillometry measures systolic, diastolic and mean blood pressure. No chest electrodes or microphones reduce operator costs. Uses a real-time clock to facilitate diary notations.

- 90207-1Q ABP Monitor w/quick Connect, English
- 90217-1Q ABP Ultralite Monitor w/quick Connect, English

BURDICK CUFFS WITH QUICK CONNECT

- 015-0067-01Q Small Adult 17-26cm
- 015-0068-02Q Standard Adult 24-32cm
- 015-0118-01Q Child 12-20cm



WELCH ALLYN AMBULATORY BLOOD PRESSURE MONITOR ACCESSORIES

- 5100-30 Microphone/Pneumatic Tubing Assembly
- 5100-31 Large Cuff & Pneumatic System Assembly
- 5100-32 Medium Cuff & Pneumatic System Assembly
- 5100-34 Large Cuff & Bladder
- 5100-35 Medium Cuff & Bladder
- 5100-36 Small Cuff & Bladder
- 5100-37 Large Cuff
- 5100-38 Medium Cuff
- 5100-41E Patient Hook-Up Kit, 25/bx
- 5100-42E Patient Diary, 50/bx
- 5100-43 Microphone Anchoring Pads, 50/bx
- 5100-44 Cuff Anchoring Pads, 50/bx
- 5100-45 Cuff Anchoring Straps, 5/bx
- 5100-46 Carrying Pouch
- 5100-47 Waist Belt
- 5100-48 Shoulder Strap
- 5100-49 Calibration T-Connector
- 5100-51 Small Bladder
- 5100-52 Medium Bladder
- 5100-53 Large Bladder
- 5100-61 Q-Trak PC Software (includes PC cable)

WELCH ALLYN ATLAS MONITOR

- 622NP-E1 Monitor with Printer



4200B-E1



4200B-E1-M

WELCH ALLYN SPOT VITAL SIGNS

Standard non-invasive blood pressure (NIBP) as well as pulse rate and MAP. Large LCD display with backlight. Extremely user-friendly and easy to use. Uses the oscillometric method to determine systolic and diastolic pressure, MAP and pulse rate. Lightweight and portable: weights just 4 pounds. Training video included.

- 4200B-E1 Device w/BP, Domestic/English
- 4200B-E1-M Device w/BP & Mobile Stand, Domestic/English
- 420TB-E1 Device w/BP & Temperature, Domestic/English
- 420TB-E1-M Device w/BP & Temperature & Mobile Stand, Domestic/English
- 42N0B-E1 Device w/BP & Pulse Oximeter, Domestic/English
- 42N0B-E1-M Device w/BP & Pulse Oximeter, Mobile Stand, Domestic/English
- 42NTB-E1 Device w/BP, Pulse Oximeter & Temperature, Domestic/English
- 42NTB-E1-M Device w/BP, Pulse Oximeter, Temperature & Mobile Stand, Domestic/English

WELCH ALLYN SPOT VITAL SIGNS ACCESSORIES

- 4200-100 Carrying Case
- 4200-145E Technical Service Manual
- 4200-60 Mobile Stand
- 4200-62 Wall Mount
- 4200-64 IV Pole Mount
- 4200-70 Anti-Theft Kit
- 4200-84 Lead Acid Battery, 6V
- 5200-01 Adult Cuff Assembly, Latex Free (Cuff, Bladder & Connector)
- 5200-02 Large Adult Cuff Assembly, Latex Free (Cuff, Bladder & Connector)
- 5200-03 Child Print Cuff Assembly, Latex Free (Cuff, Bladder & Connector)
- 5200-10 Thigh Cuff Assembly, Latex Free (Cuff, Bladder & Connector)
- D-YSPD Pedicheck Pediatric Spot Check




GE MEDICAL DINAMAP® PROCARE PATIENT MONITORS

More than 3 million times every single day critical patient care decisions are made based on blood pressure determinations using DINAMAP® technology. With the ProCare Series, you can take BP measurements that are fast, accurate and comfortable. Get reliable temperature readings in seconds and get oxygen saturation readings with your choice of leading SpO2 technologies.

- DPC120X-US DINAMAP ProCare 120 Monitor, NIBP, Heart Rate, & Printer; Typical Accessories Included At No Additional Charge; Two (2) Year Parts & Labor Warranty
- DPC220X-US DINAMAP ProCare 220 Monitor, NIBP, Heart Rate, IVAC® TurboTemp®, & Printer; Typical Accessories Included At No Additional Charge; Two (2) Year Parts & Labor Warranty
- DPC320N-US DINAMAP ProCare 320 Monitor w/Nellcor SpO2, NIBP, Heart Rate, Nellcor® Oximax® SpO2 & Printer; Typical Accessories Included At No Additional Charge; Two (2) Year Parts & Labor Warranty
- DPC420N-US DINAMAP ProCare 420 Monitor w/Nellcor SpO2, NIBP, Heart Rate, IVAC® TurboTemp®, Nellcor® Oximax® SpO2 & Printer; Typical Accessories Included At No Additional Charge; Two (2) Year Parts & Labor Warranty

Accessories


- 633178CR 12 Volt Lead Acid Battery For Use w/DPC420N-US



GE MEDICAL DINAMAP TEMPERATURE PROBES & PROBE COVERS

Reusable temperature probes. Accurate temperature reading. Probe covers reduce potential for cross-contamination between patients. Latex free.


- 088012 Oral Temperature Probe, Blue
- 088015 Electronic Probe Cover, 20/bx, 250 bx/cs
- 088016 Rectal Temperature Probe, Red
- 2008774-001 TurboTemp Oral Probe, Blue
- 2008775-001 TurboTemp Rectal Probe, Red



1026BOO




760GARCQ



768-641PEA

ADC® CARTOON CHARACTER PRODUCTS

- 1026BOO Nurse Bag, Betty Boop
- 1026GAR Nurse Bag, Garfield
- 1026HK Nurse Bag, Hello Kitty
- 1026PEA Nurse Bag, Peanuts
- 760BOOCQ Child, Aneroid, Betty Boop
- 760GARCQ Child, Aneroid, Garfield
- 760HKCQ Child, Aneroid, Hello Kitty
- 760PEACQ Child, Aneroid, Peanuts
- 760BOOCQ Adult, Aneroid, Betty Boop
- 760GARQ Adult, Aneroid, Garfield
- 760HKQ Adult, Aneroid, Hello Kitty
- 760PEAQ Child, Aneroid, Peanuts
- 760BOOXQ Large Adult, Aneroid, Betty Boop
- 760GARXQ Large Adult, Aneroid, Garfield
- 760PEAXQ Large Adult, Aneroid, Peanuts
- 768-641BOO Adult, Aneroid, Betty Boop
- 768-641GAR Adult, Aneroid, Garfield
- 768-641HK Adult, Aneroid, Hello Kitty
- 768-641PEA Adult, Aneroid, Peanuts
- 888BOO Nylon Zipper Case, Betty Boop
- 888GAR Nylon Zipper Case, Garfield
- 888HK Nylon Zipper Case, Hello Kitty
- 888PEA Nylon Zipper Case, Peanuts



ADC® ADSCOPE™ 665 DISPOSABLE STETHOSCOPE

Nurse style diaphragm only single patient use scope. Acrylesin construction throughout is super-lightweight. Non-toxic. ETO sterilizable. May be safely incinerated. Weighs about 1.5 oz for comfort, yet is remarkably durable. Ultra sensitive diaphragm for outstanding acoustic performance. Individually packaged in a sealed polybag. Available in yellow and royal blue. Latex free.

- 665RB-10 Single Patient Use Scope, Royal Blue, 10/ctn
- 665RB-100 Single Patient Use Scope, Royal Blue, 100/cs
- 665Y-10 Single Patient Use Scope, Yellow, 10/ctn
- 665Y-100 Single Patient Use Scope, Yellow, 100/cs

Accessories

- 698BK Stethoscope ID Tag
- 697BK Premium Stethoscope ID Tag

ADC PROSCOPE™ 664Y DISPOSABLE STETHOSCOPE

Diaphragm only disposable ABS chestpiece with non-chill rim for patient comfort. 22" PVC Yellow Y-tubing. Chrome-plated brass binaurals. Weighs 2.5 oz, measures 32".

- 664Y Proscope™ 664Y, Disposable Scope



WELCH ALLYN/TYCOS UNISCOPE DISPOSABLE STETHOSCOPE

17460	Adult Stethoscope, Green
17461	Adult Stethoscope, Yellow
17462	Adult Stethoscope, Dk. Blue
17462P	Pediatric Stethoscope, Lt. Blue
17463	Adult Stethoscope, Orange
17464	Adult Stethoscope, Lt. Blue
17465	Adult Stethoscope, Green, Anesthesia



OMRON DISPOSABLE NURSE STETHOSCOPE

Lightweight plastic scopes are ideal for single-patient use in contagious areas. Sound transmission is comparable to a standard nurse scope.

0300BLUE	Blue Nurse Stethoscope
0300RED	Red Nurse Stethoscope
0300YELLOW	Yellow Nurse Stethoscope
13-176H	Yellow Nurse Stethoscope w/aluminum binaurals, 50/cs
HSS-490	Lightweight Nurse Stethoscope



STARLINE SINGLE HEAD NURSE STETHOSCOPES

Diaphragm only aluminum chestpiece. Adjustable chrome-plated binaurals. Flexible 22" PVC Y tubing. White plastic eartips. Overall length 31". Latex free.

104266030	Black Stethoscope
104266031	Gray Stethoscope
104266033	Light Blue Stethoscope
104266035	Royal Blue Stethoscope
104266037	Green Stethoscope
104266039	Red Stethoscope
104266040	Pink Stethoscope
104266042	Gold Stethoscope
104266044	Lavender Stethoscope
104266045	Orchid Stethoscope
104266046	Magenta Stethoscope
104266047	Yellow Stethoscope
104266049	White Stethoscope
104266052	Neon Orange Stethoscope
104266053	Neon Pink Stethoscope
104266054	Neon Yellow Stethoscope
104266056	Burgundy Stethoscope
104266057	Dark Green Stethoscope
104266059	Teal Stethoscope
104266060	Purple Stethoscope
104266061	Navy Stethoscope

STARLINE NURSESCOPE ACCESSORIES

660-01	Chestpiece For 660 & 661
660-05	Binaurals For PROSCOPES
660-05S	Binaural Spring For 660-05 Only
660-05TBD	Y Tubing For PROSCOPES, Burgundy
660-05TBK	Y Tubing For PROSCOPES, Black
660-05TDG	Y Tubing For PROSCOPES, Dark Green
660-05TG	Y Tubing For PROSCOPES, Gray
660-05TGD	Y Tubing For PROSCOPES, Gold
660-05TGR	Y Tubing For PROSCOPES, Green
660-05TLB	Y Tubing For PROSCOPES, Light Blue
660-05TLV	Y Tubing For PROSCOPES, Lavender
660-05TM	Y Tubing For PROSCOPES, Magenta
660-05TN	Y Tubing For PROSCOPES, Navy
660-05TNO	Y Tubing For PROSCOPES, Neon Orange
660-05TNP	Y Tubing For PROSCOPES, Neon Pink
660-05TNY	Y Tubing For PROSCOPES, Neon Yellow
660-05TOD	Y Tubing For PROSCOPES, Orchid
660-05TP	Y Tubing For PROSCOPES, Pink
660-05TR	Y Tubing For PROSCOPES, Red
660-05TRB	Y Tubing For PROSCOPES, Royal Blue
660-05TTL	Y Tubing For PROSCOPES, Teal
660-05TV	Y Tubing For PROSCOPES, Purple
660-05TW	Y Tubing For PROSCOPES, White
660-05TY	Y Tubing For PROSCOPES, Yellow
660-06	Eartips For PROSCOPES, White



OMRON LIGHTWEIGHT SINGLE HEAD NURSE STETHOSCOPES

Match a favorite uniform with one of eleven brightly colored nurse stethoscopes. Silver chestpiece. These colorful stethoscopes are useful for color coding departments or identifying contagious areas.

411BLK	Nurse Stethoscope w/black tubing
411BLUE	Nurse Stethoscope w/blue tubing
411RED	Nurse Stethoscope w/red tubing



WELCH ALLYN SPECTRUM™ SINGLE-HEAD STETHOSCOPES


Excellent for taking blood pressures and general auscultation.

5079-230	38" Adult Black Stethoscope, Anesthesiology w/taper lock
----------	--

MILTEX FORD STETHOSCOPE

Shallow bell and plastic ear tips. Chrome.


1-4	Ford Stethoscope
-----	------------------



OMRON BOWLES STETHOSCOPE

Most economical professional stethoscope features a chrome-plated binaural and black PVC tubing. Standard single stem.

0407 Black Bowles Stethoscope



3M™ LITTMANN® MASTER CLASSIC II STETHOSCOPE

The excellent acoustic sensitivity of a patented tunable diaphragm is the hallmark of the low-profile Littmann Master Classic II. Award-winning design combines high performance with portable convenience in this single-sided stethoscope. Features patented Littmann Snap Tight soft-sealing eartips, patient friendly, nonchill rim and diaphragm. Three-year warranty, plus lifetime service and repair policy.

- 2141 Master Classic II Stethoscope, 27" Black Plated Chestpiece & Eartubes
- 2142G Master Classic II Stethoscope, 27" Gold-Plated Chestpiece & Eartubes
- 2143 Master Classic II Stethoscope, 27" Purple Tubing
- 2144L Master Classic II Stethoscope, 27" Black Tubing
- 2146 Master Classic II Stethoscope, 27" Burgundy Tubing
- 2147 Master Classic II Stethoscope, 27" Navy Blue Tubing
- 2630 Master Classic II Stethoscope, 27" Caribbean Blue Tubing
- 2632 Master Classic II Stethoscope, 27" Hunter Green Tubing
- 2633 Master Classic II Stethoscope, 27", Ceil Blue Tubing, ea



3M™ LITTMANN® SELECT STETHOSCOPES

A convenient tunable diaphragm in a single-sided coated alloy chestpiece providing reliable and clear acoustics.

- 2290 Select Stethoscope, 28" Black Tubing
- 2291 Select Stethoscope, 28" Caribbean Blue
- 2293 Select Stethoscope, 28" Burgundy Tubing
- 2294 Select Stethoscope, 28" Purple Tubing
- 2296 Select Stethoscope, 28" Raspberry Tubing
- 2298 Select Stethoscope, 28" Royal Blue Tubing
- 2302 Select Stethoscope, 28", Seafoam Green Tubing, ea
- 2303 Select Stethoscope, 28", Lilac Tubing, ea



STARLINE DUAL-HEAD STETHOSCOPES

Combination diaphragm/bell chestpiece for wide frequency response. Bell and diaphragm equipped with non-chill ring for patient comfort. Lightweight, durable construction. Adjustable metal binaurals. Flexible 22" PVC tubing. Plastic eartips. Overall length 31 1/2". Latex free.

- 104267030 Black Stethoscope
- 104267031 Gray Stethoscope
- 104267033 Light Blue Stethoscope
- 104267035 Royal Blue Stethoscope
- 104267037 Green Stethoscope
- 104267039 Red Stethoscope
- 104267040 Pink Stethoscope
- 104267042 Gold Stethoscope
- 104267044 Lavender Stethoscope
- 104267045 Orchid Stethoscope
- 104267046 Magenta Stethoscope
- 104267047 Yellow Stethoscope
- 104267049 White Stethoscope
- 104267052 Neon Orange Stethoscope
- 104267053 Neon Pink Stethoscope
- 104267054 Neon Yellow Stethoscope
- 104267056 Burgundy Stethoscope
- 104267057 Dark Green Stethoscope
- 104267059 Teal Stethoscope
- 104267060 Purple Stethoscope
- 104267061 Navy Stethoscope

Accessories

- 670-01 Chestpiece For 670 & 671
- 670-04 Non-Chill Ring For 670 & 671



ADC® ADSOFT™ 603 ADULT STETHOSCOPE

Combination adult size chestpiece precisely machined to exacting tolerances from surgical stainless steel for outstanding performance and rugged durability. Extra large bell (1 3/8" diameter) for unsurpassed low frequency response. Ultra sensitive diaphragm (1 7/8" diameter) for greater amplification and crisper high frequency response. Color coordinate non-chill bell and diaphragm retaining rim for patient comfort. Stainless steel binaurals with double leaf internal spring fixed at 15° angle for greater comfort. Reinforcing yoke molded into flexible 22" PVC "Y" tubing. ADSOFT™ deluxe PVC eartips. Spare diaphragm and scope ID tag included. Weights 6.5 oz. Overall length 31 1/2". Lifetime quality performance guarantee. Latex free.

- 603BD Burgundy Stethoscope
- 603BK Black Stethoscope
- 603DG Dark Green Stethoscope
- 603FG Frosted Glacier Stethoscope
- 603FL Frosted Lilac Stethoscope
- 603FS Frosted Seafoam Stethoscope
- 603G Gray Stethoscope
- 603LV Lavender Stethoscope
- 603N Navy Stethoscope
- 603R Red Stethoscope
- 603RB Royal Blue Stethoscope
- 603TL Teal Stethoscope

ADC® ADSCOPE™ 603 ADULT STETHOSCOPE ACCESSORIES

603-01BD	Burgundy Chestpiece
603-01BK	Black Chestpiece
603-01DG	Dark Green Chestpiece
603-01FG	Frosted Glacier Chestpiece
603-01FL	Frosted Lilac Chestpiece
603-01FS	Frosted Seafoam Chestpiece
603-01G	Gray Chestpiece
603-01LV	Lavender Chestpiece
603-01N	Navy Chestpiece
603-01R	Red Chestpiece
603-01RB	Royal Blue Chestpiece
603-01TL	Teal Chestpiece
603-04BD	Burgundy Non-Chill Ring For 603 & 609
603-04BK	Black Non-Chill Ring For 603 & 609
603-04DG	Dk. Green Non-Chill Ring For 603 & 609
603-04FG	Frosted Glacier Non-Chill Ring For 603 & 609
603-04FL	Frosted Lilac Non-Chill Ring For 603 & 609
603-04FS	Frosted Seafoam Non-Chill Ring For 603 & 609
603-04G	Gray Non-Chill Ring For 603 & 609
603-04LV	Lavender Non-Chill Ring For 603 & 609
603-04N	Navy Non-Chill Ring For 603 & 609
603-04R	Red Non-Chill Ring For 603 & 609
603-04RB	Royal Blue Non-Chill Ring For 603 & 609
603-04TL	Teal Non-Chill Ring For 603 & 609
603-05BD	Burgundy Binaural/Tubing
603-05BK	Black Binaural/Tubing
603-05DG	Dk. Green Binaural/Tubing
603-05FG	Frosted Glacier Binaural/Tubing
603-05FL	Frosted Lilac Binaural/Tubing
603-05FS	Frosted Seafoam Binaural/Tubing
603-05G	Gray Binaural/Tubing
603-05LV	Lavender Binaural/Tubing
603-05N	Navy Binaural/Tubing
603-05R	Red Binaural/Tubing
603-05RB	Royal Blue Binaural/Tubing
603-05TL	Teal Binaural/Tubing



ADC® ADSCOPE™ 609 ADSCOPE-LITE™ STETHOSCOPE

Combination adult chestpiece precisely machined to exacting tolerances from lightweight aluminum is remarkable durable, yet extraordinarily comfortable. Extra large 1 3/8" diameter bell. High sensitivity 1 3/4" diaphragm. Non-chill bell and diaphragm retaining rim for patient comfort. Double leaf internal spring binaurals are fixed at 15° angle for greater comfort. Custom molded 22" PVC tubing. ADSOFT™ PVC eartips and spare diaphragm included. Weights 4.5 oz. Overall length 31 1/2". Five year quality performance guarantee. Latex free.

609BD	Burgundy Stethoscope
609BK	Black Stethoscope
609DG	Dk. Green Stethoscope
609FG	Frosted Glacier Stethoscope
609FL	Frosted Lilac Stethoscope
609FS	Frosted Seafoam Stethoscope
609G	Gray Stethoscope
609LV	Lavender Stethoscope
609N	Navy Stethoscope
609R	Red Stethoscope
609RB	Royal Blue Stethoscope
609TL	Teal Stethoscope

ADC® ADSCOPE™ 609 ADSCOPE-LITE™ STETHOSCOPE ACCESSORIES

609-01BD	Burgundy Chestpiece
609-01BK	Black Chestpiece
609-01DG	Dark Green Chestpiece
609-01FG	Frosted Glacier Chestpiece
609-01FL	Frosted Lilac Chestpiece
609-01FS	Frosted Seafoam Chestpiece
609-01G	Gray Chestpiece
609-01LV	Lavender Chestpiece
609-01N	Navy Chestpiece
609-01R	Red Chestpiece
609-01RB	Royal Blue Chestpiece
609-01TL	Teal Chestpiece
609-03	Diaphragm For All
609-05BD	Burgundy Binaural/Tubing
609-05BK	Black Binaural/Tubing
609-05DG	Dark Green Binaural/Tubing
609-05FG	Frosted Glacier Binaural/Tubing
609-05FL	Frosted Lilac Binaural/Tubing
609-05FS	Frosted Seafoam Binaural/Tubing
609-05G	Gray Binaural/Tubing
609-05LV	Lavender Binaural/Tubing
609-05N	Navy Binaural/Tubing
609-05R	Red Binaural/Tubing
609-05RB	Royal Blue Binaural/Tubing
609-05TL	Teal Binaural/Tubing

ADC® ADSCOPE™ 645 STETHOSCOPE

Extraordinary design for outstanding acoustic performance and rugged durability. Chestpiece drum, retaining rings, tubing clip and binaurals generously plated with 5 microns of 18 karat gold, then polished to a brilliant shine. Complete with an accessory pouch containing two extra pairs of eartips, spare diaphragms and the three bell chestpieces. Weighs 8 oz. Overall length 30". Inspected and packaged in the U.S.A. One year warranty on finish, lifetime on performance. Latex free.

645 ADCSCOPE™ 645 18 Karat Gold Plated Stethoscope

Accessories

645-01	Chestpiece
645-01D	Chestpiece Drum
645-02	Retaining Ring, Adult
645-02S	Retaining Ring, Pediatric
645-05	Binaurals, Adult
645-05S	Binaural Spring, Pediatric
645-14A	Diaphragm Assembly
645-15A	Diaphragm Assembly
645-16	Tubing Clip Assembly



OMRON LIGHTWEIGHT DUAL HEAD STETHOSCOPES

Lightweight stethoscopes feature adjustable binaurals for a comfortable fit, one-piece molded PVC tubing for good sound transmission, and built-in rubber non-chill rings. Includes extra pair of eartips and diaphragm.

412BLK	Dual Head Stethoscope, Silver Chestpiece w/black tubing
412BLUE	Dual Head Stethoscope, Blue Chestpiece w/blue tubing
412GRAY	Dual Head Stethoscope, Silver Chestpiece w/gray tubing
412GREEN	Dual Head Stethoscope, Green Chestpiece w/green tubing
412RED	Dual Head Stethoscope, Red Chestpiece w/red tubing



ADC® ADSCOPE™ 615 PLATINUM PROFESSIONAL EDITION

Combines the acoustic response of a traditional bell and diaphragm into a convenient one-sided chestpiece. New gasketed diaphragm enables the scope to achieve multi-frequency response simply by altering pressure. Light pressure produces bell-type low frequency response. Firm pressure will produce diaphragm-type high frequency performance. New ovoid shaped chestpiece improves ergonomics. Chestpiece made from lightweight zinc alloy plated with an attractive satin finish. Non-chill diaphragm retaining rim for patient comfort. Large bore chrome-plated brass binaurals in satin finish. Double leaf internal spring fixed at 15° angle for greater comfort. Reinforcing yoke molded into flexible 22" PVC "Y" tubing. Adsoft™ deluxe PVC eartips in small and large size for the ultimate in wearing comfort and acoustic seal. Spare diaphragm included. Scope ID Tag included. Weighs 5.4 oz. Overall length 30.5". Lifetime Warranty.

615BD	Multifrequency Scope, Burgundy
615BK	Multifrequency Scope, Black
615DG	Multifrequency Scope, Dark Green
615N	Multifrequency Scope, Navy
615RB	Multifrequency Scope, Royal Blue
615ST	Multifrequency Scope, Ninja
615V	Multifrequency Scope, Purple



WELCH ALLYN ELITE® STETHOSCOPE

Stainless steel chestpiece and binaurals designed for optimal weight and durability. Dual-channel PVC tubing provides separate sound pathways to each ear. Triple-leaf binaural spring encased in polyurethane permits size adjustment while preventing breakage. Ten year warranty.

5079-122	22" Stethoscope, Black
5079-125	25" Stethoscope, Black
5079-125P	22" Stethoscope, Pediatric, Black
5079-270	25" Stethoscope, Burgundy
5079-271	25" Stethoscope, Navy
5079-284	25" Stethoscope, Forest Green



5079-135



5079-149

WELCH ALLYN SPECTRUM™ DOUBLE-HEAD STETHOSCOPES

The full-range Spectrum stethoscope is for professionals who demand the best in auscultation and comfort. Superior acoustics mean accurate auscultation of heart, lung and blood pressure sounds. Binaurals adjust to angle of ear. Stainless steel chestpiece. Portable and durable. Five year warranty.

5079-135	28" Adult Black Stethoscope
5079-137	28" Adult Navy Stethoscope
5079-139	28" Adult Burgundy Stethoscope
5079-145	28" Pediatric Black Stethoscope
5079-147	28" Pediatric Navy Stethoscope
5079-149	28" Pediatric Burgundy Stethoscope
5079-285	28" Adult Forest Green Stethoscope
5079-287	28" Pediatric Forest Green Stethoscope
5079-289	28" Adult Blue Stethoscope
5079-291	28" Pediatric Blue Stethoscope



3M™ LITTMANN® CLASSIC II S.E. STETHOSCOPE

Featuring a two-sided chestpiece, this classic stethoscope offers high acoustic sensitivity for superior performance. A patented tunable diaphragm on one side augments the traditional bell function of the opposite side. Features patented Littmann Snap Tight soft-sealing eartips plus patient friendly nonchill rim and diaphragm. Three-year warranty plus lifetime service and repair policy.

2201	S.E. Stethoscope, 28" Black Tubing
2203	S.E. Stethoscope, 28" Gray Tubing
2205	S.E. Stethoscope, 28" Navy Blue Tubing
2206	S.E. Stethoscope, 28" Caribbean Blue Tubing
2208	S.E. Stethoscope, 28" Hunter Green Tubing
2209	S.E. Stethoscope, 28" Purple Tubing
2210	S.E. Stethoscope, 28" Raspberry
2211	S.E. Stethoscope, 28" Burgundy Tubing
2812	S.E. Stethoscope, 28", Orange Tubing
2814	S.E. Stethoscope, 28", Seafoam Green Tubing
2815	S.E. Stethoscope, 28", Pink Tubing

3M™ LITTMANN® LIGHTWEIGHT II S.E. STETHOSCOPES

A general-purpose stethoscope that offers a traditional bell and diaphragm with an aluminum chestpiece to provide acoustic sensitivity.

2450	Lightweight Stethoscope, 28" Black Tubing
2451	Lightweight Stethoscope, 28" Burgundy Tubing
2452	Lightweight Stethoscope, 28" Caribbean Blue Tubing
2453	Lightweight Stethoscope, 28" Lilac Tubing
2454	Lightweight Stethoscope, 28" Ceil Blue Tubing
2455	Lightweight Stethoscope, 28" Sea Foam Green Tubing

3M™ LITTMANN® VETERINARIAN STETHOSCOPES

1392	Vet Stethoscope, Master Classic II, Single-Sided Chestpiece, 32" Navy Tubing
------	--



DRG PURETONE™ MENTOR II TRADITIONAL STETHOSCOPES

PT14B27	Mentor II Traditional Stethoscope, 27" Black Tube
PT14N27	Mentor II Traditional Stethoscope, 27" Navy Tube



PT3B27



PT4B27

DRG PURETONE™ TRADITIONAL STETHOSCOPES

The original stethoscope with the patented ENR™ (External Noise Reduction) tube technology is truly an innovation of design and function. Significantly reduces the amount of external noise while optimizing internal acoustics. 40% less weight as a result of titanium chestpiece construction; also available in stainless steel. Both with conventional diaphragm functionality. GELseal™ eartips and Anti-Slip™ rings are included. Three year warranty.

PT3B27	Traditional Stethoscope, 27" Black Tube
PT3N27	Traditional Stethoscope, 27" Navy Tube
PT4B27	Traditional Ti-Lite Stethoscope, 27" Black Tube
PT4N27	Traditional Ti-Lite Stethoscope, 27" Navy Tube



STARLINE PEDIATRIC STETHOSCOPES

Pediatric sized, combination diaphragm/bell chestpiece for wide frequency response. 1 1/4" diaphragm/bell equipped with non-chill ring for patient comfort. Lightweight aluminum construction. Adjustable chrome-plated binaurals. Flexible 22" PVC Y tubing. White plastic eartips. Overall length 31 1/2". Latex free.

104267530	Black Stethoscope
104267531	Gray Stethoscope
104267533	Light Blue Stethoscope
104267535	Royal Blue Stethoscope
104267537	Green Stethoscope
104267539	Red Stethoscope
104267540	Pink Stethoscope
104267542	Gold Stethoscope
104267544	Lavender Stethoscope
104267545	Orchid Stethoscope
104267546	Magenta Stethoscope
104267547	Yellow Stethoscope
104267549	White Stethoscope
104267552	Neon Orange Stethoscope
104267553	Neon Pink Stethoscope
104267554	Neon Yellow Stethoscope
104267556	Burgundy Stethoscope
104267557	Dark Green Stethoscope
104267559	Teal Stethoscope
104267560	Purple Stethoscope
104267561	Navy Stethoscope



ADC® PROSCOPE™ PEDIATRIC STETHOSCOPES

Pediatric sized, combination diaphragm/bell chestpiece for wide frequency response. 1 1/4" diaphragm 3/4" bell equipped with non-chill ring for patient comfort. Non-chill rim. Lightweight aluminum construction. Adjustable chrome-plated binaurals. Flexible 22" PVC Y tubing. White plastic eartips. Overall length 30". Latex free.

675BD	Burgundy Stethoscope
675BK	Black Stethoscope
675DG	Dark Green Stethoscope
675G	Gray Stethoscope
675GD	Gold Stethoscope
675GR	Green Stethoscope
675LB	Light Blue Stethoscope
675LV	Lavender Stethoscope
675M	Magenta Stethoscope
675N	Navy Stethoscope
675NO	Neon Orange Stethoscope
675NP	Neon Pink Stethoscope
675NY	Neon Yellow Stethoscope
675OD	Orchid Stethoscope
675P	Pink Stethoscope
675R	Red Stethoscope
675RB	Royal Blue Stethoscope
675TL	Teal Stethoscope
675V	Purple Stethoscope
675W	White Stethoscope
675Y	Yellow Stethoscope

ADC® PROSCOPE™ PEDIATRIC STETHOSCOPES ACCESSORIES

Accessories

- 675-01 Chestpiece For 675
- 675-02 Diaphragm Retaining Ring For 675
- 675-03 Diaphragm For 675
- 675-04 Non-Chill Ring For 675



3M™ LITTMANN® CLASSIC II PEDIATRIC & INFANT STETHOSCOPES

Two-sided models with traditional bell/diaphragm chestpieces are designed, sized and acoustically precise for children and infants. Feature patented Littmann Snap Tight soft-sealing eartips plus patient-friendly nonchill rim. Three-year warranty plus a lifetime service and repair policy.

- 2113 Pediatric Stethoscope, 28" Black Tubing
- 2113R Pediatric Stethoscope, 28" Red Tubing
- 2114 Infant Stethoscope, 28" Black Tubing
- 2114R Infant Stethoscope, 28" Red Tubing
- 2119 Pediatric Stethoscope, 28" Caribbean Blue Tubing
- 2122 Pediatric Stethoscope, 28" Raspberry Tubing
- 2124 Infant Stethoscope, 28" Caribbean Blue Tubing
- 2125 Infant Stethoscope, 28" Raspberry Tubing



ADC® ADSCOPE™ 604 PEDIATRIC STETHOSCOPE

Combination pediatric chestpiece precisely machined to exacting tolerances from surgical stainless steel for outstanding performance and rugged durability. Extra large bell (1 3/8" diameter) for unsurpassed low frequency response. Ultra sensitive diaphragm (1 3/4" diameter) for greater amplification and crisper high frequency transmission. Color coordinated non-chill bell and diaphragm retaining rim for patient comfort. Reinforcing yoke molded into flexible 22" PVC Y tubing. ADSOFT™ deluxe PVC eartips. Spare diaphragm included. Weights 5.5 oz. Overall length 31". Lifetime quality performance guarantee. Latex free.

- 604BK Black Pediatric Stethoscope
- 604G Gray Pediatric Stethoscope
- 604LB Lt. Blue Pediatric Stethoscope
- 604P Pink Pediatric Stethoscope

ADC® ADSCOPE™ 604 PEDIATRIC STETHOSCOPE ACCESSORIES

- 604-01BK Black Chestpiece
- 604-01G Gray Chestpiece
- 604-01LB Lt. Blue Chestpiece
- 604-01P Pink Chestpiece
- 604-02BD Burgundy Diaphragm Rim
- 604-02BK Black Diaphragm Rim
- 604-02DG Dark Green Diaphragm Rim
- 604-02G Gray Diaphragm Rim
- 604-02LB Lt. Blue Diaphragm Rim
- 604-02N Navy Diaphragm Rim
- 604-02P Pink Diaphragm Rim
- 604-03 Diaphragm For All
- 604-04BD Burgundy Non-Chill Ring
- 604-04BK Black Non-Chill Ring
- 604-04DG Dark Green Non-Chill Ring
- 604-04G Gray Non-Chill Ring
- 604-04LB Lt. Blue Non-Chill Ring
- 604-04N Navy Non-Chill Ring
- 604-04P Pink Non-Chill Ring



PT5B27



PT6B27




PT9B27

DRG PURETONE™ PEDIATRIC & INFANT STETHOSCOPES

Designed specifically for smaller patients. Features the superior acoustic response of the breakthrough ENR™ tube technology found on the Classic. Patented ENR™ (External Noise Reduction) tube technology significantly reduces the amount of external noise while optimizing internal acoustics. 40% less weight as a result of titanium chestpiece construction; also available in stainless steel. Both with conventional diaphragm functionality sized appropriately for use on smaller patients. GELseal™ eartips and Anti-Slip™ rings add comfort and style. Three year warranty.

- PT5B27 Pediatric Stethoscope, 27" Black Tube
- PT5N27 Pediatric Stethoscope, 27" Navy Tube
- PT6B27 Infant Stethoscope, 27" Black Tube
- PT6N27 Infant Stethoscope, 27" Navy Tube
- PT7B27 Pediatric Ti-Lite Stethoscope, 27" Black Tube
- PT7N27 Pediatric Ti-Lite Stethoscope, 27" Navy Tube
- PT9B27 Infant Ti-Lite Stethoscope, 27" Black Tube
- PT9N27 Infant Ti-Lite Stethoscope, 27" Navy Tube



STARLINE INFANT STETHOSCOPES

Infant sized, combination diaphragm/bell chestpiece for wide frequency response. 25mm diaphragm/bell equipped with non-chill ring for patient comfort. Chrome-plated brass chestpiece. Adjustable chrome-plated binaurals. Flexible 22" PVC Y tubing. White plastic eartips. Overall length 31". Latex free.

104267630	Black Stethoscope
104267631	Gray Stethoscope
104267633	Light Blue Stethoscope
104267635	Royal Blue Stethoscope
104267637	Green Stethoscope
104267639	Red Stethoscope
104267640	Pink Stethoscope
104267642	Gold Stethoscope
104267644	Lavender Stethoscope
104267645	Orchid Stethoscope
104267646	Magenta Stethoscope
104267647	Yellow Stethoscope
104267649	White Stethoscope
104267652	Neon Orange Stethoscope
104267653	Neon Pink Stethoscope
104267654	Neon Yellow Stethoscope
104267656	Burgundy Stethoscope
104267657	Dark Green Stethoscope
104267659	Teal Stethoscope
104267660	Purple Stethoscope
104267661	Navy Stethoscope



ADC® ADSCOPE™ 605 INFANT STETHOSCOPE

Combination infant chestpiece precisely machined to exacting tolerances from surgical stainless steel for outstanding performance and rugged durability. 7/8" diameter bell provides unsurpassed low frequency response. Ultra sensitive diaphragm (1 1/8" diameter) for greater amplification and crisper high frequency transmission. Non-chill bell and diaphragm retaining rim for patient comfort. Reinforcing yoke molded into flexible 22" PVC Y tubing. ADSOFT™ deluxe PVC eartips. Spare diaphragm and scope ID tag included. Weights 5 oz. Overall length 31". Lifetime quality performance guarantee. Latex free.

605BK	Black Infant Stethoscope
605G	Gray Infant Stethoscope
605LB	Lt. Blue Infant Stethoscope
605P	Pink Infant Stethoscope

Accessories

605-01BK	Black Chestpiece
605-01G	Gray Chestpiece
605-02BK	Black Diaphragm Rim
605-02G	Gray Diaphragm Rim
605-03	Diaphragm For All
605-04BK	Black Non-Chill Ring
605-04G	Gray Non-Chill Ring



ADC® PLATINUM EDITION ADSCOPE-LITE™

Unique design that combines the acoustic response of a traditional bell and diaphragm into a convenient one-sided chestpiece. Bell-type low frequency response is achieved by simply holding the chestpiece in place with light pressure. Firm pressure will produce diaphragm-type high frequency performance. A rubber gasket recessed into the chestpiece alters the diaphragm permitting the scope to achieve multi-frequency response. Chestpiece made from lightweight zinc alloy plated with an attractive satin finish. Non-chill diaphragm retaining rim for patient comfort. Large bore chrome plated brass binaurals with matching satin finish. Double leaf internal spring fixed at 15° angle for greater comfort. Reinforcing yoke molded into flexible 22" PVC Y tubing. ADSOFT™ deluxe PVC eartips in small and large size for the ultimate in wearing comfort and acoustic seal. Spare diaphragm included. Scope ID tag included. Weights 6 oz. Overall length 31". Lifetime warranty.

614BK	Black, Child Adscope-Lite™
614G	Gray, Child Adscope-Lite™
614LB	Light Blue, Child Adscope-Lite™
614P	Pink, Child Adscope-Lite™

MILTEX PINARD FETAL HEART STETHOSCOPE

Aluminum.

30-2200	Fetal Heart Stethoscope
---------	-------------------------

OMRON FETAL STETHOSCOPE

Designed for superior auscultation of the fetal heart. Stethoscope is equipped with chrome-plated brass binaurals and stem. 22" tubing can be easily shortened to facilitate use of the headrest.

13-570	Black Fetal Stethoscope
--------	-------------------------

OMRON NEONATAL STETHOSCOPE

Chrome-plated brass chestpiece with screw-on retaining ring and slip-on non-chill ring. Quality sound transmission. Chrome-plated binaurals and 22" Y-tubing. Overall length, 30".

13-310	Gray Neonatal Stethoscope
13-311B	Light Blue Neonatal Stethoscope
13-315P	Pink Neonatal Stethoscope

MILTEX DE LEE OBSTETRICAL STETHOSCOPES

Chrome. Complete with binaurals.

30-2205	Stethoscope w/padded head support, Rigid Metal Arm
30-2210	Stethoscope w/padded head support, Flexible Arm



STARLINE SPRAGUE RAPPAPORT STETHOSCOPES

Threaded chestpiece drum of chrome plated zinc permits the use of five interchangeable chestpiece fittings. 22" tubing length; overall length 30". Latex free, double tube configuration for improved acoustic separation. Proprietary valve mechanism eliminates acoustic leakage. Adjustable chrome plated brass binaurals. Complete with an accessory pouch containing the bell chestpieces, chestpiece key, 2 spare diaphragms, one pair of mushroom eartips and one pair of silicon eartips. Weighs 7.5 oz. Inspected and packaged in the U.S.A. Five year warranty.

104264130	Black Stethoscope
104264131	Gray Stethoscope
104264133	Light Blue Stethoscope
104264135	Royal Blue Stethoscope
104264137	Green Stethoscope
104264139	Red Stethoscope
104264140	Pink Stethoscope
104264142	Gold Stethoscope
104264144	Lavender Stethoscope
104264145	Orchid Stethoscope
104264146	Magenta Stethoscope
104264150	Neon Blue Stethoscope
104264151	Neon Green Stethoscope
104264152	Neon Orange Stethoscope
104264153	Neon Pink Stethoscope
104264154	Neon Yellow Stethoscope
104264155	Boysenberry Stethoscope
104264156	Burgundy Stethoscope
104264157	Dark Green Stethoscope
104264159	Teal Stethoscope
104264160	Purple Stethoscope
104264161	Navy Stethoscope

Everything you
need for your
state-of-the-art
facility.



ADC® ADSCOPE™ SPRAGUESCOPES

Threaded chestpiece drum of chrome plated zinc permits the use of five interchangeable chestpiece fittings. 22" tubing length; overall length 30" (except #640 has 16" tubing length with overall length 24"). Latex free, double tube configuration for improved acoustic separation. Proprietary valve mechanism eliminates acoustic leakage. Adjustable chrome plated brass binaurals. Complete with an accessory pouch containing the bell chestpieces, chestpiece key, 2 spare diaphragms, one pair of mushroom eartips and one pair of silicon eartips. Weighs 6 oz. Inspected and packaged in the U.S.A. Five year warranty. Latex free.

640BK	16", Black Spraguescope
641BB	Boysenberry Spraguescope
641BD	Burgundy Spraguescope
641BK	Black Spraguescope
641CA	Caribbean Blue Spraguescope
641FG	Frosted Glacier Spraguescope
641FL	Frosted Lilac Spraguescope
641FS	Frosted Seafoam Spraguescope
641G	Gray Spraguescope
641GR	Green Spraguescope
641LB	Light Blue Spraguescope
641LV	Lavender Spraguescope
641M	Magenta Spraguescope
641N	Navy Spraguescope
641NB	Neon Blue Spraguescope
641NGR	Neon Green Spraguescope
641NO	Neon Orange Spraguescope
641NP	Neon Pink Spraguescope
641NY	Neon Yellow Spraguescope
641P	Pink Spraguescope
641R	Red Spraguescope
641RB	Royal Blue Spraguescope
641TL	Teal Spraguescope
641V	Purple Spraguescope

ADC® SPRAGUE RAPPAPORT STETHOSCOPE ACCESSORIES

640-01	Chestpiece For 640 & 641
640-01D	Chestpiece Drum For 640 & 641
640-02	Retaining Ring For 640 & 641, Adult
640-02S	Retaining Ring For 640 & 641, Pediatric
640-03	Diaphragm For All Spragues, Adult
640-03S	Diaphragm For All Spragues, Pediatric
640-05	Binaurals For 640 & 641
640-05S	Binaural Spring For 640 & 641
640-05TBK	Tubing For 640, 16" Black
640-06	Tubing Caps For 640 & 641, Black
640-09	Eartips, Black, For 640 & 641
640-10	Eartips, Clear
640-11	Bell, 1 1/4", For 640, 641 & 645
640-12	Bell, 1", For 640, 641 & 645
640-13	Bell, 3/4", For 640, 641 & 645
640-14	Aluminum Disc For 640 & 641, Adult
640-14A	Diaphragm Assembly For 640 & 641
640-15	Aluminum Disc For 640 & 641, Pediatric
640-15A	Diaphragm Assembly For 640 & 641, Pediatric
640-16	Tubing Clip For 640 & 641
640-17	Accessory Kit For 640 & 641
640-18	Accessory Pouch For 640 & 641

ADC® SPRAGUE RAPPAPORT STETHOSCOPE ACCESSORIES (Continued)

640-19	Chestpiece Key For All Spragues
641-05TBB	Boysenberry 22" Tubing For 641 Stethoscope
641-05TBD	Burgundy 22" Tubing For 641 Stethoscope
641-05TBK	Black 22" Tubing For 641 Stethoscope
641-05TCA	Caribbean Blue 22" Tubing For 641 Stethoscope
641-05TDG	Dark Green 22" Tubing For 641 Stethoscope
641-05TFG	Frosted Glacier 22" Tubing For 641 Stethoscope
641-05TFL	Frosted Lilac 22" Tubing For 641 Stethoscope
641-05TFS	Frosted Seafoam 22" Tubing For 641 Stethoscope
641-05TG	Gray 22" Tubing For 641 Stethoscope
641-05TGR	Green 22" Tubing For 641 Stethoscope
641-05TLB	Light Blue 22" Tubing For 641 Stethoscope
641-05TLV	Lavender 22" Tubing For 641 Stethoscope
641-05TM	Magenta 22" Tubing For 641 Stethoscope
641-05TN	Navy 22" Tubing For 641 Stethoscope
641-05TNB	Neon Blue 22" Tubing For 641 Stethoscope
641-05TNGR	Neon Green 22" Tubing For 641 Stethoscope
641-05TNO	Neon Orange 22" Tubing For 641 Stethoscope
641-05TNP	Neon Pink 22" Tubing For 641 Stethoscope
641-05TNY	Neon Yellow 22" Tubing For 641 Stethoscope
641-05TP	Pink 22" Tubing For 641 Stethoscope
641-05TR	Red 22" Tubing For 641 Stethoscope
641-05TRB	Royal Blue 22" Tubing For 641 Stethoscope
641-05TTL	Teal 22" Tubing For 641 Stethoscope
641-05TV	Purple 22" Tubing For 641 Stethoscope

ADC® NURSE COMBO

Combines the Pocket Pal II™ Kit complete with three instruments and choice of 8 Velcro® tag colors with an ADSCOPE™ 641 Sprague Rappaport scope with 22" tubing color coordinated to the Velcro® tags and a nine piece accessory kit.

116-641BDQ	Nurse Combo, Burgundy
116-641BKQ	Nurse Combo, Black
116-641DGQ	Nurse Combo, Dark Green
116-641NQ	Nurse Combo, Navy
116-641NYQ	Nurse Combo, Neon Yellow
116-641RBQ	Nurse Combo, Royal Blue
116-641TLQ	Nurse Combo, Teal
116-641VQ	Nurse Combo, Purple



OMRON SPRAGUE RAPPAPORT-TYPE STETHOSCOPES

5-in-1 multipurpose Quintscope™ has approximately 22" long tubing (except #416) and silver chestpiece. Accessories include three different size open bells, two different size replacement diaphragms, and two pairs of eartips in a vinyl storage case.

416-22-BLK	Black Stethoscope
416-22-DB	Dark Blue Stethoscope
416-22-DKGN	Dark Green Stethoscope



3M™ LITTMANN® CARDIOLOGY III STETHOSCOPE

Outstanding acoustic performance and exceptional versatility characterize the Littmann Cardiology III stethoscope. The innovative design offers a patented tunable diaphragm on each side of the chestpiece. The large side can be used for adult patients, while the small side is especially useful for pediatric or thin patients. Features patented Littmann Snap Tight soft-sealing eartips, patient-friendly, nonchill rim and diaphragms. Five-year warranty on chestpiece and binaural, lifetime service and repair policy.

3127	Cardiology Stethoscope, 22" Black Tubing
3128	Cardiology Stethoscope, 27" Black Tubing
3129	Cardiology Stethoscope, 27" Burgundy Tubing
3130	Cardiology Stethoscope, 27" Navy Blue Tubing
3134	Cardiology Stethoscope, 27" Hunter Green Tubing
3135	Cardiology Stethoscope, 27", Plum Tubing

3M™ LITTMANN® CARDIOLOGY S.T.C. STETHOSCOPES

Outstanding acoustics are combined with a unique and easy to grasp chestpiece that allows easy movement between auscultation sites. Single-sided chestpiece provides both bell and diaphragm frequency responses. Dual lumen design helps to eliminate noise artifact. Comfortably angled, anatomically correct headset. Patient-friendly, nonchill rim. Six-year warranty.

4470	22" Black Cardiology S.T.C. Stethoscope
4471	27" Black Cardiology S.T.C. Stethoscope
4472	27" Burgundy Cardiology S.T.C. Stethoscope
4473	27" Navy Blue Cardiology S.T.C. Stethoscope
4474	27" Hunter Green Cardiology S.T.C. Stethoscope
4475	27" Plum Cardiology S.T.C. Stethoscope



3M™ LITTMANN® MASTER CARDIOLOGY STETHOSCOPE

A patented tunable diaphragm and superb craftsmanship offer the ultimate performance in a superior-single sided stethoscope meets the exacting standards of health care professionals who demand the highest acoustic response. Features patented Littmann Snap Tight soft-sealing eartips, patient-friendly, nonchill rim and diaphragm and special procedures adapter. Seven-year warranty, plus lifetime service and repair policy.

2159	Cardiology Stethoscope, 22" Black Tubing
2160	Cardiology Stethoscope, 27" Black Tubing
2161	Cardiology Stethoscope, 27" Black Plated Chestpiece & Eartubes
2163	Cardiology Stethoscope, 27" Burgundy Tubing
2164	Cardiology Stethoscope, 27" Navy Blue Tubing
2165	Cardiology Stethoscope, 27" Hunter Green Tubing
2167	Cardiology Stethoscope, 27" Plum Tubing



ADC® ADSCOPE™ 600 CARDIOLOGY STETHOSCOPE

Unique design that combines the acoustic response of a traditional bell and diaphragm into a convenient one-sided chestpiece. Bell-type low frequency response is achieved by simply holding the chestpiece in place with light pressure. Firm pressure will produce diaphragm-type frequency performance. A rubber gasket recessed into the chestpiece alters the diaphragm permitting the scope to achieve multi-frequency response. Custom molded extra thick bilumen 19" PVC tubing stays soft and flexible. Extra thick tubing walls minimize extraneous noise. Includes one pair of plastic eartips and two pairs of ADSOFT™ PVC eartips (small and large). Overall length 28". Weight 8.5 oz. Scope ID tag included. Lifetime quality performance guarantee. Latex free.

600BD	ADSCOPE™ 600 Cardiology Stethoscope, Burgundy
600BK	ADSCOPE™ 600 Cardiology Stethoscope, Black
600DG	ADSCOPE™ 600 Cardiology Stethoscope, Dark Green
600GP	ADSCOPE™ 600 Cardiology Stethoscope, Gold Plated, Black
600N	ADSCOPE™ 600 Cardiology Stethoscope, Navy

ADC® ADSCOPE™ 600 CARDIOLOGY STETHOSCOPE ACCESSORIES

600-01BD	Chestpiece For 600, Burgundy
600-01BK	Chestpiece For 600, Black
600-01DG	Chestpiece For 600, Dark Green
600-01GP	Chestpiece For 600GP, Black
600-01N	Chestpiece For 600, Navy
600-03	Diaphragm For 600
600-05GP	Binaural/Tubing For 600GP
600-06BK	Eartips, Black

ADC® ADSCOPE™ 601 CARDIOLOGY STETHOSCOPE

Threaded chestpiece drum allows rapid conversion from a pediatric diaphragm to an adult bell. Chestpiece precisely machined to exacting tolerance from surgical stainless steel for outstanding performance and rugged durability. Extra deep bell for unsurpassed low frequency response. Ultra-sensitive diaphragms for greater amplification and crisper high frequency transmission. Non-chill bell and diaphragm retaining rims for patient comfort. Large internal bore stainless steel binaurals are fixed at a 15° angle for greater comfort. Custom molded bilumen 19" PVC tubing stays soft and flexible. Latex free. Extra thick tubing walls minimize extraneous noise. Double leaf internal spring construction reduces bulk and prevents kinking when folded. Complete with one pair of PVC eartips and ADSOFT™ deluxe PVC eartips in small and large sizes for the ultimate in wearing comfort and acoustic seal. Complete with an accessory case containing the adult bell, 2 spare diaphragms, and two sets of eartips. Scope ID tag included. Weighs 9 oz (adult bell) and 10 oz (pediatric diaphragm). Overall length 28". Lifetime warranty.

601BD	Burgundy Stethoscope
601BK	Black Stethoscope
601DG	Dark Green Stethoscope
601N	Navy Stethoscope

ADC® ADSCOPE™ 601 CARDIOLOGY STETHOSCOPE ACCESSORIES

601-01BD	Chestpiece For 601, Burgundy
601-01BK	Chestpiece For 601, Black
601-01DG	Chestpiece For 601, Dark Green
601-01N	Chestpiece For 601, Navy
601-05BD	Binaural/Tubing For 600 & 601, Burgundy
601-05BK	Binaural/Tubing For 600 & 601, Black
601-05DG	Binaural/Tubing For 600 & 601, Dark Green
601-05N	Binaural/Tubing For 600 & 601, Navy
601-11BD	Bell For 601, Burgundy
601-11BK	Bell For 601, Black
601-11DG	Bell For 601, Dark Green
601-11N	Bell For 601, Navy
601-15BD	Diaphragm Assembly For 601, Burgundy
601-15BK	Diaphragm Assembly For 601, Black
601-15DG	Diaphragm Assembly For 601, Dark Green
601-15N	Diaphragm Assembly For 601, Navy
601-18	Accessory Case For 600 & 601



ADC® ADSCOPE™ 602 CARDIOLOGY STETHOSCOPE

Stainless steel construction for rugged durability. Combination chestpiece precisely machined to exacting tolerances for outstanding acoustic performance. Extra deep bell for unsurpassed low frequency response. Non-chill bell and diaphragm ring for patient comfort. Custom molded 16" bilumen PVC tubing improves acoustic separation. Internal spring construction reduces bulk. Binaurals fixed at 15° angle for greater comfort. Includes one pair of plastic eartips and two pair of PVC eartips with threaded metal inserts for the ultimate in wearing comfort and acoustic seal. Weighs 9 oz. Overall length 25". Latex free. Lifetime warranty.

602BD	Burgundy Stethoscope
602BK	Black Stethoscope
602DG	Dark Green Stethoscope

ADC® ADSCOPE™ 602 CARDIOLOGY STETHOSCOPE ACCESSORIES

- 600-17 Accessory Case For 602, 603, 604, 605, 612, 655
- 602-01BD Chestpiece For 602, Burgundy
- 602-01BK Chestpiece For 602, Black
- 602-01DG Chestpiece For 602, Dark Green
- 602-02BD Burgundy Diaphragm Rim For 600, 601, 602, 603, 609, 612
- 602-02BK Black Diaphragm Rim For 600, 601, 602, 603, 609, 612
- 602-02DG Dark Green Diaphragm Rim For 600, 601, 602, 603, 609, 612
- 602-02FG Frosted Glacier Diaphragm Rim For 600, 601, 602, 603, 609, 612
- 602-02FL Frosted Lilac Diaphragm Rim For 600, 601, 602, 603, 609, 612
- 602-02FS Frosted Seafoam Diaphragm Rim For 600, 601, 602, 603, 609, 612
- 602-02G Gray Diaphragm Rim For 602, 603, 609, 612
- 602-02LV Lavender Diaphragm Rim For 603, 609, 612
- 602-02N Navy Diaphragm Rim For 600, 601, 602, 603, 609, 612
- 602-02R Red Diaphragm Rim For 603, 609, 612
- 602-02RB Royal Blue Diaphragm Rim For 602, 603, 609, 612
- 602-02TL Teal Diaphragm Rim For 603, 609, 612
- 602-03 Diaphragm For 600GP, 601, 602, 603, 612
- 602-04BD Burgundy Non-Chill Rings For 601, 602
- 602-04BK Black Non-Chill Rings For 601, 602
- 602-04DG Dark Green Non-Chill Rings For 601, 602
- 602-04N Navy Non-Chill Rings For 601, 602
- 602-05BD Binaural/Tubing For 602, Burgundy
- 602-05BK Binaural/Tubing For 602, Black
- 602-05DG Binaural/Tubing For 602, Dark Green
- 602-07BK ADSOFT Eartips, Black, Large
- 602-07G ADSOFT Eartips, Gray, Large
- 602-08BK ADSOFT Eartips, Black, Small
- 602-08G ADSOFT Eartips, Gray, Small
- 602-09 ADSOFT Eartips, White
- 602-02V Purple Diaphragm Rim For 600, 601, 602, 603, 609, 612



PT2B27



PT8B27

DRG PURETONE™ CARDIOLOGY STETHOSCOPES

Provides the very best acoustics available through a dual lumen tube and a double-sided chestpiece. Patented dual lumen ENR™ (External Noise Reduction) tube technology significantly reduces the amount of external noise while optimizing internal acoustics. 40% less weight as a result of titanium chestpiece construction; also available in stainless steel. Both with conventional adult- and pediatric-sized diaphragm functionality and the added convenience of a traditional bell. GELseal™ eartips and Anti-Slip™ rings add to the comfort and style. Five year warranty.

- PT2B27 Cardiology Stethoscope, 27" Black Tube
- PT2N27 Cardiology Stethoscope, 27" Navy Tube
- PT8B27 Cardiology Ti-Lite Stethoscope, 27" Black Tube
- PT8N27 Cardiology Ti-Lite Stethoscope, 27" Navy Tube



DRG PURETONE™ MENTOR III CARDIOLOGY STETHOSCOPES

- PT11B27 Mentor III Cardiology Stethoscope, 27" Black Tube
- PT11N27 Mentor III Cardiology Stethoscope, 27" Navy Tube



DRG PURETONE™ STAFF CARDIOLOGY STETHOSCOPES

- PT10B27 Staff Cardiology Stethoscope, 27" Black Tube
- PT10N27 Staff Cardiology Stethoscope, 27" Navy Tube
- PT12B27 Staff Cardiology Ti-Lite Stethoscope, 27" Black Tube
- PT12N27 Staff Cardiology Ti-Lite Stethoscope, 27" Navy Tube



5079-321



5079-325

WELCH ALLYN HARVEY DELUXE DOUBLE AND TRIPLE HEAD STETHOSCOPES

The Harvey™ DLX Stethoscope is our most prestigious and durable model with dual-lumen PVC tubing, brass chestpieces and binaurals, and triple-leaf binaural springs. Optimized for listening to higher frequency and critically important sounds like murmurs, clicks, and ejection sounds. Provides clear transmission of the first and second heart sounds. Ideally sized bell bridges intercostal space between the ribs. Tightly constrained diaphragm provides superior response at large frequencies. Additional corrugated diaphragm for detecting mid-range and low frequencies.

- 5079-321 Harvey Deluxe Triple Head Stethoscope, Black
- 5079-322 Harvey Deluxe Triple Head Stethoscope, Burgundy
- 5079-323 Harvey Deluxe Triple Head Stethoscope, Navy
- 5079-324 Harvey Deluxe Triple Head Stethoscope, Forest Green
- 5079-325 Harvey Deluxe Double Head Stethoscope, Black
- 5079-326 Harvey Deluxe Double Head Stethoscope, Burgundy
- 5079-327 Harvey Deluxe Double Head Stethoscope, Navy
- 5079-328 Harvey Deluxe Double Head Stethoscope, Forest Green



WELCH ALLYN ORIGINAL HARVEY DOUBLE & TRIPLE HEAD STETHOSCOPE ACCESSORIES

- 5079-05 Corrugated Diaphragm Only
- 5079-09 Flat Pediatric Chestpiece
- 5079-102 Tubing, 2 Tubed Joined For Use w/two-piece spindle
- 5079-103 Name Clip For Use w/5079-102 & 105
- 5079-104 Flat Diaphragm Only
- 5079-105 Tubing, 2 Tubes Joined For Use w/single-piece spindle
- 5079-11 Pediatric Diaphragm Only
- 5079-20 Accessory Kit, Includes 1 Pediatric Chestpiece, 1 Flat Diaphragm, 2 Corrugated Diaphragms, 1 Pair Soft Rubber Eartips, 1 Pair Small Delrin Eartips
- 5079-315 Accessory Kit
- 5079-318 Bell Chestpiece Assembly
- 5079-32 Eartips, Rubber
- 5079-33 Eartips, Small, White, Pair
- 5079-60 Binaurals
- 5079-61 Corrugated Adult Chestpiece
- 5079-62 Flat Adult Chestpiece
- 5079-63 Tubing, 2 Separate Tubes For Use w/two-piece spindle
- 5079-65 Bell Chest Piece



5079-27




5079-28

WELCH ALLYN ORIGINAL HARVEY DOUBLE & TRIPLE HEAD STETHOSCOPES

Original Harvey Stethoscopes come with solid brass chestpieces and dual-lumen latex tubes. Features three easily rotated chestpieces (bell, corrugated, diaphragm, and flat diaphragm). Clearly detects the full range of important sounds from low-frequency gallops to the faint high-frequency blow of aortic insufficiency. Sound detection capability of the three chestpieces allows for the most complete auscultatory examination.

- 5079-27 Triple Head Stethoscope
5079-28 Double Head Stethoscope



ADC® ADSOFT™ 655 AMPLIFYING STETHOSCOPE

Diaphragm only type chestpiece fabricated from lightweight acrylic provides up to twice the amplification of conventional acoustic scopes. No electronics to malfunction, no batteries to wear out. Stainless steel binural with double leaf internal spring fixed at 15° angle for greater comfort. Reinforcing yoke molded into flexible 22" PVC Y tubing. ADSOFT™ deluxe PVC eartips in small and large size. Scope ID tag included. Weighs 5 oz. Overall length 31". Lifetime quality performance guarantee. Latex free.

- 655BD Burgundy Stethoscope
655BK Black Stethoscope
655DG Dark Green Stethoscope
655N Navy Stethoscope
655R Red Stethoscope
655RB Royal Blue Stethoscope

Accessories

- 655-01 Chestpiece For 655
655-02 Retaining Ring For 655
655-03K Diaphragm Kit For 655



ADC® ADSOFT™ ADVANCED ELECTRONIC SCOPE

Combines traditional scope appearance with state-of-the-art electronics for the ultimate in acoustic performance and ease of use. Amplifies up to 18x over traditional acoustic scopes. 8 level volume control. Three frequency modes. Bell (20-230Hz) for cardiac sounds. Diaphragm (100-800Hz) for lung sounds. Extended Range (20-800Hz). Electronic transmission of sounds to separate speakers near earpiece for improved performance. LCD with iconic display of heart and lung for more intuitive use, and volume level and battery strength. Impact muting switch protects clinician's ears from high intensity impact noises. Battery warning signals when replacement needed. Binaurals fixed at 15° angle for snug aural fit. Specially constructed Adsoft™ PVC eartips for the ultimate in wearing comfort and acoustic seal. Lightweight and compact - weighs under 7.5 oz., measures just 36". Auto shut off extends battery life. Latex Free. Uses 2 "AA" batteries. 1 year warranty.

- 657BK Advanced Electronic Scope, Black



WELCH ALLYN MEDITRON™ ELECTRONIC STETHOSCOPE SYSTEM

Welch Allyn has teamed up with Meditron AS of Norway to provide the first computer power and Auscultation combination. A unique, patented listening system that dramatically enhances sound clarity. Ideal educational tool. Records sounds and images and speeds consultant review and future analysis (based on Microsoft™ wav files). Should be used by any physical examination where additional amplification is required, including doctor's offices, EMS sites, veterinarian's office, tele-medicine site, hospitals, and medical schools. Reimbursement is the same as regular physical examination (multiple codes). 200 hours or four years of use from the battery. Excellent for use with difficult patients (weak sounds, obese patients, children). May be cleaned with a lightly moistened cotton wipe using chlorhexidine or alcohol. Does not contain any latex.

- 5079-400 Electronic Stethoscope

Accessories

- 5079-401 Accessory Package
5079-414 Headset For Additional Users
5079-430 Analyzer Software Accessory



DRG ECHO™ AMPLIFIED STETHOSCOPES

Highlighted by the benefits of both a conventional and digital stethoscope in one lightweight, portable design. Digital amplification assists in the identification of even the most subtle cardiopulmonary events, while its conventional bell and diaphragm modes provide clinicians with the security against unexpected battery failure during critical clinical interactions. 400% increase in amplification when in digital mode over typical stethoscopes, a heightened frequency response range, and batteries warranted to last up to 100 hours of use. Switches easily between three modes of operation: digital, diaphragm, and bell. Conventional acoustics accompany the bell and diaphragm modes. Patented dual lumen ENR™ (External Noise Reduction) tube technology significantly reduces the amount of external noise while optimizing internal acoustics. GELseal™ eartips and Anti-Slip™ rings are included. One year warranty.


- EC3B27 Amplified Stethoscope, 27" Black Tube



3M™ LITTMANN® TEACHING STETHOSCOPES

A stethoscope with two headsets which allow instructor and student to perform simultaneous auscultation. Patented tunable diaphragms.

- 2138 Classic II S.E. Teaching Stethoscope, Bell & Diaphragm w/2 headsets, 40" 1 Tube Black, 1 Tube Burgundy
2139 Master Classic II Teaching Stethoscope, Single-Sided Chestpiece w/2 headsets, 40" 1 Tube Black, 1 Tube Burgundy



ADC® PROSCOPE™ TEACHING STETHOSCOPE

Choice of diaphragm only (660 Series) or combination diaphragm/bell (670 Series) PROSCOPE™ chestpiece. Two sets of adjustable chrome plated brass binaurals. Two sets of flexible 22" and one straight 10" length of black PVC tubing. Two spare pair of black mushroom eartips included. Overall length 39". Inspected and packaged in the U.S.A. One year quality performance guarantee. Latex free.

- 661 PROSCOPE™ 661 Teaching Scope, Nurse Style
- 671 PROSCOPE™ 671 Teaching Scope, Dualhead

Accessories

- 661-02 Diaphragm Retaining Ring For 660, 661, 670, 671
- 661-03 Diaphragm For 660, 661, 670, 671



OMRON DUAL-HEAD TRAINING SCOPE

Two-channel stethoscope provides accurate sound transmission to both student and instructor. Dual-head chest piece converts from diaphragm to an open bell head with a single turn of the stem.

- 0405 Dual-Head Training Scope

TRIMLINE TEACHING STETHOSCOPE

- SK20358 Teaching Stethoscope w/Two Sets of Binaurals



DRG ECHO+PLUS™ DIGITAL RECORDING STETHOSCOPES

Versatility and digitally-enhanced amplification plus the diagnostic and teaching options of digital record and replay. Allows clinicians to instantly capture critical, unexpected heart sounds for future analysis or instant consultation with colleagues. Assists instructors in teaching students the fundamentals of auscultation. 400% increase in amplification when in digital mode over typical stethoscopes, a heightened frequency response range, and batteries warranted to last up to 100 hours of use. Switches easily between three modes of operation: digital, diaphragm, and bell. Allows fifteen seconds of recording time and two replay speeds. Conventional acoustics accompany the bell and diaphragm modes. Patented dual lumen ENR™ (External Noise Reduction) tube technology significantly reduces the amount of external noise while optimizing internal acoustics. GELseal™ eartips and Anti-Slip™ rings are included. One year warranty.

- EC2B27 Digital Recording Stethoscope, 27" Black Tube

DRG ECHO+PLUS™ TEACHING MODULE

Effective, yet affordable instrumentation to assist in teaching the fundamentals of auscultation as well as the identification of more complicated cardiac pathologies. The instructor can record patient sounds directly at the bedside using the DataLink and replay them through the teaching module to students who are using their own stethoscope in a simulated setting. Ten individual pushbutton sound storage locations. Nine separate listening stations and universal volume control. Self-contained design is effective, portable, and a valuable asset. Universal 110V power requirements. One year warranty.

- TM100 Teaching Module



DRG ECHO+PLUS™ WITH DATALINK™ STETHOSCOPES

The push-button record, playback, and digital amplification features allows the ability to quickly download recorded patient data. By plugging in directly to a standard PC, the DataLink feature enable .wav file outputs for patient tracking and archiving, case study slide presentations and unlimited educational applications. A simple device to more effectively communicate with your students and remain current on your patient's progress. 400% increase in amplification when in digital mode over typical stethoscopes, a heightened frequency response range, and batteries warranted to last up to 100 hours of use. Switches easily between three modes of operation: digital, diaphragm, and bell. Allows fifteen seconds of recording time and two replay speeds. Conventional acoustics accompany the bell and diaphragm modes. Patented dual lumen ENR™ (External Noise Reduction) tube technology significantly reduces the amount of external noise while optimizing internal acoustics. GELseal™ eartips and Anti-Slip™ rings are included. One year warranty.

- EC1B27 DataLink Recording Stethoscope, 27" Black Tube





3M™ DOUBLE STICK DISCS

Pre-cut, hypoallergenic discs that are ideal for securing 3M™ Littmann® stethoscope chestpieces. The discs have an 11/32" aperture and a 1 1/4" diameter.

- 2181 Double Stick Discs w/11/32" aperture, 102/pk, 6 pk/cs

Our new prepping solution
can take a few things off your mind.



New 3M™ One-Step Patient Prep (Povidone Iodine Topical Solution, 7.5% w/w)

- One-step convenience in a non-alcohol-based product
- In a 2 minute application, provides the efficacy of a 5 minute scrub followed by a paint
- Continues to kill after blood and/or saline challenge
- Delivers improved drape adhesion
- Resists wash-off by blood or saline
- Provides persistent activity on skin for up to 24 hours. To learn more, visit www.3m.com/infectionprevention.

*Do not use in the eye. © 3M 2005 All Rights Reserved

3M Medical

from the 3M Health Care family

3M



Your one-stop shop for office setups

Look to Welch Allyn for all of these high quality products:

AEDs
Defibrillators
Otoscopes
Thermometers

Audiometers
ECGs
Lighting Products
Vag Specs

Blood Pressure
Laryngoscopes
Stethoscopes
Vital Signs Monitors

Colposcopes
Ophthalmoscopes
Tympanometric Devices
Headlights

Welch Allyn's easy-to-use, intuitive products were developed with busy clinicians in mind and are designed to work together. For example, our Integrated Diagnostic System combines all the standard physical exam products—ophthalmoscope, otoscope, blood pressure gauge, thermometer and ear tips—in one easy-to-order, easy-to-install package. This flexible product can be configured to fit any doctor's need, and is a great way to start every office setup.

To learn more about Welch Allyn products, visit our website at www.welchallyn.com

© 2005 Welch Allyn MC 3197

WelchAllyn®

Advancing Frontline Care™


3M™ LITTMANN® STETHOSCOPE BINAURALS

Includes eartips, eartubes and tubing.

12237612523	Binaurals For Select & Lightweight, 28" Light Blue
1223768527	Binaurals For Classic II S.E., Classic II & Classic II Pediatric & Infant, 28" Black
12237685289	Binaurals For Classic II S.E., Classic II & Classic II Pediatric & Infant, 28" Gray
12237685305	Binaural Assembly for Littmann® Classic II Pediatric, Navy Blue Tubing, 28"
12237685321	Binaurals For Classic II S.E., Classic II & Classic II Pediatric & Infant, 28" Red
12237685339	Binaurals For Classic II S.E., Classic II & Classic II Pediatric & Infant, 28" Burgundy
7880555776	Binaurals For Classic II S.E., Classic II & Classic II Pediatric & Infant, 28" Pink
78806297442	Binaurals For Master Classic II, 27" Navy
78806297459	Binaurals For Master Classic II, 27" Burgundy
78806299273	Binaurals For Master Classic II, 27" Black
78806311375	Binaurals For Master Classic II, 27" Black, Black Plated Eartubes
78806311425	Binaurals For Master Classic II, 27" Black, 23 Karat Gold-Plated Eartubes
78807802844	Binaurals For Teaching Stethoscope, 40" Black w/2 headsets For Classic II S.E.
78807802851	Binaurals For Teaching Stethoscope, 40" Black w/2headsets For Master Classic II
78807812066	Binaurals For Select & Lightweight, 28" Black
78807812074	Binaurals For Select & Lightweight, 28" Pink
78807833831	Binaural Assembly For Littmann Master Classic II Veterinarian Stethoscope, 32" Navy Blue
78807845165	Binaurals For Classic II S.E., Classic II & Classic II Pediatric & Infant, 28" Royal Blue
78807845181	Binaurals For Select & Lightweight, 28" Burgundy
78807845207	Binaurals For Select & Lightweight, 28" Royal Blue
78807858390	Binaurals For Select & Lightweight, 28" Purple
78807858465	Binaurals For Master Classic II, 27" Purple
78807865023	Binaurals For Classic II S.E., Classic II & Classic II Pediatric & Infant, 28" Hunter Green
78807865031	Binaurals For Classic II S.E., Classic II & Classic II Pediatric & Infant, 28" Purple
78807865825	Binaurals For Master Cardiology, Cardiology III & Cardiology II S.E., 22" Black, Short
78807865833	Binaurals For Master Cardiology, Cardiology III & Cardiology II S.E., 27" Black
78807865858	Binaurals For Master Cardiology, Cardiology III & Cardiology II S.E., 27" Navy
78807865866	Binaurals For Master Cardiology, Cardiology III & Cardiology II S.E., 27" Hunter Green
78807865874	Binaurals For Master Cardiology, Cardiology III & Cardiology II S.E., 27" Hunter Green
78807866211	Binaurals For Master Cardiology, Cardiology III & Cardiology II S.E., 27" Black, Black Plated Eartubes
78807866229	Binaurals For Master Cardiology, Cardiology III & Cardiology II S.E., 27" Black, 23 Karat Gold-Plated Eartubes
78807866359	Binaurals For Classic II S.E., Classic II & Classic II Pediatric & Infant, 28" Raspberry
78807866367	Binaurals For Classic II S.E., Classic II & Classic II Pediatric & Infant, 28" Caribbean Blue
78807866508	Binaurals For Master Classic II, 27" Caribbean Blue
78807866516	Binaurals For Master Classic II, 27" Hunter Green
78807866581	Binaural Assembly For Littmann Select & Lightweight Stethoscope, 28" Raspberry
7880786659	Binaurals For Select & Lightweight, 28" Caribbean Blue
7880789464	Binaurals For Cardiology S.T.C., 27" Burgundy
7880789465	Binaurals For Cardiology S.T.C., 27" Navy Blue
7880789466	Binaurals For Cardiology S.T.C., 27" Black
7880789467	Binaurals For Cardiology S.T.C., 27" Hunter Green
7880789513	Binaurals For Cardiology S.T.C., 22" Black
78808303578	Binaural Assembly for Littmann® Master Cardiology and Cardiology III, Plum Tubing, 27"
78808303586	Binaural Assembly for Littmann® Master Classic II, Ceil Blue Tubing, 27"
78808303594	Binaural Assembly for Littman® Classic II S.E., Ceil Blue Tubing, 28"
78808303602	Binaural Assembly for Littmann® Classic II S.E., Seafoam Green Tubing, 28"
78808303610	Binaural Assembly for Littmann® Classic II Infant, Lilac Tubing, 28"
78808303628	Binaural Assembly for Littmann® Classic II S.E., Orange Tubing, 28"
78808303636	Binaural Assembly for Littmann® Select, Ceil Blue Tubing, 28"
78808303644	Binaural Assembly for Littmann® Select, Seafoam Green Tubing, 28"
78808303651	Binaural Assembly for Littmann® Select, Lilac Tubing, 28"
78808303909	Binaural Assembly for Littman® Cardiology S.T.C., Plum Tubing, 27"

**3M™ LITTMANN® STETHOSCOPE DIAPHRAGMS**

36559	Stethoscope Diaphragms For the Cardiology II, Cardiology, Classic II, Classic, Lightweight, Blood Pressure & Differential Stethoscope, Printed w/Littmann Logo, 10/bg
36560	Stethoscope Diaphragms For Classic II Pediatric Stethoscope, Printed w/Littmann Logo, 10/bg
36561	Stethoscope Diaphragms For Classic II Infant Stethoscope, Printed w/Littmann Logo, 10/bg


**3M™ LITTMANN® STETHOSCOPE EARTIPS**

37808	Snap-Tight Stethoscope Eartips, Firm, Gray, Large, 10 pr/bg
37809	Snap-Tight Stethoscope Eartips, Soft-Sealing, Black, Small, 10 pr/ctn
37810	Snap-Tight Stethoscope Eartips, Soft-Sealing, Gray, Small, 10 pr/ctn
37811	Snap-Tight Stethoscope Eartips, Soft-Sealing, Black, Large, 10 pr/ctn
37812	Snap-Tight Stethoscope Eartips, Soft-Sealing, Gray, Large, 10 pr/ctn
38742	Snap-Tight Stethoscope Eartips, Firm, Gray, Small, 10 pr/bg
51161	Soft-Sealing Eartips, Light Brown, Large, 10 pr/ctn

**3M™ LITTMANN® STETHOSCOPE IDENTIFICATION TAGS**

Fits all sizes.

2170	Stethoscope Identification Tags, Black, 5/cs
2171	Stethoscope Identification Tags, Gray, 5/cs

**3M™ LITTMANN® STETHOSCOPE NONCHILL BELL SLEEVES**

For use with traditional two-sided chestpieces on 3M™ Littmann™ stethoscopes.

36544	Stethoscope Nonchill Bell Sleeves For Classic II S.E., Classic II & Lightweight, Black, 10/bg
36545	Stethoscope Nonchill Bell Sleeves For Classic II S.E., Classic II & Lightweight, Gray, 10/bg
36546	Stethoscope Nonchill Bell Sleeves For Classic II Pediatric, Black, 5/bg
36547	Stethoscope Nonchill Bell Sleeves For Classic II Infant, Black, 5/bg
36568	Stethoscope Nonchill Bell Sleeves For Classic II Pediatric, Gray, 5/bg
36569	Stethoscope Nonchill Bell Sleeves For Classic II Infant, Gray, 5/bg
36570	Stethoscope Nonchill Bell Sleeves For Cardiology III, Cardiology II S.E. & Cardiology II, Black, 10/bg
36571	Stethoscope Nonchill Bell Sleeves For Cardiology III, Cardiology II S.E. & Cardiology II, Gray, 10/bg
51162	Stethoscope Nonchill Bell Sleeve For Littmann Lightweight II S.E. Stethoscope, Black, 10/bg
51163	Stethoscope Nonchill Bell Sleeve For Littmann Lightweight II S.E. Stethoscope, Light Brown, 10/bg

3M™ LITTMANN® STETHOSCOPE REPAIR KIT CASE

44525 Stethoscope Repair Kit Case

3M™ LITTMANN® STETHOSCOPE SNAP-ON RIMS

- 36548 Stethoscope Snap-On Rims For Cardiology III (large side), Cardiology II S.E., Classic II S.E., Select & Lightweight, Black, 10/bg
- 36549 Stethoscope Snap-On Rims For Cardiology III (large side), Cardiology II S.E., Classic II S.E., Select & Lightweight, Gray, 10/bg
- 36550 Stethoscope Snap-On Rims For Cardiology III (small side) & Classic II Pediatric, Black, 5/bg
- 36551 Stethoscope Snap-On Rims For Cardiology III (small side) & Classic II Pediatric, Gray, 5/bg
- 36552 Stethoscope Snap-On Rims For Classic II Infant, Black, 5/bg
- 36553 Stethoscope Snap-On Rims For Classic II Infant, Gray, 5/bg



3M™ LITTMANN® STETHOSCOPE SPECIAL PROCEDURE ADAPTORS

- 36563 Stethoscope Adaptors For Cardiology II S.E. & Cardiology II, Black, 3/bg
- 36564 Stethoscope Adaptors For Cardiology II S.E. & Cardiology II, Gray, 3/bg
- 36565 Stethoscope Adaptors For Master Cardiology, Black, 3/bg
- 36566 Stethoscope Adaptors For Master Cardiology, Gray, 3/bg



3M™ LITTMANN® TUNABLE DIAPHRAGM & RIM ASSEMBLY

- 36554 Diaphragm & Rim Assembly For Master Cardiology, Black Rim, 5/bg
- 36555 Diaphragm & Rim Assembly For Master Cardiology, Gray Rim, 5/bg
- 36556 Diaphragm & Rim Assembly For Cardiology III (large side), Cardiology II S.E., Master Classic II, Classic II S.E. & Select, Black Rim, 5/bg
- 36557 Diaphragm & Rim Assembly For Cardiology III (large side), Cardiology II S.E., Master Classic II, Classic II S.E. & Select, Gray Rim, 5/bg
- 36572 Diaphragm & Rim Assembly For Cardiology III (small side), Black Rim, 5/bg
- 36573 Diaphragm & Rim Assembly For Cardiology III (small side), Gray Rim, 5/bg
- 51159 Diaphragm & Rim Assembly For Littmann Lightweight II S.E. Stethoscope, Black Rim, 5/bg
- 51160 Diaphragm & Rim Assembly For Littmann Lightweight II S.E. Stethoscope, Light Brown Rim, 5/bg
- 51185 Diaphragm & Rim Assembly For Littmann Cardiology S.T.C., Black Rim, 5/bg



GET100-1



SS50-4CL-4C



SS50-4PK-4C-1

DRG STETHOSCOPE DIAPHRAGMS FOR ECHO™, ECHO+PLUS™ & PURETONE™

Durable, chemically inert and acoustically sensitive, and provide a level of protection unmatched in the industry. Add AgION™ Antimicrobial Compound to help prevent the growth of life-threatening microbes. Latex free, C-Flex, medical-grade synthetic elastomer. Chemically inert, non-allergenic, and resistant to ripping or hardening over time. Warm against patient's skin, unlike traditional anti-chill rings. Disposable diaphragms are single use (packaged clean); antimicrobial diaphragms with AgION can be used for up to one month before replacement. Standard diaphragms available in white; AgION and disposable diaphragm colors are: red, green, orange, royal blue and maroon.

- BR100-C-BK Bell Ring, Cardiology, Black
- D100-A Stethoscope Hard Diaphragm, Adult
- D100-I Stethoscope Hard Diaphragm, Infant
- D100-P Stethoscope Hard Diaphragm, Pediatric
- DR100-BK Diaphragm Ring, Black
- GET100-1 GELSeal Eartips, Add Color Code to Part #: BK-Black, BL-Blue, GR-Green, RD-Red, or YL-Yellow
- HSAD150-A Diaphragm Anti-Microbial, Adult, Specify Color: Green, Orange, Blue, Red, 50/bx
- HSAD150-I Diaphragm Anti-Microbial, Infant, Specify Color: Green, Orange, Blue, Red, 50/bx
- HSAD150-P Diaphragm Anti-Microbial, Pediatric, Specify Color: Green, Orange, Blue, Red, 50/bx
- HUM100-S Hummer Relex Hammer & Tuning Fork
- ID100-BK ID Tag
- SS50-4CL-4C SafeSeal, 4 Column 4 Color, Adult, 200/bx
- SS50-4CL-4C-I SafeSeal, 4 Color, Infant, 200/bx
- SS50-4CI-4C-P SafeSeal, 4 Color, Pediatric, 200/bx
- SS50-4PK-4C-I SafeSeal, 4 Color (Blue, Green, Orange, Red), Infant, 200/bx
- SS50-4PK-4C-P SafeSeal, 4 Color (Blue, Green, Orange, Red), Pediatric, 200/bx
- SS50-4PK-PAT SafeSeal, Patriotic (2 Red, 2 Blue), Adult, 200/bx
- BR100-CL-BK Bell Ring, Traditional, Black
- SS50-4PK-4C SafeSeal, 4 Color (Blue, Green, Orange, Red), Adult, 200/bx

OMRON STETHOSCOPE REPLACEMENT PARTS

- 0418 Binaurals For Sprague Rappaport-Type
- 0426 Chestpiece For Bowles 407
- 0431 Plastic Eartips For Lightweight Binaurals, Black
- 0431CLEAR Soft Clear Eartips For Sprague Rappaport-Types
- 0437 Plastic Eartips For Lightweight Binaurals, White
- 0463 Retaining Rings For 412 Series
- 0485 Chestpiece Sprague Rappaport-Type
- 0487 Mushroom Eartips For Lightweight Binaurals, Black
- 423 Black PVC "Y" Tubing
- 423BLUE Blue PVC "Y" Tubing
- 423GRAY Gray PVC "Y" Tubing
- 423GREN Green PVC "Y" Tubing
- 423RED Red PVC "Y" Tubing



WELCH ALLYN ELITE® STETHOSCOPE ACCESSORIES

- 5079-108 Nametag
- 5079-115 Pediatric Diaphragm Assembly, Black
- 5079-116 Adult Diaphragm Assembly, Black
- 5079-117 Pediatric Bell, Non-Chill Rim, Black
- 5079-118 Adult Bell, Non-Chill Rim, Black
- 5079-119 22" Tubing, Black
- 5079-120 25" Tubing, Black
- 5079-121 Binaural/Spring Assembly, Black
- 5079-123 Adult Diaphragm Disc Only
- 5079-124 Pediatric Diaphragm Disc Only
- 5079-126 Adult Diaphragm, Non-Chill Rim, Black
- 5079-127 Pediatric Diaphragm, Non-Chill Rim, Black
- 5079-129 Black Tubing Assembly
- 5079-184 Adult Diaphragm, Non-Chill Rim, Gray
- 5079-185 Pediatric Diaphragm, Non-Chill Rim, Gray
- 5079-232 Accessory Kit, Black, Includes a Pediatric Diaphragm, Bell Ring & 3 Additional Set of Eartips
- 5079-233 Comfort Sealing Eartips, Gray
- 5079-272 25" Tubing, Burgundy
- 5079-273 25" Tubing, Navy
- 5079-274 Binaural/Spring Assembly, Burgundy
- 5079-275 Binaural/Spring Assembly, Navy
- 5079-276 Pediatric Bell, Non-Chill Rim, Gray
- 5079-277 Adult Bell, Non-Chill Rim, Gray
- 5079-278 Pediatric Diaphragm Assembly, Gray
- 5079-279 Adult Diaphragm Assembly, Gray
- 5079-280 Elite Accessory Kit, Gray, Includes Pediatric Diaphragm, Bell Ring & 3 Additional Sets of Eartips
- 5079-305 25" Tubing, Forest Green
- 5079-307 Binaural/Spring Assembly, Forest Green
- 5079-314 Comfort Sealing Eartips, Small
- 5079-66 Hard Eartips, Large

WELCH ALLYN SPECTRUM™ ACCESSORIES & PARTS

- 5079-170 Comfort Sealing Eartips, Large
- 5079-172 Hard Eartips, Medium
- 5079-173 Hard Eartips, Large
- 5079-180 Adult Bell, Non-Chill Rim, Black
- 5079-181 Adult Bell, Non-Chill Rim, Gray
- 5079-182 Pediatric Bell, Non-Chill Rim, Black
- 5079-183 Pediatric Bell, Non-Chill Rim, Gray
- 5079-195 Binaural/Spring Assembly & Tubing, 28", Black
- 5079-196 Binaural/Spring Assembly & Tubing, 28", Gray
- 5079-197 Binaural/Spring Assembly & Tubing, 28", Navy
- 5079-199 Binaural/Spring Assembly & Tubing, 28", Burgundy
- 5079-210 Double-Head Chestpiece Assembly, Adult, Gray
- 5079-211 Double-Head Chestpiece Assembly, Adult, Black
- 5079-214 Double-Head Chestpiece Assembly, Pediatric, Black
- 5079-215 Double-Head Chestpiece Assembly, Pediatric, Gray
- 5079-313 Comfort Sealing Eartips, Small



SUMMIT DOPPLERS

Interchangeable probes. Superior sound quality. Extended battery life.

- L150ASD2 No Display Doppler w/2MHz OB probe & recorder
- L150ASD3 No Display Doppler w/3MHz OB probe & recorder
- L150ASD4 No Display Doppler w/4MHz vascular probe & recorder
- L150ASD5 No Display Doppler w/5MHz vascular probe & recorder
- L150ASD8 No Display Doppler w/8MHz vascular probe & recorder
- L150ASDW No Display Doppler w/OB waterproof probe & recorder
- L150RSD2 No Display Doppler w/2MHz OB probe & recharger
- L150RSD3 No Display Doppler w/3MHz OB probe & recharger
- L150RSD4 No Display Doppler w/4MHz vascular probe & recharger
- L150RSD5 No Display Doppler w/5MHz vascular probe & recharger
- L150RSD8 No Display Doppler w/8MHz vascular probe & recharger
- L150RSDW No Display Doppler w/OB waterproof probe & recharger
- L150SD2 No Display Doppler w/2MHz OB probe
- L150SD3 No Display Doppler w/3MHz OB probe
- L150SD4 No Display Doppler w/4MHz vascular probe
- L150SD5 No Display Doppler w/5MHz vascular probe
- L150SD8 No Display Doppler w/8MHz vascular probe
- L150SDW No Display Doppler w/OB waterproof probe
- L250ARSD2 Display Doppler w/2MHz OB probe, recharger & recorder
- L250ARSD3 Display Doppler w/3MHz OB probe, recharger & recorder
- L250ARSD4 Display Doppler w/4MHz vascular probe, recharger & recorder
- L250ARSD5 Display Doppler w/5MHz vascular probe, recharger & recorder
- L250ARSD8 Display Doppler w/8MHz vascular probe, recharger & recorder
- L250ARSDW Display Doppler w/OB waterproof probe, recharger & recorder
- L250RSD2 Display Doppler w/2MHz OB probe & recharger
- L250RSD3 Display Doppler w/3MHz OB probe & recharger
- L250RSD4 Display Doppler w/4MHz vascular probe & recharger
- L250RSD5 Display Doppler w/5MHz vascular probe & recharger
- L250RSD8 Display Doppler w/8MHz vascular probe & recharger
- L250RSDW Display Doppler w/OB waterproof probe & recharger
- L250SD2 Display Doppler w/2MHz OB probe
- L250SD3 Display Doppler w/3MHz OB probe
- L250SD4 Display Doppler w/4MHz vascular probe
- L250SD5 Display Doppler w/5MHz vascular probe
- L250SD8 Display Doppler w/8MHz vascular probe
- L250SDW Display Doppler w/OB waterproof probe
- L350RSD2 Table-Top Doppler w/2MHz OB probe & recharger
- L350RSD3 Table-Top Doppler w/3MHz OB probe & recharger
- L350RSD4 Table-Top Doppler w/4MHz Vascular probe & recharger
- L350RSD5 Table-Top Doppler w/5MHz Vascular probe & recharger
- L350RSD8 Table-Top Doppler w/8MHz Vascular probe & recharger
- L350RSDW Table-Top Doppler w/2MHz OB waterproof probe & recharger
- SD2 2MHz OB Probe Only
- SD3 3MHz OB Probe Only
- SD4 4MHz Vascular Probe Only
- SD5 5MHz Vascular Probe Only
- SD8 8MHz Vascular Probe Only
- SDW 2MHz OB Waterproof Probe Only
- L250ABI Reimbursable Vascular Doppler w/printer & 8 MHz Bi-Directional Probe
- Accessories**
- K150 ABI Kit Includes Aneroid, 2 Cuffs, 50 Forms, Chart & Video

SUMMIT DOPPLER ACCESSORIES

B150	Rechargeable Battery, 3/pk
B155	Alkaline Battery, 3/pk
C150	Transeptic Cleansing Solution
G150	Aquasonic Clean Ultrasound Gel, 60 gm
G155	Aquasonic Clean Ultrasound Gel, 250 gm
H150	Headphones
K160	ABI 2 Cuff Package w/trigger release aneroid, Includes: 1-10cm and 1-12cm Cuff
K170	ABI 4 Cuff Package w/trigger release aneroid, Includes: 2-10cm and 2-12cm Cuffs
K180	Printer Label Paper, 5rolls/pk
K185	Printer Paper, 5 rolls/pk
K200	LifeDop Hand-Held Stand w/tray, 5 Castor
K220	Pinnacle Hand-Held Stand w/tray, 5 Castor
K230	LifeDop Hand-Held Stand w/out tray, 5 Castor
K240	Pinnacle Hand-Held Stand w/out tray, 5 Castor
K250	LifeDop Universal Clamp
K260	LifeDop Carrying Case



P200



P500



P800

NICOLET IMEX POCKET-DOP II™ DOPPLER SYSTEM

Hand-held ultrasonic Doppler features interchangeable probes and built-in speaker for vascular and obstetrical applications. Supplied in soft sided case with 1 or 2 probes, recharger, walkman style headset, gel, batteries and CD ROM of pulse sounds. Includes a five-year parts and one-year labor warranty.

P200	Unit w/2 MHz OB Probe
P230	Unit w/2 & 3 MHz Probes
P250	Unit w/2 & 5 MHz Probes
P280	Unit w/2 & 8 MHz Probes
P2WP	Unit w/2 MHz Waterproof OB Probe
P300	Unit w/3 MHz OB Probe
P350	Unit w/3 & 5 MHz Probes
P380	Unit w/3 & 8 MHz Probes
P500	Unit w/5 MHz Vascular Probe
P580	Unit w/5 & 8 MHz Probes
P800	Unit w/8 MHz Vascular Probe
PAD5	Unit w/5 MHz Probe & ABI Kit
PAD8	Unit w/8 MHz Probe & ABI Kit



NICOLET IMEX STETHODOP™ DOPPLER SYSTEM

Mounts to most bell-style stethoscopes using a custom non-chill ring. Easy one-hand operation. Momentary on/off actuator operates only when depressed to save battery life. Sliding volume control to adjust to individual user comfort level. Super sensitive Doppler. Custom pager style carrying case. Includes one-year labor and parts warranty.

SD50	Imex StethoDop w/5 MHz Probe
SD52	Imex StethoDop w/Cardiology Style Stethoscope
SS10	Imex Cardiology Stethoscope
SS20	Imex Classic Stethoscope



ED20



EN2W



EN2WR

NICOLET IMEX ELITE™

A fully configurable ultrasound Doppler used to determine fetal viability and to assist in monitoring peripheral arterial and venous blood flow. Options include interchangeable OB, vascular and waterproof probes, display or non-display models and alkaline or rechargeable batteries.

ED20	Elite 200 Doppler, Digital Display w/alkaline batteries & 2 MHz OB Probe
ED20R	Elite 200R Doppler, Digital Display w/rechargeable batteries & 2 MHz OB Probe
ED23	Elite 200 Doppler, Digital Display w/alkaline batteries & 2 & 3 MHz OB Probes
ED23R	Elite 200R Doppler, Digital Display w/rechargeable batteries & 2 & 3 MHz OB Probes
ED25	Elite 200 Doppler, Digital Display w/alkaline batteries & 2 & 5 MHz Probes
ED25R	Elite 200R Doppler, Digital Display w/rechargeable batteries & 2 & 5 MHz Probes
ED28	Elite 200 Doppler, Digital Display w/alkaline batteries & 2 & 8 MHz Probes
ED28R	Elite 200R Doppler, Digital Display w/rechargeable batteries & 2 & 8 MHz Probes
ED2W	Elite 200 Doppler, Digital Display w/alkaline batteries & 2 MHz Waterproof Probe
ED2WR	Elite 200R Doppler, Digital Display w/rechargeable batteries & 2 MHz Waterproof Probe
ED30	Elite 200 Doppler, Digital Display w/alkaline batteries & 3 MHz OB Probe
ED30R	Elite 200R Doppler, Digital Display w/rechargeable batteries & 3 MHz OB Probe
ED35	Elite 200 Doppler, Digital Display w/alkaline batteries & 3 & 5 MHz Probes
ED35R	Elite 200R Doppler, Digital Display w/rechargeable batteries & 3 & 5 MHz Probes
ED38	Elite 200 Doppler, Digital Display w/alkaline batteries & 3 & 8 MHz Probes
ED38R	Elite 200R Doppler, Digital Display w/rechargeable batteries & 3 & 8 MHz Probes

NICOLET IMEX ELITE™ (Continued)

ED50	Elite 200 Doppler, Digital Display w/alkaline batteries & 5 MHz Vascular Probe
ED50R	Elite 200R Doppler, Digital Display w/rechargeable batteries & 5 MHz Vascular Probe
ED58	Elite 200 Doppler, Digital Display w/alkaline batteries & 5 & 8 MHz Probes
ED58R	Elite 200R Doppler, Digital Display w/rechargeable batteries & 5 & 8 MHz Vascular Probes
ED80	Elite 200 Doppler, Digital Display w/alkaline batteries & 8 MHz Vascular Probe
ED80R	Elite 200R Doppler, Digital Display w/rechargeable batteries & 8 MHz Probe
ED8A	Elite 200 Doppler, Digital Display w/alkaline batteries, 8 MHz Probe & ABI Kit
EN20	Elite 100 Doppler, Non-Digital Display w/alkaline batteries & 2 MHz OB Probe
EN20R	Elite 100R Doppler, Non-Digital Display w/rechargeable batteries & 2 MHz OB Probe
EN23	Elite 100 Doppler, Non-Digital Display w/alkaline batteries & 2 & 3 MHz OB Probes
EN23R	Elite 100R Doppler, Non-Digital Display w/rechargeable batteries & 2 & 3 MHz OB Probes
EN25	Elite 100 Doppler, Non-Digital Display w/alkaline batteries & 2 & 5 MHz Probes
EN25R	Elite 100R Doppler, Non-Digital Display w/rechargeable batteries & 2 & 5 MHz Probes
EN28	Elite 100 Doppler, Non-Digital Display w/alkaline batteries & 2 & 8 MHz Probes
EN28R	Elite 100R Doppler, Non-Digital Display w/rechargeable batteries & 2 & 8 MHz Probes
EN2W	Elite 100 Doppler, Non-Digital Display w/alkaline batteries & 2 MHz Waterproof Probe
EN2WR	Elite 100R Doppler, Non-Digital Display w/rechargeable batteries & 2 MHz Waterproof Probe
EN30	Elite 100 Doppler, Non-Digital Display w/alkaline batteries & 3 MHz OB Probe
EN30R	Elite 100R Doppler, Non-Digital Display w/rechargeable batteries & 3 MHz OB Probe
EN35	Elite 100 Doppler, Non-Digital Display w/alkaline batteries & 3 & 5 MHz Probes
EN35R	Elite 100R Doppler, Non-Digital Display w/rechargeable batteries & 3 & 5 MHz Probes
EN38	Elite 100 Doppler, Non-Digital Display w/alkaline batteries & 3 & 8 MHz Probes
EN38R	Elite 100R Doppler, Non-Digital Display w/rechargeable batteries & 3 & 8 MHz Probes
EN50	Elite 100 Doppler, Non-Digital Display w/alkaline batteries & 5 MHz Vascular Probe
EN50R	Elite 100R Doppler, Non-Digital Display w/rechargeable batteries & 5 MHz Vascular Probe
EN58	Elite 100 Doppler, Non-Digital Display w/alkaline batteries & 5 & 8 MHz Vascular Probes
EN58R	Elite 100R Doppler, Non-Digital Display w/rechargeable batteries & 5 & 8 MHz Vascular Probes
EN5A	Elite 100 Doppler, Non-Digital Display w/alkaline batteries, 5 MHz Probe & ABI Kit
EN80	Elite 100 Doppler, Non-Digital Display w/alkaline batteries & 8 MHz Vascular Probe
EN80R	Elite 100R Doppler, Non-Digital Display w/rechargeable batteries & 8 MHz Vascular Probe
EN8A	Elite 100 Doppler, Non-Digital Display w/alkaline batteries, 8 MHz Probe & ABI Kit
EN8AR	Elite 100R Doppler, Non-Digital Display w/rechargeable batteries, 8 MHz Probe & ABI Kit

NICOLET IMEX POCKET-DOP 3™ DOPPLER SYSTEM

A valuable aid in detecting pregnancy and monitoring the fetal heartbeat. Features interchangeable probes: the 3 MHz probe is for early fetal heartbeat detection and the 2 MHz probe is for use in later pregnancy or with larger patients. Also features auto-correlation which is a pattern recognition technique that quickly acquires and accurately tracks the complex signal of the fetal heart rate. Comes with a convenient recharging stand to be placed on a counter or mounted to a wall. Covered by a 5-year parts and 1-year labor warranty.

J200	Doppler w/2 MHz OB Probe
J230	Doppler w/2 & 3 MHz OB Probes
J300	Doppler w/3 MHz OB Probe



NICOLET IMEXDOP CT+™ DOPPLER SYSTEM

Counter top/wall mount Doppler designed for hospitals, clinics and heavy volume offices. Features interchangeable probes, convenient recharging stand and a fetal heart rate display in easy-to-read digital numbers. Supplied with recharging stand, gel, battery pack, monaural earphone and mounting kit. Includes a five-year parts, one-year labor warranty.

CT20+	Unit w/2 MHz OB Probe
CT23+	Unit w/2 & 3 MHz Probes
CT25+	Unit w/2 & 5 MHz Probes
CT28+	Unit w/2 & 8 MHz Probes
CT30+	Unit w/3 MHz OB Probe
CT35+	Unit w/3 & 5 MHz Probes
CT38+	Unit w/3 & 8 MHz Probes
CT50+	Unit w/5 MHz Vascular Probe
CT58+	Unit w/5 & 8 MHz Probes
CT80+	Unit w/8 MHz Vascular Probe
CTWP	Unit w/2 MHz Waterproof OB Probe



NICOLET IMEX VERSALAB BI-DIRECTIONAL DOPPLER SYSTEM

Continuous wave doppler transducers optimized for detecting blood flow throughout peripheral parts of the body. Optional color display feature allows quantification, documentation and qualification of blood flow using a color frequency distribution spectral analysis for accurate determination of direction and velocity of blood flow. Fast, built-in printer for quality documentation of all information obtained. A variety of industry standard calculations and indices calculated to complement the diagnostic capabilities of the VersaLab.

VLL40	Unit w/4 MHz LE Bi-Directional Probe
VLL48	Unit w/4 & 8 MHz LE Bi-Directional Probes
VLL80	Unit w/8 MHz LE Bi-Directional Probe
VLS40	Unit w/2 MHz LE Bi-Directional Probe
VLS48	Unit w/4 & 8 MHz LE Bi-Directional Probes
VLS80	Unit w/8 MHz LE Bi-Directional Probe

NICOLET IMEX VERSALAB ACCESSORIES

A220	Personal Stereo Headset
A350	Printer Pack
A355	VersaLab Printer Paper, 5 r/l/cs
A430	Soft-Sided Carrying Case
BD40	4 MHz Bi-Directional Probe
BD80	8 MHz Bi-Directional Probe
C633	Battery Pack
C644	Power Supply 110-250V
C644I	Power Supply 220V
C645	Power Cord, 110V
C645I	Power Cord 220V
DV100	VersaLab APM/APM2 Event Marker
DV400	VersaLab Toco Probe Covers, 5/pk
ST30	Tilt Stand
TV100	VersaLab APM/APM2 Toco Probe
TV180	VersaLab APM2, 1.8 MHz Doppler Probe
TV200	VersaLab APM/APM2, 2MHz Doppler Probe
V200	VersaLab IV Pole Mount



MVAPM



MVAPT

NICOLET IMEX VERSALAB ANTEPARTUM MONITOR™

Easy-to-use antepartum monitor designed for the office or the antepartum testing area of the hospital. Pulsed, directional 1.5 MHz ultrasound locks into the fetal heart rate; sensitive toco transducer tracks maternal contractions. Trace prints out on Z-fold paper for easy scanning and storage. Supplied with chart paper, belt, plastic membrane for toco, gel and patient event marker. Includes a one-year parts and labor warranty.

MVAPM	Single APM Unit w/toco & 2 MHz Doppler Transducer
MVAPT	Twin APM Unit w/toco 1.8 & 2 MHz Doppler Transducers

NICOLET IMEX ANTEPARTUM MONITOR ACCESSORIES

D200	Z-Fold Chart Paper
D300	Belts, Pre-Cut, 2/pk



FD2PUSA/EZ8



FD2PUSA/OP2HS

HUNTLEIGH FETAL DOPPLEX® II POCKET DOPPLER

Advanced pocket Doppler offers three Fetal Heart Rate (FHR) modes, LCD display, auto-correlation, active noise reduction, digital interface to Dopplex® Printa™ to produce FHR tracings. Supplied with one transducer of choice.

FD2PUSA/EZ8	Fetal Dopplex® II w/EZ8 Widebeam Probe
FD2PUSA/OP2HS	Fetal Dopplex® II w/2 MHz Probe
FD2PUSA/OP3HS	Fetal Dopplex® II w/3 MHz Probe
FD2PUSA/VP10HS	Fetal Dopplex® II w/10 MHz Probe
FD2PUSA/VP4HS	Fetal Dopplex® II w/4 MHz Probe
FD2PUSA/VP5HS	Fetal Dopplex® II w/5 MHz Probe
FD2PUSA/VP8HS	Fetal Dopplex® II w/8 MHz Probe



SD2PUSA/EZ8



SD2PUSA/OP2HS

HUNTLEIGH SUPER DOPPLEX® II DOPPLER

Sophisticated pocket Doppler with bi-directional indication of blood flow velocity, LC display, true separated audio signals, clear visual indication of forward and reverse flow, 4 level gain control, active noise reduction, connects to EKG with custom cable. Supplied with one transducer of choice.

SD2PUSA/EZ8	Super Dopplex 2 w/EZ8 Widebeam Probe
SD2PUSA/OP2HS	Super Dopplex 2 w/2 MHz Probe
SD2PUSA/OP3HS	Super Dopplex 2 w/3 MHz Probe
SD2PUSA/VP10HS	Super Dopplex 2 w/10 MHz Probe
SD2PUSA/VP4HS	Super Dopplex 2 w/4 MHz Probe
SD2PUSA/VP5HS	Super Dopplex 2 w/5 MHz Probe
SD2PUSA/VP8HS	Super Dopplex 2 w/8 MHz Probe



HUNTLEIGH FETAL DOPPLEX® I PLUS DOPPLER

Fetal dopplex with LCD display of fetal heart rate, fixed 2 MHz transducer and built-in loudspeaker. Waterproof.

FD1+-P-USA	Fetal Dopplex® I Plus
------------	-----------------------



HUNTLEIGH AUDIO DOPPLEX® DOPPLER (WITHOUT DISPLAY)

Dedicated fetal Doppler with fixed 2 MHz transducer and integral loudspeaker. Waterproof. Includes carry case, headphones, alkaline battery (9V) and gel.


D920-P-USA	Audio Dopplex® w/o display
------------	----------------------------



HUNTLEIGH TRANSDUCERS

For FD2, MD2, SD2, and D900.

- EZ8-HS EZ8 Widebeam Probe
- OP2-HS 2 MHz Probe For Fetal Heart Detection, High Sensitivity
- OP3-HS 3 MHz Probe For Fetal Heart Detection, High Sensitivity
- VP10-HS 10 MHz Probe For Very Superficial Vascular Applications, High Sensitivity
- VP4-HS 4 MHz Probe For Very Deep Vessels, High Sensitivity
- VP5-HS 5 MHz Probe Moderately Deep Vessels, High Sensitivity
- VP8-HS 8 MHz Transducer For Peripheral/Superficial Vessels, High Sensitivity



HUNTLEIGH BABY DOPPLEX® DOPPLER 3000 MARK 2

The advanced electronic design enables accurate data to be produced using autocorrelation signal processing techniques, and displayed and recorded on the thermal paper trace. The pulsed Doppler, 7 Crystal Array Ultrasound Transducer provides excellent fetal heart pick-up and wrap-around bumpers provide additional impact protection. The standard monitor also includes: ultrasound channel, contractions measurement (Toco), automatic fetal movement detection with adjustable sensitivity, patient event marker, clinical event marker and annotation facility, clear LED indication of FHR and UA, high quality audio FHR reproduction, signal quality indicator (dual for twins), integral thermal printer.

BD3000MK2 Baby Doppler® 3000 Mark 2

HUNTLEIGH TWIN BABY DOPPLEX® DOPPLER 3002

Having all of the features of the 3000 Mark 2 version, the 3002 has a unique and extremely useful feature. Ultrasound FHR traces are produced on "Double Width" paper allowing two completely separate traces to be viewed. This avoids possible confusion which can occur, if conventional twins monitors which superimpose one trace on top of another, are used. A carry case is available to transport the products from clinic to clinic.

BD3002 Baby Doppler® 3002

HUNTLEIGH BABY DOPPLEX® DOPPLER 3000 MARK 2 ACCESSORIES

- ACC-15 Printer Paper, Plain, No Grid, 10/bx
- ACC-30 Carry Case/Shoulder Bag
- ACC-31 Wall Bracket
- ACC-32 Desk Top Tray
- ACC-53 Trolley Cart



HUNTLEIGH FLEXI DOPPLEX® DOPPLER

Dedicated fetal Doppler with fixed 3 MHz transducer and integral loudspeaker. Waterproof. Includes carry case, headphones, alkaline battery (9V) and gel.

D930-P-USA Flexi Doppler®



HUNTLEIGH MINI DOPPLEX® DOPPLER (WITHOUT DISPLAY)

Non-directional pocket Doppler with integral loudspeaker. Supplied with one transducer of choice, carry case, headphones, alkaline battery (9V), gel and operating instructions.

- D900USA/EZ8 Mini Doppler w/EZ8 Widebeam Probe
- D900USA/OP2HS Mini Doppler w/2 MHz Probe
- D900USA/OP3HS Mini Doppler w/3 MHz Probe
- D900USA/VP10HS Mini Doppler w/10 MHz Probe
- D900USA/VP4HS Mini Doppler w/4 MHz Probe
- D900USA/VP5HS Mini Doppler w/5 MHz Probe
- D900USA/VP8HS Mini Doppler w/8 MHz Probe



HUNTLEIGH MULTI DOPPLEX® II DOPPLER

Sophisticated pocket Doppler with bi-directional indication of blood flow velocity, 4 level gain control, LCD display, outputs and true separated audio signals, active noise reduction, digital interface to Doppler® Printa/EKG/or Doppler® Reporter software package for printout of waveforms (RS232). Supplied with one transducer of choice.

- MD2PUSA/EZ8 Multi Doppler II Pocket w/EZ8 Widebeam Probe
- MD2PUSA/OP2HS Multi Doppler II Pocket w/2 MHz Probe
- MD2PUSA/OP3HS Multi Doppler II Pocket w/3 MHz Probe
- MD2PUSA/VP10HS Multi Doppler II Pocket w/10 MHz Probe
- MD2PUSA/VP4HS Multi Doppler II Pocket w/4 MHz Probe
- MD2PUSA/VP5HS Multi Doppler II Pocket w/5 MHz Probe
- MD2PUSA/VP8HS Multi Doppler II Pocket w/8 MHz Probe



HUNTLEIGH DOPPLER ACCESSORIES

- ACC-14 Thermal Dopplex® Printa™ Paper, Plain, No Grid
- ACC-67 Ultrasonic Gel, 12/bx
- DP100 Dopplex® Stand, Locking Device To Help Prevent The Doppler From Being Mishandled, Dropped or Stolen



PARKER LABS AQUASONIC® 100 ULTRASOUND TRANSMISSION GEL

The most widely used gel for diagnostic and therapeutic medical ultrasound. Recommended for all procedures where a viscous gel is required. Completely aqueous, will not stain clothing or damage transducers. Hypoallergenic, bacteriostatic, non-sensitizing and non-irritating. No formaldehyde. Not a spermicide.

- 01-02 Ultrasound Transmission Gel, Non-Sterile, 60g Doppler® Size Tube, 12/bx, 12 dz/cs
- 01-08 Ultrasound Transmission Gel, Non-Sterile, 0.25 Liter Bottle, 6 dz/cs
- 01-34 Ultrasound Transmission Gel, Non-Sterile, 1 Liter dispenser bottle, 6/bx, 12/cs
- 01-50 Ultrasound Transmission Gel, Non-Sterile, 5 Liter Sonicpac® w/dispenser bottle, 4/cs



PARKER LABS AQUASONIC® CLEAR ULTRASOUND GEL

Based on the same "can't be copied" formula as Aquasonic 100 Ultrasound Transmission Gel, without fragrance or color. Acoustically correct for a broad range of frequencies. Hypoallergenic, non-irritating and bacteriostatic. Water soluble. No formaldehyde. Not a spermicide.

- 03-02 Ultrasound Gel, 60g Doppler® Size Tube, 12/bx, 12 dz/cs
- 03-08 Ultrasound Gel, .25 Liter Bottle, 6 dz/cs
- 03-34 Ultrasound Gel, 1 Liter w/dispenser cap, 6/bx, 2 bx/cs
- 03-50 Ultrasound Gel, 5 Liter SONICPAC® w/disposable bottle, 4/cs
- 03-54 Ultrasound Gel, 5 Liter ECONOPAC w/2 dispensers, 4/cs



PARKER LABS STERILE AQUASONIC® 100 ULTRASOUND TRANSMISSION GEL

Recommended for diagnostic and therapeutic ultrasound applications when sterility is indicated. Overwrapped, sterilized foil pouches. Acoustically correct for the broad range of frequencies used in medical ultrasound. Non-irritating, hypoallergenic. Completely aqueous, will not stain clothing or damage equipment. Not a spermicide.

- 01-01 Ultrasound Transmission Gel, Sterile, 20g Overwrapped Sterilized Foil Pouches, 48/bx, 6 bx/cs



PARKER LABS THERMASONIC® GEL WARMER

Warms three bottles of ultrasound gel to body temperature. Durable wall mountable construction with warranty. Durable case. Low electrical demand. UL and CSA listed.

- 82-04 Three-Bottle Gel Warmer (plus 3 dispensers of ultrasound gel), 4/cs

PARKER LABS TRANSEPTIC® CLEANSING SOLUTION

Multi-purpose spray for professional use in cleansing ultrasound transducer/probe surfaces between procedures. Will not cause transducer membrane swelling. Does not contain phenol or glutaraldehyde.

- 09-25 Cleansing Solution, 250ml (8.5 fl oz) Clear Spray Bottle, 12/bx, 4 bx/cs



CHESTER LIQUASONIC® ULTRASOUND GEL

LiquaSonic Ultrasound Transmission Gel has a thicker, no-slip formula that remains viscous even when applied to salty skin. Our salt-free, formaldehyde free formula eliminates the air between the transducer and skin, offering exceptional coupling efficiency for more effective treatment/diagnosis. LiquaSonic Gel has been formulated to be acoustically correct for the broad range of frequencies used in medical ultrasound applications. Our formula is water soluble and non-greasy preventing damage to expensive transducers and clothing.


- 001205 8 oz bottle w/mustard cap, Ultrasound Gel, 12/cs
- 001222 5 Liter Collapsible Cubitainer w/empty bottle, Ultrasound Gel, 4/cs
- 001297 128 oz Multi-Purpose, Medium Viscosity Ultrasound Gel, 4/cs
- 001298 8 oz Bottle w/mustard cap, Multi-Purpose, Medium Viscosity Ultrasound Gel, 12/cs
- 001305 8 oz Bottle w/mustard cap, Clear Ultrasound Gel, 12/cs
- 001333 128 oz Clear Ultrasound Gel, 4/cs



CHESTER LIQUASONIC® ULTRASOUND LOTION

Similar in use to the ultrasound gel, the lotion is often preferred by medical professionals. The aqueous formulation is pleasing to the patient and is acoustically correct for the large spectrum of frequencies used in medical ultrasound applications. The lotion is hypoallergenic, bacteriostatic, salt free, formaldehyde free and will not harm clothing or equipment.

- 001267 128 oz Ultrasonic Lotion, 4/cs
- 001268 8 oz Bottle w/mustard cap, Ultrasonic Lotion, 12/cs



STARLINE PHYSICIAN'S PENLIGHT

Penlight has high-intensity prefocused bulb powered by two AAA batteries. Aluminum contact strip assures good contact for brighter light. Extruded case has a unique V-clip. 1-year shelf life. Useful life is 3- to 6-months with intermittent use.

- 6666S Dr. Pack Disposable Penlight, 6/pk



AARON DIAGNOSTIC LIGHT EXAM KITS

Disposable. Two AAA batteries. Cobalt blue with Wratten #47A cobalt filter. For use with fluorescein dye. One-year shelf life.


- 6303 Penlights w/cobalt filter, 3/pk



AARON FLUORO-DOT™ PENLIGHT

Battery-operated penlight designed for use with fluorescein dye. Excellent for diagnosing abrasions and fitting contact lenses. Constant on/off switch. Contains two AA alkaline batteries. Penlight has a 4-year shelf life and a 4-hour burn time.

- 6803 Fluoro Slit™ Penlight
- 6802 Fluoro Dot™ Penlight w/wratten cobalt blue filter



AARON PENLIGHTS

- 1125 Disposable Penlight
- 1135 Reusable Penlight

Accessories

- 222 Replacement Bulb For 1135
- 7156 Cobalt Filter Only



AARON PHYSICIAN'S PENLIGHT

Penlight has high-intensity prefocused bulb powered by two AAA batteries. Aluminum contact strip assures good contact for brighter light. Extruded case has a unique V-clip. 1-year shelf life. Useful life is 3- to 6-months with intermittent use.

- 6666 Dr. Pack Disposable Penlight, 6/pk
- 666R Doctors' Replaceable Battery Penlight, 6/pk
- 66RN RN Pack Disposable Penlight w/pupil gauge, 6/pk



AARON SPECIALTY LIGHTS

Individually packaged. Two AA batteries. All lights have a 4-year shelf life.

- 6800 Tech Light™ w/concentric beam
- 6801 White Slit Beam For Cataract Screening
- 6804 1mm Light (for muscle tests)



AARON WOODS LIGHT

- UV59 Hand-Held Ultra-Violet Light

Accessories

- UV59B Replaceable Black Lamp For UV59



ADC® ADLITE™ DISPOSABLE PENLIGHT

Sealed concentric beam activated by depressing pocket clip. Available with or without pupil gauge. 5" long, 1/2" diameter.

- 351 Disposable Penlight, 6/pk
- 351P Disposable Penlight w/pupil gauge, 6/pk



ADC® ADLITE II™ PENLIGHT

Reusable penlight with attractive tapered barrel. Activated by depressing pocket clip. Two AAA batteries included. Measures 5" long, 1/2" diameter.

- 354 Reuseable Penlight



ADC® ADLITE™ PLUS PENLIGHT

Sealed concentric beam, disposable diagnostic penlight. Activated by depressing pocket clip. Measures 5" long, 1/2" diameter. Sold individually or in display cylinder of 22.

- 356BK ADLITE™ Disposable Penlight, Black
- 356OR ADLITE™ Disposable Penlight, Orange
- 356R ADLITE™ Disposable Penlight, Red
- 356TL ADLITE™ Disposable Penlight, Teal
- 356V ADLITE™ Disposable Penlight, Purple
- 356W ADLITE™ Disposable Penlight, White
- 356WP ADLITE™ Disposable Penlight, White w/pupil gauge
- 358AST Disposable Penlight, Assorted Colors, 22/cylinder
- 358TVB Disposable Penlight, Black, Purple & Teal, 22/cylinder
- 358W Disposable Penlight, White, 22/cylinder
- 358WP Disposable Penlight, White w/pupil gauge, 22/cylinder



ADC® METALITE™ PENLIGHT

Reusable with brushed aluminum housing. Activated by depressing pocket clip. Two AAA batteries included. Measures 6" long, 1/2" diameter.

- 352 Penlight

ADC® POCKET PAL III™ ORGANIZER AND KIT

An indispensable uniform accessory that organizes while it protects. The organizer features washable 420D nylon pocket sewn with nylon thread. Five slide-in pockets and one covered pocket. Key ring. Two matched closures secure instruments and double as ID tags. Sewn loop permits belt wearing. Measures 6 1/4" x 4 3/4" x 1/4". Weighs 1 oz. The kit includes: Pocket Pal™ Organizer (216), 5 1/2" Lister Bandage scissors, 5" straight Mosquito forceps, 3 color pen and ADLITE™ disposable penlight. All components are latex free.

- 117BK Pocket Pal III™ Kit, Black
- 117TL Pocket Pal III™ Kit, Teal
- 117V Pocket Pal III™ Kit, Purple

ADC® TOGGELITE™ PENLIGHT

Constructed of chrome plated brass for rugged durability. Activated by a locking toggle switch. Two AA batteries included. Measures 6" long, 3/4" diameter.

- 353 Penlight

OMRON HIGH INTENSITY DIAGNOSTIC PENLIGHTS

Focused beam lens-end lamps. Each light has a clip to prevent accidental discharge of the batteries when carried in a pocket. Lights weigh less than one ounce.

- 24-100 Reusable Pen Light, 6/pk
- 24-150 Reusable Pen Light, 100/bx



WELCH ALLYN PROFESSIONAL PENLITE

Halogen light for true tissue color and consistent, long lasting illumination. Brass with black matte finish for durability.

- 76600 PenLite



21504



24031



25020

WELCH ALLYN 2.5V/3.5V HALOGEN DIAGNOSTIC OTOSCOPE

Halogen illumination with fiber optic light transmission provides brilliant, cool, distal illumination without visual obstruction or specular reflections. Wide angle viewing lens slides left or right for instrumentation under magnification.

- 24031 2.5V Halogen Otoscope w/4 polypropylene specula (sizes 2, 3, 4 & 5 mm)
- 24032 2.5V Halogen Otoscope w/insufflator bulb, tube w/tip & 4 polypropylene specula (sizes 2, 3, 4 & 5 mm)
- 24035 2.5V Halogen Otoscope w/o specula
- 25020 3.5V Halogen Otoscope Head w/4 polypropylene specula (sizes 2, 3, 4, 5 mm)
- 25021 Same as No. 25020, Including Insufflator Bulb & Tube

Accessories

- 200055-502 Lens Holder Assembly
- 21504 Insufflator Bulb & Tube w/tip



21600



21700

WELCH ALLYN 2.5V/3.5V HALOGEN OPERATING OTOSCOPE

Open system for procedures with fully rotatable 5/16" (7.9 mm) diameter magnifying lens and rotatable speculum. Halogen illumination with fiber optic light transmission.

- 21600 2.5V Otoscope w/5 polypropylene specula (sizes 2, 3, 4, 5 & 9mm)
- 21601 2.5V Otoscope w/o specula
- 21700 3.5V Otoscope w/5 polypropylene specula (sizes 2, 3, 4, 5, 9 mm)
- 21701 Same as No. 21700, w/o specula

Accessories

- 209026-501 Lens and Attachable Arm




20100



20200



20201



21501

WELCH ALLYN 2.5V/3.5V HALOGEN PNEUMATIC OTOSCOPE

Closed system, with insufflator fitting for pneumatic otoscopy. Large magnifying lens can be easily rotated to facilitate use of swabs, curettes, etc. Halogen illumination with fiber optic light transmission.

- 20100 2.5V Otoscope w/5 polypropylene specula (sizes 2, 3, 4, 5 & 9 mm)
- 20101 2.5V Otoscope w/o specula
- 20200 3.5V Otoscope w/5 polypropylene specula (sizes 2, 3, 4, 5, 9 mm)
- 20201 Same as No. 20200, w/o specula

Accessories

- 21501 Insufflator Bulb & Tube w/o tip

WELCH ALLYN 3.5V HALOGEN PNEUMATIC OTOSCOPE/THROAT ILLUMINATOR

Dual-purpose instrument allows practitioner to examine the ears, then detach otoscope head and have an instant, high-intensity throat illuminator. Wide angle viewing lens slides left or right. Halogen illumination with fiber optic light transmission.

- 25257 Illuminator Section Only




WELCH ALLYN MACROVIEW™ OTOSCOPE

30% higher magnification over standard otoscopes. Nearly a full view of the tympanic membrane. Ability to focus to fine tune the view. Greater tip retention and easier tip removal.

- 23810 3.5V Otoscope w/specula (head only)
- 23814 3.5V Otoscope (23810) w/insufflation bulb
- 23820 3.5V Otoscope w/throat illuminator
- 23824 3.5V Otoscope (23820) w/insufflation bulb

WELCH ALLYN MACROVIEW™ OTOSCOPE ACCESSORIES

- 05215-M Soft Case For 3.5V Set
- 05258-M Hard Case For PanOptic Ophthalmoscope 3.5V Set
- 05259-M Hard Case For 3.5V Set
- 05815-M Soft Case For PanOptic Ophthalmoscope 3.5V Set
- 06500-U 3.5V Halogen Lamp For Otoscope
- 23804 Insufflation Bulb For Otoscope
- 23857 3.5V Throat Illuminator Section Only For Otoscope
- 26038 Bivalve Nasal Speculum For MacroView™ Complete
- 26538 Nasal Illuminator For MacroView™ Complete



22820



22821

WELCH ALLYN 2.5V POCKETSCOPE™ OTOSCOPE/THROAT ILLUMINATOR

Halogen light for true tissue color and consistent, long lasting illumination. Fiber optics for cool light with no reflection, no obstruction. Wide angle viewing lens. Sealed system for pneumatic otoscopy. Built-in throat illuminator/penlight for added versatility.

- 21111 PocketScope Otoscope/Throat Illuminator w/o handle
- 22800 PocketScope Otoscope/Throat Illuminator w/handle & rechargeable battery
- 22811 PocketScope Otoscope/Throat Illuminator w/handle, rechargeable battery & soft case
- 22820 PocketScope Otoscope/Throat Illuminator w/AA alkaline battery handle
- 22821 PocketScope Otoscope/Throat Illuminator w/AA alkaline battery handle & soft case

Accessories

- 211005 Replacement Lens

ADC® STANDARD & VETERINARY OTOSCOPE

German 2.7V lamp for bright illumination. Bayonet locking otoscope head. Removable 4X magnifying lens. Rheostatic on/off switch. Autoclavable, polypropylene ear specula in sizes 2.5mm, 3.5mm and 4.5mm (one of each). Veterinary model includes 4mm, 5mm, and 7mm specula sizes. Knurled brass handle with a beautiful satin finish. Ergonomically designed for comfort. Compact instrument measures just under 7 1/2" and weighs approximately 7 oz. Powered by 2 C cell batteries (not included). Comes with zippered leatherette storage pouch that holds the otoscope and included speculas. Inspected and packaged in the U.S.A. from foreign components. One year quality performance guarantee.

- 5211 Standard Otoscope
- 5211V Veterinary Otoscope

ADC® STANDARD & VETERINARY OTOSCOPE ACCESSORIES

- 5211-1 Standard Otoscope Specula, 2.5mm
- 5211-2 Standard Otoscope Specula, 3.5mm
- 5211-3 Standard Otoscope Specula, 4.5mm
- 5211-4 Standard Otoscope Replacement Lamp, 2.7V
- 5211-7 Veterinary Otoscope Specula, 4.0mm
- 5211-8 Veterinary Otoscope Specula, 5.0mm
- 5211-9 Veterinary Otoscope Specula, 7.0mm
- 5211-AD Specula Adapter For 5211

MILTEX BRUENING MAGNIFIER OTOSCOPE

Chrome. 3.5x magnification with 3 specula, small, medium and large sizes. Complete with inflation bulb and tubing in a fitted case.

- 19-70 Otoscope Set
- 19-71 Set of 3 Specula: Small, Medium & Large Sizes For Otoscope
- 19-71-1 Small Speculum
- 19-71-2 Medium Speculum
- 19-71-3 Large Speculum

MILTEX SIEGEL MAGNIFIER OTOSCOPE

Chrome. 3.15x magnification with 3 specula: small, medium and large sizes. Complete with inflation bulb and tubing.

- 19-75 Rubber Bulb w/tubing For Otoscope



22002



24305

WELCH ALLYN POLYPROPYLENE REUSABLE SPECULA

Material is lightweight and durable. Can be autoclaved, boiled or cleaned with standard germicides.

- 22002 2 mm Speculum, For Use With Pneumatic, Operating & Consulting Otoscopes, Dark Green
- 22003 3 mm Speculum, For Use With Pneumatic, Operating & Consulting Otoscopes, Dark Green
- 22004 4 mm Speculum, For Use With Pneumatic, Operating & Consulting Otoscopes, Dark Green
- 22005 5 mm Speculum, For Use With Pneumatic, Operating & Consulting Otoscopes, Dark Green
- 22009 9 mm Speculum (nasal), For Use With Pneumatic, Operating & Consulting Otoscopes, Dark Green
- 22100 Set of 5 Specula (2, 3, 4, 5, 9mm), For Use With Pneumatic, Operating & Consulting Otoscopes, Dark Green
- 24302-U 2 mm Speculum, Universal Poly Reusable
- 24303-U 3 mm Speculum, Universal Poly Reusable
- 24304-U 4 mm Speculum, Universal Poly Reusable
- 24305 5 mm Speculum, For Use With Diagnostic Otoscopes, Black
- 24305-U 5 mm Speculum, Universal Poly Reusable
- 24400-U Set of 4 Speculum, Universal Poly Reusable



WELCH ALLYN KLEENSPEC® DISPOSABLE OTOSCOPE SPECULA

Make of nontoxic plastic; to be used once then discarded. No adapter required. Black.

- 52133 3 mm Specula, For Use With Pneumatic, Operating & Consulting Otoscopes, 500/bg, 10 bg/cs
- 52134 4 mm Specula, For Use With Pneumatic, Operating & Consulting Otoscopes, 500/bg, 10 bg/cs
- 52135 5 mm Specula, For Use With Pneumatic, Operating & Consulting Otoscopes, 500/bg, 10 bg/cs
- 52432-U 2.75mm Universal KleenSpec Disposable Otoscope Specula, 850/bg, 10 bg/cs
- 52434-U 4.25mm Universal KleenSpec Disposable Otoscope Specula, 850/bg, 10 bg/cs



WELCH ALLYN SOFSPEC™ OTOSCOPE SPECULA

Polypropylene base with soft, flexible tip provides improved seal for pneumatic otoscopy. Can be autoclaved, boiled or cleaned with standard germicides.

- 22023 3 mm Speculum, For Use With Pneumatic, Operating and Consulting Otoscopes
- 22025 5 mm Speculum, For Use With Pneumatic, Operating and Consulting Otoscopes
- 22027 7 mm Speculum, For Use With Pneumatic, Operating and Consulting Otoscopes
- 22120 Set of Above 3 Specula, For Use With Pneumatic, Operating and Consulting Otoscopes
- 24323 3 mm Speculum, For Use With Diagnostic Otoscopes
- 24325 5 mm Speculum, For Use With Diagnostic Otoscopes
- 24327 7 mm Speculum, For Use With Diagnostic Otoscopes
- 24420 Set of Above 3 Specula, For Use With Diagnostic Otoscopes

SPECLINC OTREC® DISPOSABLE OTOSCOPE SPECULA

Complete speculum is disposable, providing a clean tip for each use. No covers or adapters are needed. Black color for less reflection and greater clarity.

- 7225 2.5 mm Otoscope Specula, 500/pk, 10 pk/cs (5000/cs)
- 7250 Pediatric, 2.5 mm Otoscope Specula, (cross reference #52432), 1000/pk, 10 pk/cs (10,000/cs)
- 7300 Youth, 3 mm Otoscope Specula, 1000/pk, 10 pk/cs (10,000/cs)
- 7333 3 mm Otoscope Specula (cross reference #52133) , 500/pk, 10 pk/cs (5000/cs)
- 7400 Med. Adult, 4 mm Otoscope Specula (cross reference #52434), 1000/pk, 10 pk/cs (10,000/cs)
- 7444 4 mm Otoscope Specula (cross reference #51234), 500/pk, 10 pk/cs (5000/cs)
- 7500 Lg. Adult, 5 mm Otoscope Specula, 1000/pk, 10 pk/cs (10,000/cs)
- 7555 5 mm Otoscope Specula (cross reference #52135), 500/pk, 10 pk/cs (5000/cs)
- 7666 6 mm Otoscope Specula, 500/pk, 10 pk/cs (5000/cs)
- 7999 9 mm Otoscope Specula, 500/pk, 10 pk/cs (5000/cs)

SPECLINC OTREC DISPOSABLE OTOSCOPE SPECULA WALL DISPENSERS

- 7000 Wall Dispenser for Diagnostic: 7250, 7300, 7400, 7500 5/pk (includes Wall Bracket), 10pk/cs
- 7777 Wall Dispenser for Operating and Consulting: 7225, 7333, 7555, 7666 and 7999, 5/pk (includes Wall Bracket), 10pk/cs



WELCH ALLYN 2.5V/3.5V HALOGEN OPHTHALMOSCOPE

Features high-intensity halogen illumination and mirror optical system. 28 lenses with illuminated lens dial. Apertures (microspot, small spot, large spot, fixation, red-free, and slit) for all diagnostic applications.

- 11470 2.5V Halogen Ophthalmoscope, Head Only
- 11710 3.5V Halogen Ophthalmoscope, Head Only



WELCH ALLYN HALOGEN COAXIAL AUTOSTEP OPHTHALMOSCOPE

Provides 68 viewing lenses with continuous single diopter steps or automatic progression from -30 to +38. Wide range of diopter steps facilitates precise focusing. Lens values fully displayed, no mathematical calculations required. Advanced halogen coaxial optics for high intensity illumination.

- 11730 3.5-V Halogen Coaxial Autostep Ophthalmoscope



WELCH ALLYN HALOGEN COAXIAL OPHTHALMOSCOPE

23 lenses. Illuminated lens dial. Sealed optics provide better dust protection. Coaxial optical system eliminates shadows and facilitates exam through virtually any size pupil or vitreous disorder. Special contoured mirror and matching lamp filament allow axis of illumination to be practically coincidental with axis of physician's vision. Apertures include microspot, small spot, large spot, fixation, cobalt blue, and slit.

- 11720 3.5V Halogen Coaxial Ophthalmoscope, Head Only



WELCH ALLYN PANOPTIC™ OPHTHALMOSCOPE

Patented new Axial PointSource™ Optics. Illumination system converges the light to a point at the cornea, allowing easy entry into small, undilated pupils. The illumination pathway then diverges to the retina, illuminating a very wide area of the fundus. Viewing system enables the operator to view the illuminated area on virtually the same axis, creating the widest field-of-view attainable. Provides a five-times larger view of the fundus than ever before achieved with a standard ophthalmoscope in an undilated eye. Increases magnification by 26% over a standard ophthalmoscope. Patented Glare Extinguishment System prevents interference from unwanted glare and reflections. Works with all existing Welch Allyn 3.5v power sources (battery handle or wall transformer) saving additional investment. The patient eyecup establishes and maintains proper viewing distance, provides stabilization for the view, is a pivoting point for leverage in panning around the retina, occludes ambient light, and helps with orientation. The dynamic focusing wheel adjusts the focus in a continuous, smooth action for more precise control and the optimum view with a focusing range of -20 to +20 diopters. The aperture dial has micro, small, and large spot sizes, a slit aperture, and red-free filter, and cobalt blue filter with add-on corneal magnifying lens (Model 11820 only).

- 11810 PanOptic Ophthalmoscope w/o Blue Filter
- 11820 PanOptic Ophthalmoscope w/Blue Filter & Lenses

Accessories

- 03800-U PanOptic Lamp
- 05258-MBX PanOptic Diagnostic Case
- 05815-MBX PanOptic Diagnostic Soft Case
- 11870 Patient Eye Cup, 5/set
- 11875 Corneal Lens



WELCH ALLYN 2.5V HALOGEN STREAK RETINOSCOPIES

Continuously rotating external focusing sleeve is easy to grip and manipulate regardless of hand size. Magnetic fixation cards. Rubber brow rest prevents scratching of eyeglasses. Halogen light for bright, consistent, long-lasting illumination.

17200 2.5V Halogen Streak Retinoscope



18320

WELCH ALLYN 3.5V HALOGEN RETINOSCOPE SETS

- 18310 18200 Retinoscope, 11710 Ophthalmoscope, 71000 Handle, 05820 Hard Case
- 18311 18200 Retinoscope, 11710 Ophthalmoscope, 71000 Handle, 05180 Soft Case
- 18320 18200 Retinoscope, 11720 Ophthalmoscope, 71000 Handle, 05820 Hard Case
- 18321 18200 Retinoscope, 11720 Ophthalmoscope, 71000 Handle, 05180 Soft Case
- 18330 18200 Retinoscope, 11730 Ophthalmoscope, 71000 Handle, 05820 Hard Case
- 18331 18200 Retinoscope, 11730 Ophthalmoscope, 71000 Handle, 05180 Soft Case
- 18340 18200 Retinoscope, 71000 Handle, 05820 Hard Case
- 18341 18200 Retinoscope, 71000 Handle, 05180 Soft Case
- 18410 18300 Retinoscope, 11710 Ophthalmoscope, 71000 Handle, 05820 Hard Case
- 18411 18300 Retinoscope, 11710 Ophthalmoscope, 71000 Handle, 05180 Soft Case
- 18420 18300 Retinoscope, 11720 Ophthalmoscope, 71000 Handle, 05820 Hard Case
- 18421 18300 Retinoscope, 11720 Ophthalmoscope, 71000 Handle, 05180 Soft Case
- 18430 18300 Retinoscope, 11730 Ophthalmoscope, 71000 Handle, 05820 Hard Case
- 18440 18300 Retinoscope, 71000 Handle, 05820 Hard Case
- 18441 18300 Retinoscope, 71000 Handle, 05180 Soft Case



18200



18250



18300

WELCH ALLYN 3.5V RETINOSCOPIES

External focusing sleeve. Brighter halogen light. Crossed linear polarizing filter. Magnetic, age appropriate targets. One-hand operation for streak focus and 360° streak rotation. Rubber brow rest. Interchangeable streak retinoscope can be converted to a spot retinoscope by simple changing the lamp.

- 18200 3.5V Streak Retinoscope
- 18300 3.5V Spot Retinoscope

Accessories

- 18250 Complete Set of Fixation Cards



WELCH ALLYN 2.5V POCKETSCOPE™ OPHTHALMOSCOPE

Halogen light for true tissue color and consistent, long lasting illumination. Five aperture selections. 48 lenses for better resolution. Polarized filter to reduce corneal reflection. Compact size fits in pocket.

- 12800 PocketScope Ophthalmoscope w/handle & rechargeable battery
- 12811 PocketScope Ophthalmoscope w/handle, rechargeable battery & soft case
- 12820 PocketScope Ophthalmoscope w/AA alkaline battery handle
- 12821 PocketScope Ophthalmoscope w/AA alkaline battery handle & soft case
- 13000 PocketScope Ophthalmoscope w/o handle

Accessories

- 05280-1 Compact Zipper Case
- 05928-U Soft Case



WELCH ALLYN STRABISMOSCOPE™ UNIDIRECTIONAL OCCLUDER

Test for both heterophorias (latent Strabismus) and heterotropias (manifest Strabismus) at the same time. One-way mirror allows direct observation of occluded eye. View movement of both covered and uncovered eyes simultaneously. Allows for confident diagnosis of Dissociated Vertical Deviation. Rugged ABS housing withstands drops and falls. Lightweight and easy to clean. Fits all Welch Allyn 3.5V power sources. Handles must be ordered separately.

- 12400 StrabismoScope Head Only
- 12451 3.5-V Halogen StrabismoScope Set Includes StrabismoScope, Rechargeable Handle & Soft Carrying Case

ADC® COMPLETE 2.5V INSTRUMENT SET

Five different instruments in a fitted case. Features: bayonet locking mounts for each instrument head. C cell powered, satin finish chrome-plated battery handle with knurled finish. Rheostatic on/off switch. Complete in zippered case with fitted foam insert that is scuff and scratch resistant and cleans easily

- 5215 Complete Diagnostic Set

ADC® COMPLETE 2.5V INSTRUMENT SET ACCESSORIES

- 5215-M3 #3 Laryngeal Mirror
- 5215-M4 #4 Laryngeal Mirror
- 5215-NS Nasal Speculum
- 5215-TD Metal Tongue Depressor
- 5215-TH Illuminated Tongue Blade Holder
- 5215-TO Tongue Depressor Attachment For Disposable
- 5215-ZC Black Zipper Case
- 5220 Otoscope Head Only
- 5240 Ophthalmoscope Head Only
- 5260 Battery Handle
- 5312 Dermoscope



ADC® POCKET OTOSCOPE/OPHTHALMOSCOPE SET

Ophthalmoscope features five aperture selections for general and specialist use; 17 lenses from -20 to +20 diopters; illuminated dial and built-in lens dust cover. Otoscope features fiberoptics for cool, obstruction free light transmission; wide angle, removable lens and disposable polypropylene ear specula. Both units feature a 2.5V halogen light source, interchangeable metal handles and a lifetime quality performance guarantee on optics. Includes fitted case. Latex free system.

- 5110 Pocket Otoscope/Ophthalmoscope Set
- 5111 Pocket Otoscope Set
- 5112 Pocket Ophthalmoscope Set

ADC® POCKET OTOSCOPE/OPHTHALMOSCOPE SET ACCESSORIES

- 5120 Otoscope Head Only For 5111
- 5121 Insufflator For Otoscope
- 5122 Insufflator Bulb For Otoscope
- 5140 Ophthalmoscope Head Only For 5112
- 5160 Battery Handle For Pocket Instruments
- 925110-00 Hard Case Only For 5110
- 925111-00 Hard Case Only For 5111
- 925112-00 Hard Case Only For 5112
- 5182 2.75mm Disposable Specula, 850/pkg
- 5185 4.25mm Disposable Specula, 850/pkg

ADC® STANDARD 2.5V DIAGNOSTIC SET

Ophthalmoscope head with lens wheel features 25 corrective lenses from -20 to +20 diopters. Otoscope head with removable 4x magnifying lens and insufflator fitting. Bayonet head locking mounts for secure connection to the battery handle. C cell powered battery handle with rheostatic on/off switch. Auto-clavable, polypropylene ear specula in 2.5mm, 3.5mm, and 4.5mm sizes. Complete in zippered case with fitted foam insert that is scuff and scratch resistant and cleans easily. Batteries not included.

- 5210 Standard Diagnostic Set



72600



92821

WELCH ALLYN 2.5V HALOGEN POCKETSCOPE DIAGNOSTIC SETS

Complete PocketScope Set with ophthalmoscope, otoscope/throat illuminator, two rechargeable battery handles and desk charger (110-130V, 50/60 Hz) in hard case.

- 92820 PocketScope Set Includes Ophthalmoscope, Otoscope/Throat Illuminator, 2 "AA" Alkaline Battery Handles, Hard Case (w/o charger)
- 92821 Same as No. 92820 in Soft Case

Accessories

- 72600 2.5V Rechargeable Battery
- 72801 PocketScope Set Handle w/o batteries
- 72800 PocketScope Set Handle w/rechargeable battery

WELCH ALLYN 3.5V COAXIAL MACROVIEW OTOSCOPE/OPHTHALMOSCOPE SETS

- 97200-M 23820 Otoscope, 11710 Ophthalmoscope, 71000 Handle and Hard Case
- 97200-MC 11720 Ophthalmoscope, 23820 Otoscope w/hard case and Convert. Handle



05215-MBX



97120



97220

WELCH ALLYN 3.5V HALOGEN DIAGNOSTIC SETS

Portable 3.5-V halogen diagnostic sets feature ophthalmoscopes and otoscopes. Optional nasal illuminator. With 71000 handle.

- 97120 Standard Set, 11710 Ophthalmoscope, 20200 Pneumatic Otoscope w/hard case
- 97121 Standard Set, 11710 Ophthalmoscope, 20200 Otoscope w/soft case
- 97170 Standard Set, 11710 Ophthalmoscope, 21700 Otoscope w/hard case
- 97171 Standard Set, 11710 Ophthalmoscope, 21700 Otoscope w/soft case
- 97220 Coaxial Set, 11720 Ophthalmoscope, 20200 Otoscope w/hard case
- 97221 Coaxial Set, 11720 Ophthalmoscope, 20200 Otoscope w/soft case
- 97270 Coaxial Set, 11720 Ophthalmoscope, 21700 Otoscope w/hard case
- 97271 Coaxial Set, 11720 Ophthalmoscope, 21700 Otoscope w/soft case
- 97320 AutoStep Set, 11730 Ophthalmoscope, 20200 Otoscope w/hard case
- 97321 AutoStep Set, 11730 Ophthalmoscope, 20200 Otoscope w/soft case
- 97370 AutoStep Set, 11730 Ophthalmoscope, 21700 Otoscope w/hard case
- 97371 AutoStep Set, 11730 Ophthalmoscope, 21700 Otoscope w/soft case

Accessories

- 05215-MBX Soft Case



WELCH ALLYN 3.5V INSTRUMENT SETS

- 93100 11710 Ophthalmoscope, 20000 Otoscope, 28100 Tongue Blade Holder, 27000 Larynx Illuminator, 43300 Transilluminator, 71000 Handle, 05510 Case
- 93120 11710 Ophthalmoscope, 20200 Otoscope, 28100 Tongue Blade Holder, 27000 Larynx Illuminator, 43300 Transilluminator, 71000 Handle, 05510 Case
- 93150 11710 Ophthalmoscope, 25020 Otoscope, 28100 Tongue Blade Holder, 27000 Larynx Illuminator, 43300 Transilluminator, 71000 Handle, 05520 Case
- 93170 11710 Ophthalmoscope, 21700 Otoscope, 28100 Tongue Blade Holder, 27000 Larynx Illuminator, 43300 Transilluminator, 71000 Handle, 05510 Case
- 93200 11720 Ophthalmoscope, 20000 Otoscope, 28100 Tongue Blade Holder, 27000 Larynx Illuminator, 43300 Transilluminator, 71000 Handle, 05510 Case
- 93220 11720 Ophthalmoscope, 20200 Otoscope, 28100 Tongue Blade Holder, 27000 Larynx Illuminator, 43300 Transilluminator, 71000 Handle, 05510 Case
- 93250 11720 Ophthalmoscope, 25020 Otoscope, 28100 Tongue Blade Holder, 27000 Larynx Illuminator, 43300 Transilluminator, 71000 Handle, 05520 Case
- 93270 11720 Ophthalmoscope, 21700 Otoscope, 28100 Tongue Blade Holder, 27000 Larynx Illuminator, 43300 Transilluminator, 71000 Handle, 05510 Case

Accessories

05510-U Hard Case



WELCH ALLYN 719 SERIES POWER HANDLE DIAGNOSTIC SETS

Our coaxial ophthalmoscope uses a halogen HPX™ lamp that provides 30% more light output for true tissue color and consistent, long-lasting illumination. It has 28 lenses (except the Prestige Set has 68 lenses) from -25 to +40 diopters for better resolution. The coaxial optical system produces a shadow-free spot and easier entry into undilated pupils. The crossed-linear polarizing filter reduces corneal reflection. The red-free filter may be used with any aperture. It has 18 possible aperture combinations for greater versatility. The otoscope uses a halogen HPX™ lamp and fiber optics for cool light with no reflections, no obstruction. It has a wide-angle viewing lens for instrumentation under magnification. The sturdy Cycology® resin handle is half the weight of current 3.5v handles and provides twice the on-time.

- 97310-SM Prestige Plus Set, Otoscope w/Throat Illuminator & Nasal Illuminator, Coaxial-Plus Ophthalmoscope, Lithium Ion Handle, Specula & Hard Case



WELCH ALLYN COMPACSET™ PORTABLE DIAGNOSTIC SET

Integrates both a full-featured ophthalmoscope and fiber-optic otoscope into one durable, easy-to-use instrument. Small enough to fit in any pocket, weighs just 8.5 ounces. Fully protected by a durable ABS plastic housing. Halogen illumination provides brighter, whiter, longer-lasting light. Instrument heads swivel to allow access from any position.

- 92000 2.5V CompactSet w/AA batteries
- 92100 2.5V CompactSet w/rechargeable batteries

Accessories


- 72610 2.5V Rechargeable Battery



WELCH ALLYN 2.5V/3.5V DESK CHARGERS WITH RECHARGEABLE HANDLES

Compact units provide constant, automatic recharging. Can not be overcharged. Heads must be ordered separately.

- 71500 2.5V Rechargeable Handle
- 71670 Rechargeable 3.5V Handle Only



WELCH ALLYN 719 SERIES POWER HANDLE

Half the weight of current 3.5v handles. Average 110 minutes of on-time between charging. Smallest diameter of any 3.5v power handle. Made from Cycology® resin from GE plastics. Low-battery indicator. Lithium ion battery technology.

- 71900 719 Series Power Handle



76520



76710



76720



76730



79520

WELCH ALLYN 767 WALL TRANSFORMERS

Puts everything needed for fast, accurate diagnoses at arm's length. Convenient, modular package. Handle-based rheostats for quick light-intensity adjustments. Comfortable, contoured handles are easy to hold and maneuver. Modular coiled cords are guaranteed for 10 years. High-impact ABS housing is durable and easy to clean. Improved cradle design featuring OptiSense™, a reliable optical sensor that automatically turns handles on when they are removed from the cradles. New Halogen HPX™ Lamp provides 30% more light output. Optional digital clock and pulse timer. KleenSpec® Plus Dispenser with a large, 200-specula capacity and integral storage area. Optional third handle module makes expansion a snap.

- 76510 76710 Wall Transformer & 52401 KleenSpec Plus Dispenser
- 76520 76720 Wall Transformer & 52401 KleenSpec Plus Dispenser
- 76710 3.5V Wall Transformer w/o clock
- 76710-L 76710 Wall Transformer & 76752 SureTemp Thermometer w/back light
- 76710-LA 76710 Wall Transformer, 76752 SureTemp Thermometer w/back light & 767 Wall Aneroid
- 76710-LAK 76710 Wall Transformer, 76752 SureTemp Thermometer w/back light, 767 Wall Aneroid & 52401 KleenSpec Plus Dispenser
- 76710-LAP 76710 Wall Transformer, 76752 SureTemp Thermometer w/back light, 767 Wall Aneroid & 52101 KleenSpec Plus Dispenser
- 76710-LK 76710 Wall Transformer, 76752 SureTemp Thermometer w/back light & 52401 KleenSpec Plus Dispenser
- 76710-LP 76710 Wall Transformer, 76752 SureTemp Thermometer w/back light & 52101 KleenSpec Plus Dispenser
- 76710-T 76710 Wall Transformer & 76751 SureTemp Thermometer
- 76710-TA 76710 Wall Transformer, 76751 SureTemp Thermometer & 767 Wall Aneroid
- 76710-TAK 76710 Wall Transformer, 76751 SureTemp Thermometer, 767 Wall Aneroid & 52401 KleenSpec Plus Dispenser
- 76710-TAP 76710 Wall Transformer, 76751 SureTemp Thermometer, 767 Wall Aneroid & 52101 KleenSpec Plus Dispenser

WELCH ALLYN 767 WALL TRANSFORMERS (Continued)

- 76710-TK 76710 Wall Transformer, 76751 SureTemp Thermometer & 52401 KleenSpec Plus Dispenser
- 76710-TP 76710 Wall Transformer, 76751 SureTemp Thermometer & 52101 KleenSpec Plus Dispenser
- 76720 3.5V Wall Transformer w/clock
- 76720-L 76720 Wall Transformer & 76752 SureTemp Thermometer w/back light
- 76720-LA 76720 Wall Transformer, 76752 SureTemp Thermometer w/back light & 767 Wall Aneroid
- 76720-LAK 76720 Wall Transformer, 76752 SureTemp Thermometer w/back light, 767 Wall Aneroid & 52401 KleenSpec Plus Dispenser
- 76720-LAP 76720 Wall Transformer, 76752 SureTemp Thermometer w/back light, 767 Wall Aneroid & 52101 KleenSpec Plus Dispenser
- 76720-LK 76720 Wall Transformer, 76752 SureTemp Thermometer w/back light & 52401 KleenSpec Plus Dispenser
- 76720-LP 76720 Wall Transformer, 76752 SureTemp Thermometer w/back light & 52101 KleenSpec Plus Dispenser
- 76720-T 76720 Wall Transformer & 76751 SureTemp Thermometer
- 76720-TA 76720 Wall Transformer, 76751 SureTemp Thermometer & 767 Wall Aneroid
- 76720-TAK 76720 Wall Transformer, 76751 SureTemp Thermometer, 767 Wall Aneroid & 52401 KleenSpec Plus Dispenser
- 76720-TAP 76720 Wall Transformer, 76751 SureTemp Thermometer, 767 Wall Aneroid & 52101 KleenSpec Plus Dispenser
- 76720-TK 76720 Wall Transformer, 76751 SureTemp Thermometer & 52401 KleenSpec Plus Dispenser
- 76720-TP 76720 Wall Transformer, 76751 SureTemp Thermometer & 52101 KleenSpec Plus Dispenser
- 76730 Third Handle Module
- 76910 76710 Wall Transformer & 767 Wall Aneroid
- 76920 76720 Wall Transformer & 767 Wall Aneroid
- 79510 76710 Wall Transformer, 52401 KleenSpec Plus Dispenser & 767 Wall Aneroid
- 79520 76720 Wall Transformer, 52401 KleenSpec Plus Dispenser & 767 Wall Aneroid



70000



70500



70700



71000-A



71050




71000-C

WELCH ALLYN HANDLE FOR 3.5-V HALOGEN INSTRUMENTS

Rechargeable, high-performance handles provide constant 3.5V output. Use rechargeable nickel-cadmium batteries. Heads must be ordered separately.

- 70000 Standard Battery Handle "D" Size
- 70500 Standard Battery Handle "C" Size
- 70700 2.5-V Handle
- 71000-A 3.5V Handle
- 71050 Extra Bottom Section Only For 71000, 110-130V, 50/60 Hz
- 71050-C Extra Bottom Section Only For 71000-C, 110-130V, 50/60 Hz
- 71000-C 3.5-V Handle, Can Convert to C Battery Use in Emergencies



WELCH ALLYN REPLACEMENT RECHARGEABLE BATTERIES

Nickel-cadmium rechargeable batteries for Welch Allyn handles. Color coded by voltage.

- 72000 2.5V Rechargeable Batteries For Nos. 70700, 70710, 70715, 70720, 70750, 70751, 70754, 70755, 71500 & Discontinued 60700 Handles, Color-Coded Red
- 72100 2.5V Rechargeable Batteries For Discontinued Nos. 617, 715 or 715 Handles, Color-Coded Green
- 72200 3.5V Rechargeable Batteries For Nos. 71000, 71010, 71015, 71020, 71050, 71051, 71054, 71055 & 71670 Handles, Color-Coded Black

WELCH ALLYN UNIVERSAL CHARGER

Accepts all Welch Allyn 2.5V, 3.5V, and PocketScope™ (requires adapter) well rechargeable handles. Has a green LED indicating power on. Has a yellow LED under each well indicating the handles are seated and charging. Sports an updated design to fit with our old and new products. Meets all new Agency standards.

- 71630 Universal Charger w/o instrument heads, Two 71670 NiCad Handles
- 71640 Universal Charger w/o instrument heads, Two 71900 L-ion Handles
- 71641-M Universal Charger, (2) 71670 NiCad Handles, 11720 Ophthalmoscope & 23810 MacroView Otoscope
- 71641-MS Universal Charger, (2) 71900 L-ion Handles, 11720 Ophthalmoscope & 23810 MacroView Otoscope
- 71811-MP Universal Charger, (2) 71670 NiCad Handles, 11810 Ophthalmoscope & 23810 MacroView Otoscope
- 71811-MPS Universal Charger, (2) 71900 L-ion Handles, 11810 Ophthalmoscope & 23810 MacroView Otoscope


Accessories

- 71140 Universal Charger Only
- 71249 Pocket Scope Adapter Sleeve Only
- 71310 Specula Tray
- 71420 Wall Mount Kit

ADC® DIAGNOSTIC REPLACEMENT LAMPS

Manufactured to provide extra long life and reliable illumination. Unconditionally guaranteed to be free from defects and give trouble free service. Rated for 20 hours continuous use. Packaged individually in ziplock poly bags.

- 4500-1 Lamp, WA#001100
- 4501-1 Lamp, WA#00400
- 4502-1 Lamp, WA#00200
- 4504-1 Lamp, WA#00800
- 4507-1 Lamp, WA#01300
- 4513-1 Lamp, WA#03000
- 4514-1 Lamp, WA#03100
- 4516-1 Lamp, WA#03400
- 4518-1 Lamp, WA#03900
- 4522-1 Lamp, WA#04900
- 4523-1 Lamp, WA#06000
- 4530-1 Lamp, 2.5V Halogen For 5120
- 4531-1 Lamp, 2.5V Halogen For 5121
- 4532-1 Lamp, 2.5V Halogen For 4065F to 4067F




WELCH ALLYN REPLACEMENT LAMPS

- 00200-U Vacuum Replacement Lamp For 41000
- 00900-U 2.5V Vacuum Lamp for OBS17600
- 01300-U 2.5V Vacuum Replacement Lamp
- 01400-U 2.5V Vacuum Replacement Lamp For 77800
- 02000-U 7V Vacuum Replacement Lamp For 45003
- 02500-U Replacement Lamp
- 02600-U Replacement Lamp
- 02900-U 2.5V Halogen Replacement Lamp For 290S
- 03000-U Halogen Replacement Lamp
- 03100-U Halogen 3.5V Replacement Lamp
- 03300-U 2.5V Halogen Replacement Lamp
- 03400-U 2.5V Halogen Replacement Lamp
- 03600-U Vacuum Replacement Lamp For 60500
- 03700-U 3.5V Halogen Replacement Lamp For 18100
- 03900-U Halogen Replacement Lamp For 12810
- 04050-U 150-W Halogen Replacement Lamp
- 04100-U Halogen Replacement Lamp
- 04200-U Halogen Replacement Lamp
- 04400-U Halogen Replacement Lamp For 11470
- 04430-U Halogen Replacement Lamp For LS135 Light
- 04450-U Halogen Replacement Bulb
- 04800-U 2.5V Vacuum Lamp For Laryngoscope
- 04900-U 3.5V Halogen Replacement Lamp
- 06000-U Halogen Replacement Lamp
- 06200-U Halogen Replacement Lamp
- 06300-U Halogen Replacement Lamp
- 06400-U Halogen Replacement Lamp
- 07200-U Halogen Replacement Lamp For 17200
- 07600-U Halogen Replacement Lamp For 76600
- 07800-U 6V Halogen Replacement Lamp
- 08200-U Halogen Replacement Lamp For 18200
- 08300-U Halogen Replacement Lamp For 18300
- 08800-U 4.6V Halogen Replacement Lamp For 78810 Illumination System
- DA-ARC Lamp Replacement For DLS100
- DA-L400 Lamp Replacement For DLS400 & DCU2000

ADC® FIBEROPTIC LARYNGOSCOPE BATTERY HANDLES

Lightsource in handle simplifies blade cleaning, allows for use of brighter, whiter light. Satin chrome plated brass handle with a knurled finish for a positive grip. Fully compatible with all standard hook-on blades that conform to ASTM (F965) and ISO (7376-1). Color coded green to signify compatibility ONLY with other fiberoptic handles. Batteries not included.

- 4065F Laryngoscope Fiberoptic Battery Handle, Medium "C" Size
- 4066F Laryngoscope Fiberoptic Battery Handle, Small "AA" Size
- 4067F Laryngoscope Fiberoptic Battery Handle, Stubby "AA" Size



ADC® FIBEROPTIC LARYNGOSCOPE BLADES

Fiberoptic light transmission provides cool obstruction free viewing. 2.5V halogen illumination for true tissue color and consistent long lasting lamp performance. Stainless steel. Satin, non-glare finish. Lightweight design. Fully compatible with all fiberoptic hook on handles that conform to ASTM (F965) and ISO (7376-1). Color coded green to signify compatibility and differentiate from standard hook on models.

- 4070F Size 0 Macintosh Blade, Premie
- 4071F Size 1 Macintosh Blade, Infant
- 4072F Size 2 Macintosh Blade, Child
- 4073F Size 3 Macintosh Blade, Medium Adult
- 4074F Size 4 Macintosh Blade, Large Adult
- 4080F Size 0 Miller Blade, Premie
- 4081F Size 1 Miller Blade, Infant
- 4082F Size 2 Miller Blade, Child
- 4083F Size 3 Miller Blade, Medium Adult
- 4084F Size 4 Miller Blade, Large Adult

ADC® LARYNGOSCOPE SET

Choose from either a 4 blade Macintosh or 5 blade Miller set with both AA and C battery handles. Foam lined fitted hard case with carry handle. Batteries included.

- 4079 Standard Macintosh Laryngoscope Set
- 4079F Fiberoptic Macintosh Laryngoscope Set
- 4089 Standard Miller Laryngoscope Set
- 4089F Fiberoptic Miller Laryngoscope Set



ADC® STANDARD LARYNGOSCOPE BATTERY HANDLES

Satin chrome plated brass with a knurled finish for a positive grip. Fully compatible with all standard hook-on blades that conform to ASTM (F965) and ISO (7376-1). Batteries not included.

- 4065 Laryngoscope Battery Handle, Medium "C" Size
- 4066 Laryngoscope Battery Handle, Small "AA" Size
- 4067 Laryngoscope Battery Handle, Stubby "AA" Size

ADC® STANDARD LARYNGOSCOPE BLADES

Stainless steel construction with satin, non-glare finish. Brass electrical contacts enhance conductivity. Lightweight design. Fully compatible with all standard hook on handles that conform to ASTM (F965) and ISO (7376-1).

- 4070 Size 0 Macintosh Blade, Premie
- 4071 Size 1 Macintosh Blade, Infant
- 4072 Size 2 Macintosh Blade, Child
- 4073 Size 3 Macintosh Blade, Medium Adult
- 4074 Size 4 Macintosh Blade, Large Adult
- 4080 Size 0 Miller Blade, Premie
- 4080-0 Size 00 Miller Blade, Small Premie
- 4081 Size 1 Miller Blade, Infant
- 4082 Size 2 Miller Blade, Child
- 4083 Size 3 Miller Blade, Medium Adult
- 4084 Size 4 Miller Blade, Large Adult
- 4090 Size 0 Wisconsin Blade, Premie
- 4091 Size 1 Wisconsin Blade, Infant
- 40915 Size 1.5 Wisconsin Blade, Small Child
- 4092 Size 2 Wisconsin Blade, Child
- 4094 Size 4 Wisconsin Blade, Large Adult
- 4093 Size 3 Wisconsin Blade, Medium Adult



WELCH ALLYN KLEENSPEC® DISPOSABLE LARYNGOSCOPES

Designed for single use. Light source inserts in the handle. After use, remove the light source and discard the blade and handle.

- 56901 Small Laryngoscope, 10/bx
- 56903 Large Laryngoscope, 10/bx

Accessories


- 60500 Illuminator For Laryngoscope

We offer:

- Great Service
- Great Pricing
- Great Products




68040



68060



68696



69061



69211



69696



69697

WELCH ALLYN LARYNGOSCOPE BLADES

One-piece stainless steel construction.

05681-U	Fiber Optic - Miller Blade Case
05691-U	Fiber Optic - MacIntosh Blade Case
05701-U	Fiber Optic - English MacIntosh Blade Universal Case
05791-U	Fiber Optic - English MacIntosh Blade Standard Case
68040	Miller Blade, Size 0
68041	Miller Blade, Size 1
68042	Miller Blade, Size 2
68043	Miller Blade, Size 3
68044	Miller Blade, Size 4
68060	Fiber Optic - Miller Blade, Size 0
68061	Fiber Optic - Miller Blade, Size 1
68062	Fiber Optic - Miller Blade, Size 2
68063	Fiber Optic - Miller Blade, Size 3
68064	Fiber Optic - Miller Blade, Size 4
68470	Miller Blade, Size 0 (Neonatal)
68696	Fiber Optic - Miller Blade, Complete Set w/medium & small handles, 05681 Case
69041	MacIntosh Blade, Size 1
69042	MacIntosh Blade, Size 2
69043	MacIntosh Blade, Size 3
69044	MacIntosh Blade, Size 4
69061	Fiber Optic - MacIntosh Blade, Size 1
69062	Fiber Optic - MacIntosh Blade, Size 2
69063	Fiber Optic - MacIntosh Blade, Size 3
69064	Fiber Optic - MacIntosh Blade, Size 4
69211	Fiber Optic - English MacIntosh Blade, Size 1
69212	Fiber Optic - English MacIntosh Blade, Size 2
69213	Fiber Optic - English MacIntosh Blade, Size 3
69214	Fiber Optic - English MacIntosh Blade, Size 4
69242	English MacIntosh Blade, Size 2
69243	English MacIntosh Blade, Size 3
69696	Fiber Optic - MacIntosh Blade, Complete Set w/medium & small handles, 05691 Case
69697	Fiber Optic - English MacIntosh Blade, Complete Set w/medium & small handles, 05791 Case



60200



60300



60305



60400



60710



60713



60813

WELCH ALLYN LARYNGOSCOPE HANDLES


Knurled finish for sure grip.

60200	Standard Laryngoscope Handle D - Size Battery
60300	Standard Laryngoscope Handle C - Size Battery
60305	Standard Stubby Laryngoscope Handle
60400	Standard Laryngoscope Handle AA - Size Battery
60710	Standard 2.5V Rechargeable Laryngoscope Handle Only
60713	2.5V Rechargeable Laryngoscope Handle Only For Fiber Optic Blades
60715	Standard Rechargeable Battery Handle & Transformer
60813	Laryngoscope Handle For Fiber Optic Blades, Uses C - Size Batteries
60814	Laryngoscope Handle For Fiber Optic Blades, Uses AA - Size Batteries
60815	Stubby Laryngoscope Handle For Fiber Optic Blades

WELCH ALLYN 3.5V HALOGEN NASAL ILLUMINATOR

When attached to throat illuminator section of No. 20000 otoscope/throat illuminator, provides illumination of the nasal passages through a 9mm speculum.

26530	3.5-V Halogen Nasal Illuminator, Complete w/9mm polypropylene speculum
26535	Nasal Illuminator Section Only



WELCH ALLYN 3.5V HALOGEN TRANSILLUMINATORS

These instruments use halogen illumination transmitted through fiber optics for brilliant, yet cool illumination. Fits any Welch Allyn 3.5V power source. May be sterilized by steam, gas, or some germicidal solutions.

43300 Curved All Purpose Transilluminator

Accessories

41002 Cobalt Blue Filter Only

WELCH ALLYN 3.5V TONGUE BLADE HOLDER

New halogen HPX™ lamp that provides 30% more light output, for true tissue color and consistent, long lasting illumination. Fiber optic diffuser for cool distal light. Universal clamp accommodates most tongue blades.

28100 Tongue Blade Holder



WELCH ALLYN LARYNX & NASOPHARYNX ILLUMINATORS

Halogen illumination and fiber optics light transmission provide cool, intense light. Wide-angle mirror. Light aperture assists in positioning mirror. Ready for immediate use. Compatible with any Welch Allyn power source.

27010 2.5V or 3.5V Halogen Larynx Mirror Section Only

27050 2.5V or 3.5V Halogen Nasopharynx Illuminator, Complete

ADC® 2.5V DERMASCOPE ACCESSORIES

5312-1 Plain Contact Plate

5312-2 Graduated Contact Plate

5312-3 Dermoscope Head

5312-4 Lamp 2.5V Halogen

5312-ZC Black Zipper Case



WELCH ALLYN EPISCOPE™

For examining the entire skin surface and making accurate measurements of skin lesions. Bright, 3.5-V halogen illumination and 10x optics system ensure true tissue color and accurate identification of suspicious lesions. Attaches to rechargeable handle. Reticle contact plate allows easy measurement in millimeters. Removable contact plate is autoclavable to eliminate the possibility of cross-contamination. Easy to focus, large knurled focusing ring. Comfortable eyepiece with soft rubber eyewear protector. Fully portable. Soft carrying case accommodates handle and includes pouch for extra front cover. Handles must be ordered separately.

47300 EpiScope Skin Surface Microscope w/reticle contact plate

47351 EpiScope Set (includes EpiScope, Reticle Contact Plate, Rechargeable Handle, Nickel Cadmium Battery & Soft Carrying Case)

Accessories


05232-U Soft Carrying Case

47310 Reticle Contact Plate

47320 Plano Contact Plate



65002



65003



65050

WELCH ALLYN RL-150 RHINOLARYNGOSCOPE

For upper-respiratory exams. Procedure time is less than 5 minutes and easy to master. High patient acceptance rate. Simple one-step disinfection system. Free lifetime loaner guarantee policy. Short payback period. Reimbursable. Eliminates unnecessary referrals.

65050 RL-150 Rhinolaryngoscope

Accessories

65002 Disinfection Tray

65003 Hanger

65006 LX-150 Adaptor For 65051 Rhinolaryngoscope/LX-150 Light Source



AS032L




AS034L

AMSINO AMSMOOTH™ DISPOSABLE VAGINAL SPECULUM

Made of polystyrene with see through clarity. Provides smooth, one-handed operation for easier insertion, requires no lubricant except water, and allows for secure positioning without slippage. Single use eliminates cross contamination. The self-locking feature speeds examination and saves time.

- | | |
|--------|--|
| AS032L | Large Grave Style Vaginal Speculum, 10/bg, 10 bg/cs |
| AS032M | Medium Grave Style Vaginal Speculum, 10/bg, 10 bg/cs |
| AS032S | Small Grave Style Vaginal Speculum, 10/bg, 10 bg/cs |
| AS034L | Disposable Vaginal Specula, Graves Style, Individually Wrapped, Large, w/Light Source Option, 25/bx, 4bx/cs |
| AS034M | Disposable Vaginal Specula, Graves Style, Individually Wrapped, Medium, w/Light Source Option, 25/bx, 4bx/cs |
| AS034S | Disposable Vaginal Specula, Graves Style, Individually Wrapped, Small, w/Light Source Option, 25/bx, 4bx/cs |



COOPERSURGICAL SANI-SPEC® SPECULA

Clear plastic material allows 360° visibility and reduces thermal shock. Disposable. Provides protection of the doctor, the nurse and the patient from cross contamination. Plastic material avoids unpleasant temperature shock associated with metal. Smooth pre-lubricated blades with 6 position grooves.

- | | |
|-------|--|
| 82329 | Small Specula, 1"W x 4 1/2"L, 10/bg, 10 bg/cs |
| 82330 | Medium Specula, 1 1/4"W x 4 3/4"L, 10/bg, 10 bg/cs |
| 82331 | Large Specula, 1 1/3"W x 5 1/3"L, 10/bg, 10 bg/cs |
| 82340 | Small Specula, Individually Wrapped, 50/cs |
| 82341 | Medium Specula, Individually Wrapped, 50/cs |

DRG SOFTSPEC VAGINAL SPECULUM

A tampon-like, disposable device that is inserted into the vagina and gently inflated to expand the vaginal walls for guaranteed visualization of the cervix. Allows easy access for an adequate pap smear by inserting the spatula through the donut shape of the SoftSpec. Disposable so you do not have to worry about messy clean up or the spread of bacterial from patient to patient.

- | | |
|------------|---|
| VSI-25-001 | SoftSpec Vaginal Speculum Illuminated w/bulb, 25/bx |
| VSO-25-001 | SoftSpec Vaginal Speculum w/bulb, 25/bx |




STARLINE VAGINAL SPECULA

Recommended for routine gynecological exams. Disposable. Allows the examination of the vaginal mucosa at 360° angle without rotation. Wide range of adjustability. Ease of operation and manipulation with one hand. Self-holding in various opening positions. Made from hard, clear and non-toxic plastic resin of high strength.

- | | |
|----------|---|
| 16108309 | Small Vaginal Specula, 10/bg, 10 bg/cs |
| 16108310 | Medium Vaginal Specula, 10/bg, 10 bg/cs |
| 16108311 | Large Vaginal Specula, 10/bg, 10 bg/cs |

COOPER SURGICAL GALENICA GYNE-LITE™ SYSTEM ACCESSORIES

- | | |
|------|--------------------------|
| 9220 | Adaptor, 10/bg |
| 9230 | Sheaths, 100/bx, 4 bx/cs |
| 9240 | Bulb |



COOPER SURGICAL GALENICA VAGINAL SPECULA

Recommended for routine gynecological exams. Disposable. Allows the examination of the vaginal mucosa at 360° angle without rotation. Wide range of adjustability. Ease of operation and manipulation with one hand. Self-holding in various opening positions. Made from hard, clear and non-toxic plastic resin of high strength. Unit Pack (poly bag) non-sterile. Unit Pack (peeling pouch) sterile.

- | | |
|------|--|
| 8309 | Small Vaginal Specula, Non-Sterile, 10/bg, 10 bg/cs |
| 8310 | Medium Vaginal Specula, Non-Sterile, 10/bg, 10 bg/cs |
| 8311 | Large Vaginal Specula, Non-Sterile, 10/bg, 10 bg/cs |
| 8615 | Small Unit Pack Vaginal Specula, 50/cs |
| 8616 | Medium Unit Pack Vaginal Specula, 50/cs |
| 8617 | Large Unit Pack Vaginal Specula, 50/cs |

Accessories

- | | |
|------|--|
| 9310 | Smoke Evacuator For All Specula, 40/cs |
|------|--|



WELCH ALLYN KLEENSPEC® 580 DISPOSABLE VAGINAL SPECULA WITH BARRIER

Superior halogen illumination for a whiter, brighter light for truer tissue color rendition. Smooth, one-handed angulation and elevation mechanism with positive locking positions. Light handle is completely sealed when inserted into the speculum. Built-in reservoir prevents fluids from contaminating the outside of the speculum handle. Optional disposable sheath covers the exposed sections of the light handle and cord, virtually eliminating the threat of cross-contamination.

58010 Disposable Sheaths, 25/bx, 5 bx/cs



WELCH ALLYN KLEENSPEC® 580 SERIES DISPOSABLE VAGINAL SPECULA

Halogen light for true tissue color and consistent, long lasting illumination. Disposable for economy, safety and convenience. One-hand operation with angulation and elevation for convenience. Especially narrow blade width available for adolescent and geriatric patients. Individually wrapped specula.

- 58000 Small Specula (Graves Style), 25/bx, 4 bx/cs
- 58000S Small Specula w/sheath, 20/bx, 4 bx/cs
- 58000S-B Small Specula w/sheath, 20/bx
- 58001 Medium Specula (Graves Style), 25/bx, 4 bx/cs
- 58001S Medium Specula w/sheath, 20/bx, 4 bx/cs
- 58001S-B Medium Specula w/sheath, 20/bx
- 58004 Large Specula (Graves Style), 20/bx, 4 bx/cs
- 58004S Large Specula w/sheath, 15/bx, 4 bx/cs
- 58005 Small Specula w/small evacuator tube, 12/bx, 4 bx/cs
- 58006 Medium Specula w/small evacuator tube, 12/bx, 4 bx/cs



WELCH ALLYN KLEENSPEC® DISPOSABLE VAGINAL SPECULA ILLUMINATION SYSTEM

Permanent light carrier projects intense beam of light through plastic light guide. No external light source necessary.

- 78810 Complete Vaginal Specula Illumination System w/on-off switch & transformer
- Accessories**
- 73305 Transformer w/5 ft cord, 110-130V, 50/60 Hz
- 73410 Transformer w/5 ft cord, 110-130V, 50/60 Hz, UL Listed, CSA Certified
- 78800 Permanent 6V Halogen Light Handle w/on-off switch



58600



78600

WELCH ALLYN KLEENSPEC® VAGINAL SPECULA

Brilliant halogen illumination system allows diagnostic accuracy. Clear, polystyrene construction. Designed for one-hand operation. Single-use prevents cross-contamination. Expendable light carrier for safety. (Illumination system not included.)

- 58600 Small Specula, 25/bx, 4 bx/cs
- 58601 Medium Specula, 25/bx, 4 bx/cs
- 78600 Expendable Illuminator (can be resterilized 50 times before disposal), 6/bx, 6 bx/cs



SEILER 935 THREE-STEP COLPOSCOPE

- 30935-3 Model 935 Seiler 3-Step Colposcope
- 30935-3-D Model 935 Colposcope w/digital
- 30935-3-V Model 935 Colposcope w/video



SEILER 955 FIVE-STEP COLPOSCOPE


- 30955-5 Model 955 Colposcope
- 30955-5-D Model 955 Colposcope w/digital
- 30955-5-V Model 955 Colposcope w/video

Contact your representative for information on cryosurgical units.



SEILER 965Z ZOOM COLPOSCOPE

- 30965-Z Model 965 Colposcope
- 30965-Z-D Model 965 Colposcope w/digital
- 30965-Z-V Model 965 Colposcope w/video



WELCH ALLYN IMAGE CAPTURE SYSTEM FOR VIDEO COLPOSCOPE

A unique software-based documentation system to be used with the Welch Allyn video colposcope. The software can be controlled remotely from the Video Colposcope and enables the user to capture, review, store, print and export images.

- 88801 Image Capture System For Desktop Computer
- 88802 Image Capture System For Laptop Computer



37019



37023



37027

WELCH ALLYN FIBER OPTIC ANOSCOPES

Fiber optics provide cool distal ring of light, no reflections, no obstructions. Autoclavable stainless steel. 45° bevel tip (90° available on request). Set includes speculum, obturator, light handle and transformer.

- 37019 19mm x 7cm Speculum Only w/obturator
- 37023 23mm x 7cm Speculum Only w/obturator
- 37027 27mm x 7cm Speculum Only w/obturator



WELCH ALLYN KLEENSPEC® DISPOSABLE ANOSCOPE

45° bevel end. Individually wrapped. 19mm x 10cm with obturators.

- 53110 Disposable Anoscope, 25/bx, 4 bx/cs



38108



38614



38900



39614

WELCH ALLYN STANDARD ANOSCOPES

Halogen light for true tissue color and consistent, long-lasting illumination. Autoclavable. Comes with obturator.

- 38108 8mm Aperture, Anoscope Only
- 38114 14mm Aperture, Anoscope Only, Standard Speculum, 89mm Length
- 38119 19mm Aperture, Anoscope Only, Standard Speculum, 89mm Length
- 38122 22mm Aperture, Anoscope Only, Standard Speculum, 89mm Length
- 38614 14mm Aperture, Anoscope Only, Operating (Cut-Out) Speculum, 89mm Length
- 38619 19mm Aperture, Anoscope Only, Operating (Cut-Out) Speculum, 89mm Length
- 38622 22mm Aperture, Anoscope Only, Operating (Cut-Out) Speculum, 89mm Length
- 38900 20mm ID, 22mm OD, 68mm Length to Taper, Anoscope Only, Rotating Speculum
- 39614 14mm Aperture, Anoscope Only, Long Speculum, 127mm Length
- 39619 19mm Aperture, Anoscope Only, Long Speculum, 127mm Length
- 39622 22mm Aperture, Anoscope Only, Long Speculum, 127mm Length

COOPER SURGICAL GALENICA ANOSCOPES

Easy manipulation for septic anal/rectal examinations. Decreases the risks of cross-infections. 56° bevel at the insertion end. Adaptable to fiber optic lighthouse. Made from hard non-toxic plastic resin of fine medical grade. Disposable and lightweight.

8607 Anoscope, 25/bx, 4 bx/cs



32010



32410



32820



35303

WELCH ALLYN FIBER OPTIC SIGMOIDOSCOPIES

Halogen light. Fiber optics provide cool distal ring of light, no reflections, no obstructions. Autoclavable stainless steel. Hinged window with neoprene seal. Graduated tube with removable tip. Complete set includes: sigmoidoscope speculum, obturator, light handle, insufflation bulb & transformer

- 32010 11mm x 10cm Sigmoidoscope Speculum Only w/obturator
- 32020 11mm x 25cm Sigmoidoscope Speculum Only w/obturator
- 32410 15mm x 15cm Sigmoidoscope Speculum Only w/obturator
- 32420 15mm x 25cm Sigmoidoscope Speculum Only w/obturator
- 32810 19mm x 15cm Sigmoidoscope Speculum Only w/obturator
- 32820 19mm x 25cm Sigmoidoscope Speculum Only w/obturator
- 32830 19mm x 35cm Sigmoidoscope Speculum Only w/obturator
- 33220 23mm x 25cm Sigmoidoscope Speculum Only w/obturator
- 33830 19mm x 30cm Sigmoidoscope Speculum Only w/obturator
- 35303 Set Includes: 32410 Sigmoidoscope, 32820 Sigmoidoscope, 37023 Anoscope, 73210 Light Handle w/cord, 73305 Transformer, 30200 Insufflation Bulb Complete, 07800 Extra Halogen Lamps (2), 30130 Suction Tube, 05391 Case

Accessories

05391-U Hard Case


WELCH ALLYN FX-100 FLEXIBLE FIBER OPTIC SIGMOIDOSCOPE ACCESSORIES

- 33905 Water Bottle Assembly
- 33950 Biopsy Forceps

WELCH ALLYN KLEENSPEC® FIBER OPTIC DISPOSABLE SIGMOIDOSCOPE

Speculum attaches to fiber optic light head; no delicate or protruding light carriers to obstruct vision. Brilliant distal illumination is shadow free, without color distortion. For use with any Welch Allyn 3.5V or 2.5V power source and handle adapter or 6.5V handle and transformer.

36103 Complete Illumination System includes Fiber Optic Light Head, Light Handle, Transformer & Insufflation Bulb




30130



30200




34000



34100



36019



38300



53130



73210



73550

WELCH ALLYN KLEENSPEC® FIBER OPTIC DISPOSABLE SIGMOIDOSCOPE ACCESSORIES

- 30130 30cm Suction Tube
- 30140 40cm Suction Tube
- 30200 Insufflation Bulb, Bladder Net & Tubing
- 34000 Rectal Probe, 7 1/2"
- 34100 Rectal Hook, 7 1/2"
- 36019 Fiber Optic Light Head
- 38300 Incandescent Light Carrier
- 38700 Fiber Optic Light Carrier
- 53130 Disposable Specula (19 mm x 25 cm) w/obturator, 25/bx, 4 bx/cs
- 73210 Light Handle w/3 ft cord
- 73500 Halogen 3.5-V Handle Adapter
- 73550 Halogen 2.5-V Handle Adaptor



23200



23210


WELCH ALLYN AM 232™ MANUAL AUDIOMETER

Ergonomic panel design and convenient controls allow one-handed operation of this easily portable audiometer. Frequency values and HL levels are easily visible at a glance. The AM 232 also offers selection modes for steady, pulsed or frequency-modulated stimulus for increased reliability of response. +10 dB push-button prevents mistaken presentation of uncomfortable frequency tones. Provides ample storage space inside hard-shell case for standard headphones and optional accessories.

- 23200 AM 232 Audiometer w/AC Only
- 23210 AM 232 Audiometer w/AC Power & Rechargeable Battery




23220



52300

WELCH ALLYN AM 232™ MANUAL AUDIOMETER ACCESSORIES

- 23220 Response Hand Switch
- 23221 Patch Cord
- 23222 Audiocups
- 52300 Audiogram Form, 50 sheets/pad, 10 pads/cs
- 72320 Replacement Rechargeable Battery Pack Only
- 72329 Complete Battery Pack Assembly w/rechargeable batteries



23300



92680

WELCH ALLYN AUDIOSCOPE® 3 SCREENING AUDIOMETER/OTOSCOPE

Versatile instrument for screening hearing loss combines 20, 25 and 40 dB HL levels in one unit. Tones are presented in random intervals so patients can not second guess their hearing screen. This varied timing means improved results. Built-in 1000 Hz pretone, presented at 20 dB HL above the screening level, allows patients to practice listening before they're actually screened. After pretone, patient is given the easiest tones for enhanced reliability of responses. 110-130V.


- 23300 AudioScope 3 w/AudioSpecs (no transformer)
- 92600 AudioScope 3 Set (AudioScope 3, set of AudioSpecs, charging transformer, recording forms)
- 92680 Complete Set (AudioScope 3, set of AudioSpecs, charging transformer, charging stand, recording forms)



23400



71040



72300

WELCH ALLYN AUDIOSCOPE® 3 SCREENING AUDIOMETER/OTOSCOPE ACCESSORIES

- 23303 Small AudioSpec
- 23305 Medium AudioSpec
- 23307 Large AudioSpec
- 23400 Set of AudioSpecs (small, medium, large)
- 55230 Recording Forms, 100/pk, 10 pk/bx
- 71040 Charging Transformer
- 71123 Charging Stand Only
- 72300 3.5V Nickel-Cadmium Rechargeable Replacement Battery
- 94005-0000 Paper, 5.3" x 5.8", 140 shts/pk, 12 pk/cs



WELCH ALLYN EAR WASH SYSTEM PLUS

Easy to use by healthcare professionals and patient satisfaction. Chamber regulates the output pressure to prevent damage to the tympanic membrane. Handle indicator will contrast in color when water is within 10 degrees of body temperature. Chamber automatic shut-off will shut off system if water reaches 110 degrees. Multi-functional disposable tip. Device can be moved easily from location to location via a storage tray. One-year warranty.

- 29350 Ear Wash System Plus, English



23640



93690

WELCH ALLYN MICRO TYMP® 2 PORTABLE TYMPANOMETRIC INSTRUMENT

Portable, hand-held unit detects otitis media and other middle ear disorders. Measures sound returned from eardrum under positive, negative, and atmospheric pressure to determine status of middle ear. Accurately identifies many middle ear disorders in asymptomatic cases. Printer/charger keeps Micro Tym fully charged and provides hardcopy documentation. Results instantly displayed on LCD before printing. 226 Hz frequency, 85 db HL intensity.

- 23640 Micro Tym II (includes 2 sets of tips)
- 93690 Micro Tym II Set (includes printer/charger, extra paper & tips)



23630



71170

WELCH ALLYN MICRO TYMP® 2 PORTABLE TYMPANOMETRIC INSTRUMENT ACCESSORIES

- 05276-U Carrying Case
- 23630 Set of Micro Tympanometry Tips (4 sizes)
- 23639 Multi-Pack of Micro Tympanometry Probe Tips
- 71170 Printer/Charger (includes paper)
- 72900 2.4V Nickel Cadmium Rechargeable Battery



WELCH ALLYN TM 262™ AUTO TYMP™


Offers tympanometry and ipsilateral acoustic reflex testing (with optional manual audiometry). Built-in printer and display. Improved ear seal via lightweight handheld probe. Delivers precision results with each use, and stores up to eight test results.

- 26200 Auto Tympanometer w/Ipsi Reflex
- 26230 Auto Tympanometer w/audiometer




WELCH ALLYN TM 262™ AUTO TYMP™ ACCESSORIES


- 05260-U Carrying Case
- 26100 Eartips (probe), 6 sizes
- 26240 Dust Cover
- 26241 Test Cavity
- 52600 Printer Paper, 5/bx



3050



3051



3052

TECH-MED EYE CHARTS

- 3050 Snellen Eye Test Chart, 20 ft, Non-Reflective Matte Finish, 22" x 11"
- 3051 Illiterate Eye Test Chart, 20 ft, Non-Reflective Matte Finish, 22" x 11"
- 3052 Kindergarten Eye Chart 20 ft, Non-Reflective Matte Finish, 22" x 11"
- 3053 Pocket Eye Chart, Use at 14", Provides 20/800 Distance, Laminated Plastic, 6 1/2" x 3 1/2"
- 3060 Illuminated Eye Test Cabinet, Baked Enamel Finish w/handle & power cord, Includes 1 each Illuminated Eye Charts
- 3061 Illuminated Snellen Eye Test Chart, 20 ft
- 3062 Illuminated Snellen Eye Test Chart, 10 ft
- 3063 Illuminated Kindergarten Eye Test Chart, 20 ft
- 3064 Illuminated Illiterate Eye Test Chart, 20 ft
- 3065 Illuminated Illiterate Eye Test Chart, 10 ft

GRAHAM-FIELD GRAFCO® ILLITERATE PLASTIC EYE CHART

With green and red color bars. 20' distance.
1241 Illiterate Eye Chart, Plastic, 22" x 11"

GRAHAM-FIELD GRAFCO® KINDERGARTEN PLASTIC EYE CHART

Shows easy to distinguish objects. 20' distance.
1243 Kindergarten Eye Chart, Plastic, 22" x 11"


GRAHAM-FIELD GRAFCO® SNELLEN TYPE PLASTIC EYE CHART

Non-reflective, matte finish with green and red color bars. 20' distance.
1240 Snellen Eye Chart, Plastic, 22" x 11"
1264 Snellen Eye Chart, 10' Distance, Translucent for Use w/illuminated

GRAHAM-FIELD GRAFCO® ISHIHARA® TEST CHART BOOKS, FOR COLOR DEFICIENCY

Accepted by leading authorities as a simple and accurate method for discovering congenital color blindness and red-green blindness, each in two forms, complete and incomplete. Includes four special plates for tests to determine the kind and degree of defect in color vision. Makes use of the peculiarity that in red-green blindness, blue and yellow appear remarkably bright compared with red and green. Color plates encased in specialty designed album-type books for ease of handling. Detailed instructions included.

- 1254 Ishihara Test Chart Book For Color Deficiency, 14 Plate Book
- 1255 Ishihara Test Chart Book For Color Deficiency, 24 Plate Book
- 1256 Ishihara Test Chart Book For Color Deficiency, 38 Plate Book
- 1257 Ishihara Test Chart Book For Color Deficiency, 10 Plate Book
- 1263 Kindergarten Eye Chart, Translucent for Use w/illuminated



STEREO OPTICAL OPTEC® VISION TESTER ACCESSORIES

- 1000-21 Secondary Test Slide For Optec 1000
- 1000-22 Elementary Test Slide For Optec 1000
- 1000-23 Multiple Acuity Test Slide For Optec 1000
- 1000-NPS Near Point Lens For Optec 1000
- 2000-220 Purdue Job Standards Book w/recording forms (5 pads)
- 2000-221 Headrest Tissues (5000 tissues)
- 2000-226 Bulb (6) Color Corrected (120v)
- 2000-230 Plus Lens (+1.75D) For Optec 1000 & 2000
- 2000-231 Plus Lens (+2.25D) For Optec 2000
- 2000-232 Intermediate Lenses for VDT Operator Testing (set of 5)
- 2000-233 Intermediate Lens, 39 1/2", For Optec 2000
- 2000-234 Intermediate Lens, 31 1/2", For Optec 2000
- 2000-235 Intermediate Lens, 26", For Optec 2000
- 2000-236 Intermediate Lens, 22 1/2", For Optec 2000
- 2000-237 Intermediate Lens, 20", For Optec 2000
- 2000-245 Carrying Case For Optec 2000, Hard Sided
- 2000-246 Carrying Case For Optec 1000, Soft Sided
- 2000-247 Carrying Case For Optec 1000, Hard Sided
- 2000-250 Carrying Case For Optec 2000, Soft Sided
- 2000-HO2 H.O.T.V. Slide Set 20/30 (4 slides) w/Binocular Test For Optec 2000
- 2000-HO3 H.O.T.V. Slide Set 20/40 (slides) w/Binocular Test For Optec 2000
- 2000-HO4 H.O.T.V. Slide Set 20/30 (4 slides) Monocular For Optec 2000
- 2000-HO5 H.O.T.V. Slide Set 20/40 (4 slides) Monocular For Optec 2000
- 2500-226 Bulbs (6) For Optec 2500 (24v)



STEREO OPTICAL OPTEC® VISION TESTERS

Testers have homogeneous illumination making them the only vision screeners on the market that meet the ANSI standard for illumination in vision. All instruments are micro-processor controlled and continually calibrated. Choose one of the following test packages: 2000MTS - Medical Slide Package with Michigan Preschool; 2000HOTV - Medical Slide Package with H.O.T.V.; 2000ITS - Industrial Slide Package with Job Standards Book; 2000BL - International Slide Package with Job Standards Book; 2000PRS - Ophthalmic Prescreening Slide Package (FAA); 2000RHB - Driver Rehabilitation Slide Package. Accessories include: dust cover, pointer, lens cleaner, headrest tissues, recording forms, and training manual.

- 1000P System w/Peripheral Test & Accessory Kit, Excludes Video, Includes 1000S Test Package - Multiple Acuity Slide, Elementary Test Slide & Nearpoint Lens
- 5000 System w/Manual Control & Accessory Kit, w/o Peripheral Test, Choose a Standard Test Slide Package
- 5000P System w/Manual Control, Peripheral Test & Accessory Kit, Choose a Standard Test Slide Package
- 5500 System w/Remote Control, Peripheral Test & Accessory Kit, Choose a Standard Test Slide Package
- 5500P System w/Peripheral Test, Remote Control & Accessories, Choose a Standard Test Slide Package
- 6500P System w/Microprocessor Controlled Glare & Luminance Levels Peripheral Test, Remote Control, Accessory Kit (Choose a 5000 Standard Slide Pkg or a 6500 Test Pkg)

STEREO OPTICAL C.A.T.® CONTRAST SENSITIVITY SLIDE PACKAGE

Package includes: (2000-007) NEAR Point Acuity Right Eye & Left Eye Monocular & Both Eyes Simultaneously (20/100 to 20/20), (2000-003) FAR Point Acuity Right Eye & Left Eye Monocular & Both Eyes Simultaneously (20/100 to 20/20), (2000-037) 100% FAR C.A.T. Contrast Sensitivity Monocular, (2000-038) 50% Contrast Sensitivity Monocular, (2000-039) 25% Contrast Sensitivity Monocular, (2000-040) 12.5% Contrast Sensitivity Monocular, (2000-041) 6.25% Contrast Sensitivity Monocular, (2000-010) FAR Point Color Perception (Pseudo Ishihara), (2000-012) FAR Point Lateral Phoria (1 Diopter Increments), (2000-025) FAR Point Vertical Phoria (1/2 Diopter Increments), (2000-024) FAR Point Stereo Depth Perception (400 to 20 Seconds of Arc), (2000-087) FAR Point Fusion.

2000HS2 C.A.T. Contrast Sensitivity Slide Package

STEREO OPTICAL DRIVER REHABILITATION SLIDE PACKAGE

For complete vision screening needs. Ideal for school, athletic, employment physicals, driver licensing exams, D.O.T., F.A.A. and general check-ups. Package includes: (2000-147) Demonstration Highway, (2000-003) FAR Point Acuity Right Eye & Left Eye Monocular & Both Eyes Simultaneously (20/200-20/20), (2000-010) FAR Point Color Perception (Pseudo Ishihara), (2000-007) NEAR Point Acuity Right Eye & Left Eye Monocular & Both Eyes Simultaneously (20/100-20/20), (2000-012) FAR Point Lateral Phoria (1 Diopter Increments), (2000-025) FAR Point Vertical Phoria (1/2 Diopter Increments), (2000-024) FAR Point Stereo Depth Perception (400-20 Seconds of Arc), (2000-087) FAR Point Fusion, (2000-004) FAR Point Acuity Right & Left Eye Monocular Test Also Right & Left Eye Plus Lens Test (20/200-20/20), (2000-065) FAR Point Acuity Numbers Right Eye & Left Eye & Monocular Both Eyes Simultaneously (20/200-20/20), (2000-148) FAR Point Road Sign Recognition & Depth Perception (620-323-145 Seconds of Arc), (2000-149) FAR Point Color Recognition (Red, Yellow, Green, Blue).

2000RHB Driver Rehabilitation Slide Package

STEREO OPTICAL F.A.C.T.® CONTRAST SENSITIVITY SLIDE PACKAGE

Package includes: (2000-007) NEAR Point Acuity Right Eye & Left Monocular & Both Eyes Simultaneously (20/100 to 20/20), (2000-003) FAR Point Acuity Right & Left Eye Monocular & Both Eyes Simultaneously (20/200 to 20/20), (2000-171) 1.5 Cycles Per Degree FAR F.A.C.T. Contrast Sensitivity Monocular, (2000-172) 3 Cycles Per Degree Monocular, (2000-173) 6 Cycles Per Degree Monocular, (2000-174) 12 Cycles Per Degree Monocular, (2000-175) 18 Cycles Per Degree Monocular, (2000-010) FAR Point Color Perception (Pseudo Ishihara), (2000-012) FAR Point Lateral Phoria (1 Diopter Increments), (2000-025) FAR Point Vertical Phoria (1/2 Diopter Increments), (2000-024) FAR Point Stereo Depth Perception (400 to 20 Seconds of Arc), (2000-087) FAR Point Fusion.

2000FS3 F.A.C.T. Contrast Sensitivity Slide Package

STEREO OPTICAL INDUSTRIAL SLIDE PACKAGE

For complete occupational screening using accurate and easily understood targets to help determine effective vision. Includes Purdue Job Standards Booklet to easily identify pass/fail results for all occupations. Intermediate lenses available for VDT operators which is the fastest growing population in the industry today. Now test at FAR Point, INTERMEDIATE Point and NEAR Point. Package includes: (2000-147) Demonstration Highway, (2000-001) FAR Acuity Both Eyes Binocular (20/200-20/13), (2000-008) FAR Acuity Right Eye Monocular (20/200-20/13), (2000-005) FAR Acuity Left Eye Monocular (20/200-20/13), (2000-024) FAR Stereo Depth Perception (400-20 Seconds of Arc), (2000-010) FAR Color Perception (Pseudo Ishihara), (2000-025) FAR Vertical Phoria (1/2 Diopter Increments), (2000-012) FAR Lateral Phoria (1 Diopter Increments), (2000-002) NEAR Acuity Both Eyes Binocular (20/200-20/13), (2000-009) NEAR Acuity Right Eye Monocular (20/200-20/13), (2000-006) NEAR Acuity Left Eye Monocular (20/200-20/13), (2000-013) NEAR Lateral Phoria (1 Diopter Increments).

2000ITS Industrial Slide Package

It's our job to Support Yours.



**Because people
depend on you...**

*You need dependable
equipment.*


*Products of Extraordinary Design, Exceptional Quality,
Superlative Performance... with warranties to match.*

At **ADC**[®], everything we build, everything we **DO** is
designed to help make your job a little easier.

As one of the world's leading manufacturers of value priced diagnostic
instruments and accessories, **ADC**[®] delivers the products **YOU** need
to get your job done, at prices that really make sense.

**After all, what good are great products...
if you can't afford them?**

ADC[®]
American Diagnostic Corporation



Quick, find the bug.



*Also, ask about our new
OSOM® BVBLUE® Test!*

OSOM® Trichomonas Rapid Test. Traditional wet mount testing is slow, and in many cases, results are difficult to confirm. Now Genzyme, a world leader in diagnostics, introduces the first rapid test for *Trichomonas vaginalis*. A hassle-free, dipstick-simple test, using only one reagent. Accurate, two-color results in 10 minutes. Be precise. Be confident.



For more information visit our Web site at www.genzymediagnosics.com
BVBLUE is a registered trademark of Gryphus Diagnostics, LLC

genzyme
Diagnostics

STEREO OPTICAL INTERNATIONAL SLIDE PACKAGE

For complete occupational screening using accurate and easily understood targets to help determine effective vision. Includes Job Standards Booklet to easily identify pass/fail results for all occupations. Intermediate lenses available for VDT operators which is the fastest growing population in the industry today. Now test at FAR Point, INTERMEDIATE Point and NEAR Point. Package includes: (2000-255) FAR Acuity Both Eyes Binocular (20/200-20/17), (2000-256) FAR Acuity Right Eye Monocular (20/200-20/17), (2000-257) FAR Acuity Left Eye Monocular (20/200-20/17), (2000-258) FAR Stereo Depth Perception (362-9.7 Seconds of Arc), (2000-259) FAR Color Perception (Pseudo Ishihara), (2000-253) FAR Vertical Phoria (1/2 Diopter Increments), (2000-254) FAR Lateral Phoria (1 Diopter Increments), (2000-260) NEAR Acuity Both Eyes Binocular (20/200-20/17), (2000-261) NEAR Acuity Right Eye Monocular (20/200-20/17), (2000-262) NEAR Acuity Left Eye Monocular (20/200-20/17), (2000-264) NEAR Lateral Phoria (1 Diopter Increments), (2000-263) NEAR Vertical Phoria (1/2 Diopter Increments).

2000BL International Slide Package

STEREO OPTICAL MEDICAL SLIDE PACKAGE WITH HOTV

Complete vision screening needs preschool through adult. Ideal for school, athletic physicals, employment physicals and general check-ups. Package includes: (2000-026) FAR Point Acuity Binocular Test (20/70), (2000-030) FAR Point Acuity Right & Left Eye Monocular Test (20/30), (2000-031) FAR Point Acuity Right & Left Eye Monocular Test (20/30), (2000-032) FAR Point Acuity Right & Left Eye Monocular Test (20/30), (2000-004) FAR Point Acuity Right & Left Eye Monocular Test Also Right & Left Eye Plus Lens Test (20/200-20/20), (2000-020) FAR Point Muscle Balance (Combination Lateral & Vertical Phoria), (2000-019) FAR Point Color Perception (Tumbling "E"), (2000-022) NEAR Point Muscle Balance (Combination Lateral & Vertical Phoria), (2000-003) FAR Point Acuity Right & Left Eye Monocular & Both Eyes Simultaneously (20/200-20/20), (2000-012) FAR Point Lateral Phoria (1 Diopter Increments), (2000-025) FAR Point Vertical Phoria (1/2 Diopter Increments), (2000-007) NEAR Point Acuity Right & Left Eye Monocular & Both Eyes Simultaneously (20/100-20/20)

2000HOTV Medical Slide Package w/HOTV

GE MEDICAL TUFFSAT™ SENSORS

OXY-E4-H Ear Clip Sensor
 OXY-F1-H Finger Clip Sensor (one is included w/TuffSat)
 OXY-W4-H Wrap Sensor For 6 lbs & above

KEYSTONE VIEW VISION SCREENER

The VS-IV Standard is a versatile machine that provides up to 14 tests using both near (16") and far (20") point, as well as a supplemental lens for testing intermediate vision at 26". The VS-IV Universal is an advanced binocular unit that exceeds the requirements of school vision screening programs, and because the tests are completely non-language, it is particularly useful in testing preschool children and children with possible learning disabilities.

818-400 VS-IV Deluxe-Standard Tests
 818-590 VS-IV Deluxe-Universal Tests

KEYSTONE VIEW DRIVERS VISION SCREENER

The DVS-IV is designed for the rapid, yet accurate, evaluation of large groups. Vision testing is both standardized and completely confidential. The DVS-GT contains the same features as the DVS-IV plus tests for contrast sensitivity (simulates snow, fog, rain and night driving conditions) and glare recovery (the ability to adapt to bright lights such as oncoming traffic).

818-700 DVS-IV Drivers Vision Screener
 818-830 DVS-IV Drivers Vision Screener w/glare

KEYSTONE VIEW COMPUTER VISION SCREENER AND SOFTWARE

Working with nearly any PC, the CVS-I both automates the testing process, using a standardized script, and lets you control the test instrument remotely from the computer station. The Visionary Software allows you to easily record data and produce reports.

818-130 CVS-I Screener w/Visionary Software
 818-150 Visionary Software

KEYSTONE VIEW TESTING AND TRAINING TARGET SETS FOR USE IN CORRECT-EYE SCOPE

005-348 Van Orden Transparencies, 16/set
 006-166 Franzblau Tracing Patterns Forms (Seven Patterns Series)
 297-170 Red Crayons, 10/pk
 5348-1 Van Orden Glass Transparency
 5401 Leavell Exercise Pads Set (A, B, C), 3pads/pk
 5402 Leavell Exercise Pads Set (AA, BB, CC) for Stereoreader, 3pads/pk
 725351 Red Rock Games Only, 7 games/bx
 735-265 Leavell Exercise Pads CC for Stereoreader, 3 pads/pk
 735-265-1 Leavell Exercise Pads C, 3 pads/pk
 735-266 Leavell Exercise Pads AA for Stereoreader, 3 pads/pk
 735-266-1 Leavell Exercise Pads A, 3 pads/pk
 735-267 Leavell Exercise Pads BB for Stereoreader, 3 pads/pk
 735-267-1 Leavell Exercise Pads B, 3 pads/pk
 735-289 Groffman Pattern Cards (2 test cards, 1 demonstration card) Includes Record Form 735-287
 735-292 Delacato Exercise/Reading Cards, 152/pk
 735-294 Zweig-Bruno Pads #1 (Series A, B, C), 3 pads/pk
 735-295 Zweig-Bruno Pads #2 (Series 1, 2, 3), 3pads/pk
 735-296 Groffman Books #1, 3 books/pk
 735-297 Groffman Books #2, 3 books/pk
 735-426 V.O. 2 Star Pattern Pads, 100 sheets/pad, 3 pads/pk
 735-429 Universal Series Correct-Eye Graphs, 104 drawings/pad, 3 pads/pk
 735-436 V.O. Base-In Training Charts Pads, 100 sheets/pad, 3 pads/pk
 735-955 Manas Series Correct Eye Graphs, 110 drawings/pad, 3pads/pk
 736-690 Peek-A-Boo Red Coloring Books, 10/pk
 736-831 Brock Technique Transparency Set, 28/set

KEYSTONE VIEW TESTING AND TRAINING TARGET SETS FOR USE WITH TELEBINOCULAR

005-100 Visual Skills Test Set
 005-100-8 Severe Color Blind Test
 005-100-9 Mild Color Blind Test
 005-101 Driver's Vision Test Set
 005-102 IVS (OSV) Rapid Screening Test Set
 005-103 O.E.P. Test Set (Includes 005-100)
 005-104 Spache Binocular Reading Test (May also be used in Stereoscope)
 005-106 Professional Performance Test
 005-106-11 Gross Stereopsis Test
 005-106-13 Astigmatic Dial Test
 005-106-6 Macular Suppression Test
 005-107 Pilot Test Set
 005-108 (OVS) (Industrial Screening) Test Set
 005-109 Snellen Near Point Acuity Test
 005-110 Snellen Far Point Acuity Test
 005-111 Near Point Stereopsis Test
 005-112 D.C. Aciators Test Set
 005-113 Progress of Fusion Test Set
 005-115 Industrial Short Tests
 005-116 Multi-Stereo Test Set
 005-118 Mark Diplopia Test Set
 005-119 O'Brian Central Field
 005-121 AN Series Training Targets
 005-128 Broward County Test Set
 005-129 40" Intermediate, Sloan Letter Acuity
 005-130 Landolt Rings Far Point Test

KEYSTONE VIEW TESTING AND TRAINING TARGET SETS FOR USE WITH TELEBINOCULAR (Continued)

005-131	D.C. Aviators Short Test Set
005-134	Intermediate, Landolt Ring Acuity (40")
005-135	Rapid Screening Tests for Schools
005-138	Illinois School Test
005-138-4	Illinois Preschool Test
005-139	Massachusetts School Test
005-140	Michigan School Test
005-141	New York School Test
005-141-1	Tumbling E Far Point Test
005-142	Astigmatism Screening Test
005-155	Clinical Fusion Test Set
005-156	Walton Acuity/Astigmatism Test
005-162	Gesell Test Set
005-164	Anisekonia Test Set (Special Order)
005-165	Near Point Landolt Ring Acuity Test
005-166	Near Point "Tumbling E" Acuity Test
005-170	Peek-A-Boo Test Set
005-171	Sloan Letter Acuity Test - Near Point
005-178	Palmer Freed Dissociation Test (May also be used in Stereoscope)
005-179	Keystone Binocular Skills Test & Manual (May also be used in Stereoscope)
005-180	Supplemental NP Test, Sloan Type
005-184	Paramacular Suppression Test (May also be used in Stereoscope)
005-185	Universal Driver Test
005-189	Retinal Rivalry Test (May also be used in Stereoscope)
005-192	Contact Lens Card
005-193	Plus Lens Acuity Target
005-194	DB-2B Left Eye, DB-3B Right Eye, DB-1B Both Eyes Acuity Test
005-200	Illinois College of Optometry Test
005-202	Pacific College of Optometry Test
005-203	Special Far Point Peek-A-Boo Acuity Test
005-204	Physicians Visual Rating Tests
005-209	Far Point Sloan Letter Acuity Test
005-210	Sloan Intermediate Test
005-213	Special NP Peek-A-Boo Acuity Target
005-218	Combined Phoria FP
005-219	Combined Phoria NP
005-220	Sloan Digit Acuity FP
005-221	Sloan Digit Acuity NP

KEYSTONE VIEW TELEBINOCULAR AND ATTACHMENTS

A pioneer in the vision testing field. Offers unmatched testing and training versatility, including diagnosis, visual skills screening and vision training.

1103S	Visual Survey Telebinocular-Standard
728-000	Attachment: Perimeter Attachment
728-030	Pilot Plus Visual Survey Telebinocular
728-044	Attachment: Plus Lens Attachment, +1.50, Includes 5193 DB-18 Plus Lens Acuity Target
728-045	Attachment: Plus Lens Attachment, +1.75D, Includes 5193 DB-18 Plus Lens Acuity Target
728-046	Attachment: Plus Lens Attachment, +2.00D Includes 5193 DB-18 Plus Lens Acuity Target
728-047	Attachment: Plus Lens Attachment, +2.25D Includes 5193 DB-18 Plus Lens Acuity Target
728-053	Attachment: Plus Lens Attachment, +2.50D Includes 5193 DB-18 Plus Lens Acuity Target
728-056	Ophthalmic Telebinocular
728-067	Driver Screening Telebinocular
728-935	Industrial Telebinocular

KEYSTONE VIEW EQUIPMENT

004-016	Colored Circles-Transparent, 10/pk
004-112	Eccentric Circles-Transparent, 20/pk
297-337	Red/Green Filters Elastic Band, 1pr
297-449	Sherman Playing Cards (2 decks), 2 deck/pk
725-000	Correct-Eye Scope
725-016	Stereo Reader
725-019	Closed Stereoscope - Hand Held
725-025	Delacato Mask For Stereo Reader
725-026	Open (725-030) Stereoscope w/pedestal mount
725-030	Open Stereoscope - Hand Held
725-040	Closed (725-019) Stereoscope w/pedestal mount
725-090	Spectacles Red/Green Filters, 1pr
725-330	Franzblau Multi-Sensory Coordinator (includes tracing pattern forms)
725-336	Van Orden Trainer w/aux lens holder
725-338	Complete Franzblau Red Rock System
725-366	Cheiroscopic Attachment for 725-374 & 725-330
725-374	Light Box Transparency Viewer
725-395	Colored Circles-Opaque Card, 50/pk
725-396	Eccentric Circles-Opaque, 25pr/pk
725-924	Red Rock Face Plate For Use w/725-000 or 725-330
728-288	Sherman TV Training Kit

KEYSTONE VIEW PARTS & ACCESSORIES

297-390	VS & DVS Carrying Case
297-504	VS II & DSV III Lamps GE400, 10/pk
297-800	VS II & DSV III 1A, 250V Fuse, 5/pk
298-410	Lamps For Telebinocular, Transparency Viewer, Correct-Eye Scope, Van Orden Trainer, Franzblau Multi-Sensor Coordinator, 25 Watt, 2/pk
298-808	Telebinocular Large Dust Cover
298-971	DVS II Lamps, 10/pk
725-017	Pedestal Base for Hand Held Open Stereoscope
725-087	Pedestal Base For Hand Held Closed Stereoscope
725-323	Plastic Faceplate w/clip for Franzblau
725-331	Sensor Pen Assembly For Franzblau Multi-Sensor Coordinator, 1 pair
728-055	Auxiliary Lens Holders, 1 pair
728-241	Occluders, Flip Type For Telebinocular, 1 pair
728-310	Telebinocular Carrying Case
728-391	Wooden Pointers, 20/set
815-249	DVS II Hand Control
815-278	DVS II Transformer
815-499	VS II & DVS III Transformer
815-503	VS II 14-Pin Hand Control
815-503-1	VS II 16-Pin Hand Control
815-754	DVS III 14-Pin Hand Control
815-754-1	DVS III 16-Pin Hand Control
818-062	CVS-I / VS-IV / DVS-IV/GT Target LED
818-066	DVS-GT LED Glare Lamp
818-114	CVS-I / VS-IV Hand Control
818-140	DVS IV - 110V Transformer
818-145	DVS IV - 220V Transformer
818-202	DVS-GT Hand Control
818-470	CVS-I Series 12VDC Universal PS
818-720	DVS-IV Hand Control

KEYSTONE VIEW RECORD FORMS

5540	Vinyl Matrix for OVS Form #8 SPECIAL ORDER: Minimum Order 10 Packs, 3 pd/pk
727-226	VS II / VS IV Universal Target Record Form, 3 pd/pk
735-145	OVS Rapid Screening Record Form, 3 pd/pk
735-236	Orthoscope Refractions Form OTR (5" x 8"), 3 pd/pk
735-237	Professional Performance Test Form (8 1/2" x 11"), 3 pd/pk
735-269	OVS Test Record Form #8, 3 pd/pk
735-270	OVS Test Record Form #9, 3 pd/pk
735-272	Doctor's Record Form #3A Visual Skills Profile (8 1/2" x 11"), 3 pd/pk
735-272-1	Doctor's Record Form #3A Visual Skills Profile (5" x 8"), 3 pd/pk
735-273	Driver Test Set Record Form #10, 3 pd/pk
735-275	Visual Survey Short Test Form #4, 3 pd/pk
735-277	Delacato Test Summary Sheets, 3 pd/pk
735-278	Leavell Test Sheets, 3 pd/pk
735-287	Visual Tracing Test Form Groffman, 3 pd/pk
735-384	Industrial Short Test Form #12, 3 pd/pk
735-385	Professional Performance Test Form (5" x 8"), 3 pd/pk
735-386	Doctor's Form #2A VS (8 1/2" x 11"), 3 pd/pk
735-386-1	Doctor's Form #2A VS (5" x 8"), 3 pd/pk
735-388	Skills Test Form DG-P, 3 pd/pk
735-389	Skills Test Form DG-1 / DG-2, 3 pd/pk
735-390	Skills Test Form DG-3 / DG-4, 3 pd/pk
735-394	Spache Binocular Reading Test Record Form, 3 pd/pk
735-397	Vinyl Matrix Overlay for OVS Form #9 SPECIAL ORDER: Minimum Order 10 Packs, 3 pd/pk
735-399	Ready to Read Record Form #6, 3 pd/pk
735-400	Multi-Stereo Scoring Form, 3 pd/pk
735-401	Mark Diplopia Diagnostic Chart, 3 pd/pk
735-402	Visual Survey Parents Notify Form, 3 pd/pk
735-403	OVS Referral Form #66, 3 pd/pk
736-161	Stereometric Record Form M-DC, 3 pd/pk
736-386	School Survey Form #5B (2 Part Form) Use w/005-135, Follow by 005-100, 3 pd/pk
736-561	School Survey Form #5A, 3 pd/pk
736-698	Visual Skills Survey Form, Use w/005-135 / 005-100, 3 pd/pk
736-699	Peek-A-Boo Series Form (8 1/2" x 11"), 3 pd/pk
736-744	Employee Vision Record Form #9-A, 3 pd/pk
736-776	Modified Telebinocular Technique Form, 3 pd/pk
736-804	Peek-A-Book Series Form (5" x 8"), 3 pd/pk
736-832	Palmer-Freed Dissociation Form, 3 pd/pk
736-878	Pilot Vision Test Pre-Record Form, 3 pd/pk
737-017	Physician's Visual Rating Form, 3 pd/pk
737-025	DVS II / DVS III / DVS IV Record Forms, 3 pd/pk
737-025A	DVS-GT Standard Record Forms, 3 pd/pk
737-161	Walton MVS Form W-T, 3 pd/pk
737-225-1	VS II / VS IV Standard Record Form (Replaces 737-225A, 737-131, 737-225 & 737-254) English, 3 pd/pk
737-225B	VS II / VS IV Standard Target Record Form (2 Part Form) Spanish, 3 pd/pk
737-255	VS II / VS IV Universal Target Record Form (5" x 8"), 3 pd/pk
738-769	Aniseikonia Record Form R1, 3 pd/pk

KEYSTONE VIEW SERVICES AND SETS

006-100	Visual Skills Program 1 (Screen) Includes Instruction Manual, 728-056, 005-100, 735-272 & 735-272-1
006-101	Visual Skills Program 1 (Diagnose) Includes Instruction Manual, Rating Cards, 728-056, 005-100 & 735-237
006-103	Brock Technique Service, Includes 28 Large Plastic Slides & Lightbox
006-107	Combined Programs 1 & 2 Screening & Diagnostic Test Set
006-116	Hand-Eye Coordination Test Set, 3 test/set
006-129	Keystone 3-D Fusion Games Set, Includes: 20 Red, Green & Black 8 1/2" x 11" Picture Charts, Easy-to-Use Manual, Easel, Pointers, Record Forms & 1 Pair Red/Green Filters
006-130	Basic Binocular Test Set, Includes: 12 Special Targets Bound in Flip Chart Form, Non-Glare Plastic-Coated Cards, Red/Green Filters & 1 Pair of Pointers
006-169	Economy Color Book Package, Includes: 10 Books, 10 Red Crayons & 1 Pair Red/Clear Filters
6104	Complete Van Orden Program, Kit Includes: Keystone View Correct Eye Scope, The Van Orden Trainer, 16 Transparency & Opaque Slides, 1 Set of Split Cards, Red/Green Filters, 2 Pads of Van Orden "Star" Graphs w/VO-1 matrix, 2 Pads of Training Charts & Manual
6106	School Vision Screening Program: Includes 728-030 w/005-100 & 005-135, Instruction Manual, Familiarization Chart & 736-386
6193	Pediatric Vision Screening Service, Comes Equipped with the Visual Survey/School Telebinocular & Peek-A-Boo Test Set, Instruction Manual & 100 Peek-A-Boo Test Forms


**WELCH ALLYN SURESIGHT™ VISION SCREENER**

Screen with one-button touch from 14" away. Five-second, automatic test addresses compliance problems and makes testing easier. Minimal cooperation is required, making it ideal for use with young children. The lights and sounds engage the child's attention. Automatically measures and indicates abnormal readings. Addresses eye charts' sensitivity and specificity problems. Automatically screens for common vision problems, including near- and farsightedness, astigmatism, and anisometropia. Hand-held unit operates for three hours of continuous testing. Available battery-operated printer enables you to get printouts anywhere, anytime. Weighs two lbs. and features an ergonomic grip and replaceable, self-adjusting hand strap. Features optics and electronics shock-mounted inside housing to protect unit. The sealed unit helps keep out dirt. Screens babies, children, and adults. You can even test patients while they're wearing glasses or contact lenses.

14000	SureSight™ Vision Screener Set, Includes 110-130V Charging Transformer & Stand
-------	--

Accessories

05140-U	Carrying Case
14042	Thermal Printer
53600	Printer Paper, 24/cs



V50-1ANY



V52-1D

BURDICK VISION SYSTEM

The software includes: "Search by time" in full disclosure, R-R printed on all strips, Date, day, and time on all pages, Review workstation capability for Vision, Forms Review in Vision, French and Spanish language availability, Availability of a 19" flat panel monitor option for Vision and Vision Premier systems, Availability of Vision 5L recorders with software only configurations, Addition of Rxio to all domestic configurations, distribution of Rxio outside the US is prohibited.

V50-1ANY	Premier System Includes: English Language Software, Recorder NOT included
V51-1G	Premier ACCESSORY: English Language Software with PC Card, Single user license
V52-1D	Vision System with 92515 Holter Recorder, English Language Software
V53-1DG	Vision Software w/PC card interface. Vision 5L recorder and English language Single user license
V53-1G	Vision System ACCESSORY: English Language Software, PC Card Interface

BCI 3180 PULSE OXIMETER

Provides an exceptional value for everyday oximetry needs. This compact, stand-alone monitor delivers fast, reliable SPO₂, pulse rate and pulse strength measurements, complete with a plethysmographic waveform and up to 30 hours of on-screen trending. Plus, the bright liquid crystal display, easy-to-use menu and audible and visible alarms enhance clinical efficiency and quality of care for patients ranging from neonate to adult.

3180	3180 Pulse Oximeter
------	---------------------

BCI AUTOCORR™ DIGITAL PULSE OXIMETER

Uses refined noise-reducing hardware and a unique patented software algorithm. Instead of relying on a voltage threshold to determine a pulse, it analyzes a digitized signal at a rate of 120 times per second, and compares it to previous pulse data. These pulses are constantly stored and updated, and used as dynamic templates to determine the validity or invalidity for newly-encountered pulsatile signals. Patented Serial Autocorrelation (SAC) technology provides reliable monitoring.

3304	Autocorr® Pulse Oximeter
3304LS	Autocorr® Pulse Oximeter w/loud alarms
3404-000	Autocorr® Plus SpO ₂ Pulse Oximeter
3404-001	Autocorr® Plus SpO ₂ Pulse Oximeter/ECG/Respiration

BCI DIGIT™ FINGER OXIMETER

Combines monitor and sensor into one unit that provides fast, reliable SpO₂, pulse rate, and pulse strength measurements on patients from pediatric to adult. Ideal for use within emergency settings, in hospital or clinical environments, or for home use. Portable, lightweight pocket-size monitor. One-button keypad. Large LED display. Pulse-strength display bar indicates relative perfusion at the measurement site. Auto power shutdown after 8 seconds. Low battery indicator flashes when 30 minutes of battery time remains.

3420G	Finger Oximeter, Gray
3420Y	Finger Oximeter, Yellow



BCI FINGERPRINT® HAND-HELD PULSE OXIMETER

Portable pulse oximetry monitoring system with built-in printer includes: 3401 Monitor (with 3 rolls of paper); 1884 Operation Manual; 3418 Protective Rubber Boot; 3044 Oximetry Finger Sensor; (4) Alkaline AA Cell Batteries.

3401-000	FingerPrint® Hand-Held Pulse Oximeter
3403-000	FingerPrint® Sleep Hand-Held Pulse Oximeter/Summary Printer



BCI HAND-HELD PULSE OXIMETER

Portable pulse oximetry monitoring system includes: 3301 Monitor; 1818 Operation/Service Manual; (3) Alkaline C Cell Batteries.

3301	Hand-Held Pulse Oximeter, with Oximetry Finger Sensor
3301-EP	Hand-Held Pulse Oximeter, with 3078 Ear Probe



BCI HAND-HELD PULSE OXIMETER WITH ALARMS

Portable continuous or spot check pulse oximetry monitoring system with alarms and internal rechargeable battery includes: 3303 Monitor; 3311 Oximetry Cable, 5 ft.; 3044 Oximetry Finger Sensor; 1850 Clinician's Operation Manual; 1854 Home-Use instruction Book; 8210 AC Charger, 105-125V 60Hz; 3353 Protective Rubber Boot with Carrying Strap and Tilt Stand.

3303	Hand-Held Pulse Oximeter w/alarms
------	-----------------------------------

BCI MINICORR DIGITAL HAND-HELD PULSE OXIMETER

Designed to provide fast, reliable measurements on challenging patients. The device rejects spontaneous motion and other artifacts. Informs the operator of excess artifacts and improper sensor positioning and processes low perfusion signals. Designed for spot-checking or bedside monitoring. This portable unit is ideal for use in home care and clinical environments, as well as during sleep screening. Easy to use, 5-key operation. Infrared link to printer. Battery or charger. Boot is available. System includes 3402 Monitor; 3311 Oximetry Cable, 5 ft.; 3044 Oximetry Finger Sensor; 1896 Operation Manual.

3402	MiniCorr Digital Hand-Held Pulse Oximeter
------	---

BCI PULSE OXIMETER ACCESSORIES

1100	Airway Adapter, Straight w/o filter (8200, 8400), 10/pkg	3107	ECG Shielded Leads, 40", Safety Type Wires, 3-Lead (3404, 9100, 9200), 3/pkg
1101	Adult O2 Nasal Cannula + Sample Line, (8200, 8400) 10/pkg	3134	Neonate Attachment Tape, 50/pkg
1102	Pediatric O2 Nasal Cannula + Sample Line (8200, 8400), 10/pkg	3135	Infant Attachment Tape, 50/pkg
1105	Airway Adapter, Straight w/filter (8400), 10/pkg	3136	Neonate Attachment Tape, 100/pkg
1121	Adult Oral/Nasal CO2 Cannula + Sample Line (8200, 8400), 10/pkg	3137	Infant Attachment Tape, 100/pkg
1123	Adult Oral/Nasal CO2 Sample Line (8200, 8400), 10/pkg	3138	Universal "Y" Posey Wrap Cover, 10/pkg
1124	Pediatric Oral/Nasal CO2 Sample Line (8200, 8400), 10/pkg	3178	Pediatric Finger Sensor
1129	Adult Nasal CO2 Sample Line (8200, 8400), 10/pkg	3310	Seiko Printer Model #DPU-201GSU (3301, 3303, 8200)
1130	Pediatric Nasal CO2 Sample Line (8200, 8400), 10/pkg	3311	Oximetry Cable, 5 ft (3301, 3303, 3401, 3402, 8400)
1131	Infant Nasal CO2 Sample Line (8200, 8400), 10/pkg	3311L	Oximetry Extension Cable, 15 ft (3403, 9004-051)
1140	15 ft CO2 Sample Line Extension (8200, 8400), 10/pkg	3311N	Nellcor Oximetry Extension Cable 4 ft
1151	Airway Adapter, Pediatric w/o filter (8200, 8400), 10/pkg	3313	Printer Interface Cable, 1 ft (3301)
1152	Airway Adapter, Pediatric w/filter (8200, 8400), 10/pkg	3315	3301 Carrying Case w/belt clip/shoulder strap, No Logo
1178	Moisture Trap (9004), 5/pkg	3318	3301 Protective Rubber Boot w/carrying strap
1179	External CO2 Absorber (9004), 2/pkg	3322	3303 Carrying Case w/belt clip/shoulder strap
1187	Moisture Trap Adapter (1178)	3339	PC Adapter Cable (3301, 3303, 3401, 3402, 8200)
1300	Adult Disposable Oximetry Finger Sensor, 10/bx	3350	Printer Interface Cable
1301	Pediatric Disposable Oximetry Finger Sensor, 10/bx	3351	Analog Output Adapter (3304)
1302	Neonate Disposable Oximetry Finger Sensor, 10/bx	3352NC	Digital Output Adapter (Nurse Call - Closed Circuit) (3304, 3404)
1303	Infant Disposable Oximetry Finger Sensor, 10/bx	3352NO	Digital Output Adapter (Nurse Call - Open Circuit) (3304, 3404)
1310	D.O.T. Reusable Oximetry Sensor, 3/bx	3353	3303 Protective Rubber Boot w/carrying strap/tilt stand
1311	D.O.T. Sensor Replacement Tape, 24/pkg	3354	Pole Mount Bracket (3303, 8200)
1606	Oximeter/ECG Simulator (3303, 8400)	3356	Clarity® Series IV Pole Bracket (3304, 3404, 6004, 9004)
1606HH	Oximeter Simulator & 5 ft Cable (3301, 3301, 3401, 3402)	3361	Printer Interface Cable (3304)
1611	115V AC Power Supply (3402)	3362	PC Interface (3304/3404)/Graphics Printer Interface (6004)
1614	AC Charger, 115V, 60Hz (6004, 9004)	3365	Printer Interface Cable (9004)
1617	105-125V AC Power Supply (3404, 3408)	3366	PC Interface (9004)/Graphics Printer Interface (9004)
1621	Oximetry Data Management Program	3369	Clarity® Series Carrying Case w/shoulder strap (3304, 3404, 6004)
1818	3301 Operation/Service Manual	3371	Graphics Printer Interface (3404)
1818R	3301 Operation/Service Manual	3389	Interface Cable (1621)
1835	8200 Operation Manual	3406	ECG Shielded Safety Type Cable, 5-Lead (3404, 9200)
1836	8200 Service Manuals	3407	ECG Shielded Leads, 40" Safety Type Wires, 5-Lead (3404, 9200), 5/pkg
1850	3303 Clinician Operation Manual	3408	Graphics Printer
1851	3303 Service Manual	3409	Printer Paper (3402, 8400), 4/bx
1854	3303 Home-Use Instruction Book	3410A	DC Power Cable (3404, 3408)
1857	Autocorr® Oximeter Clinician Operation Manual (3304)	3414	Analog Output Cable (3404)
1858	Autocorr® Oximeter Service Manual (3304)	3415	3401 FingerPrint Carrying Case
1861	Autocorr® Oximeter Home-Use Instruction Book (3304)	3418	3401 FingerPrint Protective Rubber Boot w/carrying strap
1872E	Capnocheck® Plus/Capnocheck Sleep Operation Manual (9004)	3419	MiniCorr Oximeter Protective Rubber Boot w/carrying strap (3402)
1873	Capnocheck® Plus/Capnocheck Sleep Service Manual (9004)	3426	Digit® Lanyard (3420)
1876E	Mini-Torr Plus® NIBP Monitor Operation Manual (6004)	3427	Digit® Carrying Case (3420)
1877	Mini-Torr Plus® NIBP Monitor Service Manual (6004)	3444	Comfort-Clip® Oximetry Finger Sensor (3303, 8400)
1884	3401 FingerPrint® Operation Manual	5093	ETCO2 Calibration Gas (8200, 8400)
1884SS	FingerPrint® Sleep Oximeter Operation Manual (3403)	6007	Mini-Torr Plus® NIBP Monitor Printer Upgrade Kit (6004)
1884SSHU	FingerPrint® Sleep Oximeter Home-Use Instruction Note (3403)	6008	Mini-Torr Plus® NIBP Monitor Oximetry Upgrade Kit (6004)
1885	3401 FingerPrint® Service Manual	6012	3401 FingerPrint® Printer Paper, 4/pkg
1885SS	FingerPrint® Sleep Oximeter Service Manual (3403)	6014	Electronic Thermometry Temperature Probe Covers (6004), 1000/ctn
1886	Advisor® Vital Signs Monitor Operation Manual (9200)	6015L	Electronic Thermometry Oral Temperature Probe w/9 ft cable (6004)
1887	Advisor® Vital Signs Monitor Service Manual (9200)	6015R	Electronic Thermometry Rectal Temperature Probe w/9 ft cable (6004)
1888	Autocorr® Plus Oximeter/ECG Monitor Clinician Operation Manual (3404)	6016	Electronic Thermometry Calibration Key (6004)
1889	Autocorr® Plus Oximeter/ECG Monitor Service Manual (3404)	6110	NIBP Cuff Supply Line (Luer-Lock), 10 ft (may need Transition Supply Line)
1892	Graphics Printer Operation Manual (3408)	6110S	NIBP Cuff Supply Line (Luer-Lock), 5 ft (may need Transition Supply Line)
1893	Graphics Printer Service Manual (3408)	6111	Printer Paper (9100, 9200), 4/bx
1894	Capnocheck® II Monitor Operation Manual (8400)	6120A	Adult (25.3-34.3cm) Non-Disposable NIBP Cuff
1895	Capnocheck® II Monitor Service Manual (8400)	6120C	Child (15.8-21.3cm) Non-Disposable NIBP Cuff
1896	MiniCorr Pulse Oximeter Operation Manual (3402)	6120I	Infant (9.8-13.3cm) Non-Disposable NIBP Cuff
1897	MiniCorr Pulse Oximeter Service Manual (3402)	6120LA	Large Adult (32.1-43.4cm) Non-Disposable NIBP Cuff
1898	MiniCorr Pulse Oximeter Home Use Instruction Book (3402)	6120N	Newborn (7.7-10.5cm) Non-Disposable NIBP Cuff
1900	Digit® Operation Manual (3420)	6120SA	Small Adult (20.0-27.0cm) Non-Disposable NIBP Cuff
1904	3180 Oximeter Clinician Operation Manual	6120SC	Small Child (12.4-16.8cm) Non-Disposable NIBP Cuff
1905	3180 Oximeter Service Manual	6121	10 ft NIBP Cuff Transition Supply Line
3005	Power Cord (3180, 9100, 9200)	6121S	5 ft NIBP Cuff Transition Supply Line
3025	Infant Oximetry Wrap Sensor	6122	10 ft NIBP Cuff Supply Line (Quick Connect)
3026	Neonate Oximetry Wrap Sensor	6122S	5 ft NIBP Cuff Supply Line (Quick Connect)
3043	Universal "Y" Oximetry Sensor	6123	5 ft NIBP Cuff Neonate Transition Supply Line
3044	Standard Oximetry Finger Sensor	6124	5 ft NIBP Cuff Neonate Supply Line (Quick Connect)
3049	Adhesive Strip Tape, 40/pkg	6140	Neonate (3-6cm) Single Use NIBP Cuff, 3/pkg
3078	Ear Oximetry Probe Sensor	6141	Neonate (4-8cm) Single Use NIBP Cuff, 3/pkg
3106	ECG Shielded Safety Type Cable, 3-Lead (3404, 9100, 9200)	6142	Neonate (6-11cm) Single Use NIBP Cuff, 3/pkg
		6143	Neonate (7-13cm) Single Use NIBP Cuff, 3/pkg

BCI PULSE OXIMETER ACCESSORIES (Continued)

- 6144 Neonate (8-15cm) Single Use NIBP Cuff, 3/pkg
- 6145 Neonate Variety Pack Single Use NIBP Cuff (2 of each size), 10/pkg
- 7108 Clarity® Series/9200 Roll Stand w/mounting bracket & basket (3304, 3404, 6004, 9004)
- 8044 8 ft CO2 Sample Line (8200, 8400), 10/pkg
- 8061 Calibration Gas Flow Regulator w/gauge (8200, 8400)
- 8208 Gas Manifold Filters (8200, 8400), 10/pkg
- 8210 AC Charger, 105-125V, 60Hz (3303, 8200)
- 8211 4 ft CO2 Sample Line (8200, 8400), 10/pkg
- 8213 Printer Interface Cable, 1 ft (8200)
- 8214 Patient Attachment Kit: 1100, 8211, 8208 (8200, 8400), 10/pkg
- 8217 Capnocheck Calibration Kit (8200, 8400)
- 8223 Calibration Adapter (9004)
- 8403 Battery Eliminator 105-15VAC, 60Hz (8400)
- 8404 Battery Charger (8400)
- 8408 Battery 7.4V Lithium-Ion Rechargeable (8400)
- 8409 Pole Mount Bracket (8400)
- 8411 Infrared Printer Model #HP82240B (3402, 8400)
- 8412 Capnocheck® II Monitor Protective Rubber Boot (8400)
- 8414 Capnocheck® II Capnograph/Oximeter Carrying Case (8400, 8401)
- 9007 Capnocheck® Plus/Capnocheck Sleep Oximetry Upgrade Kit (9004)
- 9014 3 Channel Analog Output Cable (polysomnography) (9004)
- 9015 4 Channel/Distance Analog Output Cable (polysomnography) (9004)
- 9047 Capnocheck® Plus Capnograph Filter Upgrade Kit (9004)
- 9048 External 3.5 Micron Filter (9004), 2/pkg
- 9189 FIO2 Adapter T (9004), 10/pkg
- 9190 FIO2 Sensor (9004)
- 9191 FIO2 Cable (9004)
- 9193T Capnocheck® Plus Capnograph FIO2 Upgrade Kit (9004)



BURDICK OXY 100 PULSE OXIMETER

Fast, accurate SpO2 levels whenever and wherever needed. Displays SpO2 level, pulse rate and pulse strength on an easy to read LED screen. Uses 3 "C" cell batteries and automatically switches off when not in use. Compact, portable and durable. Up to 99 individual tests can be stored for printout using the optional printer.

- 92521-1A Burdick Oxy 100 w/finger sensor
- 92522-1A Burdick Oxy 200 w/finger sensor



GE MEDICAL TUFFSAT® PULSE OXIMETER

Ideal oximetry spot-check solution for physicians' offices and clinics. The hand-held oximeter fits in your pocket or is conveniently worn using the included carrying case with belt clip. Includes a Plr® (relative perfusion index measurement) for easy quantification of signal quality. Available in four colors; yellow, blue, pink and white.

- 6051-0000-160 TuffSat® Yellow w/sensor
- 6051-0000-185 TuffSat® Blue w/sensor
- 6051-0000-186 TuffSat® Pink w/sensor
- 6051-0000-187 TuffSat® White w/sensor

GE MEDICAL MASIMO PULSE OXIMETRY SENSORS & CABLES

Accurate under patient motion or movement. Accurate for low perfusion patients. Latex free.

- 2009743-001 Masimo PCO8 Cable
- 2009745-001 Masimo Finger Sensor, Adult, Reusable
- 2002799-001 Masimo Finger Sensor, Pediatric, Reusable
- 2002800-001 Masimo Finger Sensor, Adult, Reusable

GE MEDICAL NELLCOR DISPOSABLE PULSE OXIMETRY SENSORS (PROBES)

Single patient use reduces potential for cross-contamination. Sensor placement indexes printed on wrap for easy alignment. Latex free.

- 407705-001 Nellcor OxiSensor Neonatal Disposable Probes, Patient Range
- 407705-002 Nellcor OxiSensor Infant Disposable Probes, Digital Probe, Patient Range 3-20 KG, 24/cs
- 407705-003 Nellcor OxiSensor Pediatric Disposable Probes, Digital Probe, Patient Range 10-50 KG, 24/cs
- 407705-004 Nellcor OxiSensor Adult Disposable Probes, Digital Probe, Patient Range > 30 KG, 24/cs

GE MEDICAL NELLCOR REUSABLE PULSE OXIMETRY SENSORS & CABLES

Hinged shell accommodates a wide range of digits. Soft vinyl padding holds probe in place. Latex free.

- 2008773-001 Nellcor Interface Cable, OxiSmart DOC10 Cable
- 407705-006 Nellcor DuraSensor Reusable Finger Probe
- 2006644-001 12 ft Nellcor Interface Cable OxiSmart
- 4033CAX Intermediate Cables For Nellcor Sensors
- DS100A Nellcore Reusable SpO2 Fingercip Sensor

GE MEDICAL NONIN PULSE OXIMETER SENSORS & CABLES

Latex free.

- 412926-001 8" Apex S Nonin Pulse Oximetry Interface Cable
- 412926-002 12" Apex S Nonin Pulse Oximetry Interface Cable
- 418497-001 8" ApexPro Nonin Pulse Oximetry Interface Cable
- 420970-001 Reusable Adult Articulated Finger Clip, 36", Patient Range > 30 Kg
- 420970-004 Reusable Ear Finger Clip, 18", Patient Range > 40 Kg

GE MEDICAL OHMEDA PULSE OXIMETRY SENSORS & CABLES

Latex free.

- 9587-002 Ohmeda FingerClip, Reusable Finger Probe, 12 ft Integrated Cable
- 9587-011 Ohmeda EarProbe, Adult Reusable Ear Probe, 8 ft Integrated Cable
- 9587-022 Ohmeda Flex II Probe, Reusable Wrap Probe, 12 ft Integrated Cable
- 9587-208 Ohmeda Adult/Pediatric Oxy-Tip Disposable Probes, Use w/9587-209 Interconnect Cable, 25/cs
- 9587-209 Ohmeda Interconnect Cable For Oxy-Tip Disposable Probes, 12 ft
- 9587-210 Ohmeda Neonatal Oxy-Tip Disposable Probes, Use w/9587-209 Interconnect Cable, 10/cs

GE MEDICAL SENSORMEDICS PULSE OXIMETRY PROBES & CABLES

Latex free.

008731	Sensormedics Interconnect Cable
009084	OMNI-Sat Finger Sensor
009120	OMNI-Sat Ear Clip Sensor
009130	Neonatal Sensor, Semi-Disposable
009136	Adult Sensor, Semi-Disposable

MALLINCKRODT N-45 HANDHELD PULSE OXIMETER

Compatible with the complete line of OxiMax® pulse oximetry sensors, so clinicians have the convenience of Nellcor's most advanced sensor technology with a lightweight, battery-operated, handheld device. Designed with an artifact-reducing algorithm to deliver highly accurate SpO₂ and pulse rate values, even during challenging monitoring conditions. The N-45 oximeter can be used for spot-checks or continuous monitoring. It offers the security of audible and visible alarms, and features 14 or 30 hour battery life, depending on battery type. The 24-hour downloadable trend memory enhances patient management and enables overnight oximetry studies in the hospital or home. The N-45 oximeter is compatible with ProFox software.

N45	N-45 Handheld Pulse Oximeter, Includes: 9 V Lithium Battery, Cable Clip & Operator's Manual
-----	---

Accessories

N45ACADPT	N-45 110-230 V A/C Adapter
N45BOOTCORAL	N-45 Protective Boot, Coral
N45BOOTPURPLE	N-45 Protective Boot, Purple
N45DTCBL	N-45 Data Transfer Cable



MALLINCKRODT N20 SERIES HANDHELD PULSE OXIMETERS

Rugged, durable construction with protective boot. 32 to 37 hours long battery life. Easy-to-see day/night display. Three-year warranty. Comes with four C alkaline batteries, a finger clip sensor, a sensor extension cable, protective boot and moisture-resistant jacket. All models have a constant-pitch beep tone that sounds with each detected pulse. Displays oxygen saturation, pulse rate, vertical bar graph pulse amplitude indicator/battery meter, pulse search indicator, low battery indicator.

N20PA	N-20PA Handheld Pulse Oximeter, Includes: 4-C Size Batteries, 1-Finger-Clip Sensor (DS-1000A), DEC-4 Sensor Extension Cable, Integral Thermal Printer w/12 rolls paper (901624), Protective Boot (024605), Moisture-Resistant Jacket & Manual
N20PAV	Pulse Oximeter w/continuous operation, Printer & Alarms (Vet Only)
N20PE	Pulse Oximeter w/printer, 1/cs

Accessories

024603	Protective boot for N-20, N-20E & N-20V
024605	Protective Boot, N-20P Series
024607	Soft Carrying Case for N-20 Series Pulse Oximeters
029043	Dual Soft Carrying Case for Pulse Oximeters
901624	Printer Paper, 12 rolls/box




MALLINCKRODT NPB40 HANDHELD PULSE OXIMETER

Both compact and affordable. It is simple to operate and provides spot-check monitoring of oxygen saturation and pulse rate. The durable monitor is small enough to fit in your pocket. Designed for monitoring neonatal through adult patients in almost any setting, from hospital to home. Features an easy-to-read display and printing capabilities via optional Hewlett-Packard printer. Operates on AA batteries that provide approximately 19 hours of battery life.

NPB40	NPB40 Handheld Pulse Oximeter, 1/bx
NPB40V	NPB40 Handheld Pulse Oximeter (Vet Only)

Accessories

035183	Soft Carrying Case, 1/bx
035328	Protective Boot, 1/bx
902114	Printer, HP 82240B HP Infrared Printer
902115	Printer Paper, HP 82175, 6/bx




MALLINCKRODT OXIMAX®/OXISENSOR® II ADHESIVE SENSORS

ASP-3N	Oxisensor II Neonatal Sensor, Assorted, 1/bg
ASP3	Oxisensor II Adult Sensor, Assorted, 1/bg
ASP3P	Oxisensor II Pediatric Sensor, Assorted, 1/bg
D20	Oxisensor II Sensor, Pediatric, 24/cs
D25	Oxisensor II Sensor, Adult, 24/cs
D25L	Oxisensor II Adult Sensor w/36" cable, 24/cs
I20	Oxisensor II Sensor, Infant, 24/cs
MAXA	Oximax Sensor, Adult, 24/cs
MAXFAST	Oximax Sensor, Adult Forehead, 24/cs
MAXI	Oximax Sensor, Infant, 24/cs
MAXN	Oximax Sensor, Neonate/Adult, 24/cs
MAXP	Oximax Sensor, Pediatric, 24/cs
N25	Oxisensor II Sensor, Neonate/Adult, 24/cs
R15	Oxisensor II Sensor, Adult Nasal, 24/cs

Accessories

DEC4	Oximax 4 ft Extension Cable, 1/bx
DEC8	OxiMax 8 ft Sensor Extension Cable, 1/bx
DOC10	Oximax 10 ft Sensor/Patient Cable



D-YS



DS100A-1

MALLINCKRODT REUSABLE SENSORS

For spot checks or short-term monitoring. Non-sterile and easy-to-use design. Built-in shielding protects signal from electronic noise. High quality LEDs maximize tracking capabilities. Sensor head design optimizes signal and shields detector from ambient light.

A	OxiCliq Sensor, Adult, 24/cs
D-YS	Dura-Y Multisite Sensor, 1/bx
D-YSE	Dura-Y Ear clip For Use with the Dura-Y Sensor
D-YSPD	Pedi Check Pediatric Spot-Check Clip for the Dura-y Sensor, 1/bg
DS100A-1	DuraSensor Adult Finger Clip Sensor, 1/bx
I	OxiCliq Sensor, Infant, 24/cs
N	OxiCliq Sensor, Neonatal/Adult, 24/cs
OXI-A/N	Oxiband Adult/Neonatal Sensor, 1/bx
OXI-P/I	Oxiband Pediatric/Infant Sensor, 1/bx
P	OxiCliq Sensor, Pediatric, 24/cs
PDSLV	Pedi Check Replacement Sleeve
RS10	Reflectance Sensor, 6/cs
V-SAT	Veterinary Sensor
VSC-L	Veterinary Sensor Clip, Large
VSC-S	Veterinary Sensor Clip, Small

Accessories

ADH-A/N	Adhesive Wrap for Reusable Sensors, Adult/Neonatal, 100/bx
ADH-P/I	Adhesive Wrap for Reusable Sensors, Pediatric Infant, 100/bx
FOAM A/N	Foam Wrap for Reusable Sensors, Adult/Neonatal, 12/pk
FOAM P/I	Foam Wrap for Reusable Sensors, Pediatric/Infant, 12/pk
OC-3	OxiCliq Sensor Cable, 3 ft, 1/bx
POSEY	Velcro Wrap for Reusable Sensors, 12/pk



PEPPER MEDICAL PULSE OXIMETRY SENSORS

Cost effective alternative pulse oximetry sensors that will plug ("Plug and Play") directly into your customer's existing monitors. Compatible with 8 different manufactures including AvalonAire, NPB, BCI, and Ohmeda. Available in 4 sizes: adult pediatric, infant, and neonate. 2 tape styles available: Microfoam, Transpore, and Medplast. FDA approved.

2211	Disposable Pulse Oximetry Sensor, 24/bx
------	---



BCI CAPNOCHECK® II CAPNOGRAPH

Hand-held capnograph is a full-featured monitor providing ETCO₂, Inspired CO₂, Respiration Rate SpO₂, and Heart Rate measurements. Ideal for use in and out of the hospital environment. A back-light feature on the display enables easy viewing in difficult lighting. Battery operated with an optional battery bypass offers flexibility for stationary and transport applications. Includes: 8400 Monitor, 1894 Operation Manual, 3044 Oximetry Finger Sensor, 8412 Protective Rubber Boot, 8214 Patient Attachment Kit includes the following: 8208 Capnocheck Filters (10), 8211 Capnocheck Sample Lines, 4 ft (10), 1100 Airway Adapter, straight without filter (10).

8400	Capnograph/Oximeter w/battery & charger
8400E	Capnograph/Oximeter w/battery eliminator
8400NP	Capnograph/Oximeter w/o Power Supply
8401	Capnograph w/battery & charger
8401E	Capnograph w/battery eliminator
8401NP	Capnograph w/o power supply

BCI CAPNOCHECK® PLUS CAPNOGRAPH

Offers crisp, clear waveforms with large numerics for maximum visibility. High resolution vacuum fluorescent display. Simplified menu structure. Monitoring of both intubated and non-intubated patients. Flexibility to incorporate SPO₂ and FIO₂ technology. Patented moisture removal system. Four-channel analog output for polysomnography interfacing.

9004-000	Capnocheck® Plus Capnograph
9004-001	Capnocheck® Plus Capnograph/Oximeter
9004-002	Capnocheck® Plus Capnograph/FIO ₂
9004-003	Capnocheck® Plus Capnograph/Oximeter/FIO ₂

BCI CAPNOCHECK® SLEEP CAPNOGRAPH/OXIMETER

Specially designed for sleep screenings and lab-based sleep studies. Its 4-channel analog output capability lets you select ETCO₂ numeric value, ETCO₂ waveform, inspired CO₂, respiratory rate, heart rate, as well as SPO₂ numeric waveform displays. Patented Serial Autocorrelation (SAC) pulse oximetry technology provides reliable SPO₂ measurements.

9004-050	Capnocheck® Sleep Capnograph
9004-051	Capnocheck® Sleep Capnograph/Oximeter



MALLINCKRODT HANDHELD CAPNOGRAPH/PULSE OXIMETER

Delivers two highly advanced monitoring technologies in a convenient, portable monitor: Microstream® CO₂ technology and Nellcor's latest OxiMax pulse oximetry technology. Microstream technology and broad selection of breath sampling accessories ensure effective capnography monitoring for intubated and nonintubated patients alike. A low sample flow rate of 50 mL/min allows monitoring on the widest range of patients, from neonates to adults. OxiMax pulse oximetry technology provides exceptional performance, even during challenging monitoring conditions. The NPB-75 monitor is compatible with Nellcor's growing line of OxiMax sensors, including the Max-Fast™ Forehead Sensor and SoftCare™ Nonadhesive Sensors.

- NPB70 NPB70 Handheld Capnograph Includes: AC Adapter, Filterline Accessory Kit, Rechargeable NiMH Battery Pack, Operator's Manual & Quick Guide
- NPB75MAX Capnograph/Pulse Oximeter, Includes: AC Adapter, Accessory Kit, MAXPAC Oximetry Sensor Asst. Pack, DEC-4 Sensor Extension Cable, Rechargeable NiMH Battery Pack, Manual & Quick Guide

Accessories

- 047165 Carrying Case
- 047166 Extra Battery Pack Pouch
- 047184 Pole Mount Clamp
- 047203 External Battery Charger, 110-250 VAC
- 048127 Communication Adapter Kit: Printer/Computer Interface



BURDICK EK-10 SINGLE-CHANNEL ELECTROCARDIOGRAPH

Combines the cost advantages of a single-channel electrocardiograph with the increased productivity of automatic operation. Features a thermal array printer, liquid crystal display, 'soft-touch' membrane pad and 50 ECG acquisition/print on battery (optional). Compact and portable.

- 101070 EK-10 Single-Channel Electrocardiograph
- Accessories**
- 007175 1V Output Cable
 - 086027 Service Manual For EK-10
 - 086063 Operators Manual For EK-10

BURDICK SINGLE CHANNEL ECG MOUNTS

- 007036 ECG Mounts, Elite Folder, Pressure-Sensitive, Full Information, 100/bx, 5 bx/cs
- 007037 ECG Mounts, Card, Pressure-Sensitive, EK10 or Abridged Elite, 100/bx, 5 bx/cs
- 007857 ECG Mounts, Card, Pressure-Sensitive, 12-Lead, 100/bx, 5 bx/cs
- 007909 ECG Mounts, Card, Pressure-Sensitive, 3 Sec., 12-Lead, 100/bx, 5 bx/cs
- 007911 ECG Mounts, Folder, Pressure-Sensitive, 3 Sec., 12-Lead, 100/bx, 5 bx/cs



BURDICK ELITE II INTERPRETIVE ECG RECORDER

Simultaneously acquires 10 seconds of ECG data for all 12 leads. Prints interpretive reports in 4-channel format. Menu displays guide user through all functions. Stores up to 30 ECGs in memory. 11" x 13" x 3 1/4". 120/240 VAC, 50/60 Hz. Includes 10 lead accessory group and 12.5-V nickel cadmium battery. Compact and portable.

- 101158 Elite II, Single Channel Interpretative Electrocardiograph
- Accessories**
- 007697 Elite II Factory Installed Fax Option
 - 086136 Elite II Operators Manual
 - 086175 Elite II Service Manual
 - 086245 ECG Interpretation Criteria, GRI



010-1639-00



A30-1DIL01

BURDICK ATRIA™ 3000 RESTING ECG SYSTEM

Simple to use, compact and reliable. The lightweight and compact device operates either manually or automatically, printing 3-, 4-, 6- or 12-channel ECG tracings. The combination of a standard QWERTY keyboard and easy-to-navigate menu makes it easy to use. Features advanced functionality, including a stable baseline filter for greater accuracy, pacemaker algorithm, user-definable QTc calculations and Burdick's exclusive AccuPrint™ functionality, which restricts ECG printing if lead wires are not functioning properly. Variety of upgrade paths.

- A30-1DIL01 Atria™ 3000 ECG w/interpretation, extended storage 40 patients & no communication
- A30-1DKL01 Atria™ 3000 ECG w/measurements, 40 patient storage & no communication

Accessories

- 010-1639-00 Cart For The Atria 3000 ECG System
- 040-1457-00 Atria 3000 Kit Upgrade Interpretation
- 040-1458-00 Atria 3000 Kit Upgrade Extra Storage
- 040-1474-00 Atria 3000 Kit Upgrade Measurements



BURDICK ECLIPSE™ 850 ELECTROCARDIOGRAPHS

Multi-channel electrocardiograph with preview screen. Thermal printer produces crisp accurate 3, 4, 6 or 12 channel ECG reports and continuous rhythm reports on 8 1/2" paper. Acquires and prints on battery and/or AC power. Optional internal modem with fax and SCP-ECG communications. Includes: 10-Lead patient cable with 90 day warranty, power cord, one pad assurance 50 paper, cardiosens/ultra sensors & adapter clips, manual, handbook and video. Three year warranty on parts and labor.

- 92300-1AK09 Interpretive Eclipse 850i w/storage for 40 ECG
- 92300-1AL09 Non-Interpretive Eclipse 850 w/storage for 40 ECG
- 92300-1BKP09 Interpretive Eclipse 850i w/SCP-ECG & fax communications, internal modem & storage for 1 ECG
- 92300-1BLP09 Non-Interpretive Eclipse 850 w/SCP-ECG & fax communications, internal modem & storage for 1 ECG

BURDICK EQUIPMENT CARRYING CASE

- 007039 Soft Shelled Carrying Case for ELITE, EK-10 & Presto Spirometer #100231
- 007160 Soft Shell Carry Case For Use w/Extol LE & LEI
- 007732 Carrying Case For Eclipse 4
- 007771 Carry Case For Use w/Eclipse 8



BURDICK UNIVERSAL ECG

Portable and easy to use. Patient data is instantly and accurately obtained, viewing up to 12 channels. Unlimited data analysis including the ability to reanalyze data if system parameters are modified. Compatible with your favorite off-the-shelf PC, laptop, tablet PC or Compaq iPAQ.

- 92309-01 Universal ECG Non-Interpretive, AHA, English
- 92309-02 Universal ECG Interpretive, AHA, English



BURDICK ASSURANCE 50™ ECG RECORDING PAPER

Guaranteed to hold a trace for 50 years. Single and multi-channel formats. Special thermal paper resists effects of alcohol, sunlight, plastics and other common contaminants that degrades image quality on standard thermal papers.

- 716-0233-00 ECG Recording Paper, Red Grid Use w/ EK10, Elite, Elite II, Medic 5, SUB FOR 007957 & 007975, 100r/cs
- 716-0234-00 ECG Recording Paper, Z-Fold Pad Red Grid, Perforated No Header use w/ Eclipse LE, LEII, E550, E560, E600, ExTOL 200 sheets/pd, 10pd/cs
- 716-0235-00 ECG Recording Paper, Z-Fold Pad Red Grid, Perforated w/blank header use w/ Eclipse8, 800, 300 sheets/pd, 10pd/cs
- 716-0237-00 EKG Thermal Paper, 200 sheet/pad, 10 pad/cs
- 716-0239-00 ECG Recording Paper, Z-Fold Pad Red Grid, Perforated w/blank header use w/Eclipse 850, 300 sheets/pd, 10pd/cs

BURDICK HEARTLINE THERMAL PAPER

- 007868 Chemical/Thermal Heavyweight Paper, Perf, Multi-Channel w/blank header, Matrix, Red Grid/Black Trace, 200 shts/pd, 10 pd/cs
- 007940 Single Channel Thermal Paper, Grid, 45mm, Red Grid/Black Trace, 100'/roll, 100 r/cs
- 007941 Thermal Paper, Orange Grid, For Use w/EK4, EK5/5A, Ek8, E100, E200, 100r/cs
- 007944 Thermal Paper, Orange Grid, Perforated Header, For Use w/EK670 S/N 6000 & Above, EK7 S/N 7000 & Above, E310, E320, 20r/cs
- 007958 Single Channel Thermal Paper, Dot Matrix, 45mm, Red Grid/Black Trace, 150 ft/r, 80 r/cs (Sub for 007957)
- 007966 Thermal Paper, Z-Fold Pad, Red Grid, Perforated w/header, Use w/E350, 200 sheets/pd, 10pads/cs
- 007971 Thermal Paper, Red Grid, Perforated w/blank header, Use w/Eclipse 4/400, 200 sheets/pd, 40pd/cs
- 007979 Chemical/Thermal Paper, Perf, Multi-Channel, No Header, 8.5" x 5.5", Red Grid/Black Trace, 300 shts/pd, 10 pd/cs
- 007983 Chemical/Thermal Paper, Perf, Multi-Channel w/blank header, Matrix, Red Grid/Black Trace, 200 shts/pd, 10 pd/cs
- 007988 Chemical/Thermal Paper, Reinforced Perf, Multi-Channel, Blank Header, 8.5" x 5.5", Red Grid/Black Trace, 300/pd, 10 pd/cs



92400-01



B24-1AMP01

BURDICK QUEST® EXERCISE STRESS SYSTEM

Features easy, intuitive touch screen lets you focus on the patient. Performs both resting 12-lead ECGs and stress tests on one system. Accurate ECG data acquired at the patient eliminating non-cardiac electrical noise. Flexible programming means you choose the protocol. Generous storage, enhanced online data editing and several data management options. Full-range reporting offers over 15 different report configurations. Optional treadmill not included, See TM55 or TM65. Optional T-Wave Alternans helps predict Sudden Cardiac Death (SCD).

- 92400-01 TM55 Treadmill Option For Quest Exercise Stress System
- 214-1022-00 Stress Belt (Black) w/pouch and quick release buckle, 48"
- 214-1023-00 Belt Extender w/quick release buckle no pouch, 24"
- B24-1AMP01 Quest Exercise stress system, English language, with AHA accessories kit, TM55 treadmill and network communications option (120V), ea
- B24-1AP Quest System w/o treadmill. Includes standard US Thermal Paper and AHA Patient Cable, English Language (120V)
- B24-1AP01 Quest Exercise Stress System, English Language w/AHA accessories kit and TM55 Treadmill
- B24-1APT Quest System w/T-Wave Option w/o treadmill. Includes Standard US Thermal Paper and AHA Patient Cable, English Language and 120V
- B24-1APT01 Quest Exercise Stress System, English Language w/AHA accessories kit, TM55 Treadmill and T-Wave Option

Accessories

- 004115 Emergency Stop Switch For T600
- 007772 Patient Module
- 004219 Patient Data Acquisition Cable



BURDICK ELECTROCARDIOGRAPH ACCESSORIES

007001	Patient Cable, Resting, 10 Lead, Non-replaceable (E500)
007026	Cable, E600/ELITE/E500tw To Modem
007028	Interconnect Cable, E600/E500tw/ELITE To E-Scribe
007053	Junction Box For Modem On Stress Cart
007059	Interconnect Cable, Cart To Cart
007153	Plastic Cardiogram Mount Inserter
007157	Adult Limb Plate Sensor, 4 electrodes/pk
007158	Limb Strap For Sensor, Single Hole, 18", Adult
007159	Limb Strap For Sensor, 2 Holes, 20", Adult
007161	Infant Limb Plate Sensor, 4/pk
007162	Limb Strap For Sensor, Infant, Single Hole
007163	Welsh Suction Cup Sensor, Self-Retaining, Adult
007211	Patient Cable, Resting, 5 Lead, Non-replaceable (EK-5A)
007264	Alligator (Veterinary) Adapter Clip
007364	Welsh Suction Cup Sensor, Spring-Clip, Infant
007374	Stress Cable Belt, Black Buckle, 5 ft
007463	Lead Adapter For Snap End
007469	Lead Adapter For Grabber End
007517	Patient Cable, Stress, 10 Lead, Non-replaceable (ELITE, EK-10, E350, E350i, E550, E560)
007519	Patient Cable, Stress, Non-replaceable, EK-8, E100, E200, E310, E320, E600
007526	Stress Belt, 4 ft, White, Velcro
007550	Adult Limb Plate Sensor, 4mm, Set of 4
007557	How To Operate The Burdick EK-10
007558	A Guide To Good ECG Technique
007597	External Modem w/cable
007655	How To Operate The Burdick ELITE
007656	How To Operate The Burdick E350i
007679	How To Operate The Burdick E350
007690	Welsh Suction Cup Sensor 4mm, Set of 6
007704	Patient Cable, Resting, 10 Lead, Replaceable, Eclipse 4, 400, 8, 800, 850, Plus
007705	Patient Cable, Resting, 10 Lead, Replaceable, ELITE, E50, E350i, E550, E560, Eclipse LE/Lei
007706	Patient Cable, Stress, 10 Lead, Replaceable, E350, E350i, E560, Eclipse LE/Lei
007710	Patient Cable, Three Lead Wires (Replacement Kit) For 007704, 007704, Includes: one 20", one 30", one 40" Lead Wire & 10 Lead Identifiers
007711	Patient Cable, Three Lead Wires (Replacement Kit) For 007706, 007725, Includes: One 20", One 30", One 20" Lead Wire & 10 Lead Identifiers
007725	Patient Cable, Stress, 10 Lead, Replaceable, Eclipse 4, 400, 8, 800, 850, Plus
007726	Patient Cable, Three Lead Wires (Replacement Kit) For 007727, 007709, Includes: One 20", One 30", One 20" Lead Wire & 10 Lead Identifiers
007728	Patient Cable, Three Lead Wires (Replacement Kit) For 007707, 007708, Includes: One 20", One 30", One 40" Lead Wire & 10 Lead Identifiers
007751	EKS-95 Electronic Patient
007797	Lead Adapter, 0.158" (4mm) diameter to snap end
007861	Retaining Strap w/weight for 107155 Set
008003	How To Operate The Burdick Eclipse
010-1608-00	ECG Replacement Cables Interpretive
010-1609-00	ECG Replacement Cables Non-Interpretive
012-0700-00	Patient Cable, Resting, 10 Lead, Non-replaceable (E350, E350i, E550, E560)
012-0701-00	Patient Cable, Resting, Leads, Non-replaceable, (EK-8, E100, E200, E310, E320, E600)
012-0702-00	Patient Cable, Resting, 5 Lead, Non-Replaceable, EK10, Elite
012-0844-00	Nonreplaceable Lead ECG Cable, 10 Lead, AHA, Combs Located Between Lead Wires to Decrease Tangling, Use With Atria, Eclipse Premier & Eclipse 850
013-0150-00	30" Replacement Stress Leads, Green and Brown
047224	3-Lead Patient Cable
047236	5-Lead Patient Cable
047246	DC Input Cable (1V to 1mV)
047254	Patient Cable, Twist-Lock Connector For PA-101 Preamp
086061	Electrocardiography: A Better Way Workbook
097152	8 ft Power Cord- 3 Cond For EK-5A
097159	Power Cord For EK-7, EK670, DC-190, 10 ft
107067	Welsh Sensor Kit, Infant Plates/Straps & Bulbs
107068	Welsh Sensor Kit, Adult Plates/Straps, Infant Bulbs
107546	4mm Sensor Kit, Adult Plates/Straps, Infant Bulbs
147234	3-Lead Wire Kit For 047224 & 047254 Cables Only
147235	5-Lead Wire Kit For 047236 Cable Only, Snap
700-0277-00	USB to Serial Cable Adapter For Use w/Eclipse, CardioVault, Pyramid & AED Products
860062	Retaining Strap Only For 107155 Set
860403	Bulb Only For Welsh Infant Sensor (007364)
860404	Bulb Only For Welsh Adult Sensor (007163)
862278	Battery Pack For ELITE, EK-10

GE MEDICAL MAC 1200 STANDARD RESTING ECG

A 12-lead electrocardiograph with a full-size paper, alphanumeric keypad and waveform display. Optional Marquette® 12SL™ ECG analysis program, data storage, and bi-directional data transmission capability. AC line or battery operation from NiCd battery. All models listed include: Printed measurements and starter kit (Silver Mactrode® electrodes, grabber adapters, paper, patient cable, power cord).

MAC1200-AAX-CXX-XXX	MAC 1200 'Plus' Interpretive Electrocardiograph w/12SL ECG Analysis Program & Waveform Display
---------------------	--

Accessories

901168-001	MAC 1200 Office Cart
------------	----------------------

GE MEDICAL DEFIBRILLATOR SUPPLIES

009/00310/24V	Battery Pack, NICAD For 360 Defib
009/09956/USA	Charger Lead-Acid Battery USA
099/01270/000	Cable Vehicle Charger For Responder
210/34223/000	Battery For Data Event Card
414034-002	Data Event Card, 80 Minutes
850/40156/005	Adult Defibrillator Pads Responder, 5 pr/cs
850/40156/025	Adult Defibrillator Pads Responder, 25 pr/cs
850/40518/005	Pediatric Defibrillator Pads Responder, 5 pr/cs
850/40518/025	Pediatric Defibrillator Pads Responder, 25 pr/cs
900634-001	Custom Setup Card

GE MEDICAL ECG CABLES & LEADWIRES

107230	Leadwire Set, Snap, 3 LD, Gray, 100cm AHA, 3/set
107231	Leadwire Set, Grab, 3 LD, Gray, 100cm AHA, 3/set
107235	Leadwire Set, Grab, 6 LD, Gray, 100/130cm AHA, 6/set
107289	Leadwire Set, Grab, 6 LD, Gray, 100 AHA, 6/set
107290	Leadwire Set, Snap, 6 LD, Gray, 100 AHA, 6/set
107314	Leadwire Set, Grab, 3 LD, Gray, Pediatric, AHA, 3/set
107320	ECG Cable, 6 LD, AHA
107326	ECG Cable, 3 LD, AHA
107328	Leadwire Set, Snap, 3 LD, AHA, 3/set
107331	ECG Cable, 3 LD w/ESU Filter For OR, AHA
107333	ECG Cable, 6 LD w/ESU Filter For OR, AHA
107384	ECG Cable, 3 LD, Safety, AHA

GE MEDICAL ECG & DEFIB RECORDING PAPERS

Consolidate recording paper purchases. Reduce time and expense to generate multiple purchase orders.

- 9801-001 Thermal Paper, Spacelabs Equivalent 307438-001, All White, 1.97" x 100 ft, For Spacelabs ECG Systems PC Bedside, 90600 Series, 90731, 90323, 10/pk, 10 pk/cs
- 9801-002 Thermal Paper, Quinton Equivalent 15022-001, All White, Black Trace, 8.5" x 11", Queue Hole, Z-Fold, For Quinton Stress Systems Q3000, Q4000, 216 sheets/pk, 20 pk/cs
- 9801-003 Thermal Paper, Quinton Equivalent 16211-002, All White, Queue Mark, 8.5" x 11", Z-Fold, For Quinton Stress Systems Q5000 & QCATH, 200 sheets/pk, 20 pk/cs
- 9801-004 Thermal Paper, Hewlett-Packard Equivalent M1707A, Red Grid, 8.5" x 11", Z-Fold, For HP M1700A, M1701A Pagewriter XL Systems, 200 sheets/pk, 8 pk/cs
- 9801-010 Thermal Paper, Hewlett-Packard Equivalent 40453A/B, Single Channel, Black Grid, 1.97" x 100 ft, For HP Defibrillators, 10/pk, 10 pk/cs
- 9801-016 Thermal Paper, Burdick Equivalent 007961, Red Grid, Black Trace, 8.5" x 11", Z-Fold, For Burdick E500TW, E350, E-600 ECG Systems, 200 sheets/pk, 10 pk/cs
- 9801-028 Thermal Paper, General Scanning Equivalent AR-200, All White, Black Grid, 8.5" x 11", Z-Fold, Queue Marks, For General Scanning Systems AR-200, 70206, 300 sheets/pk, 10 pk/cs
- 9801-057 Thermal Paper, Hewlett-Packard Equivalent M2481A, Red Grid, Black Trace, 8.5" x 5.71", 8 pk/cs
- 9801-059 Thermal Paper, Spacelabs Equivalent 006-0196-001, Black Trace, 1.97" x 97 ft, For Spacelabs Systems Using 006-0196-001 Paper, 10/bx, 10 bx/cs

GE MEDICAL LASER PRINTER PAPERS


Quality paper stock minimizes wear on laser printers. Save time printing traces on preprinted grid paper. Red grids improve trace legibility.

- 9402-053 Laser Printer Paper, Red Grid w/blank interpretation area, 8.5" x 11", 500 sheets/pk, 5 pk/cs
- 9402-054 Laser Printer Paper, Full Page Red Grid, 8.5" x 11", 500 sheets/pk, 5 pk/cs

GE MEDICAL LEG PLATE ADAPTER

An alternative to the disposable attachment pads (2464AAO). Smooth surface against patient's skin. Easy to clean or sterilize. Latex free.

- 1594AAO 2" x 2" Leg Plate Strap Adapter



GE MEDICAL ULTRA-ARCHIVIST® THERMAL RECORDING PAPERS

100 year warranty for trace legibility. Eliminates need for special storage precautions. Heavyweight paper to minimize potential for paper jams. Clear, precise traces.

- 2007854-001 Ultra-Archivist Thermal Paper, Red Grid w/blank interpretation area, 8.5" x 11", Z-Fold, 250 sheets/pk, 8 pk/cs
- 2007854-002 Ultra-Archivist Thermal Paper, Red Grid A-4 Modified, 8.26" x 11", Z-Fold, 250 sheets/pk, 8 pk/cs
- 2007854-003 Ultra-Archivist Thermal Paper, Full Page Red Grid, 8.5" x 11", Z-Fold, 250 sheets/pk, 8 pk/cs

GE MEDICAL UNIVERSAL THERMAL RECORDING PAPERS

Bright white paper offers excellent trace legibility. Best choice for MUSE customers. Electronic record reduces need for more expensive, long-term trace retention.

- 9902-020 Universal Thermal Paper, Full Page Red Grid, 8.5" x 11", Z-Fold, Queue Hole, 300 sheets/pk, 8 pk/cs
- 9902-021 Universal Thermal Paper, Full Page Red Grid, A-4 Modified, 8.26" x 11", Z-Fold, Queue Hole, 300 sheets/pk, 8 pk/cs
- 9902-024 Universal Thermal Paper, Red Grid w/blank interpretation area, 8.5" x 11", Z-Fold, Queue Hole, 300 sheets/pk, 8 pk/cs
- 9902-024GY Universal Thermal Paper, Gray Grid w/blank interpretation area, 8.5" x 11", Z-Fold, Queue Hole, 300 sheets/pk, 8 pk/cs
- 9902-030 Universal Thermal Paper, All White, Black Queue Mark, 8.5" x 11", Z-Fold, 300 sheets/pk, 8 pk/cs
- 9902-031 Universal Thermal Paper, Green Grid w/blank interpretation area, 8.5" x 11", Z-Fold, Queue Hole, 300 sheets/pk, 8 pk/cs
- 9902-324 Universal Thermal Paper, Red Grid w/blank interpretation area, 3-Hole Punch, 8.5" x 11", Z-Fold, Queue Hole, 300 sheets/pk, 8 pk/cs
- 9902-920 Universal Thermal Paper, Full Page Red Grid, 8.5" x 11", Z-Fold, Queue Hole, 300 sheets/pk, 10 pk/cs
- 9902-924 Universal Thermal Paper, Red Grid w/blank interpretation area, 8.5" x 11", Z-Fold, Queue Hole, 300 sheets/pk, 10 pk/cs
- 9902-931 Universal Thermal Paper, Green Grid w/blank interpretation area, 8.5" x 11", Z-Fold, Queue Hole, 300 sheets/pk, 10 pk/cs

GE MEDICAL PREMIUM THERMAL RECORDING PAPERS

Exclusive chemical composition resists trace fading for 5-7 years following recommended storage procedures. Clear, precise trace. Designed for reading ease.

- 22616606 Premium Thermal Paper, Full Page Brown Grid, Queue Mark, 8.5" x 11", Z-Fold, 150 sheets/pk, 10 pk/cs
- 9402-060 Premium Thermal Paper, Full Page Red Grid, Queue Hole, 8.5" x 11", Z-Fold, 150 sheets/pk, 16 pk/cs
- 9402-061 Premium Thermal Paper, Red Grid w/blank interpretation area, Queue Hole, 8.5" x 11", Z-Fold, 150 sheets/pk, 16 pk/cs
- 22613001 Premium Thermal Paper, 55mm x 40M
- 22616802 Premium Thermal Paper, Brown Dot Grid, 3.54" x 115 ft, 10/cs
- 385-1000-015 Premium Thermal Paper, Brown Grid, 8.5" x 11", Z-Fold, 400 sheets/pk, 10 pk/cs
- 9402-018 Premium Thermal Paper, Red Grid, 4" DDW 4.25" x 150 ft, 32/cs
- 9402-020 Premium Thermal Paper, Full Page Red Grid, 8.5" x 11", Z-Fold, Queue Hole, 300 sheets/pk, 8 pk/cs
- 9402-021 Premium Thermal Paper, Full Page Red Grid, A-4 Modified, 8.26" x 11", Z-Fold, Queue Hole, 300 sheets/pk, 8 pk/cs
- 9402-023 Premium Thermal Paper, Red Grid, 4.25" x 75 ft, 2/pk, 9 pk/cs
- 9402-024 Premium Thermal Paper, Red Grid w/blank interpretation area, 8.5" x 11", Z-Fold, Queue Hole, 300 sheets/pk, 8 pk/cs
- 9402-028 Premium Thermal Paper, All White, 2" x 100 ft, 48/cs
- 9402-029 Premium Thermal Paper, All White, Queue Hole, 8.5" x 11", Z-Fold, 300 sheets/pk, 8 pk/cs
- 9402-030 Premium Thermal Paper, All White, Black Queue Mark, 8.5" x 11", Z-Fold, 300 sheets/pk, 8 pk/cs
- 9402-031 Premium Thermal Paper, Green Grid w/blank interpretation area, 8.5" x 11", Z-Fold, Queue Hole, 300 sheets/pk, 8 pk/cs
- 9402-046 Premium Thermal Paper, 2" x 100 ft, 48/cs
- 9402-050 Premium Thermal Paper, Red Grid w/blank interpretation area, 8.5" x 11", Z-Fold, Queue Hole, 500 sheets/pk, 5 pk/cs
- 9402-324 Premium Thermal Paper, Red Grid w/blank interpretation area, 3-Hole Punch, 8.5" x 11", Z-Fold, Queue Hole, 300 sheets/pk, 8 pk/cs
- 9472-103 Premium Thermal Paper, All White, 2" x 100 ft, 4/bx

GE MEDICAL CARDIOLOGY LEADWIRE ADAPTERS

Latex free. 90 day warranty.

9490-210	Set Adapter GE Wide-Mouth Universal Alligator/Banana Socket/Snap, 10/set
38401804	Set Leadwire Adapters Grabber/Banana, 10/set
9459-901	Set Adapters Grabber/Banana, 6" Wire, 10/set
9490-105	Set Adapter Grey/Black Alligator/Banana Socket, 4mm, 10/set
9490-106	Set Adapter Grey/Blue Alligator/Banana Pin, 3mm, 10/set
9490-111	Set Adapter Clip Black Alligator/Banana Socket, 4mm, 10/set
9490-112	Set Adapter Clip Black Alligator/Banana Pin, 3mm, 10/set
9490-114	Set Adapter Clip Black Alligator/Banana Socket, 4mm, 14/set
9490-115	Set Adapter Grey/Black Alligator/Banana Socket, 4mm, 14/set

GE MEDICAL CARDIOSMART, CARDIOSYS, MAC CABLES & LEADWIRES

Easy to clean, wipe with damp cloth. Built for durability. Latex free. One year warranty. Coaxial cable reduces potential noise.

22341809	10 ft ECG Cable For MAC 1100/1200/1200ST, AHA, Uses CAM 14 Leadwires
22341802	7 ft ECG Cable for CardioSmart, MAC 1000, CardioSys, AHA Uses CAM 14 Leadwires
38401742	Leadwire Set, Banana/Cardiosmart, 0.156" Banana, AHA, 10/set

**WELCH ALLYN CP 10 & CP 20 ECG UNITS**

Includes the following accessory items standard: 10-lead resting AHA patient cable; 500 disposable tab electrodes; 10 resting electrode adaptors; physician's guide to the analysis program (interpretive systems) on CD; in-servie video on CD; 200 sheets of paper; operator manual on CD; power cord.

11102-200	CP 10 ECG Non-Interpretive 12 Lead Multichannel ECG w/ECG measurements, 8.25" x 11" Paper, Alphanumeric Keyboard & AHA Patient Cable
11102-500	CP 10 Interpretive 12-Lead Multichannel ECG, 8.25" x 11" Paper, Alphanumeric Keyboard & AHA Patient Cable
11103-200	CP 20 Non-Interpretive 12-Lead Multichannel ECG w/ECG measurements, ECG Preview Display, 8.25" x 11" Paper, Alphanumeric Keyboard & AHA Patient Cable
11103-210	CP 20 Non-Interpretive 12-Lead Multichannel ECG w/ECG measurements, Communication Module, ECG Preview Display, 8.25" x 11" Paper, Alphanumeric Keyboard & AHA Patient Cable
11103-500	CP 20 Interpretive 12-Lead Multichannel ECG, ECG Preview Display, 8.25" x 11" Paper, Alphanumeric Keyboard & AHA Patient Cable
11103-510	CP 20 Interpretive 12-Lead Multichannel ECG w/communication module, ECG Preview Display, 8.25" x 11" Paper, Alphanumeric Keyboard & AHA Patient Cable

Accessories

80185-0000	Data Transfer Cable
96152-0000	CP 10 Interpretation Upgrade (must provide system serial #)

WELCH ALLYN ECG ACCESSORIES

08250-0000	CP 10/CP 20 Cart w/cable arm
08265-0000	ECG Utility Cart
08285-0000	ECG Carrying Case
25254-0000	Cable Separators, 10/set
45008-0000	Resting Tab Electrodes, 1000/bx
58548-0000	Pinch Clip Adapters, 10/set
65159-0000	Cable Arm
65160-0000	CP 10/CP 20 Cart w/o cable arm
80141-0000	Hospital Grade Power Cable
9293-032-60	10-Lead Patient Cable, Resting ECG, AHA CP 10/CP 20
9293-032-61	10-Lead Patient Cable, Resting ECG, IEC CP 10/CP 20
94016-0000	8.25" x 11" Paper For CP 10 & CP 20, 200 sheets/pk, 5 pk/cs
58581-0000	ACCESSORY: Universal Adapter Clip, 10/pk

**MIDMARK IQMARK™ DIGITAL ECG (SERIAL VERSION)**

Includes: IQmark™ Digital ECG Acquisition Module (serial version), IQmark™ Digital ECG Software Package, IQmark™ Diagnostic Workstation - Patient Management Software, 10-lead patient cable, 10-pack universal clips, 100-pack ECG electrodes, 4-pack AA batteries, mouse pad, operation manuals and quick reference guide.

4-000-0000	IQmark™ Digital ECG (Serial Version)
------------	--------------------------------------

MIDMARK IQMARK™ ACCESSORIES

1-370-0004	USB External Compact Flash Card Reader
3-370-0005	PCMCIA Compact Flash Adapter
1-136-0005	USB to Serial Port Adapter

MIDMARK IQMARK™ DIGITAL RECORDER KIT

Kit includes: IQmark™ Digital Recorder, lead wire set, recorder pouch and belt, compact flash card and operation manual.

4-000-0170	IQmark™ Digital Recorder Kit For EZ & Advanced Holter
------------	---

MIDMARK IQMARK™ DIGITAL RECORDER KIT - 3 CHANNEL

Kit includes: IQmark™ Digital Recorder, lead wire set, recorder pouch and belt, 16MB PCMCIA flash card and operation manual.

4-000-0105	IQmark™ Digital Recorder Kit - 3 Channel For 7700/8800 Holter
------------	---

MIDMARK IQMARK™ EZ HOLTER KIT WITH DIGITAL RECORDER

Kit includes: IQmark™ Digital Recorder Kit (lead wire set, recorder pouch and belt, compact flash card, operation manual), IQmark™ EZ Holter Software Package, IQmark™ Diagnostic Workstation - Patient Management Software, compact flash card reader, PCMCIA compact flash adapter, sentinel data protection key, (4) holter prep kits, mouse pad, operation manual and quick reference guide..

4-000-0150	IQmark™ EZ Holter Kit w/digital recorder
------------	--



MIDMARK IQMARK™ ADVANCED HOLTER KIT WITH RECORDER & PACEMAKER OPTION

Includes: leadwire set, recorder pouch and belt, compact flash card, operation manual, IQmark™ Advanced Holter with Pacemaker Analysis Software Package, IQmark™ Diagnostic Workstation - Patient Management Software, Compact Flash Card Reader, PCMCIA Compact Flash Adapter, Sentinel Data Protection Key for Pacemaker (parallel port key shipped standard, USB key available upon request), 4 Holter Prep Kits, mouse pad, and quick reference guide.

4-000-0158 IQmark™ Advanced Holter Kit w/recorder & Pacemaker Option



BURDICK BASELINE STRESS PREP KIT

Each kit contains: 10 Stress Electrodes, Foam, 1 Skin Scrub, 1 Instruction Sheet, 3 Alcohol Prep Pads, 1 Bacitracin Ointment, 1 Moist Wipe, 1 Surgical Prep Razor.

015-0629-00 Blue Max Stress Prep Kit, 30 kit/cs
097008 Baseline Stress Prep Kit, 30/cs



BURDICK TM55 AND TM65 TREADMILLS - DIGITALLY CONTROLLED

Thousands of clinicians worldwide rely on the TM55 and TM65 to deliver the accuracy and reliability required for exercise stress tests. Fewer moving parts reduces wear and tear and maximizes system reliability. Controlled "soft start" for patient safety. Digitally controlled AC-drive system provides high torque for accurate speed and grade, regardless of weight. Quiet Operation through entire speed range eliminates office disruption. Fully integrated into the Burdick Quest exercise stress system.

TM55HXXRBE TM55, 230V US & Canada
TM55LXXRBE TM55, 120V



BURDICK TM55 & TM 65 TREADMILLS

Proven design maximizes system reliability. Controlled soft start improves patient safety. Accurate speed and grade provides better test results. Quiet operation calms patients and eliminates office disruption. Fully integrated into the Burdick Quest® exercise stress system

000488-010 TM65 Treadmill, 120V
000488-011 TM65 Treadmill, 230V US/Canada
Accessories
032758-001 Short Handrail
032759-001 Long Handrail
037288-001 Emergency Stop

BURDICK TWAVE ALTERNANS OPTION FOR QUEST

015-0631-00 PM-3 Alternans Lead Wire Set Used For Microvolt T-Wave Alternans Test on Quest Exercise Stress Systems AHA
015-0635-00 Micro-V Alternans Sensors, 7/pk, 20 pk/bx
015-0637-00 PM-3 Patient Module & System Cable For Microvolt T-Waves Alternans Test Option
040-1219-00 Microvolt T-Wave Alternans Test S/W & H/W Upgrade Kit For Quest Systems AHA
92400-T Microvolt T-Wave Alternans Option AHA

BURDICK HEARTWORKS™ DIAGNOSTIC WORKSTATION

A suite of diagnostic devices conveniently integrated on one PC: Vision Holter, Vision 5L Recorder, ABP Monitor and software, Universal ECG, SpiroCard and Office Medic software.

D20-1ADH HeartWorks™ Diagnostic Workstation with Laptop Computer and Printer
D20-1D HeartWorks™ Diagnostic Workstation w/out Laptop and Printer

BURDICK HEARTSTRIDE EXERCISE STRESS SYSTEM

Designed to meet the unique needs of primary care physicians looking to add cardiac stress testing capabilities to their practices. With its simple, graphical user controls, plus optional network connectivity capability (including ability to interface with EMRs), integrating cardiac stress testing into the office has never been simpler. Designed to deliver value, the PC-based HeartStride allows smaller offices to leverage their PC hardware investment while increasing patient care. HeartStride is built on 50 years of diagnostic experience and reliably delivers clinical accuracy and diagnostic confidence. Plus it's backed by the Cardiac Science Corporation field service network. Features: Superior ECG trace quality; Easy, intuitive operation; Simple to learn and use; Comprehensive Burdick on-site training; In-office Practice Support Program; Fast, accurate reporting.

H20-1AEPT01 Network Configuration: English Language w/AHA Accessories, Includes: CPU, Printer, LCD Flat Panel Monitor, Laser Printer, Network Connectivity, TM55 Treadmill & Custom Cart.
H20-1AKPT01 Base Configuration: English Language w/AHA Accessories, Includes: CPU, Printer, LCD Flat Panel Monitor, Laser Printer, TM55 Treadmill & Custom Cart.
H21-1AKT01 SW Only: English Language w/AHA Accessories, Includes: TM55 Treadmill & Custom Cart.




MIDMARK IQMARK™ DIGITAL ECG (USB VERSION)

Includes: IQmark™ Digital ECG Acquisition Module (USB Version), IQmark™ Digital ECG Software Package, IQmark™ Diagnostic Workstation - Patient Management Software, 10-lead patient cable, 10-pack universal clips, 100-pack ECG electrodes, mouse pad, operation manuals and quick reference guide.

4-000-0056 IQmark™ Digital ECG (USB Version)

MIDMARK IQMARK™ EZ ADVANCED HOLTER SUPPLIES

2-100-0060 Lead-Lok Electrodes, 50/bx
 2-100-0080 Prep Kit - IQmark™ Holter, 5-Lead, 2AA's
 2-100-0090 Prep Kit - IQmark™ Holter, 7-Lead, 2AA's
 3-370-0001 Lead Wire, 5-Lead, IQmark™ Digital Recorder
 3-370-0002 Lead Wire, 7-Lead, IQmark™ Digital Recorder
 3-370-0003 Pouch w/belt for IQmark™ Digital Recorder
 1-370-0002 Compact Flash Card
 2-100-0081 Prep Kit - IQmark™ Holter, 5-Lead 2 AA's, 30/cs
 2-100-0091 Prep Kit - IQmark™ Holter, 7-Lead, 2AA's, 30/cs



MIDMARK IQMARK™ PACEMAKER OPTION FOR EXISTING ADVANCE HOLTER

IQmark™ Advanced Holter Pacemaker Analysis Software, Sentinel Data Protection Key for Pacemaker (parallel port key shipped standard, USB key available upon request), operation manual, and quick reference guide.

4-000-0140 IQmark™ Pacemaker Option For Existing Advanced Holter

MIDMARK IQMARK™ 7700/8800 HOLTER SUPPLIES

3-053-0001 Lead Wire, 7-Lead, IQmark™ Digital Recorder
 3-079-0001 Pouch w/belt for IQmark™ Digital Recorder
 3-513-0005 Junction Box For IQmark™ Digital Recorder
 1-079-0125 16MB Flash Card

MIDMARK IQMARK™ MANUALS

3-100-1006 IQmark™ Digital Recorder Manual
 3-100-1007 IQmark™ Holter 7700 Manual
 3-100-1008 IQmark™ Holter 8800 Manual
 3-100-1021 IQmark™ Digital ECG Manual
 3-100-1025 IQmark™ Digital Spirometer Manual
 3-100-1038 IQmark™ Diagnostic Workstation Manual
 3-100-1044 IQmark™ EZ Holter Manual
 3-100-1045 IQmark™ Digital Recorder Manual
 3-100-1049 IQmark™ Digital ECG PDA & SpiroPDA Manual
 3-100-1050 IQmark™ Advanced Holter Manual

MIDMARK IQMARK™ ECG SUPPLIES

2-100-0205 Disposable ECG EZ-Trodes, 500/bx
 3-047-0001 Universal ECG Clips (3mm, 4mm & Snap)
 3-100-0200 IQmark™ Patient Cable
 3-999-0011 24-pack, AA Batteries
 2-100-0206 Disposable ECG EZ-Trodes, 500/bx, 10 bx/cs

KENDALL-LTP MEDI-TRACE® CHART PAPERS

Kendall-LTP manufactures one of the largest offerings of recording chart papers in the industry. With more than 40 years of manufacturing experience, our customers can expect consistent trace quality as well as a product that is designed to meet or exceed the instrument specifications.

01091735 Fetal Monitoring Paper, Compatible to Philips 9270-0678-030, w/label sets, 30 pk/cs
 20608594 Recording Chart Paper, Compatible to Spacelabs/Burdick 007868, 10 pads/cs
 30015682 Recording Chart Paper, Compatible to Spacelabs 307438, 50 r/cs
 30540499 Recording Chart Paper, Compatible to Spacelabs 307568, 20 pad/cs
 30561354 Recording Chart Paper, Compatible to Quinton 16211-002, 20 pad/cs
 30597226 Recording Chart Paper, Compatible to GE/Corometrics 4305AAO, 40 pad/cs
 30659745 Recording Chart Paper, Compatible to Nihon Kohden FQW50-3-100, 100 pad/cs
 30718970 Recording Chart Paper, Compatible to Physio-Control 804700, 50 r/cs
 30725389 Recording Chart Paper, Compatible to Philips 40477A, 80 r/cs
 30729225 Recording Chart Paper, Compatible to Quinton 18416-001, 20 pad/cs
 30748696 Recording Chart Paper, Compatible to Philips M1910A, 40 pad/cs
 30763711 Recording Chart Paper, Compatible to Spacelabs 006-0196-00, 100 pad/cs
 30768678 Recording Chart Paper, Compatible to Philips M1707A, 8 pad/cs
 30781747 Recording Chart Paper, Compatible to GE/Marquette 9402-046, 50 r/cs
 31002176 Recording Chart Paper, Compatible to GE/Marquette 9402-020, 8 pad/cs
 31002184 Recording Chart Paper, Compatible to GE/Marquette 9402-024, 8 pad/cs
 31091427 Recording Chart Paper, Compatible to Physio-Control 805319-03, 18 r/cs

JONES-ZYLON HOLTER MONITOR KIT

Contains: patient diary, audio cassette tape, 5 electrodes, surgical prep razor, 2 AA alkaline batteries and 6" x 9" minigrip bag.

568803 Holter Monitor Kit, 50/cs



BURDICK BASELINE HOLTER CASSETTE PREP KIT

Each kit contains: 1 Hook-Up Guide, 1 Surgical Prep Razor, 1 Skin Prep Scrub, 3 Alcohol Prep Pads, 12 Stress Loop Tape, 1 60 Minute Audio Cassette, 1 9V Battery, 1 Patient Diary.

043084 2 Channel Cassette Prep Kit, w/5 Holter Electrodes, 30/cs
 043085 3 Channel Cassette Prep Kit, w/7 Holter Electrodes, 30/cs



BURDICK BASELINE HOLTER DISC PREP KIT

Each kit contains: 1 Hook-Up Guide, 1 Surgical Prep Razor, 1 Skin Prep Scrub, 3 Alcohol Prep Pads, 12 Stress Loop Tape, 1 Patient Diary, 1 Patient ID Card, 2 9V Batteries.

- 043101 2 Channel 6532 Prep Kit, w/5 Holter Electrodes, 30/cs
- 043104 3 Channel 6532 Prep Kit, w/7 Holter Electrodes, 30/cs



043253



043273

BURDICK HOLTER MONITORING SYSTEM ACCESSORIES

- 012-0141-00 Universal Test Box For 90205 & 90208
- 012-0439-00 5-Lead Bonded Cable/Lead Set For 90205 & 90208
- 012-0440-00 7-Lead Bonded Cable/Lead Set For 90205 & 90208
- 016-0475-00 Nylon Carrying Pouch For 90205 & 90208
- 016-5066-00 Nylon Adjustable Shoulder Strap For 90205 & 90208
- 043092 Test Box For 6532, 6632, 6300
- 043113 Interface Kit For 6300 Cassette Recorder
- 043138 2 Channel Safety Cable (cable only) For 6632
- 043139 2 Channel Lead Set, 5 Safety Wires For 6632
- 043140 3 Channel Lead Set, 7 Safety Wires For 6632
- 043229 Signal Averaging/Late Potential Software
- 043253 3 Channel Lead Set, 7 Safety Wires, 6732 Recorder
- 043255 2 Channel Lead Set, 5 Safety Wires, 6732 Recorder
- 043273 Compact Flash Card, 30Mb For Use w/6732 Compact Digital Recorder
- 086418 Altair PC Holter System Operating Manual
- 086419 Altair Simplicity Holter Operating Manual
- 086421 6500 DISC Recorder Operators Manual
- 086422 Altair 6300 Cassette Recorder Operators Manual
- 143207 Interface Kit For 6732 & 6632 PC Card Recorder
- 175-1401-00 3 Channel Safety Cable (cable only) For 6632
- 685-0200-00 Compact Flash Card, 64 MB For 6732
- 92514-01 Digital Holter Recorder, 3-Channel, 128 MB

BURDICK PC CARD RECORDER PREP KIT

Each kit contains: 1 Surgical Prep Razor, 3 Alcohol Prep Pads, 7 Holter Electrodes, 1 Patient ID Card, 1 Hook-Up Guide, 1 Skin Prep Scrub, 12 Stress Loop Tape, 1 Patient Diary, 4 "AA" Batteries

- 043202 PC Card Recorder Prep Kit, 30/cs



BURDICK THREE-CHANNEL COMPACT DIGITAL RECORDER PREP KIT

Each kit contains: 1 Surgical Prep Razor, 3 Alcohol Prep Pads, 7 Holter Electrodes, 1 Patient ID Card, 1 Hook-Up Guide, 1 Skin Prep Scrub, 12 Stress Loop Tape, 1 Patient Diary, 2 "AA" Batteries

- 043272 3-Channel Compact Digital Recorder Prep Kit, 30/cs

BURDICK TWO-CHANNEL COMPACT DIGITAL RECORDER PREP KIT

Each kit contains: 1 Surgical Prep Razor, 3 Alcohol Prep Pads, 5 Holter Electrodes, 1 Patient ID Card, 1 Hook-Up Guide, 1 Skin Prep Scrub, 12 Stress Loop Tape, 1 Patient Diary, 2 "AA" Batteries

- 043250 2-Channel Compact Digital Recorder Prep Kit, 30/cs

GE MEDICAL HOLTER 12-LEAD AMBULATORY CABLES

Ribbon cable design with strategic electrode snap placement. No need for separate leadwires; one piece design. Cable can be cleaned or disinfected. Lightweight wire can be used on a large cross section of patients. Latex free. 90 day warranty.

- 421825-001 Cable SEER MC Recorder, 12-Ld Ambulatory, AHA, 36", Includes Electrode Snap Connections
- 421825-002 Cable SEER MC Recorder, 12-Ld Ambulatory, IEC, 36", Includes Electrode Snap Connections

GE MEDICAL HOLTER 5 & 7 LEAD SNAP-LOCK CABLES & LEADWIRES

Locks leadwires in place to ensure a good connection. Release button allows the leadwires to be removed when depressed. Snap end of leadwire provides secure connection to snap style electrodes. Latex free. 90 day warranty.

- 409106-001 Cable 8500 Holter, Snap-lock, 5-Ld, DIN, 42"
- 409106-002 Cable 8500 Holter, Snap-lock, 7-Ld, DIN, 42"
- 409106-003 Cable SEER Holter, Snap-lock, 5-Ld, DIN, 42"
- 409106-004 Leadwire Set, 8500/SEER/SEER MC Holter, 5-Ld, 18", 5/set
- 409106-005 Leadwire Set, 8500/SEER MC Holter, 7-Ld, 18", 7/set
- 409106-006 Leadwire Set, 8500/SEER/SEER MC Holter, 5-Ld, 24", 5/set
- 409106-201 Cable 8500 Holter, Snap-lock, 5-Ld, DIN, 15"
- 419598-005 Cable SEER MC, Snap-lock, 5-Ld, DIN, 42"
- 419598-007 Cable SEER MC, Snap-lock, 7-Ld, DIN, 42"

GE MEDICAL HOLTER DOWNLOAD CABLES

Download patients data quickly and easily from SEER to Holter scanner. Durable flex cable. Latex free. 90 day warranty.

- 80784-001 20" Cable Download SEER Recorder

GE MEDICAL HOLTER MULTI-LINK 12-LEAD CABLES & LEADWIRES

Shielded leadwires with multi-link connectors. Provides 12-lead monitoring capability. Latex free. One year warranty.

- 411200-001 Multi-Link Leadwire Set, 5-Ld, Snap, AHA, 29", RA, LA, RL, LL, V, 5/set
- 411200-003 Multi-Link Leadwire Set, 5-Ld, Snap, IEC, 29", R, L, N, F, C, 5/set
- 419598-010 Cable SEER MC, Multi-Link, 12-Ld, AHA, 36"
- 419598-011 Cable SEER MC, Multi-Link, 12-Ld, IEC, 36"

GE MEDICAL HOLTER SUPPLIES

- 3801-902 9V Battery Energizer Alkaline, 12/box
- 3801-904 1.5V AA Battery Energizer Alkaline, 144/box
- 3801-905 9V Battery Energizer Alkaline, 72/box
- 400068-002 Holter Kits 5-Lead Kit Contents: 5 Vinyl Tape Electrodes, 24-Hour Cassette, 9V Battery, 2 Prep Pads & Patient Diary, 24 kits/cs
- 400068-003 Holter Kits 7-Lead Kit Contents: 7 Vinyl Tape Electrodes, 24-Hour Cassette, 9V Battery, 2 Prep Pads & Patient Diary, 24 kits/cs
- 400068-004 Holter Kits 7-Lead Kit Contents: 7 Vinyl Tape Electrodes, 48-Hour Cassette, 9V Battery, 2 Prep Pads & Patient Diary, 24 kits/cs
- 400068-009 Holter Kits 5-Lead Kit Contents: 5 Vinyl Tape Electrodes, 4 AA Batteries, 2 Prep Pads & Patient Diary, 24 kits/cs
- 400068-016 Holter Kits 5-Lead Kit Contents: 5 Foam Oval Electrodes, 4 AA Batteries, 2 Prep Pads & Patient Diary, 24 kits/cs
- 400068-021 Holter Kits 5-Lead Kit Contents: 5 Vinyl Tape Electrodes, 48-Hour Cassette, 9V Battery, 2 Prep Pads, Patient Diary & Prep Razor, 24 kits/cs
- 400068-024 Holter Kits 7-Lead Kit Contents: 7 Vinyl Tape Electrodes, 4 AA Batteries, 2 Prep Pads & Patient Diary, 24 kits/cs
- 400068-025 Holter Kits 12-Lead Kit Contents: 10 Vinyl Tape Electrodes, 4 AA Batteries, 2 Prep Pads & Patient Diary, 24 kits/cs
- 57995-002 8500 Holter Recorder Belt For Pouch 57995-005
- 57995-003 8500 Holter Recorder Shoulder Pouch w/strap
- 57995-005 8500 Holter Recorder Pouch w/o straps & belt
- 57995-901 8500 Holter Recorder Pouch w/strap & belt
- 90026-902 Holter Patient Diary, 50/pk
- 9435-001 24-Hour Capacity Cassette Tapes, 10/box
- 9435-002 48-Hour Capacity Cassette Tapes, 10/box
- 9435-003 24-Hour Capacity Cassette Tapes, 20/box
- 9435-004 48-Hour Capacity Cassette Tapes, 20/box

GE MEDICAL HOLTER TEST CABLES

Verify ECG signals on 8500 recorder. Positive fit connector to the recorder; banana socket for ECG cart leadwires. Accepts 3mm or 4mm leadwires. Latex free. 90 day warranty.

- 900248-001 8500 Holter Test Cable, 3 Channel, Fits Recorders With Mini-9 Pin Connectors
- 9401-901 8500 Holter Test Cable, 2 Channel, Fits Recorders With 3 Pin Round Connectors

GE MEDICAL HOLTER VALUE LINE CABLES & LEADWIRES

Easy insertion of leadwire into cable. Latex free. 90 day warranty.

- 2014788-001 Set Leadwire Holter, 5-Ld, DIN, 20", Value Line, 5/set
- 2014789-001 Set Leadwire Holter, 7-Ld, DIN, 20", Value Line, 7/set
- 2014790-001 Cable 8500 Holter Recorder, 5-Ld, DIN, 42", Value Line
- 2014791-001 Cable 8500 Holter Recorder, 7-Ld, DIN, 42", Value Line
- 2014792-001 Cable SEER Recorder, 5-Ld, DIN, 42", Value Line
- 2014793-001 Cable 8500 Holter Recorder, 5-Ld, DIN, 15", Value Line
- 2014794-001 Cable 8500 Holter Recorder, 7-Ld, DIN, 15", Value Line
- 2014795-001 Cable SEER Recorder, 5-Ld, DIN, 15", Value Line

WELCH ALLYN HOLTER SYSTEM 2000 ACCESSORIES

- 08239-0000 Template Matching, Software Option
- 08240-0000 Pouch For MT-100
- 08241-0000 Holter Hook-Up Kit, Includes 1.5V AA Batteries, Patient Diary, Razor, Abrading Pad, Alcohol Prep Pads, 5 Disp. Holter Electrodes
- 80188-0000 Patient Cable, Short, MT-100
- 80189-0000 Fiber Optic Transfer Cable, MT-200

KENDALL-LTP MEDI-TRACE® TELEMETRY/HOLTER POUCHES AND KITS

Kendall-LTP offers a large variety of Holter and Telemetry pouches, kits, and accessories. All of our pouches are designed for patient comfort, transmitter protection and improved monitor performance. Our kits are offered in many variations to accommodate almost any requirement.

- 30951514 Telemetry Kit, MPK-MOD5, Contains 535 Electrode, Prep Pad, (2) 1.5 Alkaline Batteries and Pouch, 50/cs
- 30951555 Telemetry Pouch, TMY-BLUE-50, Small, 50/cs
- 30951746 Holter Kit, HMPK-HOLTR1, Contains 535 & 533 Electrodes, Tape, Prep Pad, Abrasive Pad, 9V Battery, Razor and Diary, 50/cs
- 30951811 Telemetry Pouch, TMY-BLUE-2-100, Large, 100/cs
- 31027942 Telemetry Kit, MPK-MOD3, Contains 533 Electrode, Prep Pad, 8.4V Battery & Pouch, 50/cs
- 31027983 Telemetry Kit, MPK-MOD2, Contains 535 Electrode, Prep Pad, 8.4V Zinc-Air Battery & Pouch, 50/cs
- 31098448 Telemetry Kit, MPK-MOD4, Contains 535 Electrode, Prep Pad, 9V Battery & Pouch, 50/cs
- MK05140 Holter Kit, SK7828, Contains 535 Electrode, 9V Battery, Tape & Diary, 20/cs
- MK06200 Holter Kit, Digial, Contains 535 & 533 Electrodes, (2) AA Batteries, Razor, Prep Pad, Abrasive Pad & Diary, 20/cs
- MK07740 Telemetry Pouch, TP11 Large w/window, 50/cs
- MK07860 Telemetry Pouch, TP3 Non-woven, 100/cs
- MK07900 Telemetry Pouch, TP5 Non-woven w/2 straps, 100/cs
- MK07980 Telemetry Pouch, TP7S Large w/2 straps, 50/cs
- MK08000 Telemetry Pouch, TP8 w/slit & loop, 50/cs
- MK08080 Telemetry Pouch, TPV Vinyl w/clip, 50/cs



BURDICK BI-PHASIC MEDIC 6 PORTABLE DEFIBRILLATOR

Compact, streamlined design that provides easy to use, comprehensive defibrillation therapies in a compact 13 1/2 lb package. A radiant 6.4", non-glare screen makes functions easy to use in a variety of light settings. The easy to follow function keys, make system feature selections quick and easy to navigate. Includes a powerful rechargeable battery, which delivers 40+ consecutive discharges at 360j or 2 hours of continuous monitoring operation. Sync and a broad range of joule (energy) selections are available for emergency resuscitation, internal electrode procedures or other cardiac therapeutic modalities. Offers ECG capabilities via 3-lead (included) patient cable. Comprehensive pre/post treatment summary documentation which is provided on 50mm paper and stored into internal memory.

- 92450-1APV Medic 6 Defibrillator w/ext pacing & advisory
- 92450-1AV Medic 6 Defibrillator w/advisory

BURDICK DEFIBRILLATOR ACCESSORIES

015-0605-00	Patient Cable, 3-Lead, Medic 6
015-0608-00	Battery Pack For Medic 6
015-0617-00	Standard Thermal Paper, 150 ft x 45mm For Medic 6, 3 r/bx
015-0610-00	AC Power Cord For Medic 6
047262	Power Cord For DC-200, Medic 4 & Medic 5
086144	Operators Manual For DC-200
086183	Service Manual For DC-200
015-0612-00	Defibrillation/Pacer/Monitoring Pads, Adult For Medic 6, 10/bx
015-0621-00	Defibrillation/Pacer/Monitoring Pads, Pediatric For Medic 6, 10/bx
067006	External Paddle Set For Medic 5, Adult
067015	Adult Internal Paddles w/handles For DC-190
067016	Pediatric Internal Paddles w/handles For DC-190
067019	Adult Internal Paddles w/handles For Medic 5
067107	3-Lead Patient Cable For Medic 4 & Medic 5
067112	5-Lead Patient Cable For Medic 4 & Medic 5
850711	Battery Pack For Medic 5



WELCH ALLYN PIC30 DEFIBRILLATOR

973039C00E	PIC30 Manual Defibrillator
973039C0AE	PIC30 Manual Defibrillator w/shock advisory
973039CPOE	PIC30 Manual Defibrillator w/pacing
973039CPAE	PIC30 Manual Defibrillator w/pacing & shock advisory



WELCH ALLYN PIC40 DEFIBRILLATOR

971016	EtCO2 Capnography
973092E	PIC40 Defibrillator
973093E	PIC40 Defibrillator w/NIBP
973094E	PIC40 Defibrillator w/NIBP & pacing
973095E	PIC40 Defibrillator w/pacing



WELCH ALLYN PIC50 DEFIBRILLATOR

971081E	PIC50 Defibrillator, 3 Lead, Basic
971082E	PIC50 Defibrillator, 3 Lead, Deluxe
971083E	PIC50 Defibrillator, 12 Lead, Basic

WELCH ALLYN PIC™ RESCUER AED/MONITOR/DEFIBRILLATOR ACCESSORIES

001537-U	PIC Pediatric Electrode
001647-U	SuperPak Rechargeable Battery
001726-U	10-Pack Disposable Electrodes, 10 pk/bx
001739-U	Chart Plain White Paper, 40mm, 3 r/bx
001790-U	Cable: 3-Lead ECG (AHA) 6 ft
001795-U	PIC 5-Lead Patient Cable
001934-U	CO2 Sensor Kit
001950-U	Adult Airway Adapter, 5/bx
560639	PIC Tinted Lens Package
560640	PIC Clear Lens Package
601215	PIC Hanger, Black
90021-4	PIC Soft Carrying Case
900323	PIC30 Carrying Case
971029-E	PIC Paddle Tray
971091	Reliable Integral External Pacing
971107	PIC Hands Free Adapter, Grey
971108-E	PIC Deluxe External Paddle Set, Grey
980132	PIC 12V Power Adapter
999895	PIC EMS Service Plus Program (per year)
971104-PIC-E	PIC Quick Charger



3M™ DEFIB PADS

Ready-to-use, conductive gel pads for use in defibrillation. Self-contained to eliminate inadvertent gel bridging which can cause skin trauma to patient or operator. Bright orange solid gel pads are easy to see, require no cleanup, and leave no slippery residue that can interfere with CPR.

2345N	4 1/2" x 4 1/2" Defib Pads, 2/pk, 10 pk/bx, 10 bx/cs
2346N	4 1/2" x 6" Defib Pads, 2/pr, 10 pr/bx, 10 bx/cs



ZOLL AED PRO

Rugged construction and exceptional durability support the needs of Basic Life Support providers and tiered systems. Superior design and rubber over-molds are built to withstand harsh environments. ECG Monitoring with 3-lead cable and high-quality display. Manual override option satisfies Advanced Life Support professionals. Unrivaled battery compatibility and options that meet your specific needs. Internal memory data storage and IrDA or USB data transfer.

90110820499992010N	AED Pro Includes: (2) Sets STAT Padz II, (1) Charger, (2) Rechargeable Lead Acid Batteries
--------------------	--

ZOLL D 900 DEFIBRILLATOR

The Zoll D 900 offers advanced defibrillation technology in an easy-to-use unit ideal for "first response" defibrillation in hospital areas where complete ACLS-response capabilities are not required. It provides ECG monitoring with a large, easy-to-read monitor. The Zoll D 900 can be used either with Zoll Multi-Function Electrodes for "hands-off" defibrillation or with standard defibrillation paddles. A significant advantage is that pacing functions may be added to upgrade the unit at any time. The Zoll D 900 operates on both AC line and battery power and is completely compatible with other models in the Zoll cardiac resuscitation system.

8009-0006 1 Set Defib Paddles, 1 ECG Cable, Recorder Tape, 2 Guides & 1 Tube Defib Gel, Refurbished

**ZOLL PD 1200 DEFIBRILLATOR
PACEMAKER/DEFIBRILLATOR**

The Zoll PD 1200, combining defibrillation, pacing and monitoring, offers complete ACLS-response capabilities in a single, easy-to-use unit that is ideal for critical care specialty care areas. It can be configured with the Zoll "hands-off" Multi-Function Electrodes, with standard defibrillation paddles and pacing electrodes, and with internal defibrillation paddles and pacing electrodes. The Zoll PD 1200 operates on both AC line and battery power and is completely compatible with other models in the Zoll cardiac resuscitation system.

8002-0011 1 Set Standard Defib Paddles, 1 ECG Cable, 1 Pacer Output Cable, 2 Sets Electrodes (Adult, Pedi), Recorder Paper, 2 Guides, 1 Tube Gel, 1 Inservice Video, 2 Yr. Warranty, Refurbished

**ZOLL MEDICAL AED PLUS DEFIBRILLATOR**

Combines superior biphasic defibrillation with fast and simplified electrode placement. The AED Plus uses a brilliant new graphical/voice user interface to help the rescuer every step of the way. The electrode pad helps the rescuer perform CPR. A sensor on the CPR-D pads tracks the rescuer's compression depth and rate as soon as CPR begins. The electrode comes complete with rescue essentials and tamper-evident seals and it lasts four years. Features include audio and visual coaching and Rectilinear Biphasic™ Waveform. Requires only one electrode pad. Runs on inexpensive, consumer lithium batteries.

21000010102011010 AED Plus Defibrillator

ZOLL MEDICAL AED DEFIBRILLATOR ACCESSORIES

1000-0113-01 Replacement Trainer Remote Control
 8000-0800-01 Simulator/Tester For Functional Testing
 8000-0804-01 Universal Simulator Cable
 8000-0811-01 Flush Mount Wall Cabinet
 8000-0813-01 Zoll Data Review
 8000-0817-01 Surface Mount Wall Cabinet
 8000-0821 Replacement Trainer Phone Cord
 8000-0822 Replacement Trainer AC Adapter
 8008-0104-01 AED Trainer w/8 pre-programmed AED scenarios & remote control
 8900-2055 Adult Multi Function Electrode, 12/cs
 8900-2065 Pediatric Multi Function Electrode, 6/cs
 8900-0800-01 CPR-D Padz w/compression
 9310-0737 8 1/2" x 11" Wall Sign
 9310-0738 V-Shape Wall Sign
 9650-0301-01 Administration Guide
 9650-0304-01 Replacement Trainer Operator's Guide
 8000-0301 80mm Fanfold Paper Use Only w/M-Series Defibrillator, 20/cs
 8000-0802-01 Replacement Soft Case
 8000-0807-01 Type 123 Lithium Batteries, 10/pk
 8000-0808-01 Replacement PASS Cover
 8000-0809-01 Wall Mounting Bracket
 8000-0814-01 Recessed Mount Wall Cabinet
 8000-0815-01 USB IRDA Adapter
 8000-0816-01 RS-232 IRDA Adapter
 8900-0801-01 Stat-Padz II HVP Multi-Function Electrodes
 8900-0802-01 Stat-Padz II, 1 pr
 8900-0803-01 Replacement Adhesive Gels for CPR-D Padz - Training Electrode replacements, 5pair/bx
 8900-0804-01 CPR-D Padz Training Electrodes, reusable "z-design" electrode with CPR Hand Placement Indicator and 1 pair of Disposable Adhesive Gels
 8900-0805-01 Stat Padz II (reusable for Training)
 8900-0810-01 ACCESSORY: Pedi-padz II, Pediatric Multi-Function Electrodes, 1 pair
 8900-4003 HVP Multi Function Electrodes, stat-padz, 24 month shelf life, 12pr/cs
 8900-4004 Stat-padz HVP Multi-Function Electrodes, Shelf Life 24 Months, 1 pair
 8900-4005 Solid Gel Multi Function Electrodes, pro-padz, 24 month shelf life, 12pr/cs
 8900-4012 Multi-Function Electrode, Shelf-Life 9 mo, Sterile, 6 pr/cs
 9650-0300-01 Operator's Guide/Wall Poster
 9656-0156 AED Plus Brochure Literature, 25/pk
 9656-0158 AED Plus Automated External Defibrillator w/CPR feedback literature, 50/pk
 9656-0161 AED Plus Die-Cut Literature, 25/pk
 9656-0162 AED Plus Accessories & Training Products Literature, 25/pk

BURDICK FIRSTSAVE AED DEFIBRILLATOR

A consumer-grade AED. It includes many of the innovations found on the CardioVive AED at a value price. Audible alarms indicate system failure. Patented RHYTHMx® analysis software. Patented STAR® biphasic defibrillation energy waveform. Patented technology performs automated daily, weekly, and monthly self-tests on the battery, electrode presence and internal electronics. Pediatric capability. Simple, easy-to-use, one-button operation. Instructive voice prompt guides.

92531-201 FirstSave AED Defibrillator

BURDICK FIRSTSAVE AED DEFIBRILLATOR ACCESSORIES

010-1564-00 Adult Electrodes, 1 set/pk
 010-1572-00 Wire Rack
 010-1596-00 Pediatric Electrodes, 1 set/pk
 010-1649-00 Replacement Battery
 010-1651-00 Soft Carrying Case



BURDICK CARDIOVIVE AED DEFIBRILLATOR

Automatic daily self-test ensures safe, reliable operation. Pre-connected, non-polarized, disposable pads eliminates confusion by saving critical time. Easy to use: open cover, open electrode package, place electrodes on bare chest & push button as directed. Patented biphasic waveform measures impedance and adjusts to deliver appropriate defibrillation level. Advanced technology differentiates between shockable life-threatening cardiac rhythms and non-shockable arrhythmia events.

92532-201 CardioVive Defibrillator



WELCH ALLYN AED-10 AUTOMATED EXTERNAL DEFIBRILLATOR

Comes with an illuminated defib pad placement illustration on the device front panel, backlit display to provide the user with test prompts - battery status - number of shocks and energy settings, concise voice prompts and CPR timer. Includes the following standard accessory package: (1) non-rechargeable lithium battery (90+ shocks at 360 joules or 5 hours of continuous use), user's manual, quick reference card, defibrillation pads and five-year limited warranty program.

970300E AED-10 Device & Standard Accessories Only
 970302E AED-10 w/carrying case & standard accessories
 970308E AED-10 w/lead out pads
 970311E AED-10 RAPID RESPONSE KIT incl: AED-10 w/new expanded carrying case, 2 AED-10 Lithium Batteries, 2 sets of Defib Pads, 1 Pediatric Energy Reducer, 1 1st Responder Prep Kit, Quick Ref Card, In-Svc Video, and Operator Manual

Accessories

900432 AED 10 Carrying Case, 1 Year Limited Warranty

WELCH ALLYN AED AUTOMATED EXTERNAL DEFIBRILLATOR ACCESSORIES


00172-0 Electrode: Bi-Trode Limb Electrode Set (2 wrist clamp electrodes)
 00183-0E AED-20 Non-Rechargeable (Lithium) Battery
 00185-2 AED-10 Non-Rechargeable (Lithium) Battery
 00185-3 AED Defibrillation Pads, 2/pouch
 001851-U AED Trainer, Adult Pads
 001855-U JumpStart Adult Defibrillation Pads, 10 pr/bx
 001857-E Lead Out Pads For AED10 & 20
 001962-U SmartLink Software Program
 00211-9 AED-20 Parallel Printer Communications Kit
 00213-0 Cable: AED-20 ECG (AHA)
 00213-7 AED Quick Response Prep Kit
 00214-3 AED IR Communication Kit
 002155 AED Wire Mount
 002156-U AED Compact Wall Cabinet
 90021-6E Rugged AED-20 Carrying Case w/quick reference instruction placard, storage space for 1 battery & storage pockets for other supplies
 90042-1 Rugged AED-10 Soft Shell Carrying Case, Includes Quick Reference Instruction Placard & Storage Area For Supplies
 970205 AED-20 Upgrade ECG to Manual
 980136 Cardiology Datacard, 4MB "Linear"
 980143 AED-20 Trainer
 981125E PowerStick Battery Charging Station, 1 Bay
 990178 AED-20 In Service Video (VHS)



WELCH ALLYN AED-20 AUTOMATED EXTERNAL DEFIBRILLATOR

Includes a diagnostic quality backlit display which provides the user with text prompts, battery status, heart rate, event timer, number of shocks, volume control, display brightness control and energy settings. Also included are concise voice prompts and a CPR timer. Can be upgraded in the future for expanded functions. Comes complete with the standard accessory package: (1) battery (either a non-rechargeable long life lithium for infrequent use or rechargeable nickel metal hydride for frequent use/monitoring with charging station), operator's manual, quick reference card, defibrillation pads and five-year limited warranty program.

972211E AED 20 Basic w/non-rechargeable battery
 972212E AED 20 Basic w/rechargeable battery
 972213E AED 20 Basic Plus w/non-rechargeable battery
 972214E AED 20 Basic Plus w/rechargeable battery
 972215E AED 20 Deluxe w/non-rechargeable battery
 972216E AED 20 Deluxe w/rechargeable battery



3M™ RED DOT™ CLEAR PLASTIC MONITORING ELECTRODES

A tear-drop style monitoring electrode with clear plastic backing that is ideal for situations when immediate skin assessment is required without the need to remove the electrode. The gel's low chloride content minimizes skin irritation for increased patient comfort.

2235 5.1cm Dia. Monitoring Electrodes w/o abrader, 50/bg, 20 bg/cs
 2235-3 5.1cm Dia. Monitoring Electrodes w/abrader, 3/bg, 20 bg/bx, 10 bx/cs
 2235-5 5.1cm Dia. Monitoring Electrodes w/abrader, 5/bg, 20 bg/bx, 10 bx/cs
 9642 5.1cm Dia. Monitoring Electrodes w/o abrader, 50/bg, 20 bg/cs



3M™ RED DOT™ DIAPHORETIC FOAM MONITORING ELECTRODES

A rectangular-shaped electrode that offers a high-tack adhesive and foam backing designed for diaphoretic clinical situations. Patented solid gel penetrates the skin quickly allowing high quality. Low chloride content minimizes skin irritation.


2230 5.1cm x 5.5cm Monitoring Electrodes w/o abrader, 50/bg, 20bg/cs
 2270 5.1cm x 5.5cm Monitoring Electrodes w/abrader, 25/bg, 40bg/cs
 2270-3 5.1cm x 5.5cm Monitoring Electrodes w/abrader, 3/bg, 20 bg/bx, 10 bx/cs
 2270-5 5.1cm x 5.5cm Monitoring Electrodes w/abrader, 5/bg, 20 bg/bx, 10 bx/cs




3M™ RED DOT™ DIAPHORETIC SOFT CLOTH MONITORING ELECTRODES

A rectangular-shaped electrode that offers a high-tack adhesive and soft cloth backing designed for diaphoretic clinical situations. Patented solid gel penetrates the skin quickly allowing high quality. Low chloride content minimizes skin irritation.

2231	5.1cm x 5.5cm Monitoring Electrodes w/o abrader, 50/bg, 20 bg/cs
2271	5.1cm x 5.5cm Monitoring Electrodes w/abrader, 25/bg, 40 bg/cs
2271-3	5.1cm x 5.5cm Monitoring Electrodes w/abrader, 3/bg, 20 bg/bx, 10 bx/cs
2271-5	5.1cm x 5.5cm Monitoring Electrodes w/abrader, 5/bg, 20 bg/bx, 10 bx/cs



2269T



2280



2284

3M™ RED DOT™ ECG MONITORING ELECTRODES WITH PRE-ATTACHED LEAD WIRE

2269T	Neonatal, 2cm x 4cm, Pre-Wired Radiolucent Electrode w/clear tape, 3/bg, 10 bg/bx, 10 bx/cs
2280	Neonatal, 22mm x 22mm, Pre-Wired Non-Radiolucent Electrode w/soft cloth, 3/bg, 100 bg/cs
2282	Neonatal, 22mm x 22mm, Pre-Wired Radiolucent Electrode w/soft cloth, 3/bg, 100 bg/cs
2283	Infant, 33mm x 22mm, Pre-Wired Radiolucent Electrode, 3/bg, 100 bg/cs
2284	Neonatal, 17mm x 22mm, Pre-Wired Radiolucent Limb Band Electrode, 3/bg, 100 bg/cs



3M™ RED DOT™ FOAM MONITORING ELECTRODES

Tear-drop style monitoring electrodes with foam backing. Patented solid gel penetrates the skin quickly allowing high quality, uninterrupted traces. Low chloride content minimizes skin irritation.

2237	4.4cm Dia. Monitoring Electrodes w/o abrader, 50/bg, 20 bg/cs
2237-3	4.4cm Dia. Monitoring Electrodes w/o abrader, 3/bg, 20 bg/bx, 10 bx/cs
2237-5	4.4cm Dia. Monitoring Electrodes w/o abrader, 5/bg, 20 bg/bx, 10 bx/cs
2259	5.1cm Dia. Monitoring Electrodes w/abrader 25/bg, 40 bg/cs
2259-3	5.1cm Dia. Monitoring Electrodes w/abrader, 3/bg, 20 bg/bx, 10 bx/cs
9640	5.1cm Dia. Monitoring Electrodes w/o abrader, 50/bg, 20 bg/cs



3M™ RED DOT™ INFANT SOFT CLOTH MONITORING ELECTRODES

2258-3 Infant, 3.2cm Dia., Soft Cloth Electrode w/o abrader, 3/bg, 20 bg/bx, 10 bx/cs



3M™ RED DOT™ MONITORING ELECTRODES WITH FOAM TAPE & STICKY GEL

Multi-purpose monitoring electrodes with foam backing and sticky gel are ideally suited for the OR, ER and diagnostic stress testing applications. Conveniently packaged in individual strips of ten electrodes for easy use and less waste. Remains fresh for 30 days out of the bag.

2560	Monitoring Electrode w/o abrader, 10/strip, 5 strips/bg, 20 bg/cs
2570	Monitoring Electrode w/abrader, Radiolucent Stud, 50/bg, 20 bg/cs
2570-3	4cm x 3.5cm Monitoring Electrodes w/abrader, Radiolucent Stud, 3/bg, 200 bg/cs
2570-5	4cm x 3.5cm Monitoring Electrodes w/abrader, Radiolucent Stud, 5/bg, 200 bg/cs



3M™ RED DOT™ MONITORING ELECTRODES WITH MICROPORE™ TAPE BACKING

A round-style monitoring electrode with 3M™ Micropore™ Surgical Tape backing. Patented solid gel penetrates the skin quickly allowing high quality. Low chloride content minimizes skin irritation. Radiolucent to minimize the need to replace or remove electrodes in most radiological imaging applications.

2239	6cm Dia. Monitoring Electrodes w/o abrader, 50/bg, 20 bg/cs
2249	6cm Dia. Monitoring Electrodes w/abrader, 25/bg, 40 bg/cs
2249-3	6cm Dia. Monitoring Electrodes w/abrader, 3/bg, 20 bg/bx, 10 bx/cs
2248	Pediatric, 4.4cm Dia., Micropore Paper Tape Electrode, 25/bg, 40 bg/cs
2248-3	Pediatric, 4.4cm Dia., Micropore Paper Tape Electrode, 3/bg, 20 bg/bx, 10 bx/cs

3M™ RED DOT™ REPOSITIONABLE MONITORING ELECTRODES

2660-3	Repositionable Monitoring Electrodes, 1.56" x 1.25", 3/bg, 200 bg/cs
2660-5	Repositionable Monitoring Electrodes, 1.56" x 1.25", 5/bg, 200 bg/cs
2670-3	Repositionable Monitoring Electrodes, 1.56" x 1.25", 3/bg, 200 bg/cs
2670-5	Repositionable Monitoring Electrodes, 1.56" x 1.25", 5/bg, 200 bg/cs



3M™ RED DOT™ RESTING MONITORING ELECTRODES

Resting electrodes for 12-lead EKG applications.

- | | |
|------|--|
| 2330 | 2.2cm x 2.2cm Monitoring Electrodes , Tab Style, 10/strip, 10 strips/bg, 10 bg/bx, 4 bx/cs |
| 2352 | 1 3/4" x 7/8" Monitoring Electrodes , Snap Style, 10/bg, 100 bg/cs |
| 9650 | 2.2cm x 2.2cm Monitoring Electrodes , Tab Style, 100/bg, 10 bg/bx, 4 bx/cs |
| 2360 | Monitoring Electrode w/conductive adhesive, Tab Style, 10/card, 10 cards/bg, 40 bg/cs |



3M™ RED DOT™ SOFT CLOTH MONITORING ELECTRODES

A round-style monitoring electrode with soft cloth backing for general ECG monitoring. Patented solid gel penetrates the skin quickly allowing high quality. Low chloride content minimizes skin irritation.

- | | |
|--------|--|
| 2238 | 6cm Dia. Monitoring Electrodes w/o abrader, 50/bg, 20 bg/cs |
| 2238-3 | 6cm Dia. Monitoring Electrodes w/o abrader, 3/bg, 20 bg/bx, 10 bx/cs |
| 2238-5 | 6cm Dia. Monitoring Electrodes w/o abrader, 5/bg, 20 bg/bx, 10 bx/cs |
| 2255 | 6cm Dia. Monitoring Electrodes w/abrader, 25/bg, 40 bg/cs |
| 2255-3 | 6cm Dia. Monitoring Electrodes w/abrader, 3/bg, 20 bg/bx, 10 bx/cs |
| 9641 | 6cm Dia. Monitoring Electrodes w/o abrader, 50/bg, 20 bg/cs |



3M™ RED DOT™ TRACE PREP

A gentle abrasive material designed to remove the non-conductive skin layer to achieve a better trace quality.


- | | |
|------|---|
| 2236 | 3/4" x 196" Trace Prep w/dispenser, 36/cs |
|------|---|



3M™ RED DOT™ WET GEL FOAM MONITORING ELECTRODES

Tear-drop style monitoring electrodes with wet gel. Low chloride content minimizes skin irritation.

- | | |
|------|--|
| 2246 | 7.2cm Dia. Monitoring Electrodes w/abrader & Micropore Tape, 25/bg, 40 bg/cs |
| 2256 | 5.1cm Dia. Monitoring Electrodes w/abrader, 25/bg, 40 bg/cs |



BURDICK BLUEMAX ELECTRODES, HOLTER/STRESS

Offset snap.

- | | |
|-------------|---|
| 097001 | 40mm Dia. Electrodes, Vinyl, 25/pch, 40 pch/cs |
| 097004 | 50mm Dia. Electrodes, Foam, 10/pch, 100 pch/cs |
| 097010 | 68mm Dia. Electrodes, Micropore, 25/pch, 40 pch/cs |
| 097014 | 30mm Dia. Electrodes, Foam, 10/pch, 200 pch/cs |
| 015-0630-00 | Blue Max SUPATAB Resting ECG Electrode, Wet Gel, 100/pk, 24 pk/cs |
| 097002 | 50mm Dia. Electrodes, Foam, 25/pch, 40 pch/cs |



BURDICK CARDIOSENS™ DISPOSABLE ELECTRODE SENSORS

Three styles of disposable sensors: CardioSens™/Ultra II has aggressive adhesive. Fast signal stabilization with rapid defibrillator recovery and superior conductivity. CardioSens™/K has an exceptionally high conductive performance with natural Karaya adhesive. CardioSens™/Advantage has a narrow design, special formulation adhesive and unique, high conductivity tin backing.

- | | |
|--------|--|
| 047029 | CardioSens/Ultra II, 500/bx, 10 bx/cs |
| 047866 | CardioSens/K, 500/bx, 10 bx/cs |
| 047881 | CardioSens/Advantage (replaces 047057), 500/bx, 10 bx/cs |

Accessories

- | | |
|-------------|--|
| 704-0102-00 | Clear Choice Lead Adapter Pin to Universal Clip, 10/pk |
|-------------|--|



BURDICK HEARTLINE® MONITORING ELECTRODES

Designed to deliver high quality, easy-to-read tracings and long-lasting adhesion for consistent, dependable performance. Choose from give versatile electrodes, each offering a variety of performance possibilities.

007778	HeartLine Comfort Tape Wet Gel Electrodes, 48mm, 30/pouch, 600/cs
007781	HeartLine DuraFoam Wet Gel Electrodes, 48mm, 30/pouch, 600/cs
007782	HeartLine Radiolucent Universal Electrode, 30/pouch, 600/cs
007784	HeartLine SensaCloth Electrodes, 30mm, 30/pouch, 600/cs

GE MEDICAL DISPOSABLE RESTING ECG ELECTRODES

Latex free. Rapid, stable baseline. Reduces potential for cross-contamination. Secure adhesion. Silver/silver chloride construction.

9623-003P	Silver Mactrode® Plus Disposable Resting ECG Electrodes, Tab, 10/card, 10 card/pk, 10 pk/cs
2014778-001	SilverTRACE® Plug Disposable Resting ECG Electrodes, Plug B4MO Foam Rectangle, Banana Socket, 30/pch, 10 pch/bx
402207-210	BabyMAC® Pediatric Disposable Resting ECG Electrodes, Tab, 10/card, 10 card/pk, 5 pk/cs
900703-310	SilverTRACE® Disposable Resting ECG Electrodes, Adult Foam Rectangle, Wet Gel, Snap, 5/strip, 2 strip/pch, 160 pch/cs
9623-003	Silver Mactrode® Disposable Resting ECG Electrodes, Tab, 10/card, 10 card/pk, 10 pk/cs
9623-103P	Silver Mactrode® Plus Disposable Resting ECG Electrodes, Tab, 14/card, 10 card/pk, 10 pk/cs

GE MEDICAL EXERCISE/STRESS, HOLTHER, TELEMETRY SNAP ELECTRODES

Latex free foam and vinyl tape electrodes Two-piece snap connectors. Silver/silver chloride electrodes with wet or solid get.

2014768-001	SilverTRACE® Multi P21MO, Disposable Adult ECG Electrodes, Foam Oval w/tab, Solid Gel, 30/pch, 10 pch/bx
2014775-001	SilverTRACE® First P28MO, Disposable Adult ECG Electrodes, Wet Gel, 30/pch, 10 pch/bx
2014777-001	SilverTRACE® Soft P55MO, Disposable Adult ECG Electrodes, Foam Round, Wet Gel, 30/pch, 10 pch/bx
2014780-001	SilverTRACE® Soft Stress PS50MO, Disposable Adult ECG Electrodes, Foam Round, Wet Gel, 30/pch, 10 pch/bx
2014783-001	SilverTRACE® Window P50TR, Disposable Adult ECG Electrodes, Vinyl Tape Round, Solid Gel, 30/pch, 10 pch/bx
2014784-001	SilverTRACE® Window P50TR, Disposable Adult ECG Electrodes, Vinyl Tape Round, Solid Gel, 5/strip, 60 strip/bx
2014785-001	SilverTRACE® Window P50TR, Disposable Adult ECG Electrodes, Vinyl Tape Round, Solid Gel, 7/pch, 50 pch/bx
900703-205	SilverTRACE® Disposable Adult ECG Electrodes, Foam Oval, Wet Gel, 5/strip, 30 strip/bx, 4 bx/cs
900703-230	SilverTRACE® Disposable Adult ECG Electrodes, Foam Oval, Wet Gel, 30/pch, 10 pch/cs
9653-305	SilverTRACE® Disposable Adult ECG Electrodes, Clear Vinyl Tape, Solid Gel, 5/pch, 60 pch/cs
9653-507	SilverTRACE® Disposable Adult ECG Electrodes, Clear Vinyl Tape, Solid Gel, 7/pch, 50 pch/cs

GE MEDICAL MONITORING/TELEMETRY ELECTRODES

Latex free. Conductive gel.

2014779-001	SilverTRACE® Micro P55MI Disposable Adult ECG Electrodes, Paper Round, Wet Gel, 30/pch, 10 pch/bx
2014781-001	SilverTRACE® Soft Stress PS50MO Disposable Adult ECG Electrodes, Foam Round, Wet Gel, 5/strip, 120 strip/bx
2014782-001	SilverTRACE® Tissair P40CL Disposable Adult ECG Electrodes, Cloth Round, Solid Gel, 30/pch, 10 pch/bx
900690-030	SilverTRACE® Disposable Adult ECG Electrodes, Foam Round w/tab, Solid Gel, 30/pch, 10 pch/cs
900690-403	SilverTRACE® Disposable Pediatric ECG Electrodes, Cloth Rectangle, Solid Gel, 3/bg, 15 bg/pch, 10 pch/cs
900690-404	SilverTRACE® Disposable Neonatal ECG Electrodes, Cloth w/preattached leadwires, Solid Gel, 3/bg, 100 bg/cs
9431-004	SilverTRACE® Disposable Adult ECG Electrodes, Foam Round, Wet Gel, 30/pch, 10 pch/cs

GE MEDICAL QWIK CONNECT® PLUS DISPOSABLE ATTACHMENT PADS

Reduces potential for cross-contamination. Nickel-plated stud insures good conductivity. Eliminates the need for uncomfortable leg plate adapter straps. Latex free. Secure adhesive.

2464AAO	Spiral Leg Plate Attachment Pads, 1" x 3", 50/cs
---------	--

GE MEDICAL QWIK CONNECT® PLUS FETAL SPIRAL ELECTRODES

Locking mechanism secures the electrode during attachment. Electrode wires stabilized inside the guide tube until unlocked by the clinician. Latex free. Sterile.

7000AAO	Fetal Spiral Electrode, 50/cs
---------	-------------------------------

GE MEDICAL RADIOTRANSLUCENT ELECTRODES

Latex free polyethylene foam base. Carbon snap is x-ray translucent and non-magnetic. Conductive wet gel provides clear ECG traces.

2014786-001	SilverTRACE® First Carbon P28MO, Disposable Radiotranslucent ECG Electrodes, Foam Rectangle, Wet Gel, 30/pch, 10 pch/bx
2014787-001	SilverTRACE® First Carbon P28MO, Disposable Radiotranslucent ECG Electrodes, Foam Rectangle, Wet Gel, 5/strip, 60 strip/bx
320-1534-000	SilverTRACE® Disposable Adult ECG Electrodes, Cloth w/preattached leadwires, Solid Gel, 3/pch, 10 pch/pk
9408-401	SilverTRACE® Disposable Adult Radiotranslucent ECG Electrodes, Foam Oval, Wet Gel, 5/strip, 60 strip/cs

GE MEDICAL REUSABLE RESTING ECG ELECTRODES & ACCESSORIES

German silver electrodes minimize baseline drift and reduce contact resistance. GE suction electrode perfect fit for 0.156" (4mm) banana plug leadwires. Limb clamps that are latex free, easy to apply, and accepts both banana plug and snap leadwires.

5608-003	GE Gray Suction Bulb
5623-202	Metal Limb Plate Electrode
5880-004	Metal Snap Electrodes, 5/pk
5881-002	Rubber Boots For Metal Snap Electrodes, 5/pk
9033-006	Rubber Limb Strap, Latex Free
9033-014	GE Suction Electrode
9033-015	GE Suction Electrode & Gray Bulb
9033-401	Limb Clamps, 4/set
9235-001	Adhesive Tapes For Metal Snap Electrodes 100/bx




LEONHARD LANG SKINTACT® DIAGNOSTIC RESTING ECG ELECTRODES

Designed for resting 12 lead diagnostic ECG's. Easitabs are disposable eliminating the spread of cross infection associated with older style suction cups. The C-LINE technology ensures that every electrode will be fresh providing exceptional consistency in quality of recordings at very affordable prices.

- FSWA00 Diagnostic ECG Electrodes, Foam, Ultra-Tack Adhesive, Wet Gel, 500/bx
- RT-24 Diagnostic ECG Electrodes, All-Purpose, 1000/bx
- RT-41 Diagnostic ECG Electrodes, Ultra-Tack Adhesive, Oversize, 1000/bx, 4 bx/cs

Accessories


- 58100 EasiClip (Universal Alligator Clip), 10/bx
- 58727 EasiPrep (FingerTip Abrasion Strips), 500/bx



LEONHARD LANG SKINTACT® STRESS TEST & HOLTER ECG ELECTRODES

Designed to provide the highest quality tracings for stress test procedures. The offset design and Aqua-Set gel significantly reduce any motion artifacts ensuring a very clear stable trace even under the most strenuous conditions. The C-LINE technology ensures that every electrode will be fresh providing exceptional consistency in quality of recordings at very affordable prices.


- FS-40 40mm Stress Test ECG Electrode, Aqua-Wet Gel, Foam Backing, Diaphoretic, Hypoallergenic, 1200/bx
- FS-50 50mm Stress Test ECG Electrode, Aqua-Wet Gel, Snap Connector, Foam Backing, Enhanced Adhesive, 1200/bx
- FS-50C 50mm Stress Test ECG Electrode, Radio Translucent, Foam w/carbon stud, Diaphoretic, Hypoallergenic, 1200/bx
- FS-RC-4 35 x 45mm Stress Test ECG Electrode, Aqua-Wet Gel, Foam Backing, Diaphoretic, Hypoallergenic, 600/bx
- FS-TBI 50mm Stress Test ECG Electrode, Solid Gel, Lift Tab, Foam Backing, Diaphoretic, Hypoallergenic, 1500/bx
- FSVB01 36 x 53 mm Stress Test ECG Electrode, Aqua-Wet Gel, Offset Snap Connector, Foam Backing, Enhanced Adhesive, 1200/bx
- FSVB011 35 x 53 mm Stress Test ECG Electrode, Solid Gel, Offset Snap Connector, Foam Backing, Enhanced Adhesive, 1500/bx
- FS-RC-3 35 x 45mm Stress Test ECG Electrode, Aqua-Wet Gel, Foam Backing, Diaphoretic, Hypoallergenic, 3/strip, 1200/bx
- FS-RC-5 35 x 45mm Stress Test ECG Electrode, Aqua-Wet Gel, Foam Backing, Diaphoretic, Hypoallergenic, 5/strip, 1200/bx



NIKOMED DIAGNOSTIC TAB ELECTRODES

The Nikotab System with Nikomed diagnostic ECG is easy. Unique design allows for quick and easy prep. The highest quality adhesive and innovative shape prevent twisting and deliver the greatest percentage of first time recordings. Wide-body design holds Nikotab absolutely flat, securing it in place to eliminate the twist. This product is latex free.


- 0315 Economy Size Electrodes, 100/pk, 500/bx, 5000/cs
- 0515 Nikotabs E Synthetic Electrodes, Foil, 100/pk, 500/bx, 5000/cs
- 0715 Nikotabs II Synthetic Electrodes, Foil, Larger Size, 100/pk, 500/bx, 5000/cs
- 0815 Pediatric Tab Electrodes, 100/pk, 500/bx, 5000/cs



NIKOMED TRACE1™ SOLID GEL MONITORING ELECTRODES

This product is latex free.

- 1015 Solid Gel Foam Monitoring Electrode, Compact Design, Super Tack Adhesive, 32 x 38mm, 6 per strip, 60/pk, 300/bx
- 2001 Specialty Nonwoven Monitoring Electrodes, 7/8" x 7/8", 4/pk, 400/bx
- 2002 Specialty Nonwoven Monitoring Electrodes, 1 7/8" x 7/8", 4/pk, 400/bx
- 2010 Solid Gel Stress/Holter Monitoring Foam Electrodes, Eyeball Shape, 50mm, Helps Eliminate Artifacts, 30/pk, 300/bx
- 2012 Solid Gel Monitoring Electrodes, Clear Tape, Radio Translucent X-Ray Opaque, 38 x 38mm, 30/pk, 300/bx
- 2013 Solid Gel Monitoring Electrodes, Non-Woven Cloth, Helps Eliminate Artifacts, 55mm, 30/pk, 300/bx
- 2014 Solid Gel Monitoring Electrodes, Clear Tape, Multi-Purpose, 38 x 38mm, 30/pk, 300/bx
- 2015 Solid Gel Foam Monitoring Electrodes, 35 x 45mm, Helps Eliminate Artifacts, 30/pk, 300/bx
- 2016 ICU/CCU Monitoring Cloth, 55mm, 30/pk, 300/bx
- 2054 Solid Gel, Clear Tape Monitoring Electrode, Kidney Shape Design, Super Tack Adhesive, 34 x 50mm, 30/pk, 300/bx
- 5520 Solid Gel Foam Monitoring Electrode, Multi-Snap Design, Super Tack Adhesive, 32 x 45mm, 30/pk, 300/bx
- 5521 Solid Gel Foam Monitoring Electrode, Snap-Tab Design, Super Tack Adhesive, 32 x 45mm, 30/pk, 300/bx
- 5525 Solid Gel Foam Monitoring Electrode, Dynamic Off Center Design, Super Tack Adhesive, 32 x 45mm, 30/pk, 300/bx



NIKOMED TRACE1™ WET GEL MONITORING ELECTRODES

This product is latex free.

- | | |
|--------|---|
| 4500 | Wet Gel, Foam Stress/Holter Monitoring Electrode, Super Tack Adhesive, 50mm, 60/pk, 300/bx |
| 4530 | Wet Gel, Foam Direct Connect Snap Monitoring Electrode, Patient Friendly Adhesive, 32 x 45mm, 5 per strip, 60/pk, 300bx |
| 4533-5 | Wet Gel, Foam Monitoring Electrode, Super Tack Adhesive, 32 x 45mm, 5 per strip, 60/pk, 300/bx |
| 4535 | Wet Gel, Foam Direct Connect Snap Monitoring Electrode, Patient Friendly Adhesive, 32 x 45mm, 3 per strip, 60/pk, 300bx |
| 4540 | Wet Gel, Foam Direct Connect Snap Monitoring Electrode, Patient Friendly Adhesive, 32 x 45mm, 5 per strip, 60/pk, 300bx |
| 4555 | Wet Gel, Foam Monitoring Electrode, Super Tack Adhesive, 45mm, 5 per strip, 10/pk, 100/bx |
| 4560 | Wet Gel, Foam Stress Monitoring Electrode, Super Tack Adhesive, 38 x 60mm, 60/pk, 300/bx |
| 5500 | Stress/Holter Monitoring Electrodes, 50mm, Super-Tack Disphoretic Adhesive, Eyeball Shape, 30/pk, 300/bx |
| 5511 | General Purpose Monitoring Electrodes, Clear Tape, Diaphoretic Adhesive, 50mm, 30/pk, 300/bx |
| 5533-3 | General Purpose Monitoring Electrodes, 32 x 44mm, Super-Tack Diaphoretic Adhesive, 3 Per Strip, 30/pk, 300/bx |
| 5533-5 | General Purpose Monitoring Electrodes, 32 x 44mm, Super-Tack Diaphoretic Adhesive, 5 Per Strip, 30/pk, 300/bx |
| 5535 | Diagnostic Resting ECG Monitoring Electrodes, Pediatric Monitoring, 32 x 44mm, Easy-Tack Patient-Friendly Adhesive, 30/pk, 300/bx |
| 5539-3 | Thallium Stress Monitoring Electrodes, Cardiac Cath., 28 x 44mm, Super-Tack Diaphoretic Adhesive, Radio-Translucent X-Ray Opaque Snap/Sensor System, 3 Per Strip, 30/pk, 300/bx |
| 5539-5 | Thallium Stress Monitoring Electrodes, Cardiac Cath., 28 x 45mm, Super-Tack Diaphoretic Adhesive, Radio-Translucent X-Ray Opaque Snap/Sensor System, 5 Per Strip, 30/pk, 300/bx |
| 5550 | Wet Gel, Foam Monitoring Electrode, Kidney Shape Design, Super Tack Adhesive, 34 x 50mm, 30/pk, 300bx |
| 5559 | Wet Gel, Foam Monitoring Electrode, Kidney Shape Design, Radio-Translucent, Super Tack Adhesive, 34 x 50mm, 30/pk, 300bx |
| 5560 | Stress Monitoring Electrodes, 38 x 59mm, Super-Tack Disphoretic Adhesive, 30/pk, 300/bx |



CHESTER GEL WARMER

Chester Labs Gel Warmers are the perfect companion for any of the Transmission Media Gel products. Available in 2 and 3 bottle warming units, these durable warmers are made from stainless steel construction with a white powder-coat exterior. The adjustable temperature control maintains contents just above body temperature. Equipped with an illuminated on/off rocker switch and a six foot power cord, these warmers operate for just pennies a day and can be used on table tops or mounted to walls.

- | | |
|--------|--|
| 1205-2 | 2 Bottle Gel Warmer w/temperature adjustment, White Steel Case |
| 1205-3 | 3 Bottle Gel Warmer w/temperature adjustment, White Steel Case |



CHESTER LIQUASONIC® ULTRASOUND GEL

LiquaSonic Ultrasound Transmission Gel has a thicker, no-slip formula that remains viscous even when applied to salty skin. Our salt-free, formaldehyde free formula eliminates the air between the transducer and skin, offering exceptional coupling efficiency for more effective treatment/diagnosis. LiquaSonic Gel has been formulated to be acoustically correct for the broad range of frequencies used in medical ultrasound applications. Our formula is water soluble and non-greasy preventing damage to expensive transducers and clothing.

- | | |
|--------|--|
| 001205 | 8 oz bottle w/mustard cap, Ultrasound Gel, 12/cs |
| 001222 | 5 Liter Collapsible Cubitainer w/empty bottle, Ultrasound Gel, 4/cs |
| 001297 | 128 oz Multi-Purpose, Medium Viscosity Ultrasound Gel, 4/cs |
| 001298 | 8 oz Bottle w/mustard cap, Multi-Purpose, Medium Viscosity Ultrasound Gel, 12/cs |
| 001305 | 8 oz Bottle w/mustard cap, Clear Ultrasound Gel, 12/cs |
| 001333 | 128 oz Clear Ultrasound Gel, 4/cs |



CHESTER LIQUASONIC® ULTRASOUND LOTION

Similar in use to the ultrasound gel, the lotion is often preferred by medical professionals. The aqueous formulation is pleasing to the patient and is acoustically correct for the large spectrum of frequencies used in medical ultrasound applications. The lotion is hypoallergenic, bacteriostatic, salt free, formaldehyde free and will not harm clothing or equipment.

- | | |
|--------|---|
| 001267 | 128 oz Ultrasonic Lotion, 4/cs |
| 001268 | 8 oz Bottle w/mustard cap, Ultrasonic Lotion, 12/cs |



CHESTER LIQUASONIC® ULTRASOUND SCANNING GEL

This gel has been designed for diagnostic ultrasound scanning procedures requiring a less viscous formulation. Unlike other scanning materials, Ultrasound Scanning Gel possesses characteristics of a gel and fluid. Like a gel, it is applied with ease and like a fluid it has high-quality scanning abilities. The scanning gel is an acceptable alternative to mineral oil, is acoustically correct, completely aqueous and hypoallergenic.

- 001277 128 oz Scan Gel, 4/cs
- 001278 8 oz Bottle w/mustard cap, Scan Gel, 12/cs

GE MEDICAL SKIN PREP PRODUCTS & ELECTRODE ACCESSORIES

- 3704-901 Prep Razors, 100/pk
- 418195-100 GE Electrode Clips For Stress Loops, 100/pk
- 4828-004 Alcohol Pumice Prep Pads, 100/bx
- 4828-005 Sandpaper Prep Strips, 100/bx
- 4829-002 Pre-Cut Tape Strips For Stress Loops, 6/card, 100 card/pk
- 9386-001 Dry Skin Prep Pads, 200/pk

PARKER LABS REDUX® ELECTROLYTE CRÈME

For routine ECG. Cosmetic quality, pleasing to patients. Non-irritating. Highly conductive electrode cream.

- 66-04 Electrolyte Creme, 114g (4 oz) Tube, 10/bx, 5 bx/cs

PARKER LABS REDUX® ELECTROLYTE GEL

Non-abrasive, water-based gel for routine ECG and defibrillation. Highly conductive. Water soluble, non-staining, non-gritty.

- 65-04 Electrolyte Gel, 114g (4 oz) Tube, 10/bx, 5 bx/cs

PARKER LABS REDUX® ELECTROLYTE PASTE

Water-based, abrasive gel for short term monitoring, ECG and defibrillation. Highly conductive and water soluble.

- 67-05 Electrolyte Paste, 142g (5 oz) Tube, 10/bx, 5 bx/cs



PARKER LABS SIGNACREME® ELECTRODE CREAM

Highly conductive cream electrolyte for use in electromedical procedures. Bacteriostatic, non-irritating and hypoallergenic. Cosmetic quality; pleasing to patients.

- 17-05 Electrode Cream, 5 oz Dispenser, 12/bx, 6 dz/cs
- 17-20 Electrode Cream, 2 Liter (70 oz) w/2 dispensers & 1 Dispenser Pump, 4 bx/cs



PARKER LABS SIGNAGEL® ELECTRODE GEL

Highly conductive, multi-purpose electrolyte is an ideal saline electrode gel. Recommended for all electromedical procedures. Bacteriostatic, non-irritating, non-gritty. Water-soluble and non-staining.

- 15-25 Electrode Gel, 250g Tube, 6 dz/cs
- 15-60 Electrode Gel, 60g Tube, 12 dz/cs



PARKER LABS SIGNASPRAY® ELECTRODE SOLUTION & SKIN PREP

Highly conductive, multi-purpose spray electrolyte and non-gritty skin prep which provides superior conduction between skin and electrodes. For ECG and muscle stimulation. Assures clean, artifact-free ECG tracings. Bacteriostatic, hypoallergenic, non-staining. No residue build-up.


- 18-04 Signaspray Pac® Electrode Solution w/4 gallon containers, 2 Spray Dispenser Bottles, 1 Pump, 1/cs
- 18-25 Electrode Solution, 250-ml Spray Dispenser Bottle, 12/bx, 6 dz/cs
- 18-28 Electrode Solution, Gallon w/spray dispenser bottle, 4/cs



TRIAD® ULTRASOUND GEL

Conducts a broad range of frequencies, does not allow air to intrude between the transducer and skin and requires only minimal pressure to make contact. It is water soluble, non-irritating and hypoallergenic. It will not stain skin or clothing.

- 10-9208 8 oz Ultrasound Gel, 12/cs
- 10-9228 Gallon Ultrasound Gel w/dispenser, 4/cs



ADC® EKG CALIPER

A precision crafted instrument of hardened steel for accurate measurement of tracing intervals. Point protector has built-in 40mm metric scale. Measures 4 3/4"

- 394 EKG Ruler
- 395 EKG Caliper



2121



EZ-123



NIK-10

NIKOMED ECG ACCESSORIES

- 2121 Niko-Prep, 10/card, 500/cs
- 2121-E Niko-Prep, Individual, 500/cs
- EZ-123 Diagnostic Enhancement Kit (1 pk NIK-10 [choice of 0315, 0515, 0715], 1 pk NIKOtabs, 1 card 2121)
- NIK-10 Niko-Clip, Green, 10/pk
- NIK-10W Niko Clip Adapter, Fits 3-4mm, Modified, White, 10/pk



WESTERN SURGILAST® STRESS VEST

Pre-cut stretchable net material holds electrodes in place during electrocardiograms and stress tests. Does not restrict patient's circulation or breathing. No tapes or adhesives required.

- GL750 Small/Medium Stress Vest, 20/bx
- GL752 Large/X-Large Stress Vest, 20/bx



GE MEDICAL COROMETRICS® 170 SERIES FETAL MONITORS

Represents a new generation of fetal monitors designed to meet prenatal monitoring needs for the physician's office. The monitor features large LED displays allowing easy, at-a-glance confirmation of fetal heart rate and uterine activity. Offers 9-Crystal Pulsed Doppler Ultrasound technology providing a wide-beam profile for high-quality fetal heart rate detection. Includes standard Z-fold paper and an interface for the optional Model 146 Fetal Acoustic Stimulator. Available in single or twin models.

- 0171WAT-B Model 171 External Fetal Monitor w/Nautilus Ultrasound & Toco Transducers (loop style)
- 0172WAT-B Model 172 Twins External Fetal Monitor w/Nautilus Ultrasound & Toco Transducers (loop style)

GE MEDICAL FETAL CHART PAPERS

Designed for GE Corometrics monitors. Serialized pages for case documentation. Thermal array provides accurate, vivid tracings.

- 4305CAO Fetal Chart Paper, Red Grid, 30-240 BPM, Serialized Pages, Z-Fold, 160 sheets/pk, 40 pk/cs
- 4305AAO Fetal Chart Paper, Red Grid, 30-240 BPM, Serialized Pages, Z-Fold, 160/pk, 40 pk/cs
- 4305BAO Fetal Chart Paper, Red Grid, 50-210 BPM, Serialized Pages, Z-Fold, 160/pk, 40 pk/cs
- 4305DAO Fetal Chart Paper, Red Grid, 50-210 BPM, Serialized Pages, Z-Fold, 160/pk, 40 pk/cs
- 4305FAO Fetal Chart Paper, Red Grid, 50-210 BPM, Serialized Pages, Additional 1 Minute Time Marks, Z-Fold, 160/pk, 40 pk/cs
- 4306CAO Fetal Chart Paper, Green Grid, 50-210 BPM, Serialized Pages, Z-Fold, 150/pk, 40 pk/cs
- 4482AAO Fetal Chart Paper, 30-240 BPM, Serialized Pages, Pre-Loaded on Cartridges, Z-Fold, 17 cart/cs
- 4483AAO Fetal Chart Paper, Light Red Grid, 30-240 BPM, Serialized Pages, 350 sheets/rl, 6 rl/cs
- 4483BAO Fetal Chart Paper, Light Red Grid, 50-210 BPM, Serialized Pages, 350 sheets/rl, 6 rl/cs
- 4528AAO NST Patient Record Forms, For Recording From GE Corometrics Fetal Monitor Model 140, Compatible w/4482AAO recording paper, 300/cs
- 4528BAO Antepartum NST Patient Record Forms, For Recording From GE Corometrics Fetal Monitors Models 145 & 155, Compatible w/4483AAO & 4483BAO recording paper, 150/cs
- 4914AAO Chart Guard Label Packer For Model 115, 116, 118, 170, 150, 151, 120 Series Fetal Monitor Recordings, 200/cs
- 5231AAO Fetal Chart Paper, Red Grid, Serialized Pages, Z-Fold, 20 pk/cs
- 5235AAO Fetal Chart Paper, Red Grid, Serialized Pages, Z-Fold, 20 pk/cs
- 5236AAO Fetal Chart Paper, Red Grid, Serialized Pages, Z-Fold, 10 pk/cs
- 9801-012 Hewlett-Packard Fetal Chart Paper Equivalent M1910A, 2 Fetal Channels, Red Grid, Black Trace, 5.97" x 3.81", Z-Fold, 40 pk/cs
- 9801-058 Hewlett-Packard Fetal Chart Paper Equivalent M1910A, 2 Fetal Channels, Red Grid, Black Trace, 5.96" x 4.6 ft, Z-Fold, 40 pk/cs

GE MEDICAL FETAL BLOOD SAMPLING KITS

Latex free. Designed for the GE Corometrics Model 220 pH monitor. Sterile. Non-returnable and non-refundable.

- 4520AAO Fetal Blood Sampling Kit Contains: Endoscope, 6 Swabs, Lancet w/handle, 7 Heparinized Capillary Tubes, Wax Sealant, Capillary Tube Holder, 3ml Silicone Gel



MALLINCKRODT RENAISSANCE® II SPIROMETRY SYSTEM

Offers a range of sophisticated, yet user-friendly features for diagnosing respiratory function in a variety of clinical settings. Clinicians will appreciate the intuitive design and technological capabilities, including a true, real-time graphic display that enable tests to be viewed as patients perform them. Measures a broad selection of physiological parameters and provides color printout reports. It also has a 1,000-patient memory. Five-year warranty exchange program.

- P-000700-00 PB700 Portable Spirometer
- P-000710-00 Base Station
- P-000715-00 Basic Starter Kit, Includes: PB700 Spirometer, PB710 Base Station, AC Adapter, (10) FSII Flow Sensors & Syringe Adapter
- P-000720-00 Basic Starter Kit w/printer, Includes: PB700 Spirometer, PB710 Base Station, AC Adapter, (10) FSII Flow Sensors & Syringe Adapter
- P-000722-00 Kit w/printer, Includes: Basic Starter Kit & Data Flow Data Management Software, Plus Compatible Printer
- P-000725-00 Kit w/syringe, Includes: Basic Starter Kit Plus VS300 Calibration Syringe
- P-000730-00 Complete Kit, Includes: Basic Starter Kit Plus VS300 Calibration Syringe, Compatible Printer
- P-000732-00 Complete Kit, Includes Basic Starter Kit & Data Flow Data Management Software, Plus VS300 Calibration Syringe, Compatible Printer



MALLINCKRODT FLOW SENSORS

A low cost, disposable, single-patient-use flow sensor that minimizes the possibility of cross-contamination and reduces staff clean-up time.

- P-000200-00 FS200 Flow Sensor, 50/cs
- P-000201-00 FS200 Flow Sensor, 250/cs
- P-000250-00 BD250 Flow Sensor, 25/cs
- P-000251-00 BD250 Flow Sensor, 100/cs
- P-000263-00 FSII Flow Sensors, 50/cs
- P-000265-00 FSII Flow Sensors, 250/bx



MALLINCKRODT SIMPLICITY® SPIROMETRY SYSTEM

Offers a simple, affordable solution for performing basic respiratory function screening. The lightweight, handheld unit is portable and easy to use, yet performs all basic spirometry measurements. It also provides automated comparison to predicted values, pre- and post-medication comparison, interpretation of test results and comprehensive test quality control. Equipped with a 999-patient memory, and a blow and go option that allows clinicians to conduct quick, on-the-spot screening without entering patient data.

- P-000500-00 NPB-500 Spirometer
- P-000510-00 Spirometer Base Station
- P-000515-00 Starter Kit, Includes: NPB-500 Spirometer, Base, Flow Sensors, AC Adapter, Syringe Adapter & (2) Nose Clips
- P-000525-00 Starter Kit w/syringe, Includes: NPB-500 Spirometer, Base, (10) flow Sensors, Calibration Syringe, AC Adapter, Syringe Adapter, (2) Nose Clips
- P-000530-00 Complete Kit, Includes: NPB-500 Spirometer, Base, (10) Flow Sensors, Printer, Calibration Syringe, AC Adapter, Syringe Adapter, (2) Nose Clips

MALLINCKRODT SPIROMETRY ACCESSORIES

- 067100 Data Flow Data Management Software
- P-000150-01 Parallel Printer, 1/bx
- P-006548-00 Nose Clips, 25/pk
- P-010701-00 Syringe Adapter
- P-025513-SP Pressure Tube for REN II
- P-062033-00 NPB510/PB-710 Serial Cable
- P-062522-00 Renaissance Battery
- P-062596-00 Parallel Printer Cable, 1/bg
- P-495208-00 AC Adapter, 1/bx
- P-495219-00 NICAD Battery Pack
- P-000300-00 3 Liter Calibration Syringe

WELCH ALLYN PNEUMOCHECK SPIROMETER ACCESSORIES

Easy to use. PnuemoCheck will automatically calculate and display the best effort and five most needed values on its LCD screen. For your records, the printer/charger will print the data in either flow/volume or volume/time formats in as little as 20 seconds. The printer/charger recharges the PnuemoCheck.

- 56100 Thermal Paper, 5 rl/bx, 12 bx/cs
- 56120 Pediatric Mouthpieces, 100/bg, 5 bg/cs
- 56121 Adult Mouthpieces, 100/bg, 5 bg/cs
- 56130 Disposable Nose Clips, 100/bg, 4 bg/cs
- 61010 Hose for 3-L Syringe.
- 61020 3-L Calibration Syringe Only
- 61030 3-L Calibration Syringe & Hose
- 72500 7.2-V Rechargeable Batteries (2) - Replacement

WELCH ALLYN SP-2 SPIROMETER ACCESSORIES

- 26006-0001 Disposable Flow Tubes For SP-2 Only, 100/bx

WELCH ALLYN SPIROMETER FLOW TUBES

- 26001-0000 Disposable Flow Tubes, SP-150/SP-160, 10/bx
- 26001-0001 Disposable Flow Tubes, SP-150/SP-160, 100/bx

Can't find what you are looking for? Give us a call for fast, personal service.

NDD EASYONE™ DIAGNOSTIC SPIROMETRY SYSTEM

Ultrasonic flow measurement assures fast, accurate, reliable results. No moving parts, No codes to enter. Single patient use spirette has hygienic mouthpiece that reduces risk of cross infection. Alpha numeric keypad/graphic display makes navigation and data entry fast and easy. Graphic curve display allows to review curves without waiting for a printout. Meets ATS 1994 diagnostic standards. Meets NIOSH/OSHA and SSD guidelines. Customizable configuration. Automated coaching messages and QC. Automated interpretation. Two AA batteries. Stores up to 700 test sessions. Optional 8.5 x 11 or A4 printout. PC data transfer. Toll free technical support. Three-year exchange warranty. If your spirometer fails, we will replace if overnight without charge. Extended warranties are activated upon receipt of end-user information and must be received within six months of invoice date.

1000-1	EasyOne™ Diagnostic Spirometry System II Includes: #2001-2 Diagnostic Spirometry System I, #3104-4 Four Year Extended Warranty, #2040-2 EasyWare Software, #2050-1 Case of 50 Spirettes™
2000-1	EasyOne™ Frontline Spirometer Includes: Calibration Adapter & 10 Spirettes™
2001-1	EasyOne™ Diagnostic Spirometer Includes: Calibration Adapter & 10 Spirettes™
2001-2	EasyOne™ Diagnostic Spirometry System I Includes: #2001-1 Diagnostic Spirometer, #2010-4 Cradle, #2020-1Color InkJet Printer
2000-2	EasyOne™ Frontline Spirometry System Includes: #2000-1 Frontline Spirometer, #2010-4 Cradle, #2020-1Color InkJet Printer
2010-4	EasyOne™ Cradle Includes: USB Cable

NDD EASYONE™ DIAGNOSTIC SPIROMETRY SYSTEM ACCESSORIES

2020-1	Color InkJet Printer
2030-1	Calibration Adapter for Spirometer Serial #'s Below 46000
2030-2	3-Liter Calibration Syringe
2040-1	EasyWare Software for use with Parallel/Serial Cradle #2010-1
2050-1	Spirettes™, 50/cs
2050-5	Spirettes™, 200/cs
2060-2	Parallel Cable
3101-1	1 Year Warranty Extension EasyOne™ Diagnostic Spirometer (cradle included at no charge) *Extended Warranties are activated upon receipt of end-user information and must be received within 6 months of invoice date.
3101-2	2 Year Warranty Extension EasyOne™ Diagnostic Spirometer (cradle included at no charge) *Extended Warranties are activated upon receipt of end-user information and must be received within 6 months of invoice date.
3101-3	3 Year Warranty Extension EasyOne™ Diagnostic Spirometer (cradle included at no charge) *Extended Warranties are activated upon receipt of end-user information and must be received within 6 months of invoice date.
3101-4	4 Year Warranty Extension EasyOne™ Diagnostic Spirometer (cradle included at no charge) *Extended Warranties are activated upon receipt of end-user information and must be received within 6 months of invoice date.
3105-1	1 Year Warranty Extension EasyOne™ Frontline Spirometer (cradle included at no charge) *Extended Warranties are activated upon receipt of end-user information and must be received within 6 months of invoice date.
3105-2	2 Year Warranty Extension EasyOne™ Frontline Spirometer (cradle included at no charge) *Extended Warranties are activated upon receipt of end-user information and must be received within 6 months of invoice date.
3105-3	3 Year Warranty Extension EasyOne™ Frontline Spirometer (cradle included at no charge) *Extended Warranties are activated upon receipt of end-user information and must be received within 6 months of invoice date.
3105-4	4 Year Warranty Extension EasyOne™ Frontline Spirometer (cradle included at no charge) *Extended Warranties are activated upon receipt of end-user information and must be received within 6 months of invoice date.
2000-2NP	EasyOne™ Frontline Spirometry System Includes: #2000-1 Frontline Spirometer, #2010-4 Cradle, Does not include Printer
2001-2NP	EasyOne™ Diagnostic Spirometry System I Includes: #2001-1 Diagnostic Spirometer, #2010-4 Cradle, Does not include Printer
2001-2S	EasyOne™ Screen Spirometry System Includes: #2001-1 Diagnostic Spirometer, 2010-6S, EasyOne screen
2010-6S	EasyOne™ Screen, Includes: Screen Connector & #2040-2 EasyWare
2030-3	Calibration Adapter for Spirometer Serial #'s 46000 and above
2031-1	EasyOne™ Carry Case
2040-2	EasyWare Software for use with USB Cradle p/n 2010-4

**BURDICK SENSARE HANDHELD SPIROMETER**

Full functioning ATS compliant spirometer that weighs less than 13 ounces. User friendly touch screen provides quick navigation and easy interpretation. Base station allows easy battery recharging and data transfer. Print reports by connecting the spirometer or base station to a computer.

92524-01 Sensaire Handheld Spirometer, English Version, 120V

**BURDICK SPIROCARD SPIROMETER**

Features patented, pre-calibrated, disposable patient filters. Performs all the tests you would expect from a conventional device. Interface to Office Medic software for real-time data access, reporting and archival. Compatible with your favorite off-the-shelf PC, laptop, tablet PC or Compaq iPAQ.

92526-01 SpiroCard PC Card Spirometer, English

BURDICK SPIROMETRY ACCESSORIES

010-1600-00	Sensaire Rechargeable Battery Pack
010-1601-00	Sensaire User's Manual
010-1604-00	Sensaire Mouthpieces, 100/bx
010-1605-00	SpiroCard/SpiroCard Mouthpieces, 100/bx
010-1606-00	SpiroCard/SpiroCard Pressure Tubes, 25/bg
010-1607-00	OxiCard Adult Fingercip Sensors
010-1615-00	Sensaire Pressure Tube For Calibration Use Only, 10/pkg
010-1616-00	Sensaire Stylus, 10/pkg
010-1628-00	OxiCard Pediatric Fingercip Sensors
043231	Video "How To Operate Flash"
043232	Video "How To Operate SpiroSense"
043245	Video: How To Operate Spirotouch
086613	Spirometry Book: A Better Way Mouthpiece Adapter For Calibration Syringe SpiroTouch, Sensaire, SpiroCard & SpiroCard
100022	3.0L Calibration Syringe
100049	Presto™ Thermal Paper, 6 r/cs
100052	RJ-11 Cable For Use w/Presto
10005201	DB-9 Cable For Use w/SpiroSense
100160	Presto Operator Manual
100161	Plastic Mouthpieces, 100/pkg
100162	Noseclips, 5/pkg
100163	12" Calibration Hose For Use w/SpiroSense
10021718	36" Hose & Adapter For Use w/Presto & SpiroSense
100219	Flow Sensor Bracket For Use w/Presto & SpiroSense
100220	Flow Sensor Insert
100243	Cardboard Mouthpieces (Presto & SpiroSense), 100/pkg
100245	SpiroSense Operator Manual
100327	Flash Memory Card, 256K Memory (40 Tests)
100338	Flash Operator Manual
100340	Flash Battery Recharger
100351	Flash Carry Case
100362	Flash Cardboard Mouthpieces, 100/pkg
100363	Flash Flow Sensor Hose Replacement
100381	Flash Flow Sensor Replacement
100414	SpiroSense Battery Booster
100611	Flash Compatible Printer w/cable
100612	Mouthpiece, 25/pkg
100613	Flow Sensor Insert For Spirotouch
103-0208-00	Mouthpiece Adapter for Calibration Syringe SpiroTouch, Sensaire, SpiroCard & SpiroCard



BURDICK SPIROX CARD SPIROMETER/PULSE OXIMETER

Combines a spirometer and a pulse oximeter on a single PCMCIA card providing two diagnostic tests in one portable package. Performs all measurements performed on conventional devices. Features patented, disposable and biodegradable patient filters. Devices interface to Office Medic software for real-time data access, reporting and archival. Compatible with your favorite off-the-shelf PC, laptop, tablet PC or Compaq iPAQ.

92525-01 SpirOx Card PC Card Spirometer/Pulse Oximeter, English



MIDMARK IQMARK™ DIGITAL SPIROMETER (SERIAL VERSION)

Includes: IQmark™ Digital Spirometer (serial version), IQmark™ Digital Spirometer Software Package with trending, IQmark™ Diagnostic Workstation - Patient Management Software, 10-pack disposable pneumotach mouthpieces, 10-pack disposable nose clips, 4-pack AAA batteries, mouse pad, operation manuals and quick reference guide.

4-000-0001 IQmark™ Digital Spirometer (Serial Version)



MIDMARK IQMARK™ DIGITAL SPIROMETER (SERIAL VERSION) W/3-LITER CALIBRATION SYRINGE

Includes: IQmark™ Digital Spirometer (Serial Version), IQmark™ Digital Spirometer Software Package with Trending, IQmark™ Diagnostic Workstation - Patient Management Software, IQmark™ 3-Liter Calibration Syringe, 10-pack disposable pneumotach mouthpieces, 10-pack disposable nose clips, 4-pack AAA batteries, mouse pad, operation manuals and quick reference guide.

4-000-0002 IQmark™ Digital Spirometer (Serial Version) w/IQmark™ 3-liter calibration syringe



MIDMARK IQMARK™ DIGITAL SPIROMETER (USB VERSION)

Includes: IQmark™ Digital Spirometer (USB Version), IQmark™ Digital Spirometer Software Package with trending, IQmark™ Diagnostic Workstation - Patient Management Software, 10-pack disposable pneumotach mouthpieces, 10-pack disposable nose clips, mouse pad, operation manuals and quick reference guide.

4-000-0004 IQmark™ Digital Spirometer (USB Version)



MIDMARK IQMARK™ DIGITAL SPIROMETER (USB VERSION) W/3-LITER CALIBRATION SYRINGE

Includes: IQmark™ Digital Spirometer (USB Version), IQmark™ Digital Spirometer Software Package with Trending, IQmark™ Diagnostic Workstation - Patient Management Software, IQmark™ 3-Liter Calibration Syringe, 10-pack disposable pneumotach mouthpieces, 10-pack disposable nose clips, mouse pad, operation manuals and quick reference guide.

4-000-0003 IQmark™ Digital Spirometer (USB Version)

MIDMARK IQMARK™ DIGITAL SPIROPDA (PDA VERSION) W/3-LITER CALIBRATION SYRINGE

Includes: IQmark™ Digital Spirometer (PDA Version), IQmark™ SpiroPDA Digital Spirometer Software Package, IQmark™ Digital Spirometer Software Package with Trending, IQmark™ PDA/IQmark™ Sync Software Package, IQmark™ Diagnostic Workstation - Patient Management Software, IQmark™ 3-Liter Calibration Syringe, 10-pack disposable pneumotach mouthpieces, 10-pack disposable nose clips, PDA Cable (PDA brand and model number required), USB to serial port adapter, 4-pack AAA batteries, mouse pad, operation manuals and quick reference guide (Does NOT include a Pocket-PC device, i.e. an iPAQ).

4-000-0007 IQmark™ Digital SpiroPDA (PDA Version) w/IQmark™ 3-Liter Calibration Syringe

MIDMARK IQMARK™ SPIROMETRY SUPPLIES

- 1-100-0003 IQmark™ 3-Liter Calibration Syringe
- 2-100-1200 Disposable Pneumotach Mouthpieces, 25/bx
- 2-244-0001 Disposable Nose Clips, 25/bx
- 3-100-8500 Syringe Adapter
- 3-999-0005 16-pack AAA Batteries
- 2-100-1201 Disposable Pneumotach Mouthpieces, 25/bx, 8 bx/cs

MIDMARK IQMARK™ MANUALS

3-100-1025 IQmark™ Digital Spirometer Manual



PERMAXIM REDISCREEN™ STREP A TEST THROAT SWAB

97% sensitivity. Equivalent results to a traditional swab culture. Add sample, wait five minutes and read results. Kit includes reagents, controls, tubes, swabs, and test strips. Remove strip from foil pouch. Add reagents to tube. Swab throat and add to reagent solution. Place test strip into solution. Read results in five minutes.

2033 Strep A Test Throat Swab Dipstick, CLIA Waived, 25/bx



QUIDEL QUICKVUE® DIPSTICK STREP A TEST

A rapid qualitative test for the detection of Strep A antigen directly from a throat swab. Results in 5 minutes. Room temperature storage CLIA waived.

20108 Dipstick Strep A Test, 50 test/kit



QUIDEL QUICKVUE® IN-LINE® ONE-STEP GROUP A STREP KIT

Allows for the rapid detection of group A streptococcal antigen directly from patient throat swab specimens. The unique in-line extraction and rapid procedure allow testing to take place at the point of care. Two-color results give you clear readability, and the proven accuracy of QuickVue provides dependable results. The kit stores at room temperature and includes positive and negative external controls. CLIA waived.

0343 Strep A Test, 25 test kit, 12 kits/cs
 0345 Strep A Control Swabs, 1/kit, 10 kits/cs
 0347 Strep A Swab Pack, 75 swabs/pk



PERMAXIM REDISCREEN™ URINE HCG PREGNANCY CASSETTE DEVICE

99% clinical accuracy. 25 mIU/ml sensitivity. Some positive results can be seen in as soon as one minute. One step, no reagents. Remove cassette from foil pouch and lay on flat surface. Add three drops of urine to the sample well. Read results in three minutes for urine.

2036 Urine hCG Pregnancy Cassette Device, CLIA Waived, 25/bx
 2036-100 Urine hCG Pregnancy Cassette Device, 25/bx, 4bx/cs



PERMAXIM REDISCREEN™ URINE HCG PREGNANCY TEST STRIPS

99% clinical accuracy. 25 mIU/ml sensitivity. Some positive results can be seen as soon as one minute. One step, no reagents. Remove strip device from canister. Immerse the strip in urine sample. Read results in three minutes.

2032 Urine hCG Pregnancy Test Strips, CLIA Waived, 100/bx



PERMAXIM REDISCREEN™ URINE/SERUM HCG PREGNANCY CASSETTE DEVICE

99% clinical accuracy. 25 mIU/ml sensitivity. Some positive results can be seen in as soon as one minute. One step, no reagents. Remove cassette from foil pouch and lay on flat surface. Add three drops of urine or serum to the sample well. Read results in three minutes for urine and five minutes for serum.

2012 Urine/Serum hCG Pregnancy Cassette Device, CLIA Waived, 25/bx



QUIDEL RAPIDVUE® HCG KIT

Easy to interpret dipstick. Results in 3 minutes. Each test is individually foil wrapped. Room Temperature storage. CLIA waived.

0257 HCG Kit, 25 test/kit, 10 kits/cs



0317



20183

QUIDEL QUICKVUE® INFLUENZA A/B TESTS

Designed specifically for exam room use providing test results in ten minutes or less. Test accurately detects Influenza type A and B viral antigens in three easy steps from a nasal swab or nasal aspirate/nasal wash specimen. The A+B test identifies type A or type B in positive samples. The easy-to-use, easy-to-read kit can be stored at room temperature. No instrumentation is required. The 10-test kit is conveniently packaged for single-use testing and disposal. The 25-test kit offers bulk packaging that is ideal for batch testing. Both the 10-test and 25-test kits include everything to perform the test including positive and negative controls.

- 0317 Influenza A/B Test, Identifies Positive or Negative, 25 test/kit, 10 kits/cs
- 0324 Influenza A/B Test, Identifies Positive or Negative, 10 test/kit, 10 kits/cs
- 20171 Transport Swab Set, 25 swabs/set, 40 set/pk
- 20183 QuickVue® Influenza A+B Test, Dipstick Format, Identifies Type A, Type B, or Both, Two-Color Endpoint, 25 test/kit, 10 kits/cs



2080



2081

PERMAXIM REDISCREEN™ ULTRA GUAIAIC OCCULT BLOOD TEST/ACCESSORIES

- 2080 Patient Pack, comes with 100 Triple Slide Tests, 6 x 10mL Developers, Applicators, Foil Lined Envelopes, CLIA Waived, 50/pk, 2 pk/bx
- 2081 Lab Pack, comes with 100 Single Slide Tests, 2 x 10mL Developers, Applicators, CLIA Waived, 100 test/bx
- 2084 Office Pack, comes with 34 Triple Slide Tests, 2 x 10mL Developers, Applicators, 34/bx

Accessories

- 2082 8 x 10mL Developers and Applicators, CLIA Waived, 1/bx
- 2083 2 x 100 Test Tape, 2 x 10mL Developers, Applicators, CLIA Waived, 200/bx



QUIDEL QUICKVUE® IFOB TEST

- 20194 Kit-Cassette Format, 20 Test Kit, Includes: 1 Pkg Insert, 1 Procedure Card, 20 Specimen Collection & return Kits Each Containing 1 Collection Tube w/2mL FOB Buffer, 1 Collection Paper w/adhesive, 1 Pouch, 1 Absorbent Sleeve, 1 Return Mail Pouch & 1 Patient Instructions
- 20196 Specimen and Collection Kit, 10 Specimen & Return Kits Each Containing 1 Collection Tube w/2mL FOB Buffer, 1 Collection Paper w/adhesive, 1 Pouch, 1 Absorbent Sleeve, 1 Return mail Pouch & 1 Patient Instructions
- 20197 External Controls, Positive & Negative



555001




555050

HYOGUARD ASSURE™ 3 BLOOD GLUCOSE MONITORING SYSTEM

Makes testing three times faster. Uses biosensor technology. Calibration is easy - simply insert the code chip and the meter is ready to test. Testing is easy - insert a test strip and the meter turns on. Apply sample and obtain plasma-referenced result in 10 seconds. Uses a specially designed test strip. The test strip with ergonomic Guide-me-Curve™ facilitates blood application. Capillary action wicks the sample into the test chamber. Features: smaller sample - uses only 3ul blood sample, automatic on/off with insertion and removal of test strip, and wide test range 30-550 mg/dL.

- 555001 Assure 3 Meter Only
- 555006 Assure 3 Control Solution Level 1&2
- 555008 Assure 3 Check Strip
- 555050 Assure 3 Test Strips w/code chip, 50/btl
- 555100 Assure 3 Test Strips w/code chip, 100/btl



SURGILANCE™ SAFETY LANCETS

Needle is concealed before and after use. It is automatically rendered inoperative after a single use. No arming is required and the lancets are color-coded for easy recognition. The lancet is triggered by a gentle push on the test site area, giving optimum blood flow at shallow depths. High-speed penetration minimizes pain. The hidden needle and quiet operation helps increase patient comfort.

- SLB200 Blade, 1.8mm Depth, 75 to 100uL Blood Sample, Green, 100/bx
- SLB250 Blade, 2.3mm Depth, 150 to 200uL Blood Sample, Blue, 100/bx
- SLN200 Needle, 1.8mm Depth, 10 to 20uL Blood Sample, Gray, 100/bx
- SLN240 Needle, 2.2mm Depth, 20 to 40uL Blood Sample, Orange, 100/bx
- SLN300 Needle, 2.8mm Depth, 40 to 60uL Blood Sample, Pink, 100/bx

**See Laboratory Section
for a complete listing
of physician's office lab
diagnostic products.**