



N407015

N407025




N407035

PRO ADVANTAGE® SCISSORS

Expertly crafted from surgical grade stainless steel and meets the very highest quality standards and specifications. This outstanding quality and price point makes Pro Advantage instruments ideal for today's growing instrument market.

N407000	OR Scissors, 5 1/2", Straight, Sharp/Blunt
N407002	OR Scissors, 5 1/2", Straight, Sharp/Sharp
N407005	Iris Scissors, 4", Straight
N407010	Iris Scissors, 4", Curved
N407015	Iris Scissors, 4 1/2", Straight
N407020	Iris Scissors, 4 1/2", Curved
N407025	Lister Bandage Scissors, 5 1/2"
N407030	Lister Bandage Scissors, 7 1/4"
N407035	Bandage & Utility Scissors, 7 1/2", Serrated Blade
N407037	Knowles Bandage Scissors, 5"
N407040	Spencer Stitch Scissors, 3 1/2", Delicate
N407045	Littauer Jr. Stitch Scissors, 4 1/2", Light Pattern
N407050	Littauer Stitch Scissors, 5 1/2"
N407060	Esmarch Plaster Shears, 8"



ADC SCISSORS

Forged from corrosion resistant surgical stainless steel. These popular floor grade instruments are built to last.

300	Lister Bandage Scissors, 4 1/2"
3000G	Lister Bandage Scissors, Gold, 4 1/2"
3006	Lister Bandage Scissors w/clip, 4 1/2"
3007	Lister Bandage Scissors w/clip, 5 1/2"
301	Lister Bandage Scissors, 5 1/2"
3010G	Lister Bandage Scissors, Gold, 5 1/2"
302	Lister Bandage Scissors, 7 1/2"
308	Littauer Stitch Scissors, 5 1/2"
323BD	Listerette Scissors, 5 1/2", Burgundy
323BK	Listerette Scissors, 5 1/2", Black
323DG	Listerette Scissors, 5 1/2", Dark Green
323M	Listerette Scissors, 5 1/2", Magenta
323N	Listerette Scissors, 5 1/2", Navy
323RB	Listerette Scissors, 5 1/2", Royal Blue
323TL	Listerette Scissors, 5 1/2", Teal
323V	Listerette Scissors, 5 1/2", Purple
3404SB	Operating Scissors, Straight, Sharp/Blunt, 5 1/2"
3404SS	Operating Scissors, Straight, Sharp/Sharp, 5 1/2"
3410	Mayo Dissecting Scissors, Straight, 5 1/2"
3424	Iris Scissors, Straight, 4 1/2"
3425	Iris Scissors, Curved, 4 1/2"
3656	Umbilical Cord Scissors, 4"



5-1000

5-998

MILTEX ALL PURPOSE HOSPITAL SHEARS

Plastic handles. Autoclavable.

5-1000	7 1/2" Hospital Shears, Serrated
5-998	6 1/2" Hospital Shears, Serrated



5-316

5-318



5D-326

MILTEX BROPHY SCISSORS

5-316	5 1/2" Brophy Scissors, Straight, Sharp Points
5-318	5 1/2" Brophy Scissors, Curved, Sharp Points
5D-326	5 1/2" Brophy Scissors, Curved, One Serrated Blade



MILTEX EDGE™ SCISSORS SHARPENER

Quickly puts a sharp edge on dull scissors (not intended for Supercut or serrated edge scissors). Autoclavable.

5-5000	Hand Held Scissors Sharpener
--------	------------------------------




5-560

5-561

MILTEX KNOWLES BANDAGE SCISSORS

5-560	5 1/2" Bandage Scissors
5-561	Angled, 5 1/2" Bandage Scissors



MILTEX LISTER BANDAGE SCISSORS

Extra fine quality.

- 5-510 3 1/2" Bandage Scissors
- 5-512 4 1/2" Bandage Scissors
- 5-514 5 1/2" Bandage Scissors
- 5-516 7 1/4" Bandage Scissors
- 5-550 8" Bandage Scissors w/one large finger ring



MILTEX EYE SUTURE SCISSORS

Slightly curved with sharp pointed short blades.

- 18-1652 3 3/4" Eye Suture Scissors, Gradle
- 18-1653 4" Eye Suture Scissors, Ribbon Type



MILTEX KLEINERT-KUTZ TENOTOMY SCISSORS


Curved and rounded dissecting blades.

- 21-539 4 7/8" Tenotomy Scissors



MILTEX LITTAUER JUNIOR STITCH SCISSORS

- 9-102 4 1/2" Stitch Scissors, Slender Pattern




MILTEX IRIS SCISSORS

- 18-1396 3 1/2" Straight Iris Scissors, Delicate w/20mm blades
- 18-1398 3 1/2" Curved Iris Scissors, Delicate w/20mm blades
- 18-1404 4" Straight Iris Scissors, Heavy Pattern
- 18-1405 4" Curved Iris Scissors, Heavy Pattern
- 18-1415 4" Straight Iris Scissors, Standard, Ribbon Type
- 18-1416 4" Curved Iris Scissors, Standard, Ribbon Type
- 18-1417 3 3/4" Straight Iris Scissors, Miniature Pattern w/14mm blades, Ribbon Type
- 18-1418 3 3/4" Curved Iris Scissors, Miniature Pattern w/14mm blades, Ribbon Type
- 5-300 4" Straight Iris Scissors
- 5-302 4" Curved Iris Scissors
- 5-304 4 1/2" Straight Iris Scissors
- 5-304TC 4 1/2" Straight Iris Scissors, Carb-N-Sert
- 5-306 4 1/2" Curved Iris Scissors
- 5-306TC 4 1/2" Curved Iris Scissors, Carb-N-Sert
- 5-308 4 1/2" Iris Scissors, Angled On Side
- 5-SC-304 4 1/2" Straight Iris Scissors, SuperCut Blade
- 5-SC-304TC 4 1/2" Straight Iris Scissors, SuperCut Carb-N-Sert Blades
- 5-SC-306 4 1/2" Curved Iris Scissors, SuperCut Blade
- 5-SC-306TC 4 1/2" Curved Iris Scissors, SuperCut Carb-N-Sert Blades
- 5D-303 4 1/2" Curved Iris Scissors (curved sideways)




MILTEX LITTAUER STITCH SCISSORS

- 9-104 5 1/2" Stitch Scissors



MILTEX MAYO DISSECTING SCISSORS

5-120	5 1/2" Straight Dissecting Scissors, Standard Beveled Blades
5-120TC	5 1/2" Straight Dissecting Scissors, Carb-N-Sert, Standard Beveled Blades
5-122	5 1/2" Curved Dissecting Scissors, Standard Beveled Blades
5-122TC	5 1/2" Curved Dissecting Scissors, Carb-N-Sert, Standard Beveled Blades
5-124	6 3/4" Straight Dissecting Scissors, Standard Beveled Blades
5-124TC	6 3/4" Straight Dissecting Scissors, Carb-N-Sert, Standard Beveled Blades
5-126	6 3/4" Curved Dissecting Scissors, Standard Beveled Blades
5-126TC	6 3/4" Curved Dissecting Scissors, Carb-N-Sert, Standard Beveled Blades
5-128	9" Straight Dissecting Scissors, Standard Beveled Blades
5-128TC	9" Straight Dissecting Scissors, Carb-N-Sert, Standard Beveled Blades
5-130	9" Curved Dissecting Scissors, Standard Beveled Blades
5-130TC	9" Curved Dissecting Scissors, Carb-N-Sert, Standard Beveled Blades
5-136	5 1/2" Straight Dissecting Scissors, Rounded Blades
5-136TC	5 1/2" Straight Dissecting Scissors, Carb-N-Sert, Rounded Blades
5-138	5 1/2" Curved Dissecting Scissors, Rounded Blades
5-138TC	5 1/2" Curved Dissecting Scissors, Carb-N-Sert, Rounded Blades
5-140	6 3/4" Straight Dissecting Scissors, Rounded Blades
5-140TC	6 3/4" Straight Dissecting Scissors, Carb-N-Sert, Rounded Blades
5-142	6 3/4" Curved Dissecting Scissors, Rounded Blades
5-142TC	6 3/4" Curved Dissecting Scissors, Carb-N-Sert, Rounded Blades
5-SC-120	5 1/2" Straight Dissecting Scissors, SuperCut Blade
5-SC-122	5 1/2" Curved Dissecting Scissors, SuperCut Blade
5-SC-124	6 3/4" Straight Dissecting Scissors, SuperCut Blade
5-SC-126	6 3/4" Curved Dissecting Scissors, SuperCut Blade
5-SC-128	9" Straight Dissecting Scissors, SuperCut Blade
5-SC-130	9" Curved Dissecting Scissors, SuperCut Blade




MILTEX METZENBAUM SCISSORS

5-179	5 1/2" Straight Metzenbaum Scissors, Standard
5-179A-TC	5 1/2" Straight Metzenbaum Scissors, Delicate, Carb-N-Sert, Blunt Points
5-179B-TC	5 1/2" Straight Metzenbaum Scissors, Delicate, Carb-N-Sert, Sharp Points
5-179TC	5 1/2" Straight Metzenbaum Scissors, Standard, Carb-N-Sert, Blunt Points
5-180	5 1/2" Curved Metzenbaum Scissors, Standard
5-180A-TC	5 1/2" Curved Metzenbaum Scissors, Delicate, Carb-N-Sert, Blunt Points
5-180B-TC	5 1/2" Curved Metzenbaum Scissors, Delicate, Carb-N-Sert, Sharp Points
5-180TC	5 1/2" Curved Metzenbaum Scissors, Standard, Carb-N-Sert, Blunt Points
5-181	7" Straight Metzenbaum Scissors, Standard
5-181A-TC	7" Straight Metzenbaum Scissors, Delicate, Carb-N-Sert, Blunt Points
5-181B-TC	7" Straight Metzenbaum Scissors, Delicate, Carb-N-Sert, Sharp Points

Scissors

MILTEX METZENBAUM SCISSORS (Continued)

5-181TC	7" Straight Metzenbaum Scissors, Standard, Carb-N-Sert, Blunt Points
5-182	7" Curved Metzenbaum Scissors, Standard
5-182A-TC	7" Curved Metzenbaum Scissors, Delicate, Carb-N-Sert, Blunt Points
5-182B-TC	7" Curved Metzenbaum Scissors, Delicate, Carb-N-Sert, Sharp Points
5-182TC	7" Curved Metzenbaum Scissors, Standard, Carb-N-Sert, Blunt Points
5-183	8" Straight Metzenbaum Scissors, Delicate
5-184	8" Curved Metzenbaum Scissors, Standard
5-184TC	8" Curved Metzenbaum Scissors, Standard, Carb-N-Sert, Blunt Points
5-185	8" Curved Metzenbaum Scissors, Delicate
5-186	9" Straight Metzenbaum Scissors, Standard
5-186TC	9" Straight Metzenbaum Scissors, Standard, Carb-N-Sert, Blunt Points
5-187	9" Straight Metzenbaum Scissors, Delicate
5-187TC	9" Straight Metzenbaum Scissors, Delicate, Carb-N-Sert, Blunt Points
5-188	9" Curved Metzenbaum Scissors, Standard
5-188TC	9" Curved Metzenbaum Scissors, Standard, Carb-N-Sert, Blunt Points
5-189	9" Curved Metzenbaum Scissors, Delicate
5-189A-TC	10" Curved Metzenbaum Scissors, Standard, Carb-N-Sert, Blunt Points
5-189TC	9" Curved Metzenbaum Scissors, Delicate, Carb-N-Sert, Blunt Points
5-190	11" Straight Metzenbaum Scissors, Standard
5-192	11" Curved Metzenbaum Scissors, Standard
5-192TC	11" Curved Metzenbaum Scissors, Standard, Carb-N-Sert, Blunt Points
5-193	11" Curved Metzenbaum Scissors, Delicate
5-193A-TC	12" Curved Metzenbaum Scissors, Standard, Carb-N-Sert, Blunt Points
5-194	14 1/2" Straight Metzenbaum Scissors, Standard
5-196	14 1/2" Curved Metzenbaum Scissors, Standard
5-283	5" Straight Metzenbaum Scissors, Delicate
5-284	5" Curved Metzenbaum Scissors, Delicate
5-286	6" Curved Metzenbaum Scissors, Delicate
5-288	7" Curved Metzenbaum Scissors, Delicate
5-SC-180	5 1/2" Curved Metzenbaum Scissors, SuperCut Blade
5-SC-182	7" Curved Metzenbaum Scissors, SuperCut Blade
5-SC-182TC	7" Curved Metzenbaum Scissors, SuperCut, Carb-N-Sert Blades
5-SC-188	9" Curved Metzenbaum Scissors, SuperCut Blade
5-SC-283	4 1/2" Straight Metzenbaum Scissors, Delicate, SuperCut Blade
5-SC-284	4 1/2" Curved Metzenbaum Scissors, Delicate, SuperCut Blade



MILTEX MAYO OPERATING SCISSORS

Heavy pattern. Sharp points.

5-115	5 1/2" Operating Scissors, Straight
5-116	5 1/2" Operating Scissors, Curved



MILTEX NORTHBENT STITCH SCISSORS

Curved. Light pattern.

9-103	4 3/4" Stitch Scissors
-------	------------------------



MILTEX NOYES IRIS SCISSORS

18-1510	4 1/2" Iris Scissors, Straight, Sharp-Sharp Points
18-1512	4 1/2" Iris Scissors, Straight, Sharp-Blunt Points
18-1514	4 1/2" Iris Scissors, Straight, Blunt-Blunt Points
18-1516	4 1/2" Iris Scissors, Curved, Sharp-Sharp Points



MILTEX NURSE'S SCISSORS

Feature fluoride coated blades to prevent adhesive from sticking to cutting surface. 5 1/2" long. Autoclavable.

5-700	Nurses Scissors, Pink Handles & Guard w/carry case
5-702	Nurses Scissors, Blue Handles & Guard w/carry case



MILTEX BANDAGE & UTILITY SCISSORS

Non-stick fluoride coating over stainless steel blades. Neon-green plastic finger rings for easy identification. Cuts gauze, tubing, clothing, leather, light casts and wet fiberglass casting material.

5-800	8 1/2" Utility Scissors
-------	-------------------------



MILTEX O'BRIEN STITCH SCISSORS

9-110	3 3/4" Stitch Scissors, Angled, Sharp Points
-------	--



5-270



5-271



5-272




5-273

MILTEX PLASTIC SURGERY SCISSORS

5-270	4 3/4" Straight Surgery Scissors, Sharp-Sharp Points
5-271	4 3/4" Straight Surgery Scissors, Sharp-Sharp Points, Serrated Blade
5-272	4 3/4" Straight Surgery Scissors, Sharp-Blunt Points
5-273	4 3/4" Straight Surgery Scissors, Sharp-Blunt Points, Serrated Blade
5-274	4 3/4" Straight Surgery Scissors, Blunt-Blunt Points
5-275	4 3/4" Straight Surgery Scissors, Blunt-Blunt Points, Serrated Blade
5-276	4 3/4" Curved Surgery Scissors, Sharp-Sharp Points
5-277	4 3/4" Curved Surgery Scissors, Sharp-Sharp Points, Serrated Blade
5-278	4 3/4" Curved Surgery Scissors, Sharp-Blunt Points
5-279	4 3/4" Curved Surgery Scissors, Sharp-Blunt Points, Serrated Blade
5-280	4 3/4" Curved Surgery Scissors, Blunt-Blunt Points
5-281	4 3/4" Curved Surgery Scissors, Blunt-Blunt Points, Serrated Blade
5-282	4 3/4" Surgery Scissors, Angled On Side, Sharp-Sharp Points
5-282A	4 3/4" Surgery Scissors, Angled On Side, Sharp-Sharp Points, Serrated Blade



18-1460



18-1462



18-1464



18-1466



18-1470




18-1472




18-1474



18-1476



18-SC-1472



18-SC-1476

MILTEX STEVENS TENOTOMY SCISSORS

18-1460	4 1/8" Straight Stevens Tenotomy Scissors, Short Blades, Sharp Points
18-1462	4 1/8" Straight Stevens Tenotomy Scissors, Short Blades, Blunt Points
18-1464	4 1/8" Curved Stevens Tenotomy Scissors, Short Blades, Sharp Points
18-1466	4 1/8" Curved Stevens Tenotomy Scissors, Short Blades, Blunt Points
18-1470	4 1/2" Straight Stevens Tenotomy Scissors, Long Blades, Sharp Points
18-1472	4 1/2" Straight Stevens Tenotomy Scissors, Long Blades, Blunt Points
18-1474	4 1/2" Curved Stevens Tenotomy Scissors, Long Blades, Sharp Points
18-1476	4 1/2" Curved Stevens Tenotomy Scissors, Long Blades, Blunt Points
18-1477	3 3/4" Straight Stevens Tenotomy Scissors, Blunt Points, Ribbon Type
18-1478	3 3/4" Curved Stevens Tenotomy Scissors, Blunt Points Ribbon Type
18-SC-1472	4 1/2" Straight Stevens Tenotomy Scissors, SuperCut Blade
18-SC-1476	4 1/2" Curved Stevens Tenotomy Scissors, SuperCut Blade



MILTEX RAGNELL DISSECTING SCISSORS

Flat tip and curved.

5-290	5" Dissecting Scissors
5-291	7" Dissecting Scissors



5-312



5-312TC



5-314



5-314TC

MILTEX STRABISMUS SCISSORS

18-1453	4 1/4" Strabismus Scissors, Straight, Ribbon Type
18-1454	4 1/4" Strabismus Scissors, Curved, Ribbon Type
5-312	4" Strabismus Scissors, Straight
5-312TC	4 1/2" Strabismus Scissors, Straight, Carb-N-Sert
5-314	4" Strabismus Scissors, Curved
5-314TC	4 1/2" Strabismus Scissors, Curved, Carb-N-Sert




MILTEX SPENCER STITCH SCISSORS

9-100	3 1/2" Stitch Scissors
-------	------------------------

Can't find what
you're looking for?
Call your
representative.

Scissors



MILTEX STANDARD PATTERN OPERATING SCISSORS

5-12	4 1/2" Operating Scissors, Sharp-Blunt Points, Straight
5-14	5" Operating Scissors, Sharp-Blunt Points, Straight
5-16	5 1/2" Operating Scissors, Sharp-Blunt Points, Straight
5-16TC	5 1/2" Operating Scissors, Sharp-Blunt Points, Straight, Carb-N-Sert
5-17	6" Operating Scissors, Sharp-Blunt Points, Straight
5-18	6 1/2" Operating Scissors, Sharp-Blunt Points, Straight
5-19	7 1/2" Operating Scissors, Sharp-Blunt Points, Straight
5-2	4 1/2" Operating Scissors, Sharp-Sharp Points, Straight
5-22	4 1/2" Operating Scissors, Blunt-Blunt Points, Straight
5-24	5" Operating Scissors, Blunt-Blunt Points, Straight
5-26	5 1/2" Operating Scissors, Blunt-Blunt Points, Straight
5-26TC	5 1/2" Operating Scissors, Blunt-Blunt Points, Straight, Carb-N-Sert
5-27	6" Operating Scissors, Blunt-Blunt Points, Straight
5-28	6 1/2" Operating Scissors, Blunt-Blunt Points, Straight
5-32	4 1/2" Operating Scissors, Sharp-Sharp Points, Curved
5-34	5" Operating Scissors, Sharp-Sharp Points, Curved
5-36	5 1/2" Operating Scissors, Sharp-Sharp Points, Curved
5-36TC	5 1/2" Operating Scissors, Sharp-Sharp Points, Curved, Carb-N-Sert
5-37	6" Operating Scissors, Sharp-Sharp Points, Curved
5-38	6 1/2" Operating Scissors, Sharp-Sharp Points, Curved
5-4	5" Operating Scissors, Sharp-Sharp Points, Straight
5-42	4 1/2" Operating Scissors, Sharp-Blunt Points, Curved
5-44	5" Operating Scissors, Sharp-Blunt Points, Curved
5-46	5 1/2" Operating Scissors, Sharp-Blunt Points, Curved
5-46TC	5 1/2" Operating Scissors, Sharp-Blunt Points, Curved, Carb-N-Sert
5-48	6 1/2" Operating Scissors, Sharp-Blunt Points, Curved
5-49	7 1/2" Operating Scissors, Sharp-Blunt Points, Curved
5-52	4 1/2" Operating Scissors, Blunt-Blunt Points, Curved
5-54	5" Operating Scissors, Blunt-Blunt Points, Curved
5-56	5 1/2" Operating Scissors, Blunt-Blunt Points, Curved
5-56TC	5 1/2" Operating Scissors, Blunt-Blunt Points, Curved, Carb-N-Sert
5-58	6 1/2" Operating Scissors, Blunt-Blunt Points, Curved
5-6	5 1/2" Operating Scissors, Sharp-Sharp Points, Straight
5-6TC	5 1/2" Operating Scissors, Sharp-Sharp Points, Straight, Carb-N-Sert
5-7	6" Operating Scissors, Sharp-Sharp Points, Straight
5-8	6 1/2" Operating Scissors, Sharp-Sharp Points, Straight
5-SC-16	5 1/2" Operating Scissors, Sharp-Blunt Points, Straight, SuperCut Blade
5-SC-46	5 1/2" Operating Scissors, Sharp-Blunt Points, Curved, SuperCut Blade



MILTEX STITCH SCISSORS

9-108 4 1/2" Stitch Scissors, Delicate, Angular Blade



MILTEX UNIVERSAL SCISSORS

Different colors with black handles for easy identification. Designed to cut through even the toughest material. Stainless steel with one serrated blade. 7 1/2" length. The heavy plastic handles are contoured to provide maximum comfort. Autoclavable to 143°C/290°F. Will cut bandage and gauze, cast material, metal sheets and tough clothing material. Angled design allow for placement of the lower blade under the bandage for safe and easy removal. Available in high quality grade German stainless steel or in Vantage economy grade.

5-1025	High Quality Red Universal Scissors
5-1026	High Quality Orange Universal Scissors
5-1027	High Quality Blue Universal Scissors
5-1028	High Quality Purple Universal Scissors
5-1029	High Quality Green Universal Scissors
5-1030	High Quality Yellow Universal Scissors
V95-1000	Economy Black Universal Scissors
V95-1025	Economy Red Universal Scissors
V95-1026	Economy Orange Universal Scissors
V95-1027	Economy Blue Universal Scissors
V95-1028	Economy Purple Universal Scissors
V95-1029	Economy Green Universal Scissors
V95-1030	Economy Yellow Universal Scissors

MILTEX WESTCOTT TENOTOMY SCISSORS

Wide handles.

18-1487	5 1/4" Tenotomy Scissors, Right
18-1488	5 1/4" Tenotomy Scissors, Left



MILTEX WESTCOTT STITCH SCISSORS

Curved with sharp points.

18-1486	4 1/2" Stitch Scissors
---------	------------------------



8100



8100BL

TECH-MED ALL PURPOSE SHEARS

Easily cuts through the toughest materials. 7 1/2" with stainless steel blade. Autoclavable.

8100	All Purpose Shears, Black
8100BL	All Purpose Shears, Blue
8100G	All Purpose Shears, Green
8100PR	All Purpose Shears, Purple
8100R	All Purpose Shears, Red
8100W	All Purpose Shears, White



22-500TC




22-502TC


MILTEX WIRE TWISTER

Carb-N-Sert with tungsten carbide inserts.

22-500TC	6" Wire Twister
22-502TC	8" Wire Twister



9-120



9-122



9-124




9D-121



9D-123



9D-132



9D-133



9D-134



9D-135

MILTEX WIRE CUTTING SCISSORS WITH NOTCH

9-120	4 3/4" Wire Cutting Scissors w/one serrated blade, Straight
9-122	4 3/4" Wire Cutting Scissors w/one serrated blade, Curved
9-124	4 3/4" Wire Cutting Scissors w/one serrated blade, Angled to Side
9-124TC	4 3/4" Wire Cutting Scissors w/one serrated blade, Carb-N-Sert, Angled To Side
9D-121	4 3/4" Wire Cutting Scissors, Straight
9D-123	4 3/4" Wire Cutting Scissors, Curved
9D-132	4 1/4" Wire Cutting Scissors w/one serrated blade, Straight
9D-133	4 1/4" Wire Cutting Scissors, Straight
9D-134	4 1/4" Wire Cutting Scissors w/one serrated blade, Curved
9D-135	4 1/4" Wire Cutting Scissors, Curved



BUSSE GENERAL PURPOSE NURSES' SCISSORS

Sterile, 5 1/4" sharp/blunt scissors. Cuts gauze or tape easily. Each scissors is individually sterile tray packaged. Latex free.

7246	Scissors, 100/cs
------	------------------



MILTEX ESMARCH HEAVY DUTY BANDAGE & PLASTER SHEARS

5-568 8" Shears



MILTEX HERCULES HEAVY DUTY BANDAGE & PLASTER SHEARS

Serrated blade.

5-547 7 1/2" Bandage & Plaster Shears



MILTEX BANDAGE & UTILITY SCISSORS

Non-stick fluoride coating over stainless steel blades. Neon-green plastic finger rings for easy identification. Cuts gauze, tubing, clothing, leather, light casts and wet fiberglass casting material.

5-800 8 1/2" Utility Scissors



MILTEX BLEPHAROPLASTY SCISSORS

Curved flat tips with one serrated blade. Carb-N-Sert with tungsten carbide insert blades.

21-537TC 4 1/2" Blepharoplasty Scissors



18-1392




18-1394

MILTEX BONN MINIATURE IRIS SCISSORS

15mm blades.

18-1392	3 1/2" Iris Scissors, Straight
18-1394	3 1/2" Iris Scissors, Curved



MILTEX COTTLE SCISSORS


21-610	6 1/2" Cottle Scissors, Dorsal, Heavy Pattern w/rounded blades & angular shanks
21-612	4 1/2" Cottle Scissors, Bulldog Nasal, Slight Curved Blunt Blades



MILTEX FIRST AID & UTILITY SHEARS


With wire cutter and serrated blade.

5-600 8 1/2" First Aid & Utility Shears



MILTEX FOMON SCISSORS

21-600 5" Fomon Scissors, Saber-Back, Slightly Curved Blades w/semi-sharp outer edges used for dissection
21-602 5" Fomon Scissors, Lower Lateral, Strong Curve, Dissecting Blades
21-604 5 1/2" Fomon Scissors, Dorsal, Light Pattern, Angular Shanks
21-604TC 5 1/2" Fomon Scissors, Angled Shanks, Carab-N-Sert



MILTEX GOLDMAN-FOX WOUND DEBRIDEMENT SCISSORS

5-320 5" Wound Debridement Scissors, Curved Fine Tips w/one serrated blade
5-320TC 5" Wound Debridement Scissors, Curved Fine Tips w/one serrated blade, Carb-N-Sert
5-SC-320 5" Wound Debridement Scissors, Curved, SuperCut Blade
5D-322 5" Wound Debridement Scissors, Curved w/one serrated blade, Angled Shanks

MILTEX GORNEY FACELIFT SCISSORS

One serrated blade.

21-607 7 3/4" Gorney Facelift Scissors, Straight
21-608 7 3/4" Gorney Facelift Scissors, Slight Curved

MILTEX GORNEY SEPTAL SCISSORS

Double action. Angled handle. Serrated blades with 1 3/4" cut.

21-621 7 3/4" Gorney Septal Scissors

MILTEX KAHN DISSECTING SCISSORS

5-265 5 1/2" Dissecting Scissors, Straight
5-266 5 1/2" Dissecting Scissors, Curved
5-267 7" Dissecting Scissors, Straight
5-268 7" Dissecting Scissors, Curved



MILTEX KNAPP IRIS SCISSORS

18-1420 4" Iris Scissors, Straight, Sharp-Blunt Points
18-1422 4" Iris Scissors, Straight, Blunt-Blunt Points
18-1426 4" Iris Scissors, Curved, Sharp-Blunt Points
18-1428 4" Iris Scissors, Curved, Blunt-Blunt Points

MILTEX SERRATEX OPERATING SCISSORS

Standard pattern. One blade with fine serrations reduces tissue slippage.

5-70 5 1/2" Operating Scissors, Sharp-Sharp Points, Straight
5-72 5 1/2" Operating Scissors, Sharp-Blunt Points, Straight
5-73 5 1/2" Operating Scissors, Blunt-Blunt Points, Straight
5-74 5 1/2" Operating Scissors, Sharp-Sharp Points, Curved
5-76 5 1/2" Operating Scissors, Sharp-Blunt Points, Curved
5-77 5 1/2" Operating Scissors, Blunt-Blunt Points, Curved

MILTEX SERRATEX-MAYO DISSECTING SCISSORS

One fine serrated blade.

5-132 6 3/4" Dissecting Scissors, Straight
5-134 6 3/4" Dissecting Scissors, Curved



MILTEX U.S.A. GAUZE SCISSORS

Chrome.

5-564 8" Gauze Scissors

OMRON BANDAGE SCISSORS

Stainless steel. Lister.

0701 Bandage Scissors w/o clip, 5 1/2"

VANTAGE ESMARCH SHEARS

V95-568 8" Plaster Shears

VANTAGE EYE SUTURE SCISSORS

V918-1652 3 3/4" Curved Eye Suture Scissors



VANTAGE IRIS SCISSORS

- V95-304 4 1/8" Straight Iris Scissors
V95-306 4 1/8" Curved Iris Scissors




VANTAGE KNOWLES BANDAGE SCISSORS

- V95-560 5 1/2" Bandage Scissors



V95-502SS



V95-532



V95-550

VANTAGE LISTER BANDAGE SCISSORS

- V95-500SS 3 1/2" Bandage Scissors
V95-502SS 4 1/2" Bandage Scissors
V95-504SS 5 1/2" Bandage Scissors
V95-506SS 7 1/4" Bandage Scissors
V95-524 5 1/2" Gold Plated Bandage Scissors
V95-532 5 1/2" Bandage Scissors w/pocket clip
V95-550 8" Bandage Scissors w/one large finger ring

VANTAGE LITTAUER STITCH SCISSORS

- V99-102 4 1/2" Stitch Scissors
V99-104 5 1/2" Stitch Scissors



V95-120

V95-126

VANTAGE MAYO DISSECTING SCISSORS

- V95-120 5 1/2" Straight Dissecting Scissors
V95-122 5 1/2" Curved Dissecting Scissors
V95-124 6 3/4" Straight Dissecting Scissors
V95-126 6 3/4" Curved Dissecting Scissors



VANTAGE METZENBAUM SCISSORS

- V95-180 5 1/2" Metzenbaum Scissors, Curved
V95-182 7" Metzenbaum Scissors, Curved



VANTAGE OPERATING SCISSORS

- V95-12 4 1/2" Operating Scissors, Straight, Sharp/Blunt Points
V95-16 5 1/2" Operating Scissors, Straight, Sharp/Blunt Points
V95-18 6 1/2" Operating Scissors, Straight, Sharp/Blunt Points
V95-19 7 1/2" Operating Scissors, Straight, Sharp/Blunt Points
V95-2 4 1/2" Operating Scissors, Straight, Sharp/Sharp Points
V95-22 4 1/2" Operating Scissors, Straight, Blunt/Blunt Points
V95-26 5 1/2" Operating Scissors, Straight, Blunt/Blunt Points
V95-28 6 1/2" Operating Scissors, Straight, Blunt/Blunt Points
V95-29 7 1/2" Operating Scissors, Straight, Blunt/Blunt Points
V95-32 4 1/2" Operating Scissors, Curved, Sharp/Sharp Points
V95-36 5 1/2" Operating Scissors, Curved, Sharp/Sharp Points
V95-38 6 1/2" Operating Scissors, Curved, Sharp/Sharp Points
V95-42 4 1/2" Operating Scissors, Curved, Sharp/Blunt Points
V95-46 5 1/2" Operating Scissors, Curved, Sharp/Blunt Points
V95-48 6 1/2" Operating Scissors, Curved, Sharp/Blunt Points
V95-52 4 1/2" Operating Scissors, Curved, Blunt/Blunt Points
V95-56 5 1/2" Operating Scissors, Curved, Blunt/Blunt Points
V95-58 6 1/2" Operating Scissors, Curved, Blunt/Blunt Points
V95-6 5 1/2" Operating Scissors, Straight, Sharp/Sharp Points
V95-8 6 1/2" Operating Scissors, Straight, Sharp/Sharp Points
V95-9 7 1/2" Operating Scissors, Straight, Sharp/Sharp Points

VANTAGE SPENCER STITCH SCISSORS

- V99-100 3 1/2" Stitch Scissors, Delicate

VANTAGE STEVENS TENOTOMY SCISSORS

- V918-1462 4 1/8" Straight Stevens Tenotomy Scissors, Blunt
V918-1466 4 1/8" Curved Stevens Tenotomy Scissors, Blunt



V95-312

V95-314

VANTAGE STRABISMUS SCISSORS

- V95-312 4" Strabismus Scissors, Straight
V95-314 4" Strabismus Scissors, Curved

Needle Holders



PRO ADVANTAGE® NEEDLE HOLDERS


Expertly crafted from surgical grade stainless steel and meets the very highest quality standards and specifications. This outstanding quality and price point makes Pro Advantage instruments ideal for today's growing instrument market.

N407300	Webster Needle Holder, 5"
N407302	Webster Needle Holder, 5", Tungsten Carbide
N407305	Halsey Needle Holder, 5"
N407307	Halsey Needle Holder, 5", Tungsten Carbide
N407310	Baumgartner Needle Holder, 5 1/2"
N407312	Baumgartner Needle Holder, 5 1/2", Tungsten Carbide
N407315	Mayo-Hegar Needle Holder, 5"
N407317	Mayo-Hegar Needle Holder, 5", Tungsten Carbide
N407320	Mayo-Hegar Needle Holder, 6"
N407322	Mayo-Hegar Needle Holder, 6", Tungsten Carbide




MILTEX BAUMGARTNER NEEDLE HOLDERS

8-40	5 1/2" Baumgartner Needle Holder, Regular Jaws
8-40CG	5 1/2" Baumgartner Needle Holder, Ceram-A-Grip
8-40TC	5 1/2" Baumgartner Needle Holder, Carb-N-Sert T.C. Jaws




MILTEX CASTROVIEJO NEEDLE HOLDER

18-1816TC	4 1/8" Castroviejo Needle Holder, Carb-N-Sert, Round Handles w/lock, Delicate, Smooth, Straight
18-1820	5 1/2" Castroviejo Needle Holder, Straight
18-1824	5 1/2" Castroviejo Needle Holder, Curved
18-1828	5 1/2" Castroviejo Needle Holder w/lock, Straight
18-1828A-TC	5 1/2" Castroviejo Needle Holder w/lock, Carb-N-Sert, Serrated, Straight
18-1828TC	5 1/2" Castroviejo Needle Holder w/lock, Carb-N-Sert, Smooth, Straight
18-1831	5 1/2" Castroviejo Needle Holder, Extra Delicate w/lock, Straight
18-1832	5 1/2" Castroviejo Needle Holder w/lock, Curved
18-1832A-TC	5 1/2" Castroviejo Needle Holder w/lock, Carb-N-Sert, Serrated, Curved
18-1832TC	5 1/2" Castroviejo Needle Holder w/lock, Carb-N-Sert, Smooth, Curved
18-1833	5 1/2" Castroviejo Needle Holder, Extra Delicate w/lock, Curved




MILTEX CRILE-WOOD NEEDLE HOLDERS

8-50	6" Crile-Wood Needle Holder, Serrated Jaws
8-50A-TC	6" Crile-Wood Needle Holder, Smooth Jaws, Carb-N-Sert T.C. Jaws
8-50CG	6" Crile-Wood Needle Holder, Ceram-A-Grip
8-50TC	6" Crile-Wood Needle Holder, Serrated Jaws, Carb-N-Sert T.C. Jaws
8-51	7" Crile-Wood Needle Holder, Serrated Jaws
8-51CG	7" Crile-Wood Needle Holder, Ceram-A-Grip
8-51TC	7" Crile-Wood Needle Holder, Serrated Jaws, Carb-N-Sert T.C. Jaws
8-52	8" Crile-Wood Needle Holder, Serrated Jaws
8-52CG	8" Crile-Wood Needle Holder, Ceram-A-Grip
8-52TC	8" Crile-Wood Needle Holder, Serrated Jaws, Carb-N-Sert T.C. Jaws, Carb-N-Sert T.C. Jaws
8-53	9" Crile-Wood Needle Holder, Serrated Jaws
8-53TC	9" Crile-Wood Needle Holder, Serrated Jaws, Carb-N-Sert T.C. Jaws
8-55TC	10 1/2" Crile-Wood Needle Holder, Serrated, Carb-N-Sert T.C. Jaws
8-57TC	12" Crile-Wood Needle Holder, Serrated, Carb-N-Sert T.C. Jaws
8-750TC	6" Crile-Wood Needle Holder w/quick release ratchet, Serrated, Carb-N-Sert T.C. Jaws
8-751TC	7" Crile-Wood Needle Holder w/quick release ratchet, Serrated, Carb-N-Sert T.C. Jaws
8-752TC	8" Crile-Wood Needle Holder w/quick release ratchet, Serrated, Carb-N-Sert T.C. Jaws
8-753TC	9" Crile-Wood Needle Holder w/quick release ratchet, Serrated, Carb-N-Sert T.C. Jaws



MILTEX DERF NEEDLE HOLDERS

8-4	4 3/4" Derf Needle Holder, Regular Jaws
8-4CG	4 3/4", Ceram-A-Grip
8-4TC	4 3/4" Derf Needle Holder, Carb-N-Sert T.C. Jaws




8-8



8-8A-TC



8-8CG



8-8TC



8D-36

MILTEX HALSEY NEEDLE HOLDERS


8-8	5" Halsey Needle Holder, Smooth, Regular Jaws
8-8A-TC	5" Halsey Needle Holder, Serrated, Carb-N-Sert T.C. Jaws
8-8CG	5" Halsey Needle Holder, Ceram-A-Grip, Smooth
8-8TC	5" Halsey Needle Holder, Smooth, Carb-N-Sert T.C. Jaws
8D-36	5" Halsey Needle Holder, Serrated Jaws




MILTEX HEANEY NEEDLE HOLDER

Curved jaws.


8-100	8 1/2" Heaney Needle Holder
8-100TC	8 1/2" Heaney Needle Holder, Carb-N-Sert, Serrated




8-130



8-130TC



8-42CG




8-42TC




8-44



8-44CG




8-44TC




8-46




8-46CG




8-46TC




8-48CG




8-48TC



8-744TC



8-746TC




8-748TC

MILTEX MAYO-HEGAR NEEDLE HOLDERS


8-130	10 1/2" Mayo-Hegar Needle Holder
8-130TC	10 1/2" Mayo-Hegar Needle Holder, Carb-N-Sert T.C. Jaws
8-42	5" Mayo-Hegar Needle Holder, Serrated Jaws
8-42CG	5" Mayo-Hegar Needle Holder, Ceram-A-Grip
8-42TC	5" Mayo-Hegar Needle Holder, Serrated Jaws, Carb-N-Sert T.C. Jaws
8-44	6" Mayo-Hegar Needle Holder, Serrated Jaws
8-44CG	6" Mayo-Hegar Needle Holder, Ceram-A-Grip
8-44TC	6" Mayo-Hegar Needle Holder, Serrated Jaws, Carb-N-Sert T.C. Jaws
8-46	7" Mayo-Hegar Needle Holder, Serrated Jaws
8-46CG	7" Mayo-Hegar Needle Holder, Ceram-A-Grip
8-46TC	7" Mayo-Hegar Needle Holder, Serrated Jaws, Carb-N-Sert T.C. Jaws
8-48	8" Mayo-Hegar Needle Holder, Serrated Jaws
8-48CG	8" Mayo-Hegar Needle Holder, Ceram-A-Grip
8-48TC	8" Mayo-Hegar Needle Holder, Serrated Jaws, Carb-N-Sert T.C. Jaws
8-49TC	12" Mayo-Hegar Needle Holder, Carb-N-Sert T.C. Jaws
8-744TC	6" Mayo-Hegar Needle Holder w/quick release ratchet, Serrated, Carb-N-Sert T.C. Jaws
8-746TC	7" Mayo-Hegar Needle Holder w/quick release ratchet, Serrated, Carb-N-Sert T.C. Jaws
8-748TC	8" Mayo-Hegar Needle Holder w/quick release ratchet, Serrated, Carb-N-Sert T.C. Jaws

**Can't find what
you're looking for?
Call your
representative.**


Needle Holders




8-14




8-14TC




8-15



8-15CG




8-15TC



8-16CG




8-16TC




8-17TC

MILTEX OLSEN-HEGAR NEEDLE HOLDER-SCISSORS


- 8-14 4 3/4" Olsen-Hegar Needle Holder, Serrated, Extra Delicate, Regular Jaws
- 8-14A 4 3/4" Olsen-Hegar Needle Holder, Smooth, Extra Delicate, Regular Jaws
- 8-14A-TC 4 3/4" Olsen-Hegar Needle Holder, Smooth, Extra Delicate, Carb-N-Sert T.C. Jaws
- 8-14TC 4 3/4" Olsen-Hegar Needle Holder, Serrated, Extra Delicate, Carb-N-Sert T.C. Jaws
- 8-15 5 1/2" Olsen-Hegar Needle Holder, Serrated, Regular Jaws
- 8-15CG 5 1/2" Olsen-Hegar Needle Holder, Serrated, Ceram-A-Grip
- 8-15TC 5 1/2" Olsen-Hegar Needle Holder, Serrated, Carb-N-Sert T.C. Jaws
- 8-16 6 1/2" Olsen-Hegar Needle Holder, Serrated, Regular Jaws
- 8-16CG 6 1/2" Olsen-Hegar Needle Holder, Serrated, Ceram-A-Grip
- 8-16TC 6 1/2" Olsen-Hegar Needle Holder, Serrated, Carb-N-Sert T.C. Jaws
- 8-17 7 1/4" Olsen-Hegar Needle Holder, Serrated, Regular Jaws
- 8-17TC 7 1/2" Olsen-Hegar Needle Holder, Serrated, Carb-N-Sert T.C. Jaws




8-6




8-6CG




8-6TC




8-7




8-706TC



8-7A-TC



8-7CG



8-7TC

MILTEX WEBSTER NEEDLE HOLDERS

- 8-6 5" Webster Needle Holder, Smooth
- 8-6CG 5" Webster Needle Holder, Ceram-A-Grip
- 8-6TC 5" Webster Needle Holder, Smooth, Carb-N-Sert
- 8-7 4 3/4" Webster Needle Holder, Smooth, Very Delicate
- 8-706TC 5" Webster Needle Holder w/quick-release ratchet, Carb-N-Sert, Smooth
- 8-7A-TC 4 3/4" Webster Needle Holder, Carb-N-Sert, Extra Delicate, Serrated
- 8-7CG 4 3/4" Webster Needle Holder, Ceram-A-Grip
- 8-7TC 4 3/4" Webster Needle Holder, Smooth, Very Delicate, Carb-N-Sert



MILTEX BOYNTON NEEDLE HOLDER

Delicate.

- 8-26 5" Boynton Needle Holder



MILTEX BROWN NEEDLE HOLDER

Convex jaws with fine serrations.

- 8-22 5 1/4" Brown Needle Holder
- 8-23 6 3/4" Brown Needle Holder



MILTEX COLLIER NEEDLE HOLDER

Fenestrated jaws.

- 8-2 5" Collier Needle Holder

MILTEX CONVERSE NEEDLE HOLDER

Longitudinal serrations and large finger rings.

- 8-20 4 1/4" Converse Needle Holder



MILTEX COOLEY NEEDLE HOLDER

Carb-N-Sert with tungsten carbide inserts. Extra delicate. Serrated jaws 4000 teeth per square inch.

- 8-62TC 6" Cooley Needle Holder
- 8-63TC 6 3/4" Cooley Needle Holder



17-1090TC



17-1091TC



17-1092TC

MILTEX JACOBSON MICROVASCULAR NEEDLE HOLDER

Carb-N-Sert with tungsten carbide inserts. With lock, serrated jaws 16000 teeth per square inch.

- 17-1090TC 7" Jacobson Microvascular Needle Holder, Straight
- 17-1091TC 7" Jacobson Microvascular Needle Holder, Curved
- 17-1092TC 8 1/2" Jacobson Microvascular Needle Holder, Straight



VANTAGE NEEDLE HOLDERS

V918-1828	Castroviejo Needle Holder, 5 1/4" Straight w/catch
V98-15	Olsen-Hegar Needle Holder, 5 1/2"
V98-16	Olsen-Hegar Needle Holder, 6 1/2"
V98-17	Olsen-Hegar Needle Holder, 7 1/2"
V98-2	Collier Needle Holder, 5"
V98-28	Mathieu Needle Holder, 5 1/2"
V98-4	Derf Needle Holder, 4 3/4"
V98-40	Baumgartner Needle Holder, 5 1/2"
V98-42	Mayo-Hegar Needle Holder, 5"
V98-44	Mayo-Hegar Needle Holder, 6"
V98-46	Mayo-Hegar Needle Holder, 7"
V98-48	Mayo-Hegar Needle Holder, 8"
V98-50	Crile-Wood Needle Holder, 6", Delicate
V98-6	Webster Needle Holder, 5" Smooth Jaw
V98-8	Halsey Needle Holder, 5" Smooth Jaw



PRO ADVANTAGE® BROWN-ADSON FORCEPS

Manufactured to German standards with exceptional value. Lifetime guarantee against defects in material or workmanship assures years of serviceability from each instrument. Selection of the most popular patterns of surgical instruments.

N407130	4 3/4" Brown-Adson Forceps, 7x7 Teeth
---------	---------------------------------------



N407100



N407105

PRO ADVANTAGE® DRESSING FORCEPS

Manufactured to German standards with exceptional value. Lifetime guarantee against defects in material or workmanship assures years of serviceability from each instrument. Selection of the most popular patterns of surgical instruments.

N407100	5" Dressing Forceps
N407105	5 1/2" Dressing Forceps



N407120



N407125

PRO ADVANTAGE® INSTRUMENTS ADSON DRESSING FORCEPS

Manufactured to German standards with exceptional value. Lifetime guarantee against defects in material or workmanship assures years of serviceability from each instrument. Selection of the most popular patterns of surgical instruments.


N407120	4 3/4" Adson Dressing Forceps
N407125	4 3/4" Adson Tissue Forceps, 1x2 Teeth



PRO ADVANTAGE® FOERSTER FORCEPS

Manufactured to German standards with exceptional value. Lifetime guarantee against defects in material or workmanship assures years of serviceability from each instrument. Selection of the most popular patterns of surgical instruments.

N407225	9 1/2" Straight Serrated Foerster Forceps
---------	---



N407200



N407205

PRO ADVANTAGE® HALSTED MOSQUITO FORCEPS

Manufactured to German standards with exceptional value. Lifetime guarantee against defects in material or workmanship assures years of serviceability from each instrument. Selection of the most popular patterns of surgical instruments.

N407200	5" Straight Halsted Mosquito Forceps
N407205	5" Curved Halsted Mosquito Forceps



N407210



N407215

PRO ADVANTAGE® KELLY FORCEPS

Manufactured to German standards with exceptional value. Lifetime guarantee against defects in material or workmanship assures years of serviceability from each instrument. Selection of the most popular patterns of surgical instruments.

N407210	5 1/2" Straight Kelly Forceps
N407215	5 1/2" Curved Kelly Forceps



PRO ADVANTAGE® KEVORKIAN-YOUNGE BIOPSY FORCEPS

Manufactured to German standards with exceptional value. Lifetime guarantee against defects in material or workmanship assures years of serviceability from each instrument. Selection of the most popular patterns of surgical instruments.

N407430	Kevorkian-Young Biopsy Forceps
---------	--------------------------------



N407135




N407140

N407145

PRO ADVANTAGE® SPLINTER FORCEPS

Manufactured to German standards with exceptional value. Lifetime guarantee against defects in material or workmanship assures years of serviceability from each instrument. Selection of the most popular patterns of surgical instruments.

N407135	3 1/2" Splinter Forceps
N407140	4" Splinter Forceps
N407145	4 1/2" Splinter Forceps



N407110

N407115

PRO ADVANTAGE® TISSUE FORCEPS

Manufactured to German standards with exceptional value. Lifetime guarantee against defects in material or workmanship assures years of serviceability from each instrument. Selection of the most popular patterns of surgical instruments.

- N407110 5" Tissue Forceps 1x2 Teeth
- N407115 5 1/2" Tissue Forceps 1x2 Teeth



ADC SPECIALTY FORCEPS

Forged from corrosion resistant surgical stainless steel. These popular floor grade instruments are built to last.

- 315 Magill Catheter Forceps, Child, 8"
- 316 Magill Catheter Forceps, Adult, 9 3/4"
- 317 Young Tongue Seizing Forceps, Child, 6 1/2"



2-70



2-72

MILTEX MAGILL ENDOTRACHEAL CATHETER INTRODUCING FORCEPS

- 2-70 9" Catheter Introducing Forceps, Adult
- 2-72 7" Catheter Introducing Forceps, Child, Extra Delicate



19-434



20-140

MILTEX HARTMAN & HARTMAN-NOYES ALLIGATOR FORCEPS

- 19-434 Hartman Alligator Forceps, 3 1/4" Shaft, Serrated Jaws
- 20-140 Hartman-Noyes Alligator Forceps, 5 1/2" Shaft, Serrated Jaws
- 20-142 6 1/2" Shaft Alligator Forceps, Serrated Jaws
- 20-144 5 1/2" Shaft Alligator Forceps, 1 x 2 Teeth



20-124

MILTEX HARTMAN NASAL DRESSING FORCEPS

Serrated tips.

- 20-124 7" Nasal Dressing Forceps, Improved Delicate Pattern
- 20-126 7 1/4" Nasal Dressing Forceps, Regular Model

MILTEX ADSON BAYONET DRESSING FORCEPS

- 6-198 8 1/4" Dressing Forceps, Serrated
- 6-198TC 8 1/2" Dressing Forceps, Carb-N-Sert

MILTEX BERGH CILIA FORCEPS

- 18-1102 3 1/2" Cilia Forceps
- 18-1110 3 1/2" Cilia Forceps, 5mm Wide Jaws



MILTEX SWISS CILIA & SUTURE

3 3/4" forceps. Excellent for suture removal.

- 18-1107 Swiss Cilia & Suture



17-2510



6-118



6-118TC



6-118XL

MILTEX ADSON DRESSING FORCEPS

Serrated.

- 17-2510 4 3/4" Dressing Forceps, Micro, Delicate, 0.5mm Wide Tips
- 6-118 4 3/4" Dressing Forceps, Delicate
- 6-118A 4 3/4" Dressing Forceps, Delicate, Angled
- 6-118TC 4 3/4" Dressing Forceps, Carb-N-Sert
- 6-118XL 4 3/4" Dressing Forceps, Fenestrated Handle



6-308



6-310

MILTEX CARMALT SPLINTER FORCEPS

- 6-308 4 1/4" Splinter Forceps, Straight
- 6-310 4 1/4" Splinter Forceps, Curved



MILTEX CASTROVIEJO SUTURING FORCEPS

With tying platform.

- 18-950 4" Suturing Forceps, 1 x 2 Teeth 0.9mm
- 18-951 4" Suturing Forceps, 1 x 2 Teeth 0.3mm
- 18-952 4" Suturing Forceps, 1 x 2 Teeth 0.6mm
- 18-953 4" Suturing Forceps w/11mm wide handles, 1 x 2 Teeth 0.12mm
- 18-954 4" Suturing Forceps w/11mm wide handles, 1 x 2 Teeth 0.6mm
- 18-955 4" Suturing Forceps w/11mm wide handles, 1 x 2 Teeth 0.9mm



MILTEX CUSHING DRESSING FORCEPS


Serrated with scraper end.

- 6-166 7" Dressing Forceps w/Gutch Handle
- 6-166TC 7" Dressing Forceps, Carb-N-Sert
- 6-188 7" Dressing Forceps
- 6-190 7 1/4" Dressing Forceps, Bayonet

MILTEX DEPILATORY FORCEPS

Flat square tips with deep serrations.

- 33-500 5 1/2" Depilatory Forceps



MILTEX DRESSING FORCEPS

Serrated tips.

- 6-10 6" Dressing Forceps, Standard w/serrated handles
- 6-100 5" Dressing Forceps, Curved Narrow Tips
- 6-12 7" Dressing Forceps, Standard w/serrated handles
- 6-14 8" Dressing Forceps, Standard w/serrated handles
- 6-16 10" Dressing Forceps, Standard w/serrated handles
- 6-18 12" Dressing Forceps, Standard w/serrated handles
- 6-2 4" Dressing Forceps, Standard w/serrated handles
- 6-24 4 1/2" Dressing Forceps, Delicate w/fluted handles
- 6-26 5" Dressing Forceps, Delicate w/fluted handles
- 6-28 5 1/2" Dressing Forceps, Delicate w/fluted handles
- 6-28XL 5 1/2" Dressing Forceps, Delicate Lightweight w/fenestrated handles
- 6-30 6" Dressing Forceps, Delicate w/fluted handles
- 6-30TC 6" Dressing Forceps, Delicate w/fluted handles, Carb-N-Sert
- 6-4 4 1/2" Dressing Forceps, Standard w/serrated handles
- 6-6 5" Dressing Forceps, Standard w/serrated handles
- 6-8 5 1/2" Dressing Forceps, Standard w/serrated handles
- 6-8XL 5 1/2" Dressing Forceps, Lightweight w/fenestrated handles




MILTEX FEILCHENFELD SPLINTER FORCEPS

- 6-316 3" Splinter Forceps
- 6-318 4 1/2" Splinter Forceps




MILTEX GILL IRIS FORCEPS

- 18-825 3" Iris Forceps, Dressing, Fine Cross Serrated Tips, 0.6mm Wide
- 18-834 2 7/8" Iris Forceps, Tissue, Tips 0.6mm Wide, 1 x 2 Teeth



MILTEX HUDSON (EWALD) FORCEPS


- 6-112 4 3/4" Hudson Forceps, Dressing, Serrated
- 6-113 6" Hudson (EWALD) Forceps, Dressing, Serrated Tips
- 6-113TC 6" Hudson (EWALD) Forceps, Dressing, Cross Serrated Tips, Carb-N-Sert (Tungsten Carbide)
- 6-114 4 3/4" Hudson Forceps, Tissue, 1x2 Teeth
- 6-115 6" Hudson (EWALD) Forceps, Tissue, 1 x 2 Teeth
- 6-115TC 6" Hudson (EWALD) Forceps, Tissue & Suture, 1x2 Teeth, Cross Serrated Tying Platform, Carb-N-Sert (Tungsten Carbide)



MILTEX PLAIN SPLINTER FORCEPS

Fine points.

- 6-300 3 1/2" Splinter Forceps
- 6-302 4" Splinter Forceps
- 6-304 4 1/2" Splinter Forceps



MILTEX IRIS TISSUE FORCEPS

1 x 2 teeth.

- 18-785 4" Tissue Forceps, Straight, Extra Delicate, 0.5mm Wide Tips
- 18-786 4" Tissue Forceps, Straight, Standard, 0.8mm Wide Tips
- 18-787 4" Tissue Forceps, Straight, Non Magnetic, Standard, 0.8mm Wide Tips
- 18-788 4" Tissue Forceps, Half Curved, Standard, 0.8mm Wide Tips
- 18-789 4" Tissue Forceps, Half Curved, Extra Delicate, 0.5mm Wide Tips
- 18-790 4" Tissue Forceps, Full Curved, Standard, 0.8mm Wide Tips
- 18-791 4" Tissue Forceps, Full Curved, Extra Delicate, 0.5mm Wide Tips

MILTEX POTTS-SMITH DRESSING FORCEPS

Serrated.

- 6-154 7" Dressing Forceps
- 6-154TC 7" Dressing Forceps, Carb-N-Sert
- 6-156 8 1/4" Dressing Forceps
- 6-156TC 8 1/4" Dressing Forceps, Carb-N-Sert
- 6-158 9 1/2" Dressing Forceps
- 6-158TC 9 3/4" Dressing Forceps, Carb-N-Sert
- 6-159 12" Dressing Forceps



6-104



6-106


MILTEX SEMKIN FORCEPS

- 6-104 5" Semkin Forceps, Dressing, Serrated Tips
- 6-106 5" Semkin Forceps, Tissue, 1 x 2 Teeth
- 6-107 5" Semkin Forceps, Tissue, Curved, 1 x 2 Teeth
- 6-108 5" Semkin Forceps, Tissue, 2 x 3 Teeth
- 6-109 6" Semkin Forceps, Dressing, Serrated Tips
- 6-109A 6" Semkin Forceps, Tissue, 1 x 2 Teeth

MILTEX STEVENS IRIS TISSUE FORCEPS

Straight. 1 x 2 teeth. 1mm wide.

- 18-804 4" Iris Tissue Forceps



17-303



17-307

MILTEX SWISS JEWELER STYLE FORCEPS

Non-magnetic stainless steel.

- 17-301 Style 1 Jeweler Style Forceps, Fine, 4 3/4"
- 17-302 Style 2 Jeweler Style Forceps, Strong, Fine, 4 3/4"
- 17-303 Style 3 Jeweler Style Forceps, Narrow, Fine, 4 3/4"
- 17-303C Style 3C Jeweler Style Forceps, Narrow, Fine, 4 3/8"
- 17-303X Style 3 Jeweler Style Forceps, Micro-Fine, 4 3/4"
- 17-304 Style 4 Jeweler Style Forceps, Fine, 4 3/8"
- 17-305 Style 5 Jeweler Style Forceps, Super Fine, 4 3/8"
- 17-305X Style 5 Jeweler Style Forceps, Micro-Fine, 4 3/8"
- 17-306 Style 6 Jeweler Style Forceps, Pickup, 4 1/2"
- 17-307 Style 7 Jeweler Style Forceps, Fine Curved, 4 1/2"
- 17-307X Style 7 Jeweler Style Forceps, Micro-Fine, Curved, 4 1/2"



6-336



6-337

MILTEX VIRTUS SPLINTER FORCEPS

- 6-336 6" Splinter Forceps, Straight
- 6-337 6" Splinter Forceps, 45° Angle



6-312



6-314


MILTEX WALTHER SPLINTER FORCEPS

- 6-312 4 1/4" Splinter Forceps, Straight
- 6-314 4 1/4" Splinter Forceps, Curved



MILTEX DESMARRES CHALAZION FORCEPS

- 18-1200 3 1/2" Desmarres Chalazion Forceps, Small, Inside Ring
Measures 11 x 17mm
- 18-1202 3 1/2" Desmarres Chalazion Forceps, Medium, Inside Ring
Measures 12 x 23mm
- 18-1204 3 1/2" Desmarres Chalazion Forceps, Large, Inside Ring
Measures 17 x 28mm



ADC HEMOSTATIC FORCEPS

Forged from corrosion resistant surgical stainless steel. These popular floor grade instruments are built to last.

- 310 Kelly Forceps, Straight, 5 1/2"
- 3101 Crile Forceps, Straight, 5 1/2"
- 3102 Crile Forceps, Curved, 5 1/2"
- 311 Kelly Forceps, Curved, 5 1/2"
- 312 Kelly Forceps, Straight, 6 1/4"
- 313 Kelly Forceps, Curved, 6 1/4"
- 314 Halstead Mosquito Forceps, Straight, 5"
- 3141 Halstead Mosquito Forceps, Curved, 5"

MILTEX BABY OCHSNER FORCEPS

1 x 2 teeth. Extra delicate.

- 7-70 5 1/2" Baby Ochsner Forceps, Straight
- 7-72 5 1/2" Baby Ochsner Forceps, Curved

MILTEX BABY PEAN FORCEPS

Extra delicate.

- 7-110 5 1/2" Baby Pean Forceps, Straight
- 7-112 5 1/2" Baby Pean Forceps, Curved




7-24



7-26

MILTEX HARTMAN MOSQUITO FORCEPS


- 7-24 3 1/2" Hartman Mosquito Forceps, Straight
- 7-24A 3 1/2" Hartman Mosquito Forceps, Straight Jaws, Shanks Angled On Flat
- 7-26 3 1/2" Hartman Mosquito Forceps, Curved



MILTEX STERILIZER FORCEPS

Made of stainless steel. 3-prong.

- 3-10 8" Sterilizer Forceps
- 3-12 12" Sterilizer Forceps



17-2500

6-120A

6-120XL

6-121




6-122

6-123


6-123XL

MILTEX ADSON TISSUE FORCEPS


- 17-2500 4 3/4" Tissue Forceps, 1 x 2 Teeth, 0.5 mm Wide Tips, Delicate
- 6-120 4 3/4" Tissue Forceps, 1 x 2 Teeth, Delicate
- 6-120A 4 3/4" Tissue Forceps, 1 x 2 Teeth, Angled, Delicate
- 6-120XL 4 3/4" Tissue Forceps, 1 x 2 Teeth, Fenestrated Handles, Lightweight
- 6-121 4 3/4" Tissue Forceps, 1 x 2 Teeth, Cross Serrated Jaws
- 6-122 4 3/4" Tissue Forceps, 2 x 3 Teeth, Delicate
- 6-123 4 3/4" Tissue Forceps, 1 x 2 Teeth w/suturing platform
- 6-123TC 4 3/4" Tissue Forceps, 1 x 2 Teeth, Cross Serrated Typing Platform, Carb-N-Sert, Lightweight
- 6-123XL 4 3/4" Tissue Forceps, 1 x 2 Teeth, Fenestrated Handles w/platform, Lightweight
- 6-129 6" Adson Forceps, Dressing, Serrated
- 6-129TC 6" Adson Forceps, Dressing, Cross Serrated Tips, Carb-N-Sert (Tungsten Carbide)
- 6-130 6" Adson Forceps, Tissue, 1 x 2 Teeth
- 6-130TC 6" Adson Forceps, Tissue & Suture, 1x2 Teeth, Cross Serrated Typing Platform, Carb-N-Sert (Tungsten Carbide)




16-10




16-12




16-3



16-33



16-6



16-8

MILTEX ALLIS TISSUE FORCEPS

- 16-10 7 1/4" Allis Tissue Forceps, 5 x 6 Teeth
- 16-12 10" Allis Tissue Forceps, 5 x 6 Teeth
- 16-20 6" Allis Tissue Forceps, 4 x 5 Teeth, Angular Jaws
- 16-22 6" Allis Tissue Forceps, 5 x 6 Teeth, Angular Jaws
- 16-24 7 1/4" Allis Tissue Forceps, 5 x 6 Teeth, Angular Jaws
- 16-25 Allis-NC Tissue Forceps, 6" w/double row of non-traumatic teeth
- 16-26 Allis-NC Tissue Forceps, 7 1/4" w/double row of non-traumatic teeth
- 16-3 Baby Allis Tissue Forceps, 5", 4 x 5 Teeth, Extra Delicate Jaws 2.5mm Wide
- 16-30 Allis-Adair Tissue Forceps, 6 1/2", 9 x 10 Teeth
- 16-31 Boys-Allis Tissue Forceps, 6", 5 x 6 Teeth, Heavy Pattern
- 16-33 Judd-Allis Tissue Forceps, 6", 3 x 4 Teeth, Delicate
- 16-34 Judd-Allis Tissue Forceps, 7 1/2", 3 x 4 Teeth, Delicate
- 16-38 Thoms-Allis Tissue Forceps, 8", 6 x 7 Teeth
- 16-4 Baby Allis Tissue Forceps, 5 1/2", 4 x 5 Teeth, Delicate Jaws, 4mm Wide
- 16-6 6" Allis Tissue Forceps, 4 x 5 Teeth
- 16-8 6" Allis Tissue Forceps, 5 x 6 Teeth



MILTEX BONNEY TISSUE FORCEPS

- 6-148 7" Tissue Forceps, 1 x 2 Teeth, Serrated Tips
- 6-150 7" Tissue Forceps, 2 x 3 Teeth



6-124



6-124TC



6-124XL



6-126

MILTEX BROWN-ADSON TISSUE FORCEPS

7 x 7 side grasping teeth.

- 6-124 4 3/4" Tissue Forceps
- 6-124TC 4 3/4" Tissue Forceps, Cross Serrated Typing Platform, Carb-N-Sert
- 6-124XL 4 3/4" Tissue Forceps, Lightweight, Fenestrated Handles
- 6-126 4 3/4" Tissue Forceps, Angular

MILTEX MAYO RUSSIAN TISSUE FORCEPS

Fenestrated jaws and slotted grooved handles.

6-210 9" Mayo Russian Tissue Forceps



MILTEX POTTS-SMITH TISSUE FORCEPS

1 x 2 teeth and serrated tips.

6-160 7" Tissue Forceps
6-162 8 1/4" Tissue Forceps
6-164 9 1/2" Tissue Forceps
6-164A 12" Tissue Forceps




6-142



6-144

MILTEX RUSSIAN TISSUE FORCEPS

6-142 6" Tissue Forceps
6-144 8" Tissue Forceps
6-145 10" Tissue Forceps



MILTEX TISSUE FORCEPS

6-40 4 1/2" Tissue Forceps, 1 x 2 Teeth, Standard w/serrated handles
6-42 5" Tissue Forceps, 1 x 2 Teeth, Standard w/serrated handles
6-44 5 1/2" Tissue Forceps, 1 x 2 Teeth, Standard w/serrated handles
6-44XL 5 1/2" Tissue Forceps, 1 x 2 Teeth, Standard Lightweight w/fenestrated handles
6-46 6" Tissue Forceps, 1 x 2 Teeth, Standard w/serrated handles
6-47 7" Tissue Forceps, 1 x 2 Teeth, Standard w/serrated handles
6-48 8" Tissue Forceps, 1 x 2 Teeth, Standard w/serrated handles
6-50 10" Tissue Forceps, 1 x 2 Teeth, Standard w/serrated handles
6-52 12" Tissue Forceps, 1 x 2 Teeth, Standard w/serrated handles
6-60 4 1/2" Tissue Forceps, 2 x 3 Teeth, Standard w/serrated handles
6-62 5" Tissue Forceps, 2 x 3 Teeth, Standard w/serrated handles
6-64 5 1/2" Tissue Forceps, 2 x 3 Teeth, Standard w/serrated handles
6-66 6" Tissue Forceps, 2 x 3 Teeth, Standard w/serrated handles
6-72 5" Tissue Forceps, 3 x 4 Teeth, Standard w/serrated handles
6-74 5 1/2" Tissue Forceps, 3 x 4 Teeth, Standard w/serrated handles
6-76 6" Tissue Forceps, 3 x 4 Teeth, Standard w/serrated handles
6-82 5" Tissue Forceps, 4 x 5 Teeth, Standard w/serrated handles
6-84 5 1/2" Tissue Forceps, 4 x 5 Teeth, Standard w/serrated handles
6-86 6" Tissue Forceps, 4 x 5 Teeth, Standard w/serrated handles
6-90 4 1/2" Tissue Forceps, 1 x 2 Teeth, Delicate w/fluted handles
6-92 5" Tissue Forceps, 1 x 2 Teeth, Delicate w/fluted handles
6-94 5 1/2" Tissue Forceps, 1 x 2 Teeth, Delicate w/fluted handles
6-96 6" Tissue Forceps, 1 x 2 Teeth, Delicate w/fluted handles

ADC THUMB DRESSING & PLAIN SPLINTER FORCEPS

Forged from corrosion resistant surgical stainless steel. These popular floor grade instruments are built to last.

- 330 Thumb Dressing Forceps, 5"
- 3319 Adson Tissue Forceps, Serrated, 4 1/2"
- 332 Plain Splinter Forceps, 4 1/2"
- 3328 Plain Splinter Forceps, 3 1/2"



BUSSE POSI-GRIP™ FORCEPS

Plastic forceps are available sterile or non-sterile. 4". Suited for the removal of Penrose drains. Unique construction maintains tip alignment. Each is individually packaged. Latex free.

- 7190 Sterile Deluxe Plastic Posi-Grip™ Forceps, 4", 100/cs
- 7195 Non-Sterile Deluxe Plastic Posi-Grip™ Forceps, 4", 200/cs



DMI® MOSQUITO FORCEPS

- 539-5311-0000 5" Forceps, Straight, Stainless Steel
- 539-5312-0000 5 1/2" Forceps, Curved, Stainless Steel



DYNAREX PLASTIC FORCEPS, NON-STERILE

Offers superior performance. Tapered ends for precise pick-ups.

- 4917 5" Plastic Forceps, 3000/cs

MILTEX BARRAQUER CILIA & SUTURE FORCEPS

5mm smooth platform.

- 18-1111 4 1/2" Cilia & Suture Forceps



18-866



18-867



18-868

MILTEX BISHOP-HARMON FORCEPS

- 18-866 3 3/8" Bishop-Harmon Forceps, Dressing, Fine Cross Serrated Tips 0.5mm Wide
- 18-867 3 3/4" Bishop-Harmon Forceps, Micro Tissue, Delicate, 0.3mm, 1 x 2 Teeth
- 18-868 3 3/8" Bishop-Harmon Forceps, Standard, 0.7mm, 1 x 2 Teeth

MILTEX BONN MICRO TISSUE FORCEPS

1 x 2 teeth 0.12mm with tying platform.

- 18-835 2 3/4" Tissue Forceps
- 18-836 3 3/4" Tissue Forceps



MILTEX BROWN TISSUE FORCEPS

8 x 8 side grasping teeth.

- 6-138 6" Tissue Forceps
- 6-140 8" Tissue Forceps
- 6-141 10" Tissue Forceps



6-168



6-189



6-192

MILTEX CUSHING TISSUE FORCEPS

1 x 2 teeth and scraper end.

- 6-168 7" Tissue Forceps w/Gutch Handle
- 6-189 7" Tissue Forceps
- 6-192 7 1/4" Tissue Forceps, Bayonet

MILTEX DEBAKEY ATRAUMATIC TISSUE FORCEPS

- 24-550 6" Tissue Forceps, 2.5cm Wide Tips
- 24-552 7 3/4" Tissue Forceps, 2.5cm Wide Tips
- 24-554 9 1/2" Tissue Forceps, 2.5cm Wide Tips
- 24-556 12" Tissue Forceps, 2.5cm Wide Tips
- 24-570 6" Tissue Forceps, 1.5mm Wide Jaw Tips
- 24-570XL 6" Tissue Forceps, 1.5mm Wide Jaw Tips, Lightweight w/fenestrated handles
- 24-572 8" Tissue Forceps, 1.5mm Wide Jaw Tips
- 24-574 9 1/2" Tissue Forceps, 1.5mm Wide Jaw Tips



MILTEX DEBAKEY-ADSON TISSUE FORCEPS

Atraumatic tips 1.5mm wide.

- 6-125 4 3/4" Tissue Forceps

MILTEX DOUGLAS CILIA FORCEPS

Thin tips 2.5mm wide with fine horizontal serrations.

- 18-1106 3 1/2" Cilia Forceps

MILTEX FERRIS SMITH TISSUE FORCEPS

Serrated platform tips.

- 26-956 7" Tissue Forceps, 1 x 2 Teeth
- 26-958 7" Tissue Forceps, 2 x 3 Teeth



7-630



7-632

MILTEX FLETCHER SPONGE FORCEPS

- 7-630 9 1/2" Sponge Forceps, Straight
7-632 9 1/2" Sponge Forceps, Curved



18-795



18-796

MILTEX FOERSTER FORCEPS

Octagonal grip.

- 18-795 3 3/4" Foerster Forceps, Eye Dressing, Straight, Serrated Tips 0.5mm Wide
18-795A 3 3/4" Foerster Forceps, Eye Dressing, Full Curved, Serrated Tips 0.5mm Wide
18-796 3 3/4" Foerster Forceps, Iris Tissue, Straight, 1 x 2 Teeth, 0.5mm Wide
18-798 3 3/4" Foerster Forceps, Iris Tissue, Full Curved, 1 x 2 Teeth, 0.5mm Wide



7-600



7-610

MILTEX FOERSTER SPONGE FORCEPS

- 7-600 9 1/2" Sponge Forceps, Serrated, Straight
7-602 9 1/2" Sponge Forceps, Smooth, Straight
7-604 7" Sponge Forceps, Serrated, Straight
7-606 7" Sponge Forceps, Smooth, Straight
7-610 9 1/2" Sponge Forceps, Serrated
7-612 9 1/2" Sponge Forceps, Smooth, Curved
7-614 7" Sponge Forceps, Serrated, Curved
7-616 7" Sponge Forceps, Smooth, Curved

MILTEX GEMINI-MIXTER FORCEPS

Delicate with full curved jaws.

- 7-255 5 1/2" Gemini-Mixter Forceps
7-256 7" Gemini-Mixter Forceps
7-257 8" Gemini-Mixter Forceps
7-258 9" Gemini-Mixter Forceps
7-259 11" Gemini-Mixter Forceps



6-180




6-182



6-184

MILTEX GERALD FORCEPS

- 6-180 7" Gerald Forceps, Dressing, Straight, Serrated
6-182 7" Gerald Forceps, Tissue, Straight, 1 x 2 Teeth
6-184 7" Gerald Forceps, Dressing, Curved, Serrated
6-186 7" Gerald Forceps, Tissue, Curved, 1 x 2 Teeth



18-814



18-818



18-850

MILTEX GRAEFE FORCEPS

- 18-814 2 3/4" Graefe Forceps, Eye Dressing, Straight, Serrated
18-816 2 3/4" Graefe Forceps, Iris Tissue, Straight, 1 x 2 Teeth
18-818 2 3/4" Graefe Forceps, Eye Dressing, Curved, Serrated
18-820 2 3/4" Graefe Forceps, Iris Tissue, Curved, 1 x 2 Teeth
18-850 4 3/8" Graefe Forceps, Standard Jaws 4.5mm Wide w/fine teeth
18-852 4 3/8" Graefe Forceps, Standard Jaws 4.5mm Wide w/fine teeth, Non-Magnetic
18-854 4 3/8" Graefe Forceps, Standard Jaws 4.5mm Wide w/fine teeth & catch

MILTEX GRADLE CILIA FORCEPS

Precision ground jaws 1.4mm wide.

- 18-1116 4" Cilia Forceps

MILTEX GREEN FORCEPS

Delicate.

- 14-51 8 1/2" Green Forceps

MILTEX GREEN TISSUE FORCEPS


T-shape.

- 18-896 4" Tissue Forceps, 12mm Wide 19 x 20 Teeth
18-897 4" Tissue Forceps, 5mm Wide 8 x 9 Teeth

MILTEX GRUENWALD NASAL FORCEPS

4 1/2" shaft. Through-cutting.

- 20-594 Gruenwald Nasal Forceps, 2 x 10mm Bite
20-596 Gruenwald Nasal Forceps, 2.5 x 12mm Bite
20-598 Gruenwald Nasal Forceps, 3 x 12mm Bite



MILTEX HARTMAN-HERZFELD CUP-SHAPE FORCEPS

Alligator type.

- 19-429 3" Shaft Cup-Shape Forceps, 2mm Cup
- 19-430 3" Shaft Cup-Shape Forceps, 3mm Cup

MILTEX HEANEY FORCEPS

- 30-1700 8 1/4" Heaney Forceps, Heavy Pattern, Single Tooth, Curved
- 30-1710 8 1/4" Heaney Forceps, Heavy Pattern, Double Tooth, Curved
- 30-1715 9 3/4" Heaney Forceps, Extra Heavy Pattern, Double Tooth, Surved
- 30-1716 8 1/4" Heaney-Ballentine Forceps, Longitudinal Serrations, Single Tooth, Curved




MILTEX HEANEY TISSUE FORCEPS

2 x 3 teeth and cross serrated tips.

- 6-128 5 1/2" Tissue Forceps



17-2600



17-2602

MILTEX JACOBSON MICRO MOSQUITO FORCEPS

Extremely delicate.

- 17-2600 5" Jacobson Mosquito Forceps, Straight
- 17-2602 5" Jacobson Mosquito Forceps, Curved

MILTEX JANSEN DRESSING FORCEPS

Bayonet shape. Extra delicate pattern with 2mm wide serrated tips.

- 19-371 6 1/4" Dressing Forceps

MILTEX KELLY VULSELLUM FORCEPS


Curved sideways. 4 x 5 teeth. 7.7mm wide.

- 30-998 9 1/4" Vulsellum Forceps




MILTEX KERN BONE HOLDING FORCEPS

- 27-2 9 1/2" Bone Holding Forceps, Standard Size, Plain
- 27-4 9 1/2" Bone Holding Forceps, Standard Size w/ratchet
- 27-6 6" Bone Holding Forceps, Pediatric Size, Plain
- 27-7 6" Bone Holding Forceps, Pediatric Size w/ratchet



21-626



21-627

MILTEX KLEINERT-KUTZ (RUSKIN) BONE CUTTING FORCEPS

- 21-626 6" Bone Cutting Forceps, Straight
- 21-627 6" Bone Cutting Forceps, Curved

MILTEX LAHEY HEMOSTATIC FORCEPS

- 14-48 7 1/2" Hemostatic Forceps
- 14-49 9" Hemostatic Forceps
- 7-60 5 1/2" Hemostatic Forceps, Straight
- 7-62 5 1/2" Hemostatic Forceps, Curved

MILTEX LAMBERT CHALAZION FORCEPS

Inside ring diameter 8mm.

- 18-1174 3 1/2" Lambert Chalazion Forceps

MILTEX SCHAAF FOREIGN BODY FORCEPS

Grooved tips.

- 18-978 3 3/4" Foreign Body Forceps

MILTEX SINGLEY TISSUE FORCEPS

Fenestrated serrated jaws.

- 6-214 9" Singley Tissue Forceps



MILTEX SKIN STAPLE REMOVING FORCEPS

- 9-62 5 1/2" Staple Removing Forceps

MILTEX SMITHWICK CLIP APPLYING FORCEPS

Curved jaws and straight handles.

- 9-74 9" Smithwick Clip Applying Forceps

MILTEX STILLE PATTERN TISSUE FORCEPS

4 x 5 teeth.

- 6-110 5" Tissue Forceps



VANTAGE ADSON DRESSING FORCEPS

- V96-118 4 3/4" Dressing Forceps, Serrated Tips



VANTAGE ADSON TISSUE FORCEPS

- V96-120 4 3/4" Tissue Forceps, 1x2 Teeth



VANTAGE ALLIS TISSUE FORCEPS

- V916-10 7 1/4" Tissue Forceps, 5x6 Teeth
- V916-6 6" Tissue Forceps, 4x5 Teeth
- V916-8 6" Tissue Forceps, 5x6 Teeth

VANTAGE BOZEMAN DRESSING FORCEPS

- V97-622 10 1/2" Double Curved Dressing Forceps

VANTAGE BROWN-ADSON TISSUE FORCEPS


- V96-124 4 3/4" Brown-Adson Tissue Forceps, 7x7 Teeth

VANTAGE CRILE FORCEPS

- V97-42 5 1/2" Straight Crile Forceps
- V97-44 5 1/2" Curved Crile Forceps
- V97-46 6 1/4" Straight Crile Forceps
- V97-48 6 1/4" Curved Crile Forceps



V96-10



V96-4

VANTAGE DRESSING FORCEPS

- V96-1 3 1/2" Dressing Forceps
- V96-10 6" Dressing Forceps
- V96-14 8" Dressing Forceps
- V96-16 10" Dressing Forceps
- V96-18 12" Dressing Forceps
- V96-4 4 1/2" Dressing Forceps
- V96-6 5" Dressing Forceps
- V96-8 5 1/2" Dressing Forceps



VANTAGE FEILCHENFELD SPLINTER FORCEPS

- V96-316 3" Splinter Forceps
- V96-318 4 1/2" Splinter Forceps



VANTAGE FOERSTER SPONGE FORCEPS

- V97-600 9 1/2" Sponge Forceps, Straight, Serrated Jaws
- V97-604 7" Sponge Forceps, Straight, Serrated Jaws
- V97-610 9 1/2" Sponge Forceps, Curved, Serrated Jaws



V97-2



V97-4

VANTAGE HALSTEAD MOSQUITO FORCEPS

- V97-2 4 7/8" Halstead Mosquito Forceps, Straight
- V97-4 4 7/8" Halstead Mosquito Forceps, Curved



V919-434



V919-434A

VANTAGE HARTMAN ALLIGATOR FORCEPS

Serrated jaws.

- V919-434 3 1/4" Shaft Alligator Forceps
- V919-434A 3 1/4" Shaft Alligator Forceps, Extra Delicate Pattern

VANTAGE HARTMAN MOSQUITO FORCEPS

- V97-24 3 1/2" Straight Hartman Mosquito Forceps
- V97-26 3 1/2" Curved Hartman Mosquito Forceps

VANTAGE JANSEN BAYONET DRESSING FORCEPS

- V96-197 6 1/4" Dressing Forceps



V97-36



V97-38

VANTAGE KELLY FORCEPS

- V97-36 5 1/2" Kelly Forceps, Straight
- V97-38 5 1/2" Kelly Forceps, Curved




VANTAGE LISTON BONE CUTTING FORCEPS

- V927-1390 5 1/2" Straight Bone Forceps


VANTAGE LORNA NON-PERFORATING TOWEL FORCEPS

- V97-514 3 1/2" Towel Forceps
- V97-515 5 1/4" Towel Forceps




VANTAGE MAGILL CATHETER FORCEPS

- V92-70 9 3/4" Catheter Forceps, Adult
V92-72 8" Catheter Forceps, Child



VANTAGE NOYES ALLIGATOR FORCEPS

- V920-140 5 1/2" Shaft Alligator Forceps, Serrated Jaws
V920-142 6 1/2" Shaft Alligator Forceps, Serrated Jaws



VANTAGE POLYP FORCEPS

Angled shanks. Narrow jaws. Extra delicate.

- V930-996 10 1/2" Polyp Forceps

VANTAGE RANKIN FORCEPS

- V97-92 6 1/4" Rankin Forceps, Straight
V97-94 6 1/4" Rankin Forceps, Curved

VANTAGE ROCHESTER-OSHSNER FORCEPS

- V97-150 6 1/4" Straight Rochester-Ochsner Forceps
V97-152 7 1/4" Straight Rochester-Ochsner Forceps
V97-154 8" Straight Rochester-Ochsner Forceps
V97-156 9" Straight Rochester-Ochsner Forceps
V97-157 10" Straight Rochester-Ochsner Forceps
V97-160 6 1/4" Curved Rochester-Ochsner Forceps
V97-162 7 1/4" Curved Rochester-Ochsner Forceps
V97-164 8" Curved Rochester-Ochsner Forceps
V97-166 9" Curved Rochester-Ochsner Forceps
V97-167 10" Curved Rochester-Ochsner Forceps



VANTAGE ROCHESTER-PEAN FORCEPS

- V97-118 5 1/2" Straight Rochester-Pean Forceps
V97-120 6 1/4" Straight Rochester-Pean Forceps
V97-122 7 1/4" Straight Rochester-Pean Forceps
V97-124 8" Straight Rochester-Pean Forceps
V97-126 9" Straight Rochester-Pean Forceps
V97-128 10 1/2" Straight Rochester-Pean Forceps
V97-136 5 1/2" Curved Rochester-Pean Forceps
V97-138 6 1/4" Curved Rochester-Pean Forceps
V97-140 7 1/4" Curved Rochester-Pean Forceps
V97-142 8" Curved Rochester-Pean Forceps
V97-144 9" Curved Rochester-Pean Forceps
V97-146 10 1/4" Curved Rochester-Pean Forceps

VANTAGE RUSSIAN TISSUE FORCEPS

- V96-142 6" Tissue Forceps

VANTAGE SCHROEDER TENACULUM FORCEPS

- V930-965 9 1/2" Straight Tenaculum Forceps

VANTAGE SEMKIN FORCEPS

- V96-104 Semkin Dressing Forceps, Serrated
V96-106 Semkin Tissue Forceps, 1x2 Teeth



VANTAGE SPLINTER FORCEPS

- V96-300 3 1/2" Splinter Forceps, Fine Points
V96-304 4 1/2" Splinter Forceps, Fine Points
V96-329 3" Splinter Forceps w/magnifying lens in black plastic case


VANTAGE SPLINTER LIBERATORS

- V96-400 Splinter Liberators in Reversible Handle w/needle sharp tips



VANTAGE TISSUE FORCEPS

- V96-40 4 1/2" Tissue Forceps, 1x2 Teeth
V96-42 5" Tissue Forceps, 1x2 Teeth
V96-44 5 1/2" Tissue Forceps, 1x2 Teeth
V96-46 6" Tissue Forceps, 1x2 Teeth



V96-312



V96-314


VANTAGE WALTHER SPLINTER FORCEPS

V96-312 4 1/2" Straight Splinter Forceps


V96-314 4 1/2" Curved Splinter Forceps

VANTAGE WILDE EAR FORCEPS


V919-350 Wilde Ear Forceps, Serrated




7-42




7-44




7-46



7-48



7-50



7-52

MILTEK CRILE FORCEPS

7-42 5 1/2" Crile Forceps, Straight

7-44 5 1/2" Crile Forceps, Curved


7-46 6 1/4" Crile Forceps, Straight

7-48 6 1/4" Crile Forceps, Curved


7-50 5 1/2", Baby Crile Forceps, Straight, Extra Delicate

7-52 5 1/2", Baby Crile Forceps, Curved, Extra Delicate


**Fast, easy access to
the products you need.**




18-1934




7-10



7-20



7-4



7-8

MILTEK HALSTEAD MOSQUITO FORCEPS

18-1934 5" Halstead Mosquito Forceps, Straight, Non-Magnetic

18-1936 5" Halstead Mosquito Forceps, Curved, Non-Magnetic

7-10 5" Halstead Mosquito Forceps, Curved, Extra Delicate

7-14 5" Halstead Mosquito Forceps, Straight, 1 x 2 Teeth


7-16 5" Halstead Mosquito Forceps, Curved, 1 x 2 Teeth

7-2 5" Halstead Mosquito Forceps, Straight


7-20 5" Halstead Mosquito Forceps, Extra Delicate, Jaws 45° Angled On Flat

7-4 5" Halstead Mosquito Forceps, Curved

7-8 5" Halstead Mosquito Forceps, Straight, Extra Delicate



7-36




7-38

MILTEK KELLY FORCEPS

7-36 5 1/2" Kelly Forceps, Straight

7-38 5 1/2" Kelly Forceps, Curved




MILTEK KOCHER HEMOSTATIC FORCEPS

16-120 10 1/2" Hemostatic Forceps, Straight, Longitudinal Serrations

16-122 10 1/2" Hemostatic Forceps, Curved, Longitudinal Serrations

7-66 5 1/2" Hemostatic Forceps, Straight, 1 x 2 Teeth

7-68 5 1/2" Hemostatic Forceps, Curved, 1 x 2 Teeth



MILTEX PACKER MOSQUITO FORCEPS

Flat jaws.

7-11 5" Mosquito Forceps, Straight



MILTEX PENNINGTON HEMOSTATIC FORCEPS

7-220 6" Hemostatic Forceps



7-92




7-94


MILTEX RANKIN FORCEPS

7-92 6 1/4" Rankin Forceps, Straight


7-94 6 1/4" Rankin Forceps, Curved



7-170



7-172



7-180



7-182

MILTEX ROCHESTER-CARMALT FORCEPS

7-170 6 1/4" Rochester-Carmalt Forceps, Straight


7-172 6 1/4" Rochester-Carmalt Forceps, Curved

7-180 8" Rochester-Carmalt Forceps, Straight

7-182 8" Rochester-Carmalt Forceps, Curved




7-150



7-152



7-157



7-164

MILTEX ROCHESTER-OSHSNER FORCEPS

7-150 6 1/4" Straight Rochester-Ochsner Forceps

7-152 7 1/4" Straight Rochester-Ochsner Forceps

7-154 8" Straight Rochester-Ochsner Forceps

7-156 9" Straight Rochester-Ochsner Forceps

7-157 10" Straight Rochester-Ochsner Forceps


7-160 6 1/4" Curved Rochester-Ochsner Forceps

7-162 7 1/4" Curved Rochester-Ochsner Forceps


7-164 8" Curved Rochester-Ochsner Forceps

7-166 9" Curved Rochester-Ochsner Forceps


7-167 10" Curved Rochester-Ochsner Forceps



7-118



7-138



7-144

MILTEX ROCHESTER-PEAN FORCEPS

7-118 5 1/2" Straight Rochester-Pean Forceps

7-120 6 1/4" Straight Rochester-Pean Forceps

7-122 7 1/4" Straight Rochester-Pean Forceps

7-124 8" Straight Rochester-Pean Forceps

7-126 9" Straight Rochester-Pean Forceps

7-128 10 1/4" Straight Rochester-Pean Forceps

7-130 12" Straight Rochester-Pean Forceps

7-136 5 1/2" Curved Rochester-Pean Forceps

7-138 6 1/4" Curved Rochester-Pean Forceps

7-140 7 1/4" Curved Rochester-Pean Forceps

7-142 8" Curved Rochester-Pean Forceps

7-144 9" Curved Rochester-Pean Forceps

7-146 10 1/4" Curved Rochester-Pean Forceps

7-148 12" Curved Rochester-Pean Forceps

Hemostats



PRO ADVANTAGE® BACKHAUS TOWEL CLAMPS

Manufactured to German standards with exceptional value. Lifetime guarantee against defects in material or workmanship assures years of serviceability from each instrument. Selection of the most popular patterns of surgical instruments.

N407220 3 1/2" Backhaus Towel Clamp



7-504



7-506

MILTEX BACKHAUS TOWEL CLAMP

7-504 3 1/2" Towel Clamp

7-506 5 1/4" Towel Clamp



7-500



7-501

MILTEX JONES TOWEL CLAMP

7-500 3 1/2" Towel Clamp, Standard Size

7-501 2 1/4" Towel Clamp, Lightweight

AMSINO AMSURE® CLAMPS

TC05761 Towel Clamp, 2.5 CV, White, Non-Sterile, Bulk, 5000/cs

TC12431 Towel Clamp, 4.75 CV, Blue, Non-Sterile, Bulk, 2000/cs

UCC100 Umbilical Cord Clamp, Individually Packaged, Sterile, 50/cs

GALENICA DISPOSABLE EXTERNAL CLAMP

A commodity-type instrument used to avoid manipulation of contaminated or soiled items in a variety of ways. Designed for clamping tubes. Single use. Avoids deep cleaning and sterilization costs. Radiopaque and lightweight.

8901 Disposable External Clamp, 100/bx, 4 bx/cs

MEDICAL ACTION GENT-L-KARE® ALL-PURPOSE CLAMPS

Multitude of uses ranging from tube occlusion to securing gowns, drapes and towels. Polypropylene, Autoclavable.

4201-02 All Purpose Clamps w/serrated jaws, Turquoise, 100/cs

4205-02 Plastic Forceps, Turquoise, 2000/cs




MILTEX BIHRLE NEEDLE HOLDER/CLAMP

Carb-N-Sert with tungsten carbide inserts.

29-310TC 11" Bihrlle Needle Holder/Clamp

MILTEX ROEDER TOWEL CLAMP

7-510 5 1/4" Towel Clamp w/ball stops



VANTAGE BACKHAUS TOWEL CLAMP

V97-504 3 1/2" Towel Clamp

V97-506 5 1/4" Towel Clamp

VANTAGE JONES TOWEL CLAMP

V97-500 3 1/2" Towel Clamp



V97-585



V97-590

VANTAGE TUBE OCCLUDING CLAMP

V97-584 St. Vincents Occluding Clamp, 6 1/2", Curved, Smooth Jaws

V97-585 5 3/4" Occluding Clamp, Straight, Serrated Jaws w/guard

V97-588 4 1/2" Occluding Clamp, Straight, Smooth Jaws

V97-590 7" Occluding Clamp, Straight, Smooth Jaws

V97-591 7 1/2" Occluding Clamp, Straight, Serrated Jaws w/guard



MILTEX ALM RETRACTOR

11-12 2 3/4" Alm Retractor, Spread 2 1/4", Prongs 1/8" Deep

11-13 3 7/8" Alm Retractor, Spread 3", Prongs 1/2" Deep

Redefining Excellence

Miltex[®] Surgical Instruments



Retractors

MILTEX BECKMAN-WEITLANER RETRACTORS

Hinged blades. 3 x 4 teeth.

11-630-BL	5 1/2" Beckman-Weitlaner Retractor, Blunt
11-630-SH	5 1/2" Beckman-Weitlaner Retractor, Sharp
11-634-BL	6 1/2" Beckman-Weitlaner Retractor, Blunt
11-634-SH	6 1/2" Beckman-Weitlaner Retractor, Sharp
11-642-BL	9" Beckman-Weitlaner Retractor, Blunt
11-642-SH	9" Beckman-Weitlaner Retractor, Sharp

MILTEX COTTLE RETRACTORS

21-108	Thumb Retractor, 1 1/2" Working Length, Sharp Double Hook 8mm Wide
21-150	5 1/2" Cottle Retractor, 4 Blunt Prongs, 2 Center Prongs Recessed
21-165	8" Knife Guide & Retractor, Double End, Ball Tip Prongs 10mm Wide, 5mm Wide Slotted Hook

MILTEX DAVIS RETRACTOR

11-10	4 1/2" Davis Retractor, 3 x 3 Blunt Prongs
-------	--

MILTEX DEAVER RETRACTORS

11-200	1" x 9" Deaver Retractor
11-202	1" x 12" Deaver Retractor
11-204	1" x 13" Deaver Retractor
11-206	1 1/2" x 12" Deaver Retractor
11-208	2" x 12" Deaver Retractor
11-210	3" x 12" Deaver Retractor
11-218	5/8" x 8 1/2" Baby Deaver Retractor
11-219	3/4" x 7" Baby Deaver Retractor
11-220	7/8" x 8" Baby Deaver Retractor
11-221	1" x 12" Deaver Retractor w/hollow grip handle
11-222	1" x 13" Deaver Retractor w/hollow grip handle
11-223	1 1/2" x 12" Deaver Retractor w/hollow grip handle
11-224	2" x 12" Deaver Retractor w/hollow grip handle
11-224A	3" x 12" Deaver Retractor w/hollow grip handle
11-225	4" x 12" Deaver Retractor w/hollow grip handle



18-110



18-112

MILTEX DESMARRES RETRACTOR

18-110	5 1/2" Desmarres Retractor, Size 1, 11mm
18-112	5 1/2" Desmarres Retractor, Size 2, 13mm
18-114	5 1/2" Desmarres Retractor, Size 3, 15mm
18-116	5 1/2" Desmarres Retractor, Size 4, 17mm



21-368



21-369

MILTEX FOMON RASP


Double end.

21-368	8 1/4" Fomon Rasp, 4-Sided, Fine & Coarse Teeth, Both In Flat & Convex
21-369	8 1/4" Fomon Rasp, Coarse Teeth & Convex
21-370	8 1/4" Fomon Rasp, Fine Teeth & Convex

MILTEX FOMON RETRACTOR

Two ball end prongs 11mm wide.

21-148	6 1/2" Fomon Retractor
--------	------------------------



11-365



11-366

MILTEX GELPI RETRACTORS

11-365	5 1/2" Pediatric Gelpi Retractor w/grip lock
11-366	6 3/4" Gelpi Retractor w/grip lock
11-368	6 3/4" Gelpi Retractor w/ball stops & grip lock

MILTEX GOELET RETRACTOR

Double end.

11-120	7 1/2" Goelet Retractor
--------	-------------------------

MILTEX HARRINGTON RETRACTORS

11-314	12" Harrington Retractor, Blade 2 1/2" Wide
11-316	12" Harrington Retractor, Blade 1 1/2" Wide
11-318	9 1/2" Harrington Retractor, Blade 1" Wide



11-182



11-184

MILTEX HEANEY RETRACTORS

11-180	10 1/2" Heaney Retractor, Blade 1" x 1 3/4"
11-182	10 1/2" Heaney Retractor, Blade 1 1/8" x 2 7/8"
11-184	10 1/2" Heaney Retractor, Blade 1 1/8" x 3 1/4"
11-185	11 1/2" Heaney-Simon Retractor, Blade 1" x 4 1/2"

MILTEX HIBBS RETRACTOR

11-590	9 1/4" Hibbs Retractor, Blade 1" x 3"D
--------	--

MILTEX ISRAEL RETRACTORS

11-350	9" Israel Retractor, Blunt, 4 Prongs 1 1/2" x 1 5/8"
11-352	9" Israel Retractor, Blunt, 6 Prongs 2 1/4" x 1 3/4"

MILTEX JANSEN RETRACTOR

11-2	4 1/4" Jansen Retractor, Spread 3", 3 x 3 Blunt Prongs 10mm Deep, Standard Model
11-4	4 1/4" Jansen Retractor,, Spread 3", 4 x 4 Blunt Prongs 20mm Deep

MILTEX KELLY RETRACTORS


11-266	9 3/4" Kelly Retractor w/loop handle, 2" x 1 1/2"
11-267	9 3/4" Kelly Retractor w/loop handle, 2 1/2" x 2"
11-268	9 3/4" Kelly Retractor w/loop handle, 3" x 2 1/2"
11-269	9 3/4" Kelly Retractor w/loop handle, 3 1/2" x 2 3/4"
11-270	9 1/2" Kelly Retractor w/hollow grip handle, 2" x 1 1/2"
11-272	10" Kelly Retractor w/hollow grip handle, 2 1/2" x 2"
11-274	10" Kelly Retractor w/hollow grip handle, 3" x 2 1/2"
11-276	10" Kelly Retractor w/hollow grip handle, 3 1/2" x 2 3/4"

MILTEX KNAPP RETRACTORS

- 18-100 5 1/4" Knapp Retractor, 4 Sharp Prongs
18-102 5 1/4" Knapp Retractor, 4 Blunt Prongs

MILTEX LAHEY RETRACTOR

- 11-20 7 3/4" Lahey Retractor



MILTEX NECK RAKE RETRACTORS

- 11-30 6" Neck Rake Retractor, 1 Sharp Prong, Flexible
11-32 6" Neck Rake Retractor, 2 Sharp Prongs, Flexible
11-34 6" Neck Rake Retractor, 3 Sharp Prongs, Flexible
11-36 6" Neck Rake Retractor, 1 Blunt Prong, Flexible
11-38 6" Neck Rake Retractor, 2 Blunt Prongs, Flexible
11-40 6" Neck Rake Retractor, 3 Blunt Prongs, Flexible
11-50 6" Neck Rake Retractor, 1 Sharp Prong, Rigid
11-52 6" Neck Rake Retractor, 2 Sharp Prongs, Rigid
11-54 6" Neck Rake Retractor, 3 Sharp Prongs, Rigid
11-56 6" Neck Rake Retractor, 1 Blunt Prong, Rigid
11-58 6" Neck Rake Retractor, 2 Blunt Prongs, Rigid
11-60 6" Neck Rake Retractor, 3 Blunt Prongs, Rigid
11D-27 5" Neck Rake Retractor, 3 Sharp Prongs, Flexible



MILTEX RAGNELL RETRACTOR

Double end 3 x 8mm and 5 x 15mm blades. Delicate.

- 11-73 6" Ragnell Retractor

MILTEX RIBBON RETRACTORS

Malleable.

- 11-130 3/4" x 13" Ribbon Retractor
11-132 1" x 13" Ribbon Retractor
11-134 1 1/4" x 13" Ribbon Retractor
11-136 1 1/2" x 13" Ribbon Retractor
11-140 2" x 13" Ribbon Retractor
11-144 5/8" x 8" Baby Ribbon Retractor
11-146 1" x 8" Baby Ribbon Retractor

MILTEX RICHARDSON RETRACTORS

- 11-246 9 1/2" Richardson Retractor w/loop handle, 1" x 3/4"
11-247 9 1/2" Richardson Retractor w/loop handle, 1 1/4" x 1"
11-248 9 1/2" Richardson Retractor w/loop handle, 1 1/2" x 1 1/2"
11-249 9 1/2" Richardson Retractor w/loop handle, 2" x 3/4"
11-250 9 1/2" Richardson Retractor w/hollow grip handle, 1" x 3/4"
11-252 9 1/2" Richardson Retractor w/hollow grip handle, 1 1/4" x 1"
11-254 9 1/2" Richardson Retractor w/hollow grip handle, 1 1/2" x 1"
11-256 9 1/2" Richardson Retractor w/hollow grip handle, 2" x 3/4"
11-260 Richardson-Eastman Retractor, Double End, 2/set (10" & 11")
11-260-A Richardson-Eastman Retractor, Double End, 11"
11-260-B Richardson-Eastman Retractor, Double End, 10"



MILTEX RIGBY APPENDECTOMY RETRACTOR

- 11-374 6 3/4" Appendectomy Retractor w/grip lock



MILTEX SENN RETRACTORS

- 11-74 6 3/8" Senn Retractor, 3 Sharp Prongs
11-76 6 3/8" Senn Retractor, 3 Blunt Prongs

MILTEX VOLKMAN FINGER RETRACTORS

- 11-42 4 1/2" Finger Retractor, 2 Prongs, Sharp
11-43 4 1/2" Finger Retractor, 3 Prongs, Sharp
11-44 4 1/2" Finger Retractor, 4 Prongs, Sharp
11-45 4 1/2" Finger Retractor, 6 Prongs, Sharp
11-46 4 1/2" Finger Retractor, 2 Prongs, Blunt
11-47 4 1/2" Finger Retractor, 3 Prongs, Blunt
11-48 4 1/2" Finger Retractor, 4 Prongs, Blunt
11-49 4 1/2" Finger Retractor, 6 Prongs, Blunt




MILTEX U.S. ARMY RETRACTOR

Double end.


- 11-122 8 1/2" U.S. Army Retractor, 2/set

MILTEX VOLKMAN RETRACTORS


11-500	8 1/2" Volkman Retractor, Hook, Sharp, Hollow Handle
11-502	8 1/2" Volkman Retractor, Hook, Sharp, Ring Handle
11-504	8 1/2" Volkman Retractor, 2 Prongs, Sharp, Hollow Handle
11-506	8 1/2" Volkman Retractor, 3 Prongs, Sharp, Hollow Handle
11-508	8 1/2" Volkman Retractor, 4 Prongs, Sharp, Hollow Handle
11-510	8 1/2" Volkman Retractor, 6 Prongs, Sharp, Hollow Handle
11-512	8 1/2" Volkman Retractor, 2 Prongs, Blunt, Hollow Handle
11-514	8 1/2" Volkman Retractor, 3 Prongs, Blunt, Hollow Handle
11-516	8 1/2" Volkman Retractor, 4 Prongs, Blunt
11-530	8" Volkman Retractor, 2 Prongs, Sharp, Ring Handle
11-532	8" Volkman Retractor, 3 Prongs, Sharp, Ring Handle
11-534	8" Volkman Retractor, 4 Prongs, Sharp, Ring Handle
11-536	8" Volkman Retractor, 6 Prongs, Sharp, Ring Handle
11-538	8" Volkman Retractor, 2 Prongs, Blunt, Ring Handle
11-540	8" Volkman Retractor, 3 Prongs, Blunt, Ring Handle
11-542	8" Volkman Retractor, 4 Prongs, Blunt, Ring Handle
11-544	8" Volkman Retractor, 6 Prongs, Blunt, Ring Handle
11-546	8 3/4" Volkman Retractor, 2 Prongs, Sharp, USA Pattern Ring Handle
11-547	8 3/4" Volkman Retractor, 3 Prongs, Sharp, USA Pattern Ring Handle
11-548	8 3/4" Volkman Retractor, 4 Prongs, Sharp, USA Pattern Ring Handle
11-549	8 3/4" Volkman Retractor, 6 Prongs, Sharp, USA Pattern Ring Handle
11-551	8 3/4" Volkman Retractor, 2 Prongs, Blunt, USA Pattern Ring Handle
11-552	8 3/4" Volkman Retractor, 3 Prongs, Blunt, USA Pattern Ring Handle
11-553	8 3/4" Volkman Retractor, 4 Prongs, Blunt, USA Pattern Ring Handle
11-554	8 3/4" Volkman Retractor, 6 Prongs, Blunt, USA Pattern Ring Handle
27-230	9" Volkman Retractor, Hook, Sharp, 15mm Deep, Ring Handle
27-232	9" Volkman Retractor, Hook, Blunt, 15mm Deep, Ring Handle



11-600




11-608



11-610

MILTEX WEITLANER RETRACTORS

11-600	4" Weitlaner Retractor, 2 x 3 Teeth, Sharp
11-602	4" Weitlaner Retractor, 2 x 3 Teeth, Blunt
11-608	5 1/2" Weitlaner Retractor, 3 x 4 Teeth, Sharp
11-610	5 1/2" Weitlaner Retractor, 3 x 4 Teeth, Blunt
11-611	5 1/2" Weitlaner Retractor, Solid Blades, Blunt
11-614	6 1/2" Weitlaner Retractor, 3 x 4 Teeth, Sharp
11-616	6 1/2" Weitlaner Retractor, 3 x 4 Teeth, Blunt
11-617	8" Weitlaner Retractor, 3 x 4 Teeth, Sharp
11-618	8" Weitlaner Retractor, 3 x 4 Teeth, Blunt
11-620	9 1/2" Weitlaner Retractor, 3 x 4 Teeth, Sharp
11-622	9 1/2" Weitlaner Retractor, 3 x 4 Teeth, Blunt



MILTEX BARON SUCTION TUBES

7.5 cm tube with finger cut-off.

19-580	3 Fr Suction Tubes
19-582	5 Fr Suction Tubes
19-584	7 Fr Suction Tubes

MILTEX ANDREWS TONGUE DEPRESSOR

Serrated blade 21mm wide.

22-10	Tongue Depressor
-------	------------------

MILTEX BOSWORTH TONGUE DEPRESSOR

Fenestrated blade 25mm wide.

22-6	Tongue Depressor
------	------------------



22-2



22-4

MILTEX WIEDER TONGUE DEPRESSOR

Serrated.

22-2	Small Blade Tongue Depressor, 28mm Wide
22-4	Large Blade Tongue Depressor, 36mm Wide

OLON CARE™ TONGUE DEPRESSORS

52100	Solon Care™ 6" X 11/16" Non Sterile Tongue Depressor, 100/bx, 12 bx/cs
-------	--

MILTEX TONGUE DEPRESSOR

Double end, 20mm and 25mm wide. Cross serrated.

22-16	6 3/4" Tongue Depressor
-------	-------------------------



19-314



19-316


MILTEX BILLEAU FLEXIBLE EAR LOOPS

6 1/2".

19-314	Small Ear Loop, Size No. 1
19-316	Medium Ear Loop, Size No. 2
19-318	Large Ear Loop, Size No. 3

MILTEX BUCK (SPENCER) EAR PROBE

19-210 5 1/2" Ear Probe, Sterling Shaft

**MILTEX BUCK EAR CURETTES**

6 1/2"

19-260	00 Ear Curette, Sharp, Straight Size
19-262	0 Ear Curette, Sharp, Straight Size
19-264	1 Ear Curette, Sharp, Straight Size
19-266	2 Ear Curette, Sharp, Straight Size
19-268	3 Ear Curette, Sharp, Straight Size
19-270	00 Ear Curette, Sharp, Angled Size
19-272	0 Ear Curette, Sharp, Angled Size
19-274	1 Ear Curette, Sharp, Angled Size
19-276	2 Ear Curette, Sharp, Angled Size
19-278	3 Ear Curette, Sharp, Angled Size
19-280	00 Ear Curette, Blunt, Straight Size
19-282	0 Ear Curette, Blunt, Straight Size
19-284	1 Ear Curette, Blunt, Straight Size
19-286	2 Ear Curette, Blunt, Straight Size
19-288	3 Ear Curette, Blunt, Straight Size
19-290	00 Ear Curette, Blunt, Angled Size
19-292	0 Ear Curette, Blunt, Angled Size
19-294	1 Ear Curette, Blunt, Angled Size
19-296	2 Ear Curette, Blunt, Angled Size
19-298	3 Ear Curette, Blunt, Angled Size

**MILTEX DISPOSABLE EAR CURETTES**


Each size is color coded for easy identification. The handle is designed with a scored area for finger placement. Single use packaged in an individual clear sleeve.

19-320	Ring Tip Ear Curette, Pink, 50/bx
19-321	Loop Tip Ear Curette, White, 50/bx
19-322	Small Spoon Tip Ear Curette, Gray, 50/bx
19-323	Large Spoon Tip Ear Curette, Peach, 50/bx

MILTEX EAR SUCTION TUBE


Two interchangeable tips 20 ga and 22 ga.

19-1820	Ear Suction Tube
---------	------------------

**MILTEX EAR SYRINGES**

Chrome-plated, all brass syringes with precision ground plunger, and shaft graduated in 1/2 ounce. Complete with 1 tapered and 1 short bulbous tip. Easily disassembled for sterilization.

19-380	2 oz Ear Syringe w/o shield
19-381	2 oz Ear Syringe w/shield
19-382	3 oz Ear Syringe w/o shield
19-383	3 oz Ear Syringe w/shield
19-384	4 oz Ear Syringe w/o shield
19-385	4 oz Ear Syringe w/shield
19-390	Ear Syringe Short Bulbous Tip, Chrome
19-394	Ear Syringe Taper Tip, Chrome
19-396	Ear Syringe Shield, Chrome

**MILTEX FRAZIER FERGUSON SUCTION TUBES**

With finger cut-off.

19-570	6 Fr Suction Tubes, Angled
19-572	7 Fr Suction Tubes, Angled
19-573	8 Fr Suction Tubes, Angled
19-574	9 Fr Suction Tubes, Angled
19-575	10 Fr Suction Tubes, Angled
19-576	11 Fr Suction Tubes, Angled
19-577	12 Fr Suction Tubes, Angled
20-460	9 Fr Suction Tubes, Straight

MILTEX HARTMAN EAR FORCEPS

Delicate serrated jaws 2.5 x 6mm.

19-410 5" Hartman Ear Forceps

MILTEX HARTMAN EAR SPECULA

Round.

- 19-40 Set of 3 Ear Specula Sizes 3, 4 & 5mm
 19-48 Set of 4 Ear Specula Sizes 4, 5, 6 & 7mm
 19-48-A Single Size 3mm Ear Specula
 19-48-B Single Size 4mm Ear Specula
 19-48-C Single Size 5mm Ear Specula
 19-48-D Single Size 6mm Ear Specula
 19-48-E Single Size 7mm Ear Specula

MILTEX HOUSE SUCTION-IRRIGATION TUBES

9cm working length.

- 19-562 1.8mm Suction Tube & 1.8mm Irrigating Tube
 19-563 2.5mm Suction Tube & 2mm Irrigating Tube
 19-564 3mm Suction Tube & 2mm Irrigating Tube
 19-565 4.5mm Suction Tube & 2mm Irrigating Tube

**MILTEX JANSEN (GRUENWALD) EAR FORCEPS**

6-197 6 1/2" Ear Forceps, Bayonet Shape



20-28



20-32

MILTEX KILLIAN SEPTUM SPECULUM

6" with set screw.

- 20-28 2" Long Blades Septum Speculum
 20-29 2 1/2" Long Blades Septum Speculum
 20-30 3" Long Blades Septum Speculum
 20-32 3 1/2" Long Blades Septum Speculum



19-190



19-192

MILTEX LATHBURY APPLICATORS

- 19-190 5 1/2" Lathbury Applicator, Triangular Tip
 19-192 5 1/2" Lathbury Applicator, Cross Serrated Tip

**MILTEX LUCAE EAR FORCEPS**

19-370 5 1/2" Serrated Ear Forceps, Bayonet Shape



19-2081-B



19-2090-B




19-2100



19-2110


MILTEX MICRO ALLIGATOR EAR FORCEPS

- 19-2081-B Ear Forceps, Serrated Jaws 0.6mm Wide, Ebony Finish
 19-2090-B Ear Forceps, Oval Cup Jaws 0.5mm Wide, Ebony Finish
 19-2100 Ear Forceps, 3 1/4" Shaft, Serrated Jaws 0.8mm Wide
 19-2100-B Ear Forceps, Serrated Jaws 0.8mm Wide, Ebony Finish
 19-2110 Ear Forceps, Oval Cup Jaws 0.8mm Wide
 19-2110-B Ear Forceps, Oval Cup Jaws 0.8mm Wide, Ebony Finish

**MILTEX NOYES ALLIGATOR FORCEPS**

Spring handle.

19-404 4" Noyes Alligator Forceps

**MILTEX STRUEMPEL EAR FORCEPS**

Alligator type. Fenestrated oval spoon jaws.

19-438 3" Shaft Struempel Ear Forceps

**MILTEX WEINGARTNER EAR FORCEPS**

Alligator type. Concave serrated jaws.

19-437 3" Shaft Weingartner Ear Forceps




19-350




20-120

MILTEX WILDE EAR FORCEPS

- 19-350 5" Ear Forceps, Angular w/serrated tips
 19-356 5 1/4" Ear Forceps, Cross Action, Angular w/serrated tips
 20-120 6 1/2" Ear Forceps, Angular

**MILTEX WULLSTEIN EAR FORCEPS**

- 19-2120 3" Shaft Ear Forceps, Oval Cup, Jaws 1mm Wide
 19-2130 3" Shaft Ear Forceps, Round Cup, Jaws 1mm Diameter

**MILTEX WULLSTEIN EAR SCISSORS**


Straight and very delicate.

- 19-2060 3" Shaft Ear Scissors

**ADC TUNING FORKS**

Non magnetic, corrosion resistant aluminum alloy construction weighs 1/3 of comparable steel tuning forks. Produced from 3/8" x 1" bar stock for superior performance and consistent frequency accuracy. Extra long 2" handle of turned smooth aluminum to facilitate bone conduction tests.

- 500128 Tuning Fork w/fixed weight, 128cps Frequency
 500256 Tuning Fork w/fixed weight, 256cps Frequency
 500512 Tuning Fork w/o weight, 512cps Frequency
 501024 Tuning Fork w/o weight, 1024cps Frequency
 502048 Tuning Fork w/o weight, 2048cps Frequency
 504096 Tuning Fork w/o weight, 4096cps Frequency



19-100

19-102

19-104

19-106

19-108



19-110



19-112



19-120

MILTEX ALLOY TUNING FORKS

Aluminum alloy.

- 19-100 C-64 Vibrations Tuning Fork
 19-102 C-128 Vibrations Tuning Fork
 19-104 C-256 Vibrations Tuning Fork
 19-106 C-512 Vibrations Tuning Fork
 19-108 C-1024 Vibrations Tuning Fork
 19-110 C-2048 Vibrations Tuning Fork
 19-112 C-4096 Vibrations Tuning Fork
 19-120 Set of 5 Tuning Forks, C-128, C-256, C-512, C-1024 & C-2048 in Fitted Soft Case
 19-120A Fitted Zippered Case Only For Tuning Fork Set

MILTEX LEFORT URETHRAL SOUNDS

11" with standard 1-56 thread for filiform (not included).

29-10-10	10 Fr Urethral Sound
29-10-12	12 Fr Urethral Sound
29-10-14	14 Fr Urethral Sound
29-10-16	16 Fr Urethral Sound
29-10-18	18 Fr Urethral Sound
29-10-20	20 Fr Urethral Sound
29-10-22	22 Fr Urethral Sound
29-10-24	24 Fr Urethral Sound
29-10-26	26 Fr Urethral Sound
29-10-28	28 Fr Urethral Sound
29-10-30	30 Fr Urethral Sound
29-10-8	8 Fr Urethral Sound

MILTEX BELLUCCI MICRO EAR SCISSORS

4mm blades and ebony finish.

19-2150-B	2 3/4" Shaft Ear Scissors
-----------	---------------------------

**MILTEX BOUCHERON EAR SPECULA**

Chrome. Round.

19-2	Set of 4 Ear Specula: Sizes 1, 2, 3 & 4
19-2-1	Single Size 1 Ear Specula, 4.5mm
19-2-2	Single Size 2 Ear Specula, 5.5mm
19-2-3	Single Size 3 Ear Specula, 6.5mm
19-2-4	Single Size 4, 7.5mm

MILTEX CUT SNARE WIRE FOR NASAL SNARES

4" long with formed ends.

20-186-4	Size 4 Snare Wire
20-186-5	Size 5 Snare Wire

MILTEX DEAN TONSIL KNIFE

Straight. Double edge.

22-810	7 1/2" Tonsil Knife
--------	---------------------

MILTEX DENHARDT MOUTH GAG

With grip lock.

2-154	5" Denhardt Mouth Gag
-------	-----------------------



18-1455



19-211



19-2114



19-2115



20-170




20-171



20-172



5-269TC



8-3TC



9-115



9-116

MILTEX ENT INSTRUMENTS

Highly specialized for the ENT physician. Manufactured in Germany to the highest standards. Lifetime guarantee.

11-625	Paparella-Weitlaner Retractor, 2x2 Semi-Sharp Prongs, 5.1"
18-1455	Stevens Tenotomy Scissors Curved, 5 1/2"
18-1456	Stevens Tenotomy Scissors Curved, 6 1/4"
18-1997	Kerrison-Costen Ronguer, 4mm x 9mm Bite, Angled 90°, Curved Shaft, 5 3/4"
19-1341	Lempert Elevator 7.2" Curved, Heavy Pattern
19-211	Rosen Pick Curved 6 1/2" (3mm curve from tip)
19-2114	House Alligator & Crimper Forceps, 5.1"
19-2115	Bellucci Ear Forceps, 8mm Serrated Jaws, 5 1/2"
19-2116	Goodhill Ear Forceps, 5mm Serrated Jaws, 5 1/2"
19-2117	House Forceps, 6mm Serrated Jaws, Side Opening to the Right, 5.2"
19-2118	House Forceps, 9mm Oval Cups, 7.1cm Shaft, Angled 15° Right, 5.1"
19-2119	House Forceps, 9mm Oval Cups, 7.1cm Shaft, Angled 15° Left, 5 1/2"
19-212	Rosen Pick Curved 6 1/2" (1.5mm curve from tip)
19-2155	Bellucci Scissors, 7.5cm Shaft, Curved Left, 5.4"
19-2156	Bellucci Scissors, 7.5cm Shaft, Curved Right, 5.4"
19-2157	House-Bellucci Alligator Scissors, Straight, 7.4cm Shaft, 5 1/2"
19-2158	House-Bellucci Alligator Scissors, Curved Left, 7.4cm Shaft, 5 1/2"
19-2159	House-Bellucci Alligator Scissors, Curved Right, 7.4cm Shaft, 5 1/2"
19-2160	Bellucci Scissors, Straight, 7.7cm Shaft, 5.4"
21-205-11	Sheehan Osteotome 6 1/4, 11mm Hexagon Handle
21-205-3	Sheehan Osteotome 6 1/4, 3mm Hexagon Handle
21-205-5	Sheehan Osteotome 6 1/4, 5mm Hexagon Handle
21-205-7	Sheehan Osteotome 6 1/4, 7mm Hexagon Handle
21-205-9	Sheehan Osteotome 6 1/4, 9mm Hexagon Handle
21-61	Boies Nasal Fracture Elevator, 9mm Wide Tip, 7 1/2"
21-65	Goldman Displacer Septum Elevator, 6mm Wide Tip, 7 1/3"
21-661	Ferris-Smith Fragment Forceps, 6mm x 8mm Serrated Jaws, Fenestrated, Shaft 10.5cm, 7 1/3"
21-98	Single Skin Hook, Fine, Short Working End, 4 1/3"
5-263	Kaye Fine Dissecting Scissors, Serrated, 4 1/4"
5-269TC	Kaye Dissecting Scissors, Curved, Tungsten Carbide, 4 1/2"
5-400	Utility Scissors, Curved, Blunt Tips, Ribbon Handle Style, 4.1"
6-127	Griffiths-Brown Tissue Forceps, Delicate, 1.5mm x 5mm Brown-Type Teeth, 4 1/4"
8-3TC	Par Needle Holder, Straight, Tungsten-Carbide Inserts, Delicate 1.2cm Serrated Jaws, Ring Handle w/ratchet lock, 4 1/2"
8-5TC	Par Needle Holder, Straight, Tungsten-Carbide Inserts, Delicate 1.5cm Narrow Serrated Jaws, Ring Handle w/ratchet lock, 5"
9-115	Reeh Stitch Scissors, Sharp Pointed Tips, Small Hook On 1 Blade, Ribbon Style Handles, 3 3/4"
9-116	MilTEX Stitch Scissors, Curved, Sharp Pointed Tips, Ribbon Style Handles, 3 2/3"


MILTEX ENT INSTRUMENTS (Continued)

- 9-117 Miltex Stitch Removal Scissors, 4 3/4"
- 19-531 Rosen Suction Tube, 6cm Working Length, 15GA, 3 3/4"
- 19-532 Rosen Suction Tube, 6cm Working Length, 16GA, 3 3/4"
- 19-533 Rosen Suction Tube, 6cm Working Length, 17GA, 3 3/4"
- 19-534 Rosen Suction Tube, 6cm Working Length, 18GA, 3 3/4"
- 19-535 Rosen Suction Tube, 6cm Working Length, 19GA, 3 3/4"
- 19-536 Rosen Suction Tube, 6cm Working Length, 20GA, 3 3/4"
- 19-537 Rosen Suction Tube, 6cm Working Length, 22GA, 3 3/4"
- 19-538 Rosen Suction Tube, 6cm Working Length, 24GA, 3 3/4"
- 19-539 Rosen Suction Tube, 6cm Working Length, 26GA, 3 3/4"
- 19-540 House Cutoff Adaptor, 4mm Tubing End, Thumb Plate w/cutoff (bypass hole), For Use w/19-530 & 19-539, 2" Stainless Steel
- 19-541 House Cutoff Adaptor, 4mm Tubing End, Thumb Plate w/cutoff (bypass hole), For Use w/19-530 & 19-539, 2" Chrome
- 19-550 William House Suction-Irrigation Tube 2 x 2.5FR, 8.5cm Working Length, 5.1"
- 19-551 William House Suction-Irrigation Tube 2.5 x 4FR, 7.0cm Working Length, 4 1/2"
- 19-552 William House Suction-Irrigation Tube 2.5 x 4FR, 9.5cm Working Length, 5 1/2"
- 19-553 William House Suction-Irrigation Tube 4 x 5FR, 7cm Working Length, 4 1/2"
- 19-554 William House Suction-Irrigation Tube 5 x 7FR, 9.5cm Working Length, 5 1/2"
- 19-555 William House Suction-Irrigation Tube 4 x 5FR, 9.5cm Working Length, 5 1/2"
- 19-557 William House Suction-Irrigation Tube 7 x 10FR, 9.5cm Working Length, 5 1/2"
- 19-558 William House Suction-Irrigation Tube 8 x 12FR, 9.5cm Working Length, 5 1/2"
- 19-559 William House Suction-Irrigation Tube 12 x 14FR, 9.5cm Working Length, 5 1/2"
- 19-585 Schuknecht Suction Tube, 7.8cm Working Length, 20GA, 4 2/3"
- 19-586 Schuknecht Suction Tube, 7.8cm Working Length, 24GA, 4 2/3"
- 19-60 Shsea Speculum Holder, 11 2/3"
- 2-115 Dingman Mouth Gag, Complete Set of 3 Tongue Blades
- 2-116 Tongue Blade, Small, 24mm x 52mm Blade, 5 1/2"
- 2-117 Tongue Blade, Medium, 28mm x 55mm Blade, 5 1/2"
- 2-118 Tongue Blade, Small, 30mm x 68mm Blade, 5 1/2"
- 20-170 Nasal Scissors, Straight, Blades 1.1cm, Shaft 11cm, 7.1"
- 20-171 Nasal Scissors, Left, Blades 1.1cm, Shaft 11cm, 7.1"
- 20-172 Nasal Scissors, Right, Blades 1.1cm, Shaft 11cm, 7.1"
- 20-175 Becker Septum Scissors, Serrated, 7.1"
- 20-240 Ball End Elevator, Double Ended, 7 1/3"
- 20-426 Luer-Lok® Tonsil Needle, Straight, 23GA, 4.2"
- 20-427 Luer-Lok® Tonsil Needle, Curved, 23GA, 3 3/4"
- 20-470 Suction Nasal Forceps, 4mm Straight Jaws, Straight Shaft
- 20-471 Suction Nasal Forceps, 4mm Angled Jaws, Straight Shaft

**MILTEX FARRIOR EAR SPECULA**

Oval and angled end. Nonglare stainless.

- 19-50-1 Single Size 1 Ear Specula, 3 x 4mm
- 19-50-2 Single Size 2 Ear Specula, 4 x 5mm
- 19-50-3 Single Size 3 Ear Specula, 4.5 x 5.5mm
- 19-50-4 Single Size 4 Ear Specula, 5 x 6mm
- 19-50-5 Single Size 5 Ear Specula, 5.5 x 6.5mm
- 19-50-6 Single Size 6 Ear Specula, 6 x 7mm
- 19-50-7 Single Size 7 Ear Specula, 6.5 x 7.5mm
- 19-50-8 Single Size 8 Ear Specula, 7 x 8mm
- 19-50-9 Single Size 9 Ear Specula, 7.5 x 8.5mm

**MILTEX GRUBER EAR SPECULA**

- 19-10 Set of 4 Ear Specula Sizes 1, 2, 3 & 4 w/screw top case, Light Short Pattern, Chrome
- 19-12 Set of 4 Ear Specula Sizes 1, 2, 3 & 4 w/o case, Light Short Pattern, Chrome
- 19-14 Single Size 1 Ear Specula, 2x3mm, Light Short Pattern, Chrome
- 19-16 Single Size 2 Ear Specula, 3x4mm, Light Short Pattern, Chrome
- 19-18 Single Size 3 Ear Specula, 4x5mm, Light Short Pattern, Chrome
- 19-20 Single Size 4 Ear Specula, 5x6mm, Light Short Pattern, Chrome
- 19-21 Set of 4 Ear Specula Sizes 1, 2, 3 & 4 Oval, Standard Adult Sizes, Long Pattern, Chrome
- 19-22 Single Size 1 Ear Specula, 4x5mm Oval, Standard Adult Sizes, Long Pattern, Chrome
- 19-23 Single Size 2 Ear Specula, 5x6mm Oval, Standard Adult Sizes, Long Pattern, Chrome
- 19-24 Single Size 3 Ear Specula, 6x7mm Oval, Standard Adult Sizes, Long Pattern, Chrome
- 19-25 Single Size 4 Ear Specula, 7x8mm Oval, Standard Adult Sizes, Long Pattern, Chrome

Accessories

- 19-10A Case Only For Set of 4 Specula (Child & Infant Sizes), Light Short Pattern, Chrome

MILTEX HARTMAN TONSIL PUNCH

9mm round cut with basket.

- 22-1052 4 1/2" Shaft Tonsil Punch

MILTEX JOSEPH NASAL KNIVES

6". Double edge.

- 21-14 Large Blade Nasal Knife, Straight, Blunt Tip
- 21-16 Small Blade Nasal Knife, Curved, Blunt Tip
- 21-18 Large Blade Nasal Knife, Curved, Blunt Tip
- 21-2 Small Blade Nasal Knife, Straight, Pointed Tips
- 21-4 Large Blade Nasal Knife, Straight Pointed Tips

**MILTEX JOSEPH NASAL SCISSORS**

Sharp points.

- 21-630 5 3/4" Joseph Nasal Scissors, Straight
- 21-632 5 3/4" Joseph Nasal Scissors, Curved

**MILTEX KNIGHT NASAL FORCEPS**


Standard pattern. Serrated tips.

- 20-148 6 3/4" Knight Nasal Forceps

MILTEX KNIGHT SEPTUM FORCEPS

Cup jaws 15 x 5mm. 3 1/4" shaft.

20-522 6 3/4" Knight Septum Forceps



MILTEX TIECH (HALLE) INFANT NASAL SPECULUM

20-12 5 3/4" Infant Nasal Speculum

MILTEX TONSIL SNARE WIRE

Pre-cut 4" with bent ends.

22-1060-6 Size 6 Tonsil Snare Wire


22-1060-7 Size 7 Tonsil Snare Wire

22-1060-8 Size 8 Tonsil Snare Wire

22-1060-9 Size 9 Tonsil Snare Wire



20-2



20-3



20-4



20-7

MILTEX VIENNA NASAL SPECULA

20-2 5 3/4" Nasal Specula, Large

20-3 5 3/4" Nasal Specula, Medium

20-4 5 3/4" Nasal Specula Small

20-6 Lightweight, 5 3/4" Nasal Specula, Adult Size

20-7 Lightweight, 5 3/4" Nasal Specula, Child Size

VANTAGE VIENNA NASAL SPECULUM

V920-2 5 1/2" Nasal Specula, Large

V920-3 5 1/2" Nasal Specula, Medium

V920-4 5 1/2" Nasal Specula, Small

MILTEX MURDOCK EYE SPECULUM

Blades 15mm wide. Chrome.

18-50 2 1/4" Murdock Eye Speculum



18-132



18-132A

MILTEX TONOMETERS

Autoclavable Schiotz tonometers made with jeweled hammers to minimize friction. Scales are included, permitting reading from above. Scale markings etched for accuracy. Supplied in case with weights and test block.

18-132 Tonometer w/3 weights (5.5, 7.5 & 10 gr)

18-132A Same as No. 18-132, w/Allen plunger retractor

18-139 Sterilizing Case For Tonometers, 6 1/2" x 3 1/2" x 1 1/4", Fitted To Firmly Hold Tonometer, Weights & Test Block

BD VISIMARK™ MARKING PADS

581057 Marking Pad (Gentian Violet), 10/bx

BD VISITEC™ NUCLEUS CHOPPERS

585181 Angled 45° Choppers, 6mm From Tip (P. Koch), 3/bx

BD VISITEC™ ANTERIOR CHAMBER MAINTAINERS

585061 Anterior Chamber Maintainers, .90 x 3.5mm, OD - 2mm, Self Retaining, 3/bx

585149 Anterior Chamber Maintainers, 1.25 to .90mm Tapered, OD - 2mm, Oval, 5/bx

585163 Anterior Chamber Maintainers, .90 x 4mm, OD - 1.6mm, Micro Etched, 10/bx

BD VISITEC™ CAPSULE POLISHERS

581604 .80 x 22mm Capsule Polisher, Angled 35°, .30mm Top Port (Kratz), 10/bx

581605 .60 x 22mm Capsule Polisher, Angled 35°, .30mm Top Port (Kratz), 10/bx

581614 .80 x 22mm Capsule Polisher, Angled 30°, .50mm End Port, 10/bx

581645 .60 x 22mm Capsule Polisher, Angled 60°, .33mm End Port (Simcoe), 10/bx

585005 .40 x 22mm Capsule Polisher, Flex Tip, Angled 45°, 10mm Sleeve, 10/bx

585041 .60 x 22mm Capsule Polisher, Flex Tip, Angled 45°, 10mm Sleeve, 10/bx

585095 .50 x 24mm Capsule Polisher, Angled 30°, .30mm Top Port, Olive Tip (Pearce), 10/bx

585143 .60 x 25mm Capsule Polisher, Angled 30°, .33mm End Port, Olive Tip, 10/bx

BD VISITEC™ CORTEX EXTRACTORS

585075 .60 x 16mm Extractor, .30mm Top Port, 10/bx

585106 .55 x 12.5mm Extractor, .30mm End Port Angled 10° (Charleux), 10/bx

585193 .60 x 16mm Extractor, TW, .40mm Top Port, 10/bx

585227 .60 x 23.5 (Bmm Extractor, TW, .40mm Top Port (Blumenthal), 10/bx

BD VISITEC™ DELMAINATOR & VISCODISSECTOR

585092 .90 x 25mm Cannula w/.30 x 5mm Tip Extension, Flattened & Angled 90°, 1mm From End (Glaser), 5/bx

BD VISITEC™ EYE DRAINS & SHIELDS

581037	Shield Corneal Light 8.5mm, 10/pk
581098	40 x 14mm Eye Drain, 10mm Opening (Chayet), 20/bx
581706	Mini-Shield 4mm & 8mm w/straight edges (Banaji), 20/bx
581707	13mm Sponge (Salz), 20/bx
581713	11mm Sponge (Salz), 20/bx

BD VISITEC™ EYE PROTECTION PRODUCTS

581019	Eye Shield, 8.5 x 6.5cm (25 left/25 right), 50/bx
581040	Universal Eye Shield, 7.5 x 6.5cm, 50/bx
581062	Corneal Light Shield, 8mm, 20/bx
581705	ProOphta™ S Eye Patch, Non-Sterile, 50/bx
581708	Universal Pediatric Eye Shield, 5.1 x 4.1cm, 25/bx
583113	UV Sun Glasses w/case, Grey, Non-Sterile, 10/bx

BD VISITEC™ FRONTALIS SUSPENSION SET

585192	.90mm x 6.3cm Suspension Set, Silicon Tube Length 40cm, Tube Sleeve 10mm, ID - 1mm; OD - 2mm (Seiff), 3/bx
--------	--

BD VISITEC™ GLAUCOMA REVISION PICKS

585133	.40 x 16mm Picks, Angled 60°, 7mm From End (Lederer), 5/bx
585135	.40 x 16mm Picks, Angled 30°, 13.5mm From End (Lederer), 5/bx

BD VISITEC™ IOL MANIPULATORS

585022	Manipulator Fits .20mm Holes & Larger, Blunt Tip, Angled 45° (Jaffe/Maltzman), 3/bx
585038	Manipulator, .25mm Diameter Tip, Angled 45° (Lester), 3/bx
585062	Manipulator Fits .20mm Holes & Larger, Blunt Tip, Straight (Sinskey), 3/bx
585063	Manipulator Fits .20mm Holes & Larger, Blunt Tip, Angled 45° (Sinskey), 3/bx
585069	Manipulator, .50mm Button Tip, Angled 45°, 3/bx

BD VISITEC™ IRIS RETRACTORS

585209	Micro Iris Hooks, 8mm L x .08mm D (Engles), 5/pk, 5 pk/bx
--------	---

BD VISITEC™ IRRIGATING CYSTOTOMES

581594	.60 x 16mm Cystotomes, Straight, 10/bx
581599	.60 x 16mm Cystotomes, Formed, 10/bx
581602	.50 x 16mm Cystotomes, Formed, Side Cutting, 10/bx
581610	.50 x 16mm Cystotomes, Formed (Pearce), 10/bx
581617	.80mm/.80mm Cystotomes, .50mm End Aspiration Port, 10/bx
581618	.50 x 16mm Cystotomes, Straight, 10/bx
581630	.30 x 16mm Cystotomes, Formed, 10/bx
585004	.50 x 16mm Cystotomes, Formed, 10/bx
585008	.40 x 16mm Cystotomes, Formed, 10/bx
585030	.50 x 16mm Cystotomes, Curved (Joyce), 10/bx
585034	.50 x 16mm Cystotomes, Tip Angled 90° (Reverse), 10/bx
585039	.40 x 16mm Cystotomes, Curved (Joyce), 10/bx
585042	.40 x 16mm Cystotomes, Formed (Blumenthal), 10/bx
585049	.50 x 16mm Cystotomes, Tip Angled 45°, 10/bx
585093	.40 x 12.5mm Cystotomes, Angled 45° (Simcoe), 10/bx
585113	.50 x 16mm Cystotomes, Formed (Blumenthal), 10/bx

BD VISITEC™ IRRIGATING VECTIS

581619	.50mm Vectis, 5mm Wide, One .25mm End Port (Pearce), 5/bx
585032	.40mm Vectis, 3mm Wide, One .20mm End Port (Pearce), 5/bx

BD VISITEC™ LACRIMAL DILATOR & PROBE

585124	25mm Medium Taper, Blunt, Wilder Style Dilator Mated w/a Size 0 Bowman Lacrimal Probe .70mm x 45mm (Khan), 3/bx
--------	---

BD VISITEC™ LACRIMAL INTUBATION & DCR SETS

585010	.60mm x 17.5cm DCR Set, .30mm ID, .64mm OD, Silicon Tube Length 35cm, 3/bx
585011	.90mm x 4.5cm DCR Set, Straight, Silicon Tube Length 40cm (O'Donoghue), 3/bx
585012	.90mm x 4.5cm DCR Set, Angled 45°, 15mm From End, Silicon Tube Length 40cm (O'Donoghue), 3/bx
585013	.60mm x 11cm DCR Set, .30mm ID, .64mm OD, Silicon Tube Length 35cm w/retrieval device, 3/bx
585126	.40mm x 11cm DCR Set, .30mm ID, .64mm OD, Silicon Tube Length 35cm, Olive Shaped Tips, 3/bx
585151	.40mm x 11cm DCR Set, .30mm ID, .64mm OD, Silicon Tube Length 35cm, Olive Shaped Tips w/retrieval device, 3/bx

BD VISITEC™ LACRIMAL TREPHINE

585031	.80mm x 38mm Trephine, Ultra Sharp (Sisler), 3/bx
--------	---

BD VISITEC™ LASEK INSTRUMENTS

585284	LASEK Dissector, 5/bx
585285	LASEK Edge Lifter, 5/bx
585286	LASEK Microhoe Yee, 5/bx
585287	LASEK Cannula Yee 25ga, 10/bx

BD VISITEC™ LASIK RE-POSITIONERS

585233	Flap Roller, 14mm Sleeve, 5/bx
585254	Cap/Flap Elevator (Kritzinger/Updegraff), 5/bx

BD VISITEC™ LASIK SPEARS

581709	Expanded Spears, 6.85cm, 2/pk, 20 pk/bx
581751	PVA Spears, 7.11cm, 10/pk, 10 pk/bx

BD VISITEC™ MEMBRANE DISSECTORS

585114	.90 x 25mm Cannula w/.50 x 5mm Tip Extension, Sharp Tip Angled 50° (Glaser), 10/bx
585139	90mm Brush, Sleeve Extends 3mm Beyond End (Grizzard), 5/bx
585212	.90 x 25mm Cannula w/.30 x 5mm Tip Extension, Sharp Tip Angled 50° (Eaton), 10/bx

BD VISITEC™ NUCLEUS ROTATOR

585077	Blunt Rotator, Forked Tip, Angled 45°, 9mm From End (Jaffe/Bechert), 3/bx
--------	---

BD VISITEC™ NUCLEUS SPATULA

585141	.80mm Wide Spatula, Angled 45°, 10mm From End (Pearce), 5/bx
--------	--

BD VISITEC™ SHARP CYSTOTOMES

585020	.70 x 27mm Cystotomes, Straight, Vertical Sharp Cutting Edge, 5/bx
585053	.55 x 28mm Cystotomes, Straight, Curved Sharp Cutting Edge, 3/bx
585055	.55 x 16mm Cystotomes, Formed, Curved Sharp Cutting Edge, 3/bx

BD VISITEC™ SPATULAS

585189	Excimer Spatula (Lindstrom/Vorkas), 5/bx
585213	PRK Spatula, 5/bx

BD VISITEC™ SPECULUM

585094	40 x 26.5mm Speculum Overall Size, 13 x 6mm Retraction Size (Kratz/Barraquer), 3/bx
--------	---

BD VISITEC™ SUBRETINAL FLUID INSTRUMENTS

- 585021 Cannula, .90mm, Sleeve Extends 6mm Beyond End, Straight, Overall Length 35mm (Grizzard), 5/bx
- 585035 Cannula, .90mm, Sleeve Extends 6mm Beyond End, Curved Tip, Overall Length 35mm (Grizzard), 5/bx
- 585186 Cannula, .90mm, Sleeve Extends 3mm Beyond End, Straight, Overall Length 35mm (Grizzard), 5/bx

BD VISITEC™ SURGICAL GLIDES

- 581033 35mm L x 5mm W Surgical Glides, 2.5mm Radius, 6 mil Thick (sheets), 50/bx
- 581049 30mm L Surgical Glides, 6 mil Thick, Phaco 10mm L x 2.6mm W, Tapered to 1.75mm, IOL 20mm L x 5mm W, 2.5mm Radius (Fichman), 50/bx

BD VISITEC™ VITRECTOMY UNIT & ACCESSORIES

- 585088 Replacement Tubing 1mm ID, 2mm OD, May Be Used w/585101, 585198 & 585215, 10/bx
- 585100 Vitrectomy Unit Includes: Motor Handpiece, Control & Battery Power Unit, Power Cord, Screwdriver, Storage Case, 9V Batteries, Non-Sterile)
- 585101 Guillotine Cutter, .90mm Cutting Probe w/1.4mm removable irrigation sleeve, Includes Clear Aspirating Tubing, Luer-Lok™ Syringe 20mL, White Irrigation Tubing, Single Fluid Administration Set, 5/bx
- 585127 Vitrectomy Motor Unit, Non-Sterile
- 585128 Vitrectomy Power Cord, Non-Sterile
- 585198 Guillotine Cutter.60mm, .90mm Cutting Probe .90mm Tapers to .60mm, .5mm Aspiration Port Includes Clear Aspirating Tubing, Luer-Lok™ Syringe 3mL, Stopcock (Josephberg), 3/bx
- 585204 Vitreous Biopsy Kit Includes 585198 & 585231
- 585215 Vitreous Biopsy Probe, .75mm MVR Blade, Cutting Probe .90mm Tapers to .60mm, Includes Clear Aspirating Tubing, Luer-Lok™ Syringe 3mL, Stopcock (Josephberg), 3/bx



MILTEX BARRAQUER EYE SPECULUM

- 18-35 1 3/4" Eye Speculum, Large Blades 15mm Wide
- 18-36 1 3/4" Eye Speculum, Small Blades 10mm Wide



MILTEX CASTROVIEJO CORNEAL SCISSORS


- 18-1558 4" Corneal Scissors, Right Sub-Miniature For Micro Surgery, 6mm Blades
- 18-1559 4" Corneal Scissors, Left, Sub-Miniature For Micro Surgery, 6mm Blades
- 18-1560 4 1/2" Corneal Scissors, Right, Inner Blade 1mm Longer, Blunt Tips
- 18-1562 4 1/2" Corneal Scissors, Left, Inner Blade 1mm Longer, Blunt Tips
- 18-1570 3 3/4" Corneal Scissors, Keratoplasty, Angled Blades, 11mm, Sharp Points
- 18-1572 3 3/4" Corneal Scissors, Angled Blades, 11mm, Blunt Points
- 18-1576 3 3/4" Corneal Scissors, Curved, Sharp Points
- 18-1578 3 3/4" Corneal Scissors, Curved, Blunt Points

MILTEX CASTROVIEJO EYE SPECULUM


- 18-10 4" Castroviejo Eye Speculum, Large Blades 18 x 5mm O.D.
- 18-12 3 3/4" Castroviejo Eye Speculum, Medium Blades 15 x 5mm O.D.

MILTEX COOK EYE SPECULUM

- 18-63 1 7/8" Cook Eye Speculum, Infant Size, 9mm Blades
- 18-64 1 7/8" Cook Eye Speculum, Child Size, 11mm Blades



18-779



18-782



18-784

MILTEX EYE DRESSING FORCEPS

Serrated tips.

- 18-779 4" Eye Dressing Forceps, Straight, Delicate, 0.5mm Wide Tips
- 18-780 4" Eye Dressing Forceps, Straight, Standard, 0.8mm Wide Tips
- 18-781 4" Eye Dressing Forceps, Half Curved, Delicate, 0.5mm Wide Tips
- 18-782 4" Eye Dressing Forceps, Half Curved, Standard, 0.8mm Wide Tips
- 18-783 4" Eye Dressing Forceps, Full Curved, Delicate, 0.5mm Wide Tips
- 18-784 4" Eye Dressing Forceps, Full Curved, Standard, 0.8mm Wide Tips

MILTEX EYE MAGNET

A powerful magnetic cylinder with one round and one pointed interchangeable tip. Velvet lined steel case.

- 18-399 Eye Magnet



MILTEX EYE MAGNET & FOREIGN BODY LOOP

Magnet probe and Perlon loop are reversible in pencil-shaped handle. Plastic case.

- 18-398 Eye Magnet/Loop

MILTEX EYE SCISSORS

With probe points.

- 18-1430 4" Eye Scissors, Straight
- 18-1432 4" Eye Scissors, Curved
- 5-310 4 1/2" Eye Scissors, Angled On Side

MILTEX GILL CORNEAL KNIFE

Slight curved 3.5mm wide blade.

- 18-330 4 3/4" Gill Corneal Knife

MILTEX KNAPP CATARACT KNIFE

- 18-154 1 x 20mm Knapp Cataract Knife

MILTEX TOOKE CORNEAL KNIFE

2.5 x 15mm dissecting blade.

- 18-350 4 1/2" Tooke Corneal Knife

MILTEX WILLIAMS EYE SPECULUM

- 18-22 3 1/4" Williams Eye Speculum, Large Blades 15 x 5mm
18-24 2 3/4" Williams Eye Speculum, Small Blades 11 x 5mm




MILTEX FOX DERMAL CURETTES

5 1/2", round.

- 33-2 6 mm Dermal Curette
33-3 5 mm Dermal Curette
33-4 4 mm Dermal Curette
33-5 3 mm Dermal Curette
33-5A 2 mm Dermal Curette
33-5B 1 mm Dermal Curette

MILTEX 4" KEYES DERMAL PUNCHES


- 33-20 2 mm Dermal Punch
33-22 3 mm Dermal Punch
33-24 4 mm Dermal Punch
33-25 5 mm Dermal Punch
33-26 6 mm Dermal Punch
33-28 8 mm Dermal Punch
33-30 Keyes Dermal Punch Set w/3 1/2" Handle, Screw-In Keyes Knives
2,3,4,5,6 & 7 mm Mounted on Storage & Sterilization Rack
Handle For Dermal Punch
33-30-1 2 mm Replacement Head For Dermal Punch
33-30-2 3 mm Replacement Head For Dermal Punch
33-30-3 4 mm Replacement Head For Dermal Punch
33-30-4 5 mm Replacement Head For Dermal Punch
33-30-5 6 mm Replacement Head For Dermal Punch
33-30-6 7 mm Replacement Head For Dermal Punch



MILTEX GREG COMEDONE EXTRACTOR

One fenestrated cup and chuck with interchangeable No. 65 lancet blade, protected by cap.

- 33-215 5" Comedone Extractor



MILTEX PIFFARD DERMAL CURETTES

- 33-10 Size 3 Dermal Curette, 6 1/4" Oval
33-11 Size 4 Dermal Curette, 6 1/4" Oval
33-11A Size 5 Dermal Curette, 6 1/4" Oval
33-16-0 Size 0 Dermal Curette, 5 1/2" Oval Narrow Handles
33-16-00 Size 00 Dermal Curette, 5 1/2" Oval Narrow Handles
33-16-1 Size 1 Dermal Curette, 5 1/2" Oval Narrow Handles
33-16-2 Size 2 Dermal Curette, 5 1/2" Oval Narrow Handles
33-16-3 Size 3 Dermal Curette, 5 1/2" Oval Narrow Handles
33-6 Size 1 Dermal Curette, 6 1/4" Oval
33-8 Size 2 Dermal Curette, 6 1/4" Oval



MILTEX COMEDONE EXTRACTORS

- 33-203 4" Comedone Extractor, Double End, Flat & Wire End Expressors
33-214 6" Comedone Extractor, One Fenestrated Cup & One Delicate Curved Lancet w/protective cap

MILTEX DERMAL CURETTES

6 1/4". Oval spoon.

- 33-12 Size 0 Dermal Curette
33-13 Size 1 Dermal Curette
33-14 Size 2 Dermal Curette
33-15 Size 3 Dermal Curette



MILTEX SAALFIELD COMEDONE EXTRACTOR

With acne lancet. Professional pattern.

- 33-204 5 3/4" Comedone Extractor
33-206 4 1/4" Comedone Extractor



MILTEX SCHAMBERG COMEDONE EXTRACTOR


3 3/4".

- 33-200 Square Loop Ends Comedone Extractor
33-201 Standard Pattern Comedone Extractor w/crimped small loop
33-202 Standard Pattern Comedone Extractor w/extra fine loops, Small Loop Crimped

MILTEX DISPOSABLE DERMAL CURETTES

Features a unique symmetrically ground beveled edge. Manufactured of surgical grade stainless steel providing a clean, precise edge every time guaranteed. Individually packaged in sterile pouches.


- 33-52 2mm Dermal Curette, 50/bx
33-53 3mm Dermal Curette, 50/bx
33-54 4mm Dermal Curette, 50/bx
33-55 5mm Dermal Curette, 50/bx
33-57 7mm Dermal Curette, 50/bx



MILTEX STERILE DISPOSABLE BIOPSY PUNCHES

Seamless, stainless steel cutting edge with ribbed handle for maximum control.

33-31	2 mm Biopsy Punch, 50/bx
33-31-P/25	2 mm Biopsy Punch w/plunger, 25/bx
33-31A	1.5 mm Biopsy Punch, 50/bx
33-31A-P/25	1.5 mm Biopsy Punch w/plunger, 25/bx
33-31AA	1mm Biopsy Punch, 50/bx
33-31AA-P/25	1 mm Biopsy Punch w/plunger, 25/bx
33-31B	2.5mm Biopsy Punch, 50/bx
33-32	3 mm Biopsy Punch, 50/bx
33-32-P/25	3mm Biopsy Punch, 25/bx
33-33	3.5 mm Biopsy Punch, 50/bx
33-34	4 mm Biopsy Punch, 50/bx
33-34-P/25	4mm Biopsy Punch, 25/bx
33-35	5 mm Biopsy Punch, 50/bx
33-36	6 mm Biopsy Punch, 50/bx
33-37	8 mm Biopsy Punch, 50/bx
33-38	Assorted Biopsy Punch, 50/bx



PRO ADVANTAGE® NAIL NIPPERS & SPLITTERS

Expertly crafted from surgical grade stainless steel and meets the very highest quality standards and specifications. This outstanding quality and price point makes Pro Advantage instruments ideal for today's growing instrument market.

N407525	4 1/2" Nail Nippers, Stainless, Straight Jaws, Double Spring
N407540	5" Nail Nippers, Stainless, Straight Jaws, Double Spring
N407550	6 Nail Nippers, Stainless, Straight Jaws, Double Spring
N407560	5" Nail Splitter, Anvil Type



AARON NAIL DRILL & BITS

0007 Bit, Replaceable Drill Tip, 10/bx



MILTEX SWISS CILIA & SUTURE

3 3/4" forceps. Excellent for suture removal.

18-1107 Swiss Cilia & Suture



MILTEX NAIL SPLITTER

40-230 5" Nail Splitter, Anvil Type



MILTEX PLATYPUS NAIL PULLING FORCEPS

Stainless.

40-234 5 1/2" Nail Pulling Forceps, Standard Wide Jaws
40-235 5 1/2" Nail Pulling Forceps, Narrow Jaws




40-220



40-245-CH



40-245-SS



40-250-SS



40-251-SS



40-252

MILTEX TISSUE & CUTICLE NIPPERS

40-220	4" Cuticle Nipper, Ring Handle, Chrome, Straight Jaws
40-245-CH	4" Cuticle Nipper, Chrome, 5mm, Convex Jaws
40-245-SS	4" Cuticle Nipper, Stainless, 5mm, Convex Jaws
40-250-CH	4 1/2" Cuticle Nipper, Chrome, 10mm, Convex Jaws
40-250-SS	4 1/2" Cuticle Nipper, Stainless, 10mm, Convex Jaws
40-251-CH	5" Cuticle Nipper, Chrome, 10mm, Convex Jaws
40-251-SS	5" Cuticle Nipper, Stainless, 10mm, Convex Jaws
40-252	5" Cuticle Nipper, Heavy Pattern, Stainless, Convex Jaws



MILTEX CUTICLE SCISSORS

Curved blades. Chrome.

- 40-435 3 1/2" Cuticle Scissors, Standard
40-445 3 1/2" Cuticle Scissors, Extra Delicate



MILTEX INGROWN NAIL SPLITTING SCISSORS

Straight with one serrated blade.

- 40-260 6" Nail Splitting Scissors



MILTEX NAIL SCISSORS

Curved blades. Chrome.

- 40-425 3 1/2" Nail Scissors



MILTEX FINGER NAIL DRILL

- 33-230 4" Finger Nail Drill, Standard Pattern, Complete w/3 drill points stored in handle
33-232 4 1/2" Finger Nail Drill w/chuck head, Complete w/3 drill points stored in handle
33-234 Set of 3 Drill Points: Small, Medium & Large, Chrome



DMI® NAIL NIPPERS

- 539-5406-0000 Toe Nail Nipper, 5"

GRAHAM FIELD NAIL NIPPER

- 1790 5 1/2" Nail Nipper
1792 4 1/2" Nail Nipper, Heavy Duty



PRO ADVANTAGE® TISCHLER BIOPSY PUNCH

Manufactured to German standards with exceptional value. Lifetime guarantee against defects in material or workmanship assures years of serviceability from each instrument. Selection of the most popular patterns of surgical instruments.

- N407435 Baby-Tischler Biopsy Punch
N407440 Tischler Biopsy Punch 3x7mm



BUSSE LAMINARIA CERVICAL DILATOR

Latex free.

- 206 +3mm (Small) Cervical Dilator, 12/dispenser, 4 dispensers/cs
207 +4mm (Medium) Cervical Dilator, 12/dispenser, 4 dispensers/cs
208 +5mm (Large) Cervical Dilator, 12/dispenser, 4 dispensers/cs
209 +6mm (Extra Large) Cervical Dilator, 12/dispenser, 4 dispensers/cs

BUSSE POSI-GRIP™ UMBILICAL CORD CLAMP

Designed to seal the umbilical cord, preventing loss of blood as the cord dries beneath it. The clamps locks permanently and can be easily placed with one hand, even if gloved and slippery. Latex free.

- 384 Sterile Umbilical Cord Clamp, 64/bx, 2 bx/cs
385 Non-Sterile Umbilical Cord Clamp, 250/cs



COOPERSURGICAL SANI-SPEC® SPECULA

Clear plastic material allows 360° visibility and reduces thermal shock. Disposable. Provides protection of the doctor, the nurse and the patient from cross contamination. Plastic material avoids unpleasant temperature shock associated with metal. Smooth pre-lubricated blades with 6 position grooves.

- 82329 Small Specula, 1"W x 4 1/2"L, 10/bg, 10 bg/cs
82330 Medium Specula, 1 1/4"W x 4 3/4"L, 10/bg, 10 bg/cs
82331 Large Specula, 1 1/3"W x 5 1/3"L, 10/bg, 10 bg/cs
82340 Small Specula, Individually Wrapped, 50/cs
82341 Medium Specula, Individually Wrapped, 50/cs

COOPERSURGICAL SANI-SCOPE®

Clear plastic material allows 360° visibility and reduces thermal shock. Disposable. Allows protection for the doctor, the nurse and the patient from cross contamination. Single use obturator.

82420 9/10" Dia x 3 3/8" L Sani-Scope®, 25/bg, 4 bg/cs

**GALENICA VAGINAL SPECULA**

Recommended for routine gynecological exams. Disposable. Allows the examination of the vaginal mucosa at 360° angle without rotation. Wide range of adjustability. Ease of operation and manipulation with one hand. Self-holding in various opening positions. Made from hard, clear and non-toxic plastic resin of high strength. Unit Pack (poly bag) non-sterile. Unit Pack (peeling pouch) sterile.

8309 Small Vaginal Specula, Non-Sterile, 10/bg, 10 bg/cs
 8310 Medium Vaginal Specula, Non-Sterile, 10/bg, 10 bg/cs
 8311 Large Vaginal Specula, Non-Sterile, 10/bg, 10 bg/cs
 8615 Small Unit Pack Vaginal Specula, 50/cs
 8616 Medium Unit Pack Vaginal Specula, 50/cs
 8617 Large Unit Pack Vaginal Specula, 50/cs

Accessories

9310 Smoke Evacuator For All Specula, 40/cs

MILTEX BARRETT UTERINE TENACULUM FORCEPS

30-985 7 1/4" Uterine Tenaculum Forceps, Straight

MILTEX BOZEMAN SPONGE FORCEPS

10 1/2".

7-622 10 1/2" Double Curved Sponge Forceps

**MILTEX BOZEMAN UTERINE DRESSING FORCEPS**

7-620 10 1/2" Straight Dressing Forceps
 7-624 Double Curved Dressing Forceps w/1 large finger ring

MILTEX BRAUN EPISIOTOMY SCISSORS

Angled to side. Guarded lower blade.

30-2190 5 1/2" Episiotomy Scissors
 30-2195 8 1/2" Episiotomy Scissors

MILTEX AUWARD WEIGHTED VAGINAL SPECULA

9".

30-186 2 lbs Weighted Vaginal Specula, Straight Blade 1 1/2" x 3 1/4", Slightly Angled 75°
 30-187 2.5 lbs Weighted Vaginal Specula, Straight Blade 1 1/2" x 3 1/4", Slightly Angled 75°
 30-188 3 lbs Weighted Vaginal Specula, Straight Blade 1 1/2" x 3 1/4", Slightly Angled 75°
 30-193 2.5 lbs Weighted Vaginal Specula, Blade 4" x 1 3/8", Tapering To 1 3/4" At Tip, Blade At Acute Angle 45°
 30-195 2.5 lbs Weighted Vaginal Specula, Steiner-Auward, Slightly Curved Blade 5 1/2" x 1 1/4" At 90° Angle


MILTEX BUMM UTERINE CURETTES

12". Extra large sizes on malleable shanks.

30-1270-1 Size 1 Uterine Curette, Sharp
 30-1270-2 Size 2 Uterine Curette, Sharp
 30-1270-3 Size 3 Uterine Curette, Sharp
 30-1270-4 Size 4 Uterine Curette, Sharp
 30-1275-1 Size 1 Uterine Curette, Blunt
 30-1275-4 Size 4 Uterine Curette, Blunt

MILTEX COLLIN VAGINAL SPECULUM

30-120 Small Vaginal Specula, 1 1/4" x 3 3/4"
 30-125 Medium Vaginal Specula, 1 3/8" x 4"
 30-130 Large Vaginal Specula, 1 5/8" x 4 1/2"

**MILTEX CUP JAW UTERINE BIOPSY FORCEPS**

Plain not interfitted cups. 6mm bite.

30-1480 9" Shaft Uterine Biopsy Forceps

MILTEX DEVILBISS VAGINAL SPECULA

Two sets of screws permit parallel or bivalve opening of the blades for examination or operative use.

30-140 Small Vaginal Specula, 3 1/2" x 7/8"
 30-141 Medium Vaginal Specula, 4" x 1 1/4"
 30-142 Large Vaginal Specula, 5" x 1 3/8"

**MILTEX DUPLAY UTERINE TENACULUM FORCEPS**

Double curved.


30-975 11" Uterine Tenaculum Forceps



MILTEX EMMETT UTERINE TENACULUM HOOKS

9". Sharp points.

- 30-951 Style 1 Uterine Tenaculum Hook, Half Curved, Small
- 30-952 Style 2 Uterine Tenaculum Hook, Acute Angle
- 30-953 Style 3 Uterine Tenaculum Hook, Right Angle
- 30-954 Style 4 Uterine Tenaculum Hook, Half-Curved, Large
- 30-955 Style 5 Uterine Tenaculum Hook, Double Angle



MILTEX ENDOMETRIAL BIOPSY CANNULA

Offered in either a SoftFlex™ or FirmFlex™ style for easy one-step sampling. The SoftFlex™ is ideal for standard or routine patient sampling applications where a more flexible cannula is desired by the clinician. The FirmFlex™ is a slightly more rigid cannula ideal for more controlled insertion and sampling in patient cases where atrophic endometrium or post-menopausal patients with stenosis of the cervix is a concern. Both the SoftFlex™ and the FirmFlex™ patterns are sterilized via ETO, packaged individually for single-use applications, and provided in convenient quantities of 25 per box. The cannula are made of a clean, latex-free and flexible polypropylene outer sleeve with inner plunger. The cannula sleeve, with plunger in closed position, is inserted through the cervix and into the uterine cavity. Once properly positioned, the plunger is withdrawn and the design generates sufficient vacuum for maximum endometrial sample yield to be captured within the cannula sleeve without the need for external suction devices. Detailed directions for use are packaged with cannula and should be reviewed thoroughly and carefully prior to each use or demonstration.

- 30-3010 SoftFlex™ Endometrial Biopsy Cannula, Sterile, Single Use, 25/bx
- 30-3012 FirmFlex™ Endometrial Biopsy Cannula, Sterile, Single Use, 25/bx

MILTEX ENDOMETRIAL SAMPLING SET

A disposable plastic curette and syringe system which is very similar in function and appearance to the stainless steel Randall endometrial suction biopsy curette pattern. Includes a disposable curette and a 10ml Twist-and-Lock syringe. The curette is made of clear, latex-free polypropylene and has a semi-rigid design with a slightly curved and tapered Randall-type cutting edge tip which optimizes access to the uterine contours and efficient sample collection. The semi-rigid design of the curette allows for controlled insertion in female patients with stenosis of the cervix, and the syringe system creates 21" Hg of negative pressure for optimized suction and sample capture in situations where the endometrium is atrophic and more difficult to capture sufficient sample yield.

- 30-3015 Endometrial Sampling Set, Includes: Disposable Curette w/a Twist-and-Lock Syringe, Sterile, Single Use, 10/bx



MILTEX EPPENDORFER UTERINE BIOPSY FORCEPS

- 30-1485 8 1/2" Shaft Uterine Biopsy Forceps



30-1492



30-1493

MILTEX EPPENDORFER-KRAUSE BIOPSY FORCEPS

Shaft interchangeable on pistol grip spring handle. Oval 4.5 mm x 6.5 mm bite. Complete with pistol grip handle. Rotating punch.

- 30-1492 9" Shaft Biopsy Forceps
- 30-1493 Biopsy Forceps Shaft Only



30-10




30-10B

MILTEX GRAVES VAGINAL SPECULA

- 30-10 Small Vaginal Specula, 3/4" x 3", Original Model
- 30-10B Small Vaginal Specula, 3/4" x 3", Ebony Finish
- 30-15 Medium Vaginal Specula, 1 3/8" x 4", Original Model
- 30-15B Medium Vaginal Specula, 1 3/8" x 4", Ebony Finish
- 30-20 Large Vaginal Specula, 1 1/2" x 4 1/2", Original Model
- 30-20B Large Vaginal Specula, 1 1/2" x 4 1/2", Ebony Finish
- 30-25 Medium Vaginal Specula, 1 3/8" x 4", Wide Angle Blades, Improved
- 30-30 Medium Vaginal Specula, 1 3/8" x 4", Wide Angle Blades, Open Side

**Please ask before
ordering a product based
on the pictured item.**

**MILTEX HANK UTERINE DILATORS**


11". Double ended.

30-500	Uterine Dilators, Set of 6, sizes: 11-12 Fr to 21-22 Fr in a Khaki Roll
30-505-1112	11-12 Fr Uterine Dilator
30-505-1314	13-14 Fr Uterine Dilator
30-505-1516	15-16 Fr Uterine Dilator
30-505-1718	17-18 Fr Uterine Dilator
30-505-1920	19-20 Fr Uterine Dilator
30-505-2122	21-22 Fr Uterine Dilator
30-505-2324	23-24 Fr Uterine Dilator
30-505-2526	25-26 Fr Uterine Dilator
30-505-2728	27-28 Fr Uterine Dilator
30-505-2930	29-30 Fr Uterine Dilator
30-505-78	7-8 Fr Uterine Dilator
30-505-910	9-10 Fr Uterine Dilator

MILTEX HEANEY UTERINE BIOPSY CURETTE

Serrated loop.

30-1385	9" Uterine Biopsy Curette
---------	---------------------------

**MILTEX HEGAR UTERINE DILATORS**

30-515	Uterine Dilator, 7", Single End, Set of 14: Sizes 4 to 17 mm In Khaki Roll
30-515A	Khaki Roll Only For Set of 14 Uterine Dilator
30-520	Uterine Dilator, 7", Single End, Set of 26: Sizes 1 to 26 mm In Khaki Roll
30-520A	Khaki Roll Only For Set of 26 Uterine Dilator
30-530-1	1 mm Uterine Dilator, 7", Single End
30-530-10	10 mm Uterine Dilator, 7", Single End
30-530-105	10.5 mm Uterine Dilator, 7", Single End
30-530-11	11 mm Uterine Dilator, 7", Single End
30-530-115	11.5 mm Uterine Dilator, 7", Single End
30-530-12	12 mm Uterine Dilator, 7", Single End
30-530-125	12.5 mm Uterine Dilator, 7", Single End
30-530-13	13 mm Uterine Dilator, 7", Single End
30-530-135	13.5 mm Uterine Dilator, 7", Single End
30-530-14	14 mm Uterine Dilator, 7", Single End
30-530-145	14.5 mm Uterine Dilator, 7", Single End
30-530-15	1.5 mm Uterine Dilator, 7", Single End
30-530-150	15 mm Uterine Dilator, 7", Single End
30-530-155	15.5 mm Uterine Dilator, 7", Single End
30-530-16	16 mm Uterine Dilator, 7", Single End
30-530-165	16.5 mm Uterine Dilator, 7", Single End
30-530-17	17 mm Uterine Dilator, 7", Single End
30-530-175	17.5 mm Uterine Dilator, 7", Single End
30-530-18	18 mm Uterine Dilator, 7", Single End
30-530-185	18.5 mm Uterine Dilator, 7", Single End
30-530-19	19 mm Uterine Dilator, 7", Single End
30-530-2	2 mm Uterine Dilator, 7", Single End
30-530-20	20 mm Uterine Dilator, 7", Single End
30-530-21	21 mm Uterine Dilator, 7", Single End
30-530-22	22 mm Uterine Dilator, 7", Single End
30-530-23	23 mm Uterine Dilator, 7", Single End
30-530-24	24 mm Uterine Dilator, 7", Single End
30-530-25	2.5 mm Uterine Dilator, 7", Single End
30-530-250	25 mm Uterine Dilator, 7", Single End
30-530-26	26 mm Uterine Dilator, 7", Single End
30-530-3	3 mm Uterine Dilator, 7", Single End
30-530-35	3.5 mm Uterine Dilator, 7", Single End
30-530-4	4 mm Uterine Dilator, 7", Single End
30-530-45	4.5 mm Uterine Dilator, 7", Single End
30-530-5	5 mm Uterine Dilator, 7", Single End
30-530-55	5.5 mm Uterine Dilator, 7", Single End
30-530-6	6 mm Uterine Dilator, 7", Single End
30-530-65	6.5 mm Uterine Dilator, 7", Single End
30-530-7	7 mm Uterine Dilator, 7", Single End
30-530-75	7.5 mm Uterine Dilator, 7", Single End
30-530-8	8 mm Uterine Dilator, 7", Single End
30-530-85	8.5 mm Uterine Dilator, 7", Single End
30-530-9	9 mm Uterine Dilator, 7", Single End
30-530-95	9.5 mm Uterine Dilator, 7", Single End
30-550	Uterine Dilators, 7 1/2", Double End, Set of 8: 3-5 mm & 17-18 mm In Khaki Roll
30-550A	Khaki Roll Only For Set of 8 Uterine Dilator
30-555-1112	11-12 mm Uterine Dilator, 7 1/2", Double End
30-555-1314	13-14 mm Uterine Dilator, 7 1/2", Double End
30-555-1516	15-16 mm Uterine Dilator, 7 1/2", Double End
30-555-1718	17-18 mm Uterine Dilator, 7 1/2", Double End
30-555-34	3-4 mm Uterine Dilator, 7 1/2", Double End
30-555-56	5-6 mm Uterine Dilator, 7 1/2", Double End
30-555-78	7-8 mm Uterine Dilator, 7 1/2", Double End
30-555-910	9-10 mm Uterine Dilator, 7 1/2", Double End

Can't find what
you're looking for?
Call your
representative.

MILTEX KIELLAND OBSTETRICAL FORCEPS

30-2350 15 1/2" Obstetrical Forceps

MILTEX KIELLAND-LUIKART OBSTETRICAL FORCEPS

30-2355 15 1/2" Obstetrical Forceps, Solid Blades



30-1350

30-1351

MILTEX KOGAN ENDOCERVICAL SPECULUM

30-1350 9 1/2" Endocervical Speculum w/gauge & set screw, Offset Ring Handles, Fenestrated Jaws 7 x 25mm

30-1351 11" Endocervical Speculum, Improved Pattern w/ratchet & gauge, Extra Delicate Fenestrated Jaws 3mm Wide At Tip x 25mm Long


MILTEX SIMPSON OBSTETRICAL FORCEPS

30-2300 14" Obstetrical Forceps, Long Model

30-2305 12" Obstetrical Forceps, Short Model

MILTEX SIMPSON-LUIKART OBSTETRICAL FORCEPS

30-2400 14" Obstetrical Forceps, Solid Blades

**MILTEX SIMS RECTAL SPECULUM**

Fenestrated blades. Stainless steel.

28-2 Sims Rectal Speculum

MILTEX SIMS VAGINAL SPECULUM

30-174 Single End Vaginal Specula, Small 1 1/16" x 2 1/2"

30-176 Single End Vaginal Specula, Medium 1 1/4" x 3"

30-178 Single End Vaginal Specula, Large 1 1/2" x 3 1/2"

30-210 Double End Vaginal Specula, Small 1" x 2 1/2" & 1 1/4" x 3"

30-215 Double End Vaginal Specula, Medium 1 1/8" x 3" & 1 1/4" x 3 1/2"

30-220 Double End Vaginal Specula, Large 1 3/8" x 3 1/2" & 1 1/2" x 4"

**MILTEX TOWNSEND CERVICAL BIOPSY PUNCH FORCEPS**

4.2 x 2.3 x 1mm bite with locking mechanism.

30-1445 7 1/2" Shaft Cervical Biopsy Punch Forceps


MILTEX UMBILICAL SCISSORS

30-2575 4 1/8" Umbilical Scissors, American Pattern

**PRECISION DYNAMICS SECURLINE™ UMBILICAL CORD CLAMP CLIPPER**

Designed to safely cut the clamp off umbilical cord. Single-use eliminates concerns of cross-contamination. Each clipper is individually packaged in a label-sealed glassine pouch.

3540-00-PDO Umbilical Cord Clamp Clipper, Clean, 6/bx

**VANTAGE GRAVES VAGINAL SPECULUM**

V930-10 Small Vaginal Speculum


V930-15 Medium Vaginal Speculum

V930-15A Screw Nuts For Speculum

V930-15X Center Nut For Speculum

V930-20 Large Vaginal Speculum

V930-22 X-Large Vaginal Speculum

**MILTEX KEVORKIAN-YOUNGE BIOPSY FORCEPS**

Shaft interchangeable on pistol grip spring handle. 3.5 mm x 8 mm bite, teeth on lower jaw. Complete with pistol grip handle. Rotating punch.

- 30-1490 9" Shaft Biopsy Forceps
30-1491 Biopsy Forceps Shaft Only

**MILTEX KEVORKIAN-YOUNGE UTERINE BIOPSY FORCEPS**

3.5 mm x 8 mm bite, teeth on lower jaw.

- 30-1482 9 1/2" Shaft Uterine Biopsy Forceps

**MILTEX LATERAL VAGINAL RETRACTOR**

Thumb ratchet retracts up to 2 1/4" for extra wide viewing area. Fits inside any vaginal speculum.

- 30-1355 8 1/4" Lateral Vaginal Retractor

**MILTEX LEEP DRESSING FORCEPS**

For holding and removing tissue.

- 31-14 8" Dressing Forceps, Serrated w/coated tips
31-16 10" Dressing Forceps, Serrated w/coated tips

**MILTEX LEEP EMMETT UTERINE TENACULUM HOOKS**

Used during excision, these hooks have sharp points and are 9" long. All tips are exposed for better tissue grasping.

- 31-951 Style 1 Tenaculum Hook, Half Curved, Small
31-952 Style 2 Tenaculum Hook, Acute Angle
31-953 Style 3 Tenaculum Hook, Right Angle
31-954 Style 4 Tenaculum Hook, Half Curved, Large

**MILTEX LEEP IRIS HOOK**

Used for holding tissue during excision. It is smaller than the Emmett Hook and is slightly angled.

- 31-949 10" Iris Hook w/exposed tip & delicate hook

**MILTEX LEEP LATERAL VAGINAL RETRACTORS**

Has an extra-long thumb ratchet and retracts up to a full 2 1/4" for an extra wide viewing area while providing excellent stability to the outer vaginal wall.

- 31-1355 8 1/4" Lateral Vaginal Retractor

**MILTEX LEEP PEDERSON SPECULUM**

Has narrow blades for patients with a smaller vaginal opening. The integrated smoke evacuation tube provides excellent smoke reduction.

- 31-551 Medium Pederson Speculum
31-601 Large Pederson Speculum
31-611 X-Large Pederson Speculum

**MILTEX LEEP REDUNDANT WALL VAGINAL RETRACTOR**

Provides the widest lateral opening for the vaginal canal.

- 31-70 7 3/4" Medium Redundant Wall Vaginal Retractor
31-71 7 3/4" Large Redundant Wall Vaginal Retractor

**MILTEX LEEP SCHROEDER-BRAUM TENACULUM FORCEPS**

Maximum stabilization and safety with fully coated tips.

- 31-965 9 3/4" Tenaculum Forceps, Straight w/coated tips

**MILTEX LEEP TISSUE FORCEPS**

1x2 teeth to grasp tissue and hold for removal.

- 31-48 8" Tissue Forceps w/exposed tips
31-49 8" Tissue Forceps w/coated tips
31-50 10" Tissue Forceps w/exposed tips
31-51 10" Tissue Forceps w/coated tips



MILTEX LEEP VIEW-MORE SPECULUM


The wider 1 1/2" yoke area on this speculum provides a better view of the cervix and allows for easier access with other instruments. The integrated smoke evacuation tube provides excellent smoke reduction.

31-66	Medium View-More Speculum
31-67	Large View-More Speculum

MILTEX LEKSELL RONGEUR

Standard curve.

25-510	9 1/2" Leksell Rongeur, 3mm Bite
25-511	9 1/2" Leksell Rongeur, 8mm Bite



MILTEX NOVAK ENDOMETRIAL SUCTION BIOPSY CURETTES

9-3/4", Luer hub.

30-1390	4 mm Dia. Biopsy Curette, 12 Fr
30-1391	2 mm Dia. Biopsy Curette, Extra Delicate
30-1392	3 mm Dia. Biopsy Curette, Delicate



MILTEX PRATT UTERINE DILATORS

11-1/2". Solid. Double ended.

30-560-1315	13-15 Fr Pratt Uterine Dilator
30-560-1719	17-19 Fr Pratt Uterine Dilator
30-560-2123	21-23 Fr Pratt Uterine Dilator
30-560-2527	25-27 Fr Pratt Uterine Dilator
30-560-2931	29-31 Fr Pratt Uterine Dilator
30-560-3335	33-35 Fr Pratt Uterine Dilator
30-560-3739	37-39 Fr Pratt Uterine Dilator
30-560-4143	41-43 Fr Pratt Uterine Dilator
30-561	Set of 8 Pratt Uterine Dilator, sizes: 13-15 Fr to 41-43 Fr in a Khaki Roll

Accessories

30-561A	Khaki Roll Only For Set of 8
---------	------------------------------

MILTEX RANDALL ENDOMETRIAL SUCTION BIOPSY CURETTE

Luer hub. 4mm diameter.

30-1395	9 3/4" Endometrial Suction Biopsy Curette
---------	---



MILTEX SCHROEDER UTERINE VULSELLUM FORCEPS

30-1000	10" Uterine Vulsellum Forceps, Straight
30-1005	10" Uterine Vulsellum Forceps, Curved On Flat
30-1010	10" Uterine Vulsellum Forceps, Curved Sideways



MILTEX SCHROEDER-BRAUN UTERINE TENACULUM FORCEPS

30-965	9 3/4" Uterine Tenaculum Forceps, Straight
30-965ATR	10" Braun Tenaculum Forceps, Straight/Squared Jaw, Non-Overlapping Atraumatic Tips
30-966ATR	10" Schroeder Tenaculum Forceps, Rounded Jaw, Non-Overlapping Atraumatic Tips


MILTEX SCHUBERT UTERINE BIOPSY FORCEPS

30-1420	8 1/2" Uterine Biospy Forceps, Straight
30-1425	8 1/2" Uterine Biospy Forceps, Curved
30-1430	11" Uterine Biospy Forceps, Angled Shanks
30-1435	11" Uterine Biospy Forceps, S-Shape, Triangular Bite

MILTEX SIMS UTERINE BIOPSY CURETTE

Extra small sharp loop. 2.5 x 5mm outside measurement.

30-1380	10" Uterine Biopsy Curette
---------	----------------------------



MILTEX SIMS UTERINE CURETTES


11". Sharp blades on malleable shanks.

30-1205-0	Size 0 Sims Uterine Curette
30-1205-00	Size 00 Sims Uterine Curette
30-1205-1	Size 1 Sims Uterine Curette
30-1205-2	Size 2 Sims Uterine Curette
30-1205-3	Size 3 Sims Uterine Curette
30-1205-4	Size 4 Sims Uterine Curette
30-1205-5	Size 5 Sims Uterine Curette
30-1205-6	Size 6 Sims Uterine Curette

**MILTEX SIMS UTERINE SOUND**

13". Malleable. Silver plated.

30-650	Sims Uterine Sound Graduated in Centimeters
30-650SS	13" Sims Uterine Sound
30-655	Sims Uterine Sound Graduated in Inches

**MILTEX THOMAS UTERINE CURETTES**

11". Blunt blades on malleable shanks.


30-1225-1	Size 1 Thomas Uterine Curette
30-1225-2	Size 2 Thomas Uterine Curette
30-1225-3	Size 3 Thomas Uterine Curette
30-1225-4	Size 4 Thomas Uterine Curette
30-1225-5	Size 5 Thomas Uterine Curette
30-1225-6	Size 6 Thomas Uterine Curette

VANTAGE PEDERSON VAGINAL SPECULUM

V930-50	Small Vaginal Speculum
V930-55	Medium Vaginal Speculum
V930-60	Large Vaginal Speculum
V930-61	X-Large Vaginal Speculum

**MILTEX TISCHLER BIOPSY FORCEPS**


30-1494	3 x 7mm Bite Biopsy Forceps, Complete w/pistol grip handle
30-1495	Shaft Only For Biopsy Forceps

**MILTEX TISCHLER CERVICAL BIOPSY PUNCH FORCEPS**

30-1442	9 3/4" Shaft (3 x 7 mm bite) Biopsy Punch Forceps
30-1443	Burke Modified Baby Tischler Biopsy Punch Forceps (4.5 x 3 x 1.5 mm bite)



30-1486



30-1487

MILTEX WITTNER UTERINE BIOPSY FORCEPS

8 1/2" shaft with teeth on lower jaw.

30-1486	Standard Model Uterine Biopsy Forceps w/angled tip
30-1487	Straight Tip Uterine Biopsy Forceps

MILTEX-KERRISON CERVICAL RONGEURS

7" shaft. Improved handle.

26-388	2mm Cervical Rongeur, 40° Forward Angle
26-389	3mm Cervical Rongeur, 40° Forward Angle
26-390	5mm Cervical Rongeur, 40° Forward Angle
26-392	5mm Up-Bite Cervical Rongeur
26-392A	2mm Up-Bite Cervical Rongeur
26-392B	3mm Up-Bite Cervical Rongeur
26-393	5mm Down-Bite Cervical Rongeur
26-393A	2mm Down-Bite Cervical Rongeur
26-393B	3mm Down-Bite Cervical Rongeur

OWEN MUMFORD AMIELLE CARE VAGINAL DILATORS

Following radiotherapy to the vagina, cervix or pelvis, the woman may experience some internal scarring (adhesions) to the treated area, which can make the vagina narrower, drier and less elastic. Preventing the formation of scar tissue within the vagina helps ensure that it remains supple so that sexual intercourse and vaginal examinations can be performed without discomfort or pain. Designed in conjunction with the UK National Forum of Gynecological Oncology Nurses, Amielle Care vaginal dilators help to reduce the risk of scarring and offer an effective way for women to take an active role in their treatment.

SM2150	Vaginal Dilator
--------	-----------------

OWEN MUMFORD AMIELLE COMFORT VAGINAL DILATORS

Vaginismus is a condition where the vaginal muscles go into spasm, causing the vagina to feel tight and often making sexual intercourse impossible. This may be due to a number of reasons including fear and pain of penetration, gynecological surgery or trauma associated with childbirth. Dyspareunia can also inhibit sexual intercourse as it may cause pain and discomfort upon penetration. This is often associated with vaginal dryness or infection. Designed in conjunction with psychosexual therapists, Amielle Comfort offers an effective way for women to take an active role in the treatment of these and other conditions.

SM2100	Vaginal Dilator
--------	-----------------

RIESTER VAGINAL SPECULA

The gynecologist can select from two vaginal specula models and different sizes. Regardless which quality instrument is chosen, all specula excel through their particular reliability and durability. The speculight provides for the right illumination and is the optimal addition to all vaginal specula. The graeve vaginal speculum has an external polished finish and the cusco vaginal speculum has an external matt finish. No reflections thanks to matt inner surface. Made of stainless steel type 18/8. Can be sterilized up to 134°C in an autoclave. Each speculum is supplied in a cardboard box. The speculight has a battery handle C, chrome-plated. Uses two C type batteries. 85cm cable with clip and 2.5V standard illumination (vacuum).

5201	Vaginal Specula Graeve, 2.95 x .78"
5202	Vaginal Specula Graeve, 3.74 x 1.37"
5203	Vaginal Specula Graeve, 4.52 x 1.37"
5230	Vaginal Specula Cusco, 2.95 x .66"
5231	Vaginal Specula Cusco, 2.95 x .25"
5232	Vaginal Specula Cusco, 3.34 x 1.37"
5233	Vaginal Specula Cusco, 3.93 x 1.37"
5240	Lighting Attachment For Vaginal Specula, Speculight®

DYNAREX DISPOSABLE VAGINAL SPECULA - GRAVES STYLE

4914	Vaginal Specula, Disposable, Graves Style w/light source option, Small, 25/bx, 4 bx/cs
4915	Vaginal Specula, Disposable, Graves Style w/light source option, Medium, 25/bx, 4 bx/cs
4916	Vaginal Specula, Disposable, Graves Style w/light source option, Large, 25/bx, 4 bx/cs

PRO ADVANTAGE® KNIFE HANDLE

N407500 #3 Knife Handle

MILTEX COTTLE SEPTUM ELEVATOR

Double end. Graduated in cm.


21-60 9" Cottle Septum Elevator



TECH-MED FINGER RING CUTTER

Emergency tool designed to quickly remove ring from finger. Chrome plated handle with stainless steel blade.

8120	Finger Ring Cutter, Standard
8122	Finger Ring Cutter, Heavy Duty



MILTEX DIX FOREIGN BODY NEEDLE & SPUD


18-374	Spud, 4 3/4", 1.7 x 15mm Blade
18-400	Needle & Spud, Protected In Sliding Handle
18-402	Needle & Spud, Protected In Reversible Screw Handle
18-402A	Needle w/screw socket
18-402B	Spud w/screw socket
18-404	Needle & Gouge, Protected In Reversible Screw Handle



MILTEX EYE MAGNET & FOREIGN BODY LOOP WITH MAGNIFIER

Adjustable magnifier attached. Reversible plastic handle, plated metal parts. Plastic case.

18-398M Eye Magnet/Loop w/magnifier



MILTEX GERZOG Mallet

8 oz head lead filled with stainless jacket.

19-790 7 1/2" Gerzog Mallet



MILTEX KNAPP EYE SPECULUM

Blades 15 x 6mm.

18-4 3 1/4" Knapp Eye Speculum



MILTEX JACKSON TRACHEA TUBES WITH ROTATING LOCK

Improved model with rotating lock. Stainless steel.

23-1139SS	Size 00 Trachea Tube
23-1140SS	Size 0 Trachea Tube
23-1142SS	Size 1 Trachea Tube
23-1144SS	Size 2 Trachea Tube
23-1146SS	Size 3 Trachea Tube
23-1148SS	Size 4 Trachea Tube
23-1150SS	Size 5 Trachea Tube
23-1152SS	Size 6 Trachea Tube
23-1154SS	Size 7 Trachea Tube
23-1156SS	Size 8 Trachea Tube
23-1158SS	Size 9 Trachea Tube
23-1159SS	Size 10 Trachea Tube



MILTEX MEAD Mallet

19-796 7 1/2" Mead Mallet, Head 8 oz w/2 replaceable nylon caps

Accessories

19-797 Nylon Cap Replacement For 19-796

MILTEX POOLE SUCTION TUBE

10-310	23 Fr Poole Suction Tube, Curved
10-312	30 Fr Poole Suction Tube, Straight

General Surgery

MILTEX RANEY SCALP CLIP APPLYING FORCEPS

- 26-590 6 1/4" Scalp Clip Applying Forceps
26-591 Reusable Scale Clips, Black Finish, 1 dz/pkg



MILTEX SPLINTER LIBERATOR

Trocar point in reversible handle.

- 6-400 5" Splinter Liberator

MILTEX SPRATT (BRUN) MASTOID CURETTES

Oval cups. 6 1/2".

- 19-700 Size 5/0 Mastoid Curette
19-702 Size 4/0 Mastoid Curette
19-704 Size 3/0 Mastoid Curette
19-706 Size 2/0 Mastoid Curette
19-708 Size 0 Mastoid Curette
19-710 Size 1 Mastoid Curette
19-712 Size 2 Mastoid Curette
19-714 Size 3 Mastoid Curette
19-716 Size 4 Mastoid Curette
19-718 Size 5 Mastoid Curette
19-720 Size 6 Mastoid Curette

MILTEX SPRING WIRE RETRACTOR

- 11-19 3 1/4" Wire Retractor w/blades 10mm W x 13mm D

MILTEX WILDER LACRIMAL DILATORS

- 18-694 4" Wilder Lacrimal Dilator, Fine Taper
18-696 4" Wilder Lacrimal Dilator, Medium Taper
18-698 4" Wilder Lacrimal Dilator, Heavy Taper

MILTEX SWISS PATTERN OSTEOTOMES

5".

- 27-599 2mm Osteotome
27-600 4mm Osteotome
27-601 6mm Osteotome
27-601A 8mm Osteotome
27-602 10mm Osteotome
27-602A 12mm Osteotome
27-603 14mm Osteotome
27-604 18mm Osteotome
27-605 20mm Osteotome
27-610 Set of 9 Osteotomes 2mm To 20mm w/case

Accessories

- 27-609 Case Only, 5 1/2" x 6 1/2" x 1 1/2", Fitted To Hold 9 Osteotomes 2mm To 20mm



PERSONNA DERMABLADE® SHAVE BIOPSY INSTRUMENT


Designed specifically for shave biopsies. Special sure-grip side teeth facilitate a three digit grip to provide stability, control and direction, prevent slipping, and make it easier to pick up. It is individually packaged to guarantee sterility and protection. Plus, it utilizes Personna's MicroCoated® surgical blade to provide smooth, clean excisions with less tissue trauma.

- 72-0001 Blade, Sterile, 50/ctn, 8 ctn/cs

MILTEX WILLIAMS LACRIMAL PROBES

5" sterling.

- 18-722 Size 0000-000 Lacrimal Probe
18-724 Size 00-0 Lacrimal Probe
18-726 Size 1-2 Lacrimal Probe
18-728 Size 3-4 Lacrimal Probe
18-730 Size 5-6 Lacrimal Probe
18-732 Size 7-8 Lacrimal Probe



PRECISION DYNAMICS SECURLINE™ AMNIOTIC MEMBRANE PERFORATOR

Designed to snap-not poke or prick membrane. Overall length of 10 1/2" provides ample reach for practitioner. Perforator is sterile. 100 individual pouches are placed in a convenient dispenser box that contains a level window that allows one to check product levels at a glance.

- 3550-00-PDG Amniotic Membrane Perforator, Sterile, 100/bx

COVIDIEN/US SURGICAL SINGLE USE LAPAROSCOPIC HAND INSTRUMENTS

Curved stainless steel jaws with atraumatic serrations measure approximately 19mm in length, with a maximum jaw opening of 21mm.

- 173019 Retractor, Laparoscopic, Single Use, 5mm, 3/pk, 4 pk/cs

COVIDIEN/US SURGICAL SURGIWANT™ II 5MM W/O CAUTERY

Without cautery suction/irrigation device consists of a 5mm diameter shaft and a 34cm shaft length with (9) orifices for pool suction and irrigation. Also has a sealed portal to accommodate up to a 2.4mm laser probe. Available with non-removable spatula tip or L-hook top. Can accommodate the gravity feed system, electronic pump system and the nitrogen pump system.

- 178081 Suction/Irrigation Device w/o cautery, 5mm, 3/bx

COVIDIEN/US SURGICAL VERSAPORT™ REPOSABLE TITANIUM CANNULAS

Can be used with the Versaport™ V2 and Visiport™ obturators. Attaches and detaches from a single use seal system.

- 177791 Cannula, Repposable Titanium, 12mm

COVIDIEN/US SURGICAL VERSASTEP™ REPOSABLE CANNULAS

Composed of non conductive material for use with the radially expandable sleeve to minimize capacitive coupling, and can be re-sterilized.

- 177798 Cannula, Repposable, 12mm
177800 Cannula, Repposable, 11mm



MILTEX BIOPBLADE

- 33-100 Biopblade, Sterile, 50/bx

MILTEX SIMPSON UTERINE SOUND

Graduated in inches. Malleable, Stainless steel.

30-655SS 13" Simpson Uterine Sound

MILTEX RING INSTRUMENT HOLDERS


Latch has snap-lock feature to facilitate loading and removing of instruments.

3-70 8" Ring Instrument Holder
 3-71 10" Ring Instrument Holder
 3-72 12" Ring Instrument Holder
 3-74 16" Ring Instrument Holder

MILTEX STAINLESS WIRE BRUSH


For cleaning rasps and files. Autoclavable.

21-379 6 1/4" Wire Brush

**MILTEX TIP IT® INSTRUMENT GUARDS**

Protect instrument tips and edges from damage during handling and sterilization. PVC medical grade plastic for use in medical applications. Latex-free material that is non-allergenic. Large holes (vented styles) to reduce vacuum and moisture build up for easier removal from instrument. Steam and gas permeable to allow for effective sterilization of instruments. Four styles in nine colors and sizes for a variety of choices.

3-2501V Size 1, 1/16" x 3/4" Solid Vented Instrument Guard, White, Non-Sterile, 50/pkg
 3-2502V Size 2, 5/64" x 3/4" Solid Vented Instrument Guard, Blue, Non-Sterile, 50/pkg
 3-2503V Size 3, 7/64" x 3/4" Solid Vented Instrument Guard, Green, Non-Sterile, 50/pkg
 3-2504V Size 4, 1/8" x 1" Solid Vented Instrument Guard, Red, Non-Sterile, 50/pkg
 3-2505V Size 5, 3/16" x 1" Solid Vented Instrument Guard, Yellow, Non-Sterile, 50/pkg
 3-2506V Size 6, 13/32" x 3/4" Solid Vented Instrument Guard, Maroon, Non-Sterile, 50/pkg
 3-2507V Size 7, 1/16" x 3/8" x 1" Solid Vented Instrument Guard, Orange, Non-Sterile, 50/pkg
 3-2508V Size 8, 5/64" x 5/8" x 1" Solid Vented Instrument Guard, Brown, Non-Sterile, 50/pkg
 3-2509V Size 9, 1/8" x 1" x 1" Solid Vented Instrument Guard, Black, Non-Sterile, 50/pkg
 3-2510V Size 10, Assorted Solid Vented Instrument Guards (equal number of sizes plus 5), Non-Sterile, 50/pkg

**MILTEX TIP-IT™ INSTRUMENT GUARDS**

Used to protect instrument tips and edges from damage during instrument handling and sterilization. Non-sterile. Disposable. Steam and gas sterilizable. Radiopaque. Medical grade plastic. Latex free. Packaging: 50/package.

3-2501 Size 1 Instrument Guard, 1/16" x 3/4", White
 3-2502 Size 2 Instrument Guard, 5/64" x 3/4", Blue
 3-2503 Size 3 Instrument Guard, 7/64" x 3/4", Green
 3-2504 Size 4 Instrument Guard, 1/8" x 1", Red
 3-2505 Size 5 Instrument Guard, 3/16" x 1", Yellow
 3-2506 Size 6 Instrument Guard, 13/32" x 3/4", Maroon
 3-2507 Size 7 Instrument Guard, 1/16" x 3/8" x 1", Orange
 3-2508 Size 8 Instrument Guard, 5/64" x 5/8" x 1", Brown
 3-2509 Size 9 Instrument Guard, 1/8" x 1" x 1", Black
 3-2510 Assorted Instrument Guard

We offer:

- Great Service
- Great Pricing
- Great Products



MILTEX MAGNIFICATION LOUPE

Offers the clinician a unique design featuring the "Set and Forget™" locking mechanism. This feature simplifies the fitting process by allowing the clinician to set his or her interpupillary distance in one easy step. The "Set and Forget™" feature will automatically retain the individual setting and require to further adjustment, no matter how often the loupes are taken on and off.

- 1-5000 Magnifying Loupe, Red Frame, 2.5x Magnification, 10" x 13" Working Distance
- 1-5001 Magnifying Loupe, Red Frame, 2.5x Magnification, 12" x 15" Working Distance
- 1-5002 Magnifying Loupe, Red Frame, 2.5x Magnification, 14" x 17" Working Distance
- 1-5003 Magnifying Loupe, Red Frame, 2.5x Magnification, 16" x 20" Working Distance
- 1-5004 Magnifying Loupe, Red Frame, 3.0x Magnification, 10" x 13" Working Distance
- 1-5005 Magnifying Loupe, Red Frame, 3.0x Magnification, 13" x 16" Working Distance
- 1-5006 Magnifying Loupe, Red Frame, 3.0x Magnification, 16" x 21" Working Distance
- 1-5007 Magnifying Loupe, Red Frame, 3.0x Magnification, 18" x 22" Working Distance
- 1-5008 Magnifying Loupe, Red Frame, 3.5x Magnification, 13" x 16" Working Distance
- 1-5009 Magnifying Loupe, Red Frame, 3.5x Magnification, 16" x 20" Working Distance
- 1-5010 Magnifying Loupe, Red Frame, 3.5x Magnification, 18" x 22" Working Distance
- 1-5011 Magnifying Loupe, Black Frame, 2.5x Magnification, 10" x 13" Working Distance
- 1-5012 Magnifying Loupe, Black Frame, 2.5x Magnification, 12" x 15" Working Distance
- 1-5013 Magnifying Loupe, Black Frame, 2.5x Magnification, 14" x 17" Working Distance
- 1-5014 Magnifying Loupe, Black Frame, 2.5x Magnification, 16" x 20" Working Distance
- 1-5015 Magnifying Loupe, Black Frame, 3.0x Magnification, 10" x 13" Working Distance
- 1-5016 Magnifying Loupe, Black Frame, 3.0x Magnification, 13" x 16" Working Distance
- 1-5017 Magnifying Loupe, Black Frame, 3.0x Magnification, 16" x 21" Working Distance
- 1-5018 Magnifying Loupe, Black Frame, 3.0x Magnification, 18" x 22" Working Distance
- 1-5019 Magnifying Loupe, Black Frame, 3.5x Magnification, 13" x 16" Working Distance
- 1-5020 Magnifying Loupe, Black Frame, 3.5x Magnification, 16" x 20" Working Distance
- 1-5021 Magnifying Loupe, Black Frame, 3.5x Magnification, 18" x 22" Working Distance

MILTEX EKG CALIPER

- 1-3015 Economy Line Caliper w/protective black plastic cap, Graduated in mm, Nickel Plated

ADC NEUROLOGICAL HAMMERS

Four of the most popular hammers for diagnosis of neurological function.

- 3691 Buck Hammer, 7 1/4", Chrome Plated Handle w/2 sided rubber head, Handle Conceals "screw-in" Brush, Needle Contained Within The Head
- 3693 Taylor Hammer, 7 1/2", Chrome Handle w/triangular rubber head, Orange
- 3693BK Taylor Hammer, 7 1/2", Chrome Handle w/triangular rubber head, Black
- 3693DG Taylor Hammer, 7 1/2", Chrome Handle w/triangular rubber head, Dark Green
- 3693RB Taylor Hammer, 7 1/2", Chrome Handle w/triangular rubber head, Royal Blue
- 3695 Wartenberg Pinwheel, 7 1/2", Stainless Steel Handle w/textured grip, Rotating Spur
- 3697 Babinski Hammer, 8 1/2", Octagonal Stainless Steel Handle w/concealed needle, Rubber Head




1-200




1-202



1-210



1-212



1-220



1-222

MILTEX HAMMERS

- 1-200 Taylor Percussion Hammer, 7", Solid Handle, Stainless Steel
- 1-202 Taylor Percussion Hammer, 7 1/2", Solid Handle, Chrome
- 1-208 Taylor Percussion Hammer, 9", Hoop Handle, Chrome
- 1-210 Dejerine Percussion Hammer, 9 1/2", Chrome
- 1-211 Tromner Percussion Hammer, 9 3/4" w/large & small rubber tips, Hollow Tapered Handle w/sharp point, Chrome
- 1-212 Combo-5 Neuro Diagnostic Set Includes Buck Hammer, C-128 Tuning Fork, Pin Wheel, Brush & Needle
- 1-213 Combo-5 Neuro Diagnostic Set Includes Buck Hammer, C-256 Tuning Fork, Pin Wheel, Brush & Needle
- 1-214 Combo-4 Neuro Diagnostic Set Includes Buck Hammer, Pin Wheel, Brush & Needle
- 1-220 Buck Neurological Hammer, 7 3/4" w/brush & needle (screw into handle), Chrome
- 1-221-A Needle Only For Buck Hammer
- 1-221-SM Rubber End For Buck Hammer, Small
- 1-222 Rabiner Neurological Hammer, 9" w/brush & needle (screw into handle), Chrome




MILTEX WARTENBERG PINWHEEL

For neurological examination. Stainless steel.

- 1-226 Pinwheel

RIESTER PERCUSSION HAMMERS

5030	Taylor Stainless Steel 7.87"
5035	Buck w/needle & brush 7.08"
5038	Dejerine w/needle 8.26"
5040	Dejerine w/o needle 8.26"
5045	Troemner .40 lb 9.84"
5046	Troemner .50 lb 9.84"
5052	Babinsky 2 Threads 9.44"
5055	Berliner 7.87"
5058	Queens Plastic Handle 13.77"

**MILTEX TONGUE SEIZING FORCEPS**

2-50	Collin Tongue Seizing Forceps, 6 3/4" Jaws, 25mm Wide
2-60	Young Tongue Seizing Forceps, 6 1/2" w/soft rubber jaws

Accessories

2-62	Rubber Jaws Only For 2-60
------	---------------------------

**MILTEX BRUNS PLASTER SHEARS**

Serrated blade.

27-3070	9 1/4" Plaster Shears
---------	-----------------------

**MILTEX STILLE PLASTER SHEARS**

Serrated blades.

27-3060	11" Plaster Shears, Medium
27-3062	15" Plaster Shears, Standard

MILTEX U.S. PATTERN CAST KNIFE

Wood handle. 1 1/2" blade.

27-3002	5 1/2" Cast Knife
---------	-------------------

MILTEX WOLFE-BOEHLER CAST BREAKER

Jaws 15mm wide with deep serrations.

27-3080	10" Cast Breaker
---------	------------------

MILTEX BARR FISTULA PROBE

28-104	11 1/2" Overall Fistula Probe w/8" sterling shaft, Eye At Tip
--------	---

**MILTEX BODENHAMMER RECTAL SPECULUM**

Blades 3 1/2" long x 3/4" wide at distal end.

28-4	6" Rectal Speculum w/set screw
------	--------------------------------

MILTEX FANSLER (IVES) RECTAL SPECULUM

3 1/4" long x 1" largest diameter with plastic obturator (boilable). Chrome-plated.


28-12	Fansler (Ives) Rectal Speculum
28-14	Plastic Obturator Only

MILTEX HIRSCHMAN PROCTOSCOPE

28-48	5" Long x 15/16" O.D. Proctoscope
-------	-----------------------------------

**MILTEX HIRSCHMANN ANOSCOPES**

28-20	Small Anoscope, 2 1/4" Long, 9/16" I.D.
28-22	Medium Anoscope, 2 5/8" Long, 11/16" I.D.
28-24	Large Anoscope, 3 1/2" Long, 7/8" I.D.

**MILTEX IMPROVED BIOPSY FORCEPS WITH ROTATING SHAFTS**

In pistol-grip handle. Shafts easily disassemble from handle. Constructed without small assembly parts. All shafts and handles are interchangeable. Jaws 3.5 mm x 10 mm.

28-380	25 cm (10") Straight Complete Biopsy Forceps w/pistol grip handle
28-382	35 cm (14") Straight Complete Biopsy Forceps w/pistol grip handle
28-386	35 cm (14") Angled Complete Biopsy Forceps w/pistol grip handle
28-400	Pistol Grip Spring Handle Only
28-410	25 cm (10") Straight Shaft Biopsy Forceps Only
28-412	35 cm (14") Straight Shaft Biopsy Forceps Only
28-416	35 cm (14") Angled Shaft Biopsy Forceps Only



28-150



28-152



28-153



28-151

MILTEX MCGIVNEY HEMORRHOIDAL LIGATOR

Offset handle for better vision and comfort. 7" working length. Can be disassembled entirely for cleaning. Supplied with drum loading cone and 100 latex O-rings.


28-150 Hemorrhoidal Ligator, 7", w/drum loading cone & O-rings

28-152 Hemorrhoid Grasping Forceps, 7 1/2", Angled Shanks


28-153 Loading Cone For 28-150

Accessories

28-151 Latex O-Rings For 28-150, 100/pkg



28-300



28-303

MILTEX YEOMAN BIOPSY FORCEPS

28-294 9" Shaft Biopsy Forceps, 3.5 x 10mm Bite w/basket

28-296 13" Shaft Biopsy Forceps, 3.5 x 10mm Bite w/basket

28-298 15" Shaft Biopsy Forceps, 3.5 x 10mm Bite w/basket

28-300 10" Rotating Shaft Biopsy Forceps w/ring handle, Straight 4 x 8mm Bite w/basket

28-302 14" Rotating Shaft Biopsy Forceps w/ring handle, Straight 4 x 8mm Bite w/basket

28-303 14" Rotating Shaft Biopsy Forceps w/teeth & ring handle, Straight 4 x 8mm Bite w/basket


28-304 16" Rotating Shaft Biopsy Forceps w/ring handle, Straight 4 x 8mm Bite w/basket

28-306 10" Rotating Shaft Only, Straight 4 x 8mm Bite w/basket

28-308 14" Rotating Shaft Only, Straight 4 x 8mm Bite w/basket

28-309 14" Rotating Shaft Only w/teeth, Straight 4 x 8mm Bite w/basket

28-310 16" Rotating Shaft Only, Straight 4 x 8mm Bite w/basket



MILTEX VON GRAEFE HOOKS

18-450 5 1/2" Von Graefe Hook, Small, 8mm Deep

18-452 5 1/2" Von Graefe Hook, Medium, 10mm Deep

18-454 5 1/2" Von Graefe Hook, Large, 11mm Deep



MILTEX WIENER HOOKS

18-418 5" Wiener Hook, Sharp

18-420 5" Wiener Hook, Blunt



29-805


29-850

MILTEX NO SCALPEL VASECTOMY INSTRUMENTS

29-805 No Scalpel Vasectomy Ring Clamp, 3.5mm ID, 5 1/2"

29-850 No Scalpel Vasectomy Hemostat, Sharp Point, Smooth Jaws, Curved, 5 1/2"

29-855 No Scalpel Vasectomy Instrument Set



ADC FINGER RING CUTTER

Chrome handle. Safety lever slides protective plate between finger and ring. Large thumb screw operates a razor sharp serrated saw.

380 Finger Ring Cutter

381 Replacement Blade

BD BEAVER™ RING CUTTER & PARTS

374100 Ring Cutter

374103 Ring Cutter Thumbpiece

374105 Ring Cutter Hub

374106 Ring Cutter Washer

374150 Ring Cutter Saw Blade


374151 Ring Cutter Saw Blade, High Speed Steel

BUSSE STERILE FLOOR GRADE SINGLE USE INSTRUMENTS

Individually packaged in sterile pouches.

773 Nurses Scissors, 5 1/2", Sharp/Blunt, 50/cs

775 Kelly Hemostat, 5 1/2", Straight, 50/cs



MILTEX BANE MASTOID RONGEUR

19-808 7" Bane Mastoid Rongeur, Curved

MILTEX BARRAQUER BLADE BREAKER & HOLDER

18-170 4 1/2" Blade Breaker & Holder w/ratchet lock

MILTEX BEYER RONGEUR

Double action. Slightly curved jaws 3.5mm wide.

19-844 7" Beyer Rongeur

**MILTEX BONE WIRE GUIDE**

U.S. pattern.

27-110 8 3/4" Bone Wire Guide

MILTEX BOWMAN LACRIMAL PROBES


5" sterling.

- 18-740 Size 0000-000 Lacrimal Probe
- 18-742 Size 00-0 Director w/probe tip & tongue tie
- 18-744 Size 1-2 Lacrimal Probe
- 18-746 Size 3-4 Lacrimal Probe
- 18-748 Size 5-6 Lacrimal Probe
- 18-750 Size 7-8 Lacrimal Probe


MILTEX BOWMAN NEEDLE

With stop. Straight blade 3mm.

18-308 4 3/4" Bowman Needle

**MILTEX BRINKERHOFF RECTAL SPECULUM**

- 28-30 Small Rectal Speculum, 3 1/2"L, 1" Tapering to 1/2"
- 28-32 Medium Rectal Speculum, 4"L, 1" Tapering to 5/8"
- 28-34 Large Rectal Speculum, 4 3/4"L, 1 1/4" Tapering to 3/4"

**MILTEX BROWN APPLICATORS**


Copper.

- 19-180 6" Brown Applicator, Triangular Tip
- 19-182 6" Brown Applicator, Cross Serrated Tip


MILTEX BRUN BONE CURETTES

6 3/4". Round cups.


- 26-1650 Size 1 Bone Curette, 7mm
- 26-1652 Size 2 Bone Curette, 8mm
- 26-1654 Size 3 Bone Curette, 9mm
- 26-1656 Size 4 Bone Curette, 11mm
- 26-1658 Size 5 Bone Curette, 12mm
- 26-1660 Size 6 Bone Curette, 14mm




27-159TC




27-161TC




27-163TC



27-164TC



27-165TC



27-176TC

MILTEX C-N-S PIN & WIRE PULLER

Carb-N-Sert instruments with tungsten carbide insert blades.

- 27-159TC 7" Double Action Pin & Wire Puller, Cuts Up To 2mm (5/64") Soft Wire or 1.5mm (1/16") Hard Wire
- 27-161TC 15" Double Action Pin Cutter, Heavy Duty, Cuts Up To 4.8mm (3/16") Diameter
- 27-163TC 7" Double Action Pin & Wire Puller Flush Cutting Blades, Cuts Up To 2mm (5/64") Soft Wire or 1.6mm (1/16") Hard Wire
- 27-164TC 7" Double Action Pin & Wire Puller w/4mm wide jaws
- 27-165TC 7" Double Action Pin & Wire Puller w/6mm wide jaws
- 27-176TC 5 1/2" Pin Puller, Serrated Top Jaw Measures 1cm in Length, Grooved Cross-Serrated Opposing Jaw

**MILTEX CANNULATED PIN & WIRE CUTTER**

End cutting up to 3.2mm(5/64") soft wire or 1.6mm(1/16") hard wire.

27-158 7 1/2" Cannulated Pin & Wire Cutter

MILTEX CASTROVIEJO BLADE BREAKER & HOLDER

18-171 5 1/4" Blade Breaker & Holder w/push-pull lock

**MILTEX CASTROVIEJO CALIPER**

A combination measuring and marking caliper, graduated 0-20mm.

18-650 3 1/4" Castroviejo Caliper

MILTEX ALVIS FOREIGN BODY SPUD

Oval curette 1 x 3mm with fine point.

18-384 4 3/4" Alvis Foreign Body Spud

General Surgery

MILTEX CLAYMAN INTRAOCULAR LENS FORCEPS

15mm jaws with smooth tips.

18-1081 4 1/2" Intraocular Lens Forceps

MILTEX COLLIN DYNAMOMETER

1-228 Adult Size Dynamometer, 5" x 2 1/4"



MILTEX COTTLE CARTILAGE CRUSHER

16mm wide anvil with removing forceps.

21-300 Cartilage Crusher

MILTEX COTTLE CHISEL

Curved 5mm with depth graduations in cm.

21-212 7 1/4" Cottle Chisel



MILTEX COTTLE HOOKS

21-144 Double Hook, 5 1/2", 2 Sharp Prongs 12mm Wide

21-90 Skin Hook, 5 1/2", Small, Shallow Curve

21-92 Skin Hook, 5 1/2", Small, Deep Curve

21-94 Skin Hook, 5 1/2", Large, Deep Curve

MILTEX COTTLE OSTEOTOME

7". Rounded corners.

21-206-12 12mm Wide Cottle Osteotome

21-206-4 4mm Wide Cottle Osteotome

21-206-7 7mm Wide Cottle Osteotome

21-206-9 9mm Wide Cottle Osteotome

MILTEX COTTLE RASP

Backward cutting.

21-338 8" Cottle Rasp



MILTEX COTTLE SKIN ELEVATOR

Slightly curved blade 9mm wide.

21-72 7 1/2" Cottle Skin Elevator

MILTEX CUSHING "LITTLE JOKER" PERIOSTEAL ELEVATOR

26-1403 6" Periosteal Elevator

MILTEX CUSHING DURA HOOK

26-1080 6 1/2" Dura Hook, Blunt

MILTEX DANDY HOOK

26-1150 9" Dandy Hook, Straight

MILTEX DAVIEL LENS SPOON

18-472 5" Daviel Lens Spoon, 3mm Wide

MILTEX DAVIS FOREIGN BODY SPUD

18-372 4 3/4" Davis Foreign Body Spud, 0.8mm Blade



MILTEX DAY HOOKS

19-220 6 1/2" Hook & Spoon, Small

19-222 6 1/2" Hook & Spoon, Large

MILTEX DESCHAMPS LIGATURE CARRIER

9-140 8" Ligature Carrier, Blunt, Left

9-142 8" Ligature Carrier, Blunt, Right

9-146 8" Ligature Carrier, Sharp, Left

9-148 8" Ligature Carrier, Sharp, Right

MILTEX DIAMOND PIN & WIRE CUTTER

Cuts up to 2mm (5/64") diameter hard wire. Stainless.

27-160 6 1/2" Diamond Pin & Wire Cutter

MILTEX DOYEN-JANSEN MOUTH GAG


With grip lock.

2-150 5 1/2" Doyen-Jansen Mouth Gag

MILTEX ELLIS FOREIGN BODY SPUD

1.2mm blade. Curved.


18-380 4 1/2" Ellis Foreign Body Spud

**MILTEX EXCAVATOR CURETTES**

40-57/0	1mm Dia. Excavator Curette, Single End 5" w/o hole
40-57/1	1.5mm Dia. Excavator Curette, Single End 5" w/o hole
40-57/2	2mm Dia. Excavator Curette, Single End 5" w/o hole
40-57/3	2.5mm Dia. Excavator Curette, Single End 5" w/o hole
40-58/0	1mm Dia. Excavator Curette, Single End 5" w/hole
40-58/1	1.5mm Dia. Excavator Curette, Single End 5" w/hole
40-58/1-2	Small 1.5mm & Medium 2mm Dia. Excavator Curette, Double End 5 1/2"
40-58/1-3	Small 1.5mm & Large 2.5mm Dia. Excavator Curette, Double End 5 1/2"
40-58/2	2mm Dia. Excavator Curette, Single End 5" w/hole
40-58/3	2.5mm Dia. Excavator Curette, Single End 5" w/hole

MILTEX FANSLER OPERATING SPECULUM

28-50	2 3/8" Long x 1 3/8" O.D. Operating Speculum, Slotted Tube
-------	--



19-200



19-202

MILTEX FARRELL APPLICATORS

19-200	6 1/2" Farrell Applicator, Triangular Tip
19-202	6 1/2" Farrell Applicator, Cross Serrated Tip

MILTEX FERRIS-SMITH PITUITARY CUP JAW RONGEURS

7" shank.

26-461	3 x 10mm Jaws, Straight Pituitary Rongeur
26-462	4 x 10mm Jaws, Straight Pituitary Rongeur
26-466	3 x 10mm Jaws, Curved-Up Pituitary Rongeur
26-474	6 x 12mm Jaws, Straight Pituitary Rongeur



33-140



33-142

MILTEX FINGER RING CUTTER

Cutter is chrome. Replacement blade is stainless steel.

33-140	Ring Cutter, 6 1/2" w/circular saw blade
--------	--

Accessories

33-142	Replacement Blade
--------	-------------------

**MILTEX FLEXIBLE STAINLESS RULER**

Graduated in 1/32" fractions and mm.

18-660	6" x 1/2" Stainless Ruler
--------	---------------------------

MILTEX FOMON PERIOSTEAL ELEVATOR

Slightly curved blade 4.5mm wide. Standard pattern.

21-68	6 1/4" Fomon Periosteal Elevator
-------	----------------------------------

**MILTEX FRAZIER DURA HOOKS**

26-1084	5" Dura Hooks, Sharp Point
26-1086	5" Dura Hooks, Blunt Point

MILTEX FRAZIER DURA SEPARATOR

Angled blade 1/2" long, from 1/8" to 1/16" wide.

26-1054	6 1/2" Dura Separator
---------	-----------------------

**MILTEX FREER ELEVATOR**

Double end sharp and blunt blades 5mm wide.

18-1968	7" Elevator
---------	-------------

MILTEX FREER SEPTUM KNIFE

20-200	Septum Knife, 6 1/4", "D" Model Blade
20-202	Submucous Knife, 7 1/2", Angled Blade

**MILTEX FREER SKIN HOOKS**

21-102	6" Skin Hooks, Sharp Double Prongs 2.5mm Wide
21-103	6" Skin Hooks, Blunt Double Prongs 2.5mm Wide

MILTEX FRENCH PATTERN LACRIMAL PROBES

5" bayonet shape.

18-710	Size 00 Lacrimal Probe
18-712	Size 0 Lacrimal Probe
18-714	Size 1 Lacrimal Probe
18-716	Size 2 Lacrimal Probe
18-718	Size 3 Lacrimal Probe
18-720	Size 4 Lacrimal Probe



26-119



26-132

MILTEX GIGLI SAWS & HANDLES

Standard twisted wire type.

- 26-119 T-Shaped Saw Handle, Turnable
- 26-120 T-Shape Saw Handle, Rigid w/snap-lock
- 26-122 T-Shape Saw Handle, Turnable w/snap-lock
- 26-124 Loop Grip Saw Handle, Rigid w/snap-lock
- 26-130 12" Gigli Saw
- 26-132 20" Gigli Saw
- 26-134 30" Gigli Saw
- 26-138 Saw Guide, 12", Flexible

MILTEX GILLIES (CONVERSE) SKIN HOOKS

- 21-88 7" Skin Hook, Large
- 21-89 7" Skin Hook, Small

MILTEX GRAHAM HOOK

- 26-1156 7" Graham Hook, Blunt, 6mm Deep

MILTEX GREEN HOOK

- 18-442 5" Green Hook, 9mm Deep

MILTEX GREEN SPATULA

One end 4 x 14mm, other end 5 x 16mm.

- 18-570 6" Green Spatula



MILTEX GROOVED DIRECTORS

- 10-69 4 1/2" Director w/tongue tie
- 10-70 5" Director w/tongue tie
- 10-72 5 1/2" Director w/tongue tie
- 10-74 6" Director w/tongue tie
- 10-76 8" Director w/tongue tie
- 10-79 4 1/2" Director w/probe tip & tongue tie
- 10-80 5" Director w/probe tip & tongue tie
- 10-82 5 1/2" Director w/probe tip & tongue tie
- 10-84 6" Director w/probe tip & tongue tie



MILTEX GROSS HOOK & SPOON

- 19-228 5" Gross Hook & Spoon



MILTEX GUTHRIE (IRIS) HOOK

Very delicate, 1.5 mm wide.

- 18-416 5" Iris Hook, 2 Sharp Prongs



MILTEX HARTMAN MASTOID RONGEUR

- 19-804 5 3/4" Hartman Mastoid Rongeur, Strong Curve

MILTEX HIBBS CHISELS

9".

- 27-418 1/4" Chisel, Straight
- 27-420 3/8" Chisel, Straight
- 27-422 1/2" Chisel, Straight
- 27-424 5/8" Chisel, Straight
- 27-426 3/4" Chisel, Straight
- 27-426C 3/4" Chisel, Curved
- 27-427 7/8" Chisel, Straight
- 27-428 1" Chisel, Straight
- 27-432 1 1/2" Chisel, Straight

MILTEX HIBBS GOUGES

9".

- 27-458 1/4" Gouge, Straight
- 27-458C 1/4" Gouge, Curved
- 27-460 3/8" Gouge, Straight
- 27-460C 3/8" Gouge, Curved
- 27-462 1/2" Gouge, Straight
- 27-464C 5/8" Gouge, Curved
- 27-466 3/4" Gouge, Straight
- 27-466C 3/4" Gouge, Curved
- 27-468 1" Gouge, Straight

MILTEX HIBBS OSTEOTOMES

9".

- 27-438 1/4" Osteotome, Straight
- 27-438C 1/4" Osteotome, Curved
- 27-440 3/8" Osteotome, Straight
- 27-440C 3/8" Osteotome, Curved
- 27-442 1/2" Osteotome, Straight
- 27-442C 1/2" Osteotome, Curved
- 27-444 5/8" Osteotome, Straight
- 27-444C 5/8" Osteotome, Curved
- 27-446 3/4" Osteotome, Straight
- 27-446C 3/4" Osteotome, Curved
- 27-447 7/8" Osteotome, Straight
- 27-447C 7/8" Osteotome, Curved
- 27-448 1" Osteotome, Straight
- 27-448C 1" Osteotome, Curved
- 27-449 1 1/8" Osteotome, Straight
- 27-449C 1 1/8" Osteotome, Curved
- 27-450 1 1/4" Osteotome, Straight
- 27-450C 1 1/4" Osteotome, Curved
- 27-452 1 1/2" Osteotome, Straight
- 27-452C 1 1/2" Osteotome, Curved

MILTEX HIBBS SPRATT SPINAL FUSION CURETTES

9". Oval cup.

26-1627	Size 3/0 Spinal Fusion Curette
26-1628	Size 2/0 Spinal Fusion Curette
26-1629	Size 0 Spinal Fusion Curette
26-1630	Size 1 Spinal Fusion Curette
26-1632	Size 2 Spinal Fusion Curette
26-1634	Size 3 Spinal Fusion Curette
26-1636	Size 4 Spinal Fusion Curette
26-1638	Size 5 Spinal Fusion Curette
26-1640	Size 6 Spinal Fusion Curette

**MILTEX HOKE OSTEOTOMES**

5 1/2" straight.

27-331	1/8" Hoke Osteotome
27-332	3/16" Hoke Osteotome
27-333	1/4" Hoke Osteotome
27-334	5/16" Hoke Osteotome
27-335	3/8" Hoke Osteotome
27-336	1/2" Hoke Osteotome

MILTEX I.U.D. EXTRACTOR HOOK

30-940 10 1/4" I.U.D. Extractor Hook

MILTEX INSTRUMENT RACK

3-084107 Series 7, 7 Instrument Single Rack

MILTEX IOWA TRUMPET PUDENDAL NEEDLE GUIDE

30-2632 5 1/2" Pudendal Needle Guide

**MILTEX JACKSON TRACHEA TUBES WITH WING LOCK**

Original model.

23-1199SS	00 Regular Length Trachea Tube
23-1200SS	0 Regular Length Trachea Tube
23-1202SS	1 Regular Length Trachea Tube
23-1204SS	2 Regular Length Trachea Tube
23-1206SS	3 Regular Length Trachea Tube
23-1208SS	4 Regular Length Trachea Tube
23-1210SS	5 Regular Length Trachea Tube
23-1212SS	6 Regular Length Trachea Tube
23-1214SS	7 Regular Length Trachea Tube
23-1216SS	8 Regular Length Trachea Tube
23-1218SS	9 Regular Length Trachea Tube
23-1219SS	10 Regular Length Trachea Tube
23-1243SS	2 Short Length Trachea Tube
23-1244SS	3 Short Length Trachea Tube
23-1245SS	4 Short Length Trachea Tube
23-1246SS	5 Short Length Trachea Tube
23-1247SS	6 Short Length Trachea Tube
23-1248SS	7 Short Length Trachea Tube
23-1249SS	8 Short Length Trachea Tube
23-1250SS	9 Short Length Trachea Tube
23-1251SS	10 Short Length Trachea Tube

**MILTEX JAEGER LID PLATE**

Blades 20mm and 23mm.

18-120 4" Jaeger Lid Plate

**MILTEX JAMESON CALIPER**

Graduated in inches and mm with thin tips. Chrome.

18-658 3 3/4" Jameson Caliper

MILTEX JAMESON HOOK

18-446 5" Jameson Hook, 10mm Deep

MILTEX JENNINGS MOUTH GAG

2-160	6" Jennings Mouth Gag, Adult
2-162	4 3/4" Jennings Mouth Gag, Child
2-164	3 3/4" Jennings Mouth Gag, Infant

**MILTEX JOSEPH BUTTON END KNIFE**

21-20	6" Button End Knife, Straight Blade 4 x 30mm
21-40	6 1/4" Button End Knife, Angled Blade 4 x 12mm



21-153



21-154

MILTEX JOSEPH HOOKS

21-153	Single Hook, 6 1/4", One Sharp Prong
21-154	Double Hook, 6 1/4", 2 Sharp Prongs 2mm
21-156	Double Hook, 6 1/4", 2 Sharp Prongs 5mm
21-158	Double Hook, 6 1/4", 2 Sharp Prongs 7mm
21-160	Double Hook, 6 1/4", 2 Sharp Prongs 10mm

MILTEX JOSEPH NASAL SAW

21-280	7 1/2" Joseph Nasal Saw, Straight
21-282	7 1/2" Joseph Nasal Saw, Bayonet, Right
21-284	7 1/2" Joseph Nasal Saw, Bayonet, Left
21-288	Joseph Nasal Saw, Angled Right
21-290	Joseph Nasal Saw, Angled Left

MILTEX JOSEPH PERIOSTEAL ELEVATOR

Blade 6.5mm wide. Strong curve.

21-64 7" Joseph Periosteal Elevator

MILTEX JOSEPH RASP

Fine cross serrations.

21-330 6 1/4" Joseph Rasp



MILTEX KEVORKIAN-YOUNGE ENDOCERVICAL BIOPSY CURETTES

Loop size is 3 mm x 12 mm.

- 30-1382 12" Biopsy Curette, w/basket
30-1383 12" Biopsy Curette, w/o basket

MILTEX KEY PERIOSTEAL ELEVATORS

- 27-820 7" Periosteal Elevator, 1/4" Wide
27-821 7 1/2" Periosteal Elevator, 3/8" Wide
27-822 7 1/2" Periosteal Elevator, 1/2" Wide
27-824 8" Periosteal Elevator, 3/4" Wide
27-826 8 1/4" Periosteal Elevator, 1" Wide



MILTEX KLEINERT-KUTZ HOOKS

- 21-85 5" Kleinert-Kutz Hook, Small, 3mm Dia.
21-86 5" Kleinert-Kutz Hook, Medium, 5mm Dia.
21-87 5" Kleinert-Kutz Hook, Large, 7mm Dia.



27-490



27-504

MILTEX LAMBOTTE OSTEOTOMES

9"

- 27-488 1/4" Osteotome, Straight
27-490 1/2" Osteotome, Straight
27-492 3/4" Osteotome, Straight
27-494 1" Osteotome, Straight
27-496 1 1/4" Osteotome, Straight
27-498 1 1/2" Osteotome, Straight
27-499 1/4" Osteotome, Curved
27-500 1/2" Osteotome, Curved
27-502 3/4" Osteotome, Curved
27-504 1" Osteotome, Curved
27-506 1 1/4" Osteotome, Curved
27-508 1 1/2" Osteotome, Curved




MILTEX LANGENBECK PERIOSTEAL ELEVATORS

- 26-1410 7 1/2" Periosteal Elevator, Narrow Model, 7mm Wide Blunt Blade
26-1412 7 1/2" Periosteal Elevator, Wide Model, 17mm Wide Sharp Blade



MILTEX SKIN HOOK

- 21-152 4 3/4" Skin Hook, Sharp Prong 2mm Dia.



MILTEX STAINLESS STEELWARE AND SURGICAL INSTRUMENT STRINGERS

The Miltex® Stainless Steel Ware product line offers a wide selection of bowls, basins, cups, jars, instrument trays, and covers. The product is manufactured using 304 stainless steel and measures 0.6mm in thickness. The stainless steel Miltex® Surgical Instrument Stringers are available in various lengths and the most popular shapes; closed and open with lock - side and center. All closed styles offer a secure locking mechanism and smooth edges and are designed for effortless opening and closing.

- 3-75 Stringer w/lock, 2-1/2" x 6"
3-76 Stringer w/lock, 2-1/2" x 8"
3-77 Stringer w/lock, 2-1/2" x 10"
3-78 Stringer w/lock, 2-1/2" x 12"
3-79 Stringer w/lock, 2-1/2" x 16"
3-80 Stringer w/lock, 2-1/2" x 18"

MILTEX STAINLESS STEELWARE AND SURGICAL INSTRUMENT STRINGERS
(Continued)

3-81	Stringer w/lock, 5" x 8"
3-82	Stringer w/lock, 5" x 10"
3-83	Stringer w/lock, 5" x 12"
3-84	Stringer w/lock, 5" x 16"
3-85	Stringer w/lock, 5" x 18"
3-86	Stringer, Open, 2-1/2" x 6"
3-87	Stringer, Open, 2-1/2" x 8"
3-88	Stringer, Open, 2-1/2" x 10"
3-89	Stringer, Open, 2-1/2" x 12"
3-90	Center Lock Stringer, 2-1/2" x 6"
3-900	Sponge Bowl, 4-5/16" x 2-13/64", 12oz
3-902	Sponge Bowl, 5-5/8" x 29/32", 22oz
3-903	Sponge Bowl, 6-27/64 x 21/4", 30oz
3-904	Sponge Bowl, 6-27/64" x 2-3/4", 36oz
3-905	Sponge Bowl, 7-9/16" x 3-3/16", 47oz
3-906	Sponge Bowl, 8-3/16" x 3-5/16", 64oz
3-908	Solution Basin, 12-33/64" x 3-35/64", 4-1/2qt
3-909	Solution Basin, 13-25/32" x 4-17/32", 6qt
3-91	Center Lock Stringer, 2-1/2" x 8"
3-910	Solution Basin, 13-25/32" x 4-59/64", 7qt
3-912	Iodine Cup, 3-5/16" x 2", 6oz
3-913	Iodine Cup, 5-1/32" x 2-23/64", 14oz
3-914	Medicine Cup, Graduated, 2-5/32" x 2-5/32", 2oz
3-915	Forceps Jar, 1-19/64" x 3-5/32"
3-916	Forceps Jar, 2-3/16" x 4-21/64"
3-917	Forceps Jar, 2-3/16" x 7-5/64"
3-918	Emesis Basin, 7-7/8" x 3-3/4" x 1-19/64", 10oz
3-919	Emesis Basin, 8-5/32" x 3-55/64" x 1-17/32", 12oz
3-92	Center Lock Stringer, 2-1/2" x 10"
3-920	Emesis Basin, 9-23/32" x 4-51/64" x 1-11/16", 20oz
3-921	Size 10, Mayo Tray, Perforated, 10" x 6-1/2" x 23/32"
3-922	Size 13, Mayo Tray, Perforated, 14" x 10" x 5/8"
3-925	Size 19, Mayo Tray, Perforated, 19-11/64" x 12-23/32" x 33/64"
3-926	Size 10, Mayo Tray, Non-Perforated, 10" x 6-1/2" x 23/32"
3-927	Size 13, Mayo Tray, Non-Perforated, 14" x 10" x 5/8"
3-93	Center Lock Stringer, 2-1/2" x 12"
3-930	Size 19, Mayo Tray, Non-Perforated, 19-11/64" x 12-23/32" x 33/64"
3-934	Instrument Tray, 10" x 6-1/2" x 23/32", Rolled Edge
3-935	Instrument Tray, 14" x 10" x 5/8", Rolled Edge
3-936	Instrument Tray, 21-11/32" x 16-7/16" x 7/8", Rolled Edge
3-937	Cath Tray, 8-25/32" x 4-61/64" x 1-49/64"
3-939	Flat Cath Cover, 8-31/32" x 5-5/32" x 7/16"
3-94	Center Lock Stringer, 2-1/2" x 16"
3-940	Dome Cath Cover, 8-31/32" x 5-5/32" x 7/16"
3-941	Cath Tray, Long, 8-35/64" x 3-5/16" x 19/32"
3-943	Cath Tray, Long, 15-3/4" x 2-3/4" x 2"
3-945	Instrument Tray, Solid, 12-23/32" x 10-27/64" x 2-1/2"
3-946	Instrument Tray, Solid, 12-23/32" x 10-27/64" x 3-3/4"
3-947	Flat Cover for 3-945 and 3-946, 12-29/32" x 10-35/64" x 7/16"
3-948	Instrument Tray, Solid, 10-25/64" x 6-49/64" x 1-27/32"
3-949	Flat Cover for 3-948, 10-15/32" x 6-61/64" x 7/16"
3-95	Center Lock Stringer, 2-1/2" x 18"
3-953	Dressing Jar w/cover, 4" x 2-5/8", 1/2qt
3-954	Dressing Jar w/cover 5" x 6-3/8", 2qt
3-955	Dressing Jar w/cover, 5-7/8" x 7-5/64", 3qt
3-96	Closed U Stringer, 2-1/2" x 6"
3-97	Closed U Stringer, 2-1/2" x 8"
3-98	Closed U Stringer, 2-1/2" x 10"
3-99	Closed U Stringer, 2-1/2" x 12"

MILTEX STEVENS HOOK

18-440	5" Stevens Hook, Right Angle 5mm Deep
--------	---------------------------------------

**MILTEX STILLE-LUER RONGEUR**

25-454	8 3/4" Stille-Luer Rongeur, Straight Jaws 9 x 15mm
25-456	8 3/4" Stille-Luer Rongeur, Curved Jaws 9 x 15mm
25-458	9" Stille-Luer Rongeur, Duckbill, Angled On Side Jaws 6 x 13mm

**MILTEX SURGIGUIDE TEMPLATE**

For outlining lesions. 110 x 75mm. Flexible stainless steel.

21-278	Surgiguide Template
--------	---------------------

**MILTEX TENACULUM HOOK**

30-950	8" Tenaculum Hook, Extra Delicate, Handle Graduated in cm
34-250	6 1/2" Tenaculum Hook, Dissecting, Sharp Point

MILTEX THORPE CALIPER

Graduated in inches and mm, permits measurements in deep areas.

18-657	4 1/2" Thorpe Caliper
--------	-----------------------

**MILTEX TYRELL (IRIS) HOOK**

Very delicate.

18-412	5" Iris Hook, Sharp
18-414	5" Iris Hook, Blunt

MILTEX U.S. ARMY PATTERN CHISEL

6 1/2".

27-260	6mm Wide Chisel
27-262	12mm Wide Chisel
27-264	18mm Wide Chisel

MILTEX U.S. ARMY PATTERN GOUGE

6 1/2".


27-280	6mm Wide Gouge
27-282	12mm Wide Gouge

MILTEX U.S. ARMY PATTERN OSTEOTOME

6 1/2".

27-270	6mm Wide Osteotome
27-272	12mm Wide Osteotome
27-274	18mm Wide Osteotome

General Surgery

**MILTEX UEBE APPLICATORS**

- 19-170 7" Uebe Applicator, Triangular Tip
 19-172 7" Uebe Applicator, Cross Serrated Tip

**MILTEX UNNA COMEDONE EXTRACTOR**

Spoon end.

- 33-210 5 3/4" Comedone Extractor

**MILTEX VERUCCA CURETTE**

Cup 4mm diameter.

- 40-54 5" Verucca Curette

OMRON TAYLOR PERCUSSION HAMMER

Large 7 1/2" chrome steel handle and triangular rubber head.

- 0500 Taylor Percussion Hammer

COVIDIEN/US SURGICAL VALLEYLAB LIGASURE LAP VESSEL SEALING INSTRUMENT

For sealing and fine dissection of structures in small surgical spaces. Maryland style jaw provides fine dissection and sealing capabilities. Minimizes thermal spread to approximately 2mm. 5mm single use instrument.

- LS1000 Single Use Laparoscopic Device, 32cm, (12.6"), 5mm Diameter, 6/cs

VANTAGE COMEDOME EXTRACTOR

- V933-200 Extractor, Open Ends
 V933-201 Schamberg Extractor, Small End Crimped

VANTAGE FINGER RING CUTTER


- V933-140 6 1/4" Finger Ring Cutter

Accessories

- V933-142 Blades

VANTAGE WARTENBERG PINWHEEL

- V91-226 Wartenberg Pinwheel

**MILTEX VAN BUREN URETHRAL SOUNDS**

11".

- 29-2 Set of 6 Urethral Sounds In Khaki Roll, Sizes 14 Fr To 24 Fr
 29-4 Set of 12 Urethral Sounds In Khaki Roll, Sizes 8 Fr To 30 Fr
 29-6-10 10 Fr Urethral Sound
 29-6-12 12 Fr Urethral Sound
 29-6-14 14 Fr Urethral Sound
 29-6-16 16 Fr Urethral Sound
 29-6-18 18 Fr Urethral Sound
 29-6-20 20 Fr Urethral Sound
 29-6-22 22 Fr Urethral Sound
 29-6-24 24 Fr Urethral Sound
 29-6-26 26 Fr Urethral Sound
 29-6-28 28 Fr Urethral Sound
 29-6-30 30 Fr Urethral Sound
 29-6-32 32 Fr Urethral Sound
 29-6-34 34 Fr Urethral Sound
 29-6-36 36 Fr Urethral Sound
 29-6-38 38 Fr Urethral Sound
 29-6-40 40 Fr Urethral Sound
 29-6-8 8 Fr Urethral Sound

Accessories

- 29-2A Khaki Roll Only For 6 Sounds
 29-4A Khaki Roll For 12 Sounds

**MILTEX WALTHER FEMALE DILATOR-CATHETERS**

5 1/4".

- 29-33-12 12 Fr Female Dilator-Catheter
 29-33-14 14 Fr Female Dilator-Catheter
 29-33-16 16 Fr Female Dilator-Catheter
 29-33-18 18 Fr Female Dilator-Catheter
 29-33-20 20 Fr Female Dilator-Catheter
 29-33-22 22 Fr Female Dilator-Catheter
 29-33-24 24 Fr Female Dilator-Catheter
 29-33-26 26 Fr Female Dilator-Catheter
 29-33-28 28 Fr Female Dilator-Catheter
 29-33-30 30 Fr Female Dilator-Catheter
 29-33-32 32 Fr Female Dilator-Catheter
 29-33-34 34 Fr Female Dilator-Catheter
 29-33-36 36 Fr Female Dilator-Catheter
 29-33-38 38 Fr Female Dilator-Catheter

MILTEX WALTHER URETHRAL SOUNDS

11". Shaft tapered downward toward handle.


- 29-18-16 16 Fr Urethral Sound
 29-18-18 18 Fr Urethral Sound
 29-18-20 20 Fr Urethral Sound
 29-18-22 22 Fr Urethral Sound
 29-18-24 24 Fr Urethral Sound
 29-18-26 26 Fr Urethral Sound
 29-18-28 28 Fr Urethral Sound
 29-18-30 30 Fr Urethral Sound
 29-18-32 32 Fr Urethral Sound
 29-18-34 34 Fr Urethral Sound
 29-18-36 36 Fr Urethral Sound

MILTEX POST MORTEM HAMMER

34-220 9 1/2" Post Mortem Hammer w/hook handle



35-400



35-401



35-402

MILTEX LAPAROSCOPIC BIOPSY PUNCHES & FORCEPS

Lightweight one-piece design. Ergonomic standard free-moving handle with 360° finger rotation knob which allows for single-handed operation. Rotation lock/brake feature engages when the handle is in fully closed position. 5mm OD and 35 cm shaft. Luer lock flush port. Extra screw-on and lanyard caps. Monopolar HF hook-up port and insulated sheath and handle. Ideal for assisting with reliable biopsy specimen collection.

35-400 Biopsy Forceps w/tooth on each jaw, Single Action

35-401 Biopsy Punch w/lower serrated jaw, Single Action

35-402 Biopsy Forceps, Oval Fenestrated Jaws, Single Action



35-200



35-201




35-202



35-203



35-204



35-205

MILTEX LAPAROSCOPIC DISSECTING FORCEPS

Lightweight one-piece design. Ergonomic standard free-moving handle with 360° finger rotation knob which allows for single-handed operation. Rotation lock/brake feature engages when the handle is in fully closed position. 5mm OD and 35 cm shaft. Luer lock flush port. Extra screw-on and lanyard caps. Monopolar HF port and insulated sheath and handle. Ideal for precise dissection and separation of varying types of tissue. Dual action jaws.

35-200 Tapered Maryland Dissector

35-201 Maryland Dissector

35-202 Micro Dissecting Forceps

35-203 Bullet-Nose Dissecting Forceps

35-204 Tapered Dissecting Forceps

35-205 Mixer Spreader-Dissector



35-100



35-101



35-102



35-103



35-104



35-105



35-106




35-107



35-108



35-109



35-110

MILTEX LAPAROSCOPIC GRASPING FORCEPS

Lightweight one-piece design. Insulated sheath for added safety. Ergonomic ratchet handle with 360° finger rotation knob which allows for single-handed operation. Rotation lock/brake feature engages when the handle is in fully closed position. 5mm OD and 35 cm shaft. Luer lock flush port. Extra screw-on and lanyard caps. #35-110 available with standard non-ratchet handle, monopolar HF port and fully insulated sheath and handle. Ideal for precise manipulation and firm grasp of varying types of tissue. Dual action jaws.

35-100 Dual Cup Grasping Forceps

35-101 Retraction Grasping Forceps, 2x3 Blunt Teeth

35-102 Atraumatic Allis Grasping Forceps

35-103 Atraumatic Allis Grasping Forceps, Long Jaw 20mm Length

35-104 Allis Grasping Forceps

35-105 Babcock (Atraumatic) Grasping Forceps

35-106 Bowel Grasping Forceps, 20mm Fenestrated

35-107 Retrograde Grasping Forceps

35-108 DeBakey Forceps, Straight, 45mm Fenestrated

35-109 Grasping Forceps

35-110 Grasping Forceps, Non-Ratchet Handle, Mono HF Port



35-500



35-501



35-502

MILTEX LAPAROSCOPIC NEEDLE HOLDERS - TUNGSTEN CARBIDE

Lightweight one-piece design. Tungsten carbide jaws provide firm grip on needles. In-line ring handle design places surgeon's hand in line with the needle to help facilitate natural suturing motion. German stainless steel with gold-plated finger rings. 5mm OD and 35 cm shaft. Luer lock flush port. Extra screw-on and lanyard caps.

- 35-500 In-Line Tungsten Carbide Needle Holder, Left Curved Single Action Jaw
- 35-501 In-Line Tungsten Carbide Needle Holder, Right Curved Single Action Jaw
- 35-502 In-Line Tungsten Carbide Needle Holder, Straight Single Action Jaw



35-300



35-301



35-302



35-303


MILTEX LAPAROSCOPIC SCISSORS

Lightweight one-piece design. Ergonomic standard free-moving handle with 360° finger rotation knob which allows for single-handed operation. SuperCut patterns have a razor-sharp upper blade and finely-serrated lower blade. Rotation lock/brake feature engages when the handle is in fully closed position. 5mm OD and 35 cm shaft. Luer lock flush port. Extra screw-on and lanyard caps. Monopolar hook-up port and insulated sheath and handle. #35-303 is a standard non-insulated stainless steel pattern without rotation and monopolar features. Ideal for reliable cutting of varying types of tissue.

- 35-300 Supercut Metzenbaum Scissors, Curved Dual Action Jaws
- 35-301 Hook Scissors, Single Action Jaw
- 35-302 Supercut Scissors, Straight Single Action Jaw
- 35-303 Supercut Scissors, Stainless Steel, Non-Insulated, Straight Single Action Jaw



35-503



35-504



35-505

MILTEX LAPAROSCOPIC SPECIALITY FORCEPS

Lightweight one-piece design. German stainless steel with spring-loaded handle. Claw forceps available in 5mm & 10mm OD and are ideal for cholecystectomy and appendectomy procedures. Spoon forceps available in 10mm OD with wide and large cup-shaped jaws which are ideal for retrieval of stones, fragments and tissue. Luer lock flush port. Extra screw-on and lanyard caps. 35 cm shaft. Ideal for strong grasp of varying types of tissue and reliable extraction.

- 35-503 Claw Forceps, 5mm OD
- 35-504 Claw Forceps, 10mm OD
- 35-505 Spoon Forceps, 10mm OD

MILTEX LAPAROSCOPES

Rod-lens system and technology. 5mm and 10mm OD patterns available. Stainless steel with all possible leakage paths laser welded and/or solder-sealed. Distal and proximal ends made from scratch resistant sapphire glass. Triple tube (lumen) design. Comes equipped with special lens polishing paste and three standard light-source adapters: ACMI/British Standard, Richard Wolf & Storz/Olympus. High quality German-crafted instruments are fully autoclavable at 132°C.

- 35-600 Laparoscope, 5mm, 0 Degree, 30cm Working Length
- 35-601 Laparoscope, 5mm, 30 Degree, 30cm Working Length
- 35-603 Laparoscope, 10mm, 0 Degree, 33cm Working Length